

# **Technical Data**

## NNN Modified Medium (Twin Pack)

**M681** 

NNN Modified Medium (Twin Pack) is used for cultivation of Leishmaniae and Trypanosomes.

## Composition\*\*

I	
Ingredients	Gms / Litre
Part A	-
Meat extract	3.000
Peptone	5.000
Sodium chloride	8.000
Agar	15.000
Final pH ( at 25°C)	$7.3 \pm 0.2$
Part B	-
Sodium chloride	8.000
Potassium chloride	0.200
Calcium chloride	0.200
Monopotassium dihydrogen phosphate	0.300
Dextrose	2.500
Final pH ( at 25°C)	$7.0 \pm 0.2$

<sup>\*\*</sup>Formula adjusted, standardized to suit performance parameters

#### **Directions**

Part A: Suspend 31 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C and aseptically add 10% of sterile defibrinated rabbit or human blood after inactivation at 56°C for 30mins. Mix well and dispense in 5 ml amounts in test tubes or 25 ml amounts in flasks. Allow tubed media to cool in slanted position.

Part B: Suspend 11.2 grams of Part B in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool and add approximately 2 ml in tubes or 10-15 ml in flasks over solidified Part A medium.

## **Principle And Interpretation**

The protozoan family *Trypanosomatidae* includes members from the genera *Leishmania* and *Trypanosoma*, which are flagellates that inhabit the blood and tissues of humans.

NNN Medium was developed by Novy, McNeal (1) and modified by Nicolle (2). NNN Modified Medium is a modification of the original medium and consists of two phases, blood agar (Part A) and Lockes solution (Part B) (3). This modified medium is commonly used for diagnostic work (4, 5).

This medium consists of a blood agar base and an overlay medium. The blood agar base is a highly nutritious medium that supports the growth of fastidious organisms like *Leishmania* and *Trypanosoma*. The specimens are inoculated into the liquid phase of the diphasic medium and incubated. This favours the development of organisms in the insect vector. The amastigotes transform to promastigotes in about 24 hours (5).

#### **Quality Control**

#### Appearance

Part A: Cream to tan homogeneous free flowing powder Part B: White to cream homogeneous free flowing powder

#### **Gelling**

Firm, comparable with 1.5% Agar gel

#### **Colour and Clarity of Prepared medium**

Basal medium :Light amber clear to slightly opalescent gel. After addition of sterile defibrinated rabbit or human blood : Red coloured opaque gel Part B : Colourless clear liquid

#### Reaction

Reaction of 3.1% w/v aqueous solution (Part A) at  $25^{\circ}$ C. pH :  $7.3\pm0.2$  Reaction of 1.12% w/v aqueous solution (Part B) at  $25^{\circ}$ C. pH :  $7.0\pm0.2$ 

HiMedia Laboratories Technical Data

#### **Cultural Response**

M681: Cultural characteristics observed after an incubation at 21-26°C for 48-72 hours, with added sterile defibinated rabbit or human blood.

Organism Growth

**Cultural Response** 

Leishmania donovani luxuriant Trepanosoma cruzi luxuriant

## Storage and Shelf Life

Store below 30°C in tightly closed container and the prepared medium at 2 - 8°C. Use before expiry date on the label.

### Reference

- 1. Novy F. G. and McNeal W. J., 1904, J. Inf. Diseases B, 1:1.
- 2. Nicolle A (1908) Comptes rendus de l Academie des Sciences (Paris) 146:842.
- 3. Cruickshank R., Duguid J. P., Marmion B. P., Swain R. H. A. (Eds) 1975, Medical Microbiology, 12th Edition, Vol. II, Churchill Livingstone
- 4. Taylor A. R., Baker J. R., (Eds.), 1978, Methods of Cultivating Parasites in vitro, Academic Press, London, pp 55-88
- 5. Collee J. G., Fraser A. G., Marmion B. P., Simmons A., (Eds.), Mackie and McCartney, Practical Medical Microbiology, 1996, 14th Edition, Churchill Livingstone

Revision: 02 / 2016

CE

#### Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia<sup>™</sup> publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia<sup>™</sup> Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.