

WORLD CLASS QUALITY

HiCrome™

Single Streak Rapid Differentiation Series

Single streak
24hr
Results



Let colors

Confirm your identity

HIMEDIA®

For Life is Precious

HiMediaLaboratories™

himedialabs.com



Welcome home to **HIMEDIA**

5th generation Chemically Defined Microbiology Media



5th Gen Chemically Defined Microbiology Media
Free from TSE / BSE / GMO

4th Gen Vegetable based dehydrated culture media
Free from TSE / BSE risk

3rd Gen Animal based dehydrated culture media

2nd Gen Dehydrated culture media prepared using raw materials

1st Gen Media preparation in lab using meat and other ingredients

Index

Product	Code	Page	Product	Code	Page
HiCrome™ ECC Agar	M1293	6	HiCrome™ Clostridial Agar Base	M2026	74
HiCrome™ ECC Selective Agar Base / Modified	M1294 / M2056	8	HiCrome™ Lactobacillus Selective Agar Base	M2065	76
HiCrome™ Coliform Agar w/ SLS HiCrome™ Coliform Broth w/ SLS	M1300 / M1826	10	HiCrome™ Enterococci Agar / Broth	M1414 / M1376	78
HiCrome™ Rapid ECC Broth	M2011	12	HiCrome™ Enterococcus faecium Agar Base	M1580	80
HiCrome™ M-Coliform Differential Agar Base	M1951	14	HiCrome™ Strep B Selective Agar Base / Modified	M1840 / M1966	82
HiCrome™ M-Coliconfirm Broth Base	M2064	16	HiCrome™ VRE Agar Base / Modified	M1830 / M1925	84
M-E. coli Broth	M1426	18	HiCrome™ Listeria Agar Base / Modified	M1417F / M1417	86
HiCrome™ M-Lauryl Sulphate Agar	M1569	20	L. Mono Differential Agar Base	M1540	88
HiCrome™ Chromogenic Coliform Agar (CCA)	M1991I	22	HiCrome™ L. mono Rapid Differential Agar Base	M1924	90
HiCrome™ E. coli Agar	M1295 / M1295I	24	HiCrome™ L. mono Differential Agar Base	M2009	92
HiCrome™ Coliform Agar Modified	M1832	26	HiCrome™ Aureus Agar Base	M1468	94
HiCrome™ M-TEC Agar / HiCrome M-TEC Broth	M1571 / M1713	28	HiCrome™ Staph Agar Base, Modified	M1837	96
HiCrome™ MacConkey Sorbitol Agar Base	M1340	30	HiCrome™ Staph Selective Agar	M1931	98
HiCrome™ EC 0157: H7 Agar, Modified	M1574A	32	HiCrome™ MeReSa Agar Base	M1674	100
HiCrome™ EC 0157:H7 Selective Agar Base, Modified	M1575A	34	HiCrome™ MRSA Agar Base, Modified	M1953	102
HiCrome™ Enrichment Broth Base for EC 0157:H7	M1598	36	HiCrome Rapid MRSA Agar Base	M1974	104
HiCrome™ M-Modified ECO157:H7 Selective Agar Base	M1862	38	HiCrome™ Bacillus Agar	M1651	106
HiCrome™ Enterobacter sakazakii Agar/Modified	M1577 / M1641	40	M-CP Agar Base	M1354	108
HiCrome™ Cronobacter Isolation Agar (CCI Agar)	M2062I	42	HiCrome™ Nickels and Leesment Medium	M1712	110
Salmonella Differential Agar / Modified (Twin pack) (RajHans Medium)	M1078 / M1082	44	HiCrome™ Bifidobacterium Agar	M1960	112
HiCrome™ RajHans Medium/Modified (Salmonella Agar/Modified)	M1633 / M1634	46	HiCrome™ Acinetobacter Agar Base	M1938	114
HiCrome™ Salmonella Agar /HiCrome™ Improved Salmonella Agar	M1296 / M1466	48	HiCrome™ OGYE Agar Base	M1467	116
HiCrome™ Selective Salmonella Agar Base	M1842	50	HiCrome™ Candida Differential Agar/Base, Modified	M1297A / M1456A	118
HiCrome™ MM Agar	M1393	52	HiCrome™ Mueller Hinton Agar (For Antifungal testing)	M2067	120
HiCrome™ MM Agar Modified	M1816	54	HiCrome™ Candida Differential Agar Base	M1297AR	122
HiCrome™ Klebsiella Selective Agar Base	M1573	56	HiCrome™ Malassezia Agar (Twin Pack)	M1985	124
HiCrome™ ESBL Agar Base	M1829	58	HiCrome™ PA Broth	M1633	128
HiCrome™ KPC Agar Base	M1831	60	HiColiform™ Broth, Modified	M1850	130
HiCrome™ Vibrio Agar	M1682	62	Rapid HiColiform™ Agar / Broth	M1465 / M1453	132
HiCrome™ UTI Agar, Modified / HiCrome™ UTI Selective Agar	M1418 / M1505	64	HiCrome™ ECD Agar w/ MUG	M1488	134
HiCrome™ UTI Agar	M1353 / M1353R	66	HiFluoro™ Pseudomonas Agar Base	M1469	136
HiCrome™ Universal Differential Medium	M1600	68	HiCrome™ EC Broth w/ RUG	M2073	138
HiCrome™ Mueller Hinton Agar	M2010	70	Supplement		140-149
HiCrome™ Yersinia Agar Base	M2025	72	HiDtect		150-153

28 Chromogenic Media are also available in the HiCromeVeg Category where animal based nutrients have been substituted with Vegetable based nutrients



Corresponding lists of Chromogenic media containing Animal peptone and HiVeg™ peptones

HiCrome™ Animal Peptone Based Media

M1651 - HiCrome™ Bacillus Agar
M1297A- HiCrome™ Candida Differential Agar
M1456A- HiCrome™ Candida Differential Agar Base, Modified
M1991I - HiCrome™ Chromogenic Coliform Agar (CCA Agar)
M1300 - HiCrome™ Coliform Agar w/ SLS
M1293 - HiCrome™ ECC Agar
M1294 - HiCrome™ ECC Selective Agar Base
M1488 - HiCrome™ ECD Agar w/MUG
M1295 - HiCrome™ <i>E. coli</i> Agar
M1575A - HiCrome™ EC 0157:H7 Selective Agar Base
M1577- HiCrome™ <i>Enterobacter sakazakii</i> Agar
M1641 - HiCrome™ <i>Enterobacter sakazakii</i> Agar, Modified
M1376 - HiCrome™ Enterococci Broth
M1580- HiCrome™ <i>Enterococcus faecium</i> Agar Base
M1466 - HiCrome™ Improved <i>Salmonella</i> Agar
M1573 - HiCrome™ <i>Klebsiella</i> Selective Agar Base
M1393 - HiCrome™ MM Agar
M1674 - HiCrome™ MeReSa Agar Base
M1712 - HiCrome™ Nickels and Leesment Agar Base
M1296 - HiCrome™ <i>Salmonella</i> Agar
M1353 - HiCrome™ UTI Agar
M1418 - HiCrome™ UTI Agar, Modified
M1505 - HiCrome™ UTI Selective Agar
M1682 - HiCrome™ <i>Vibrio</i> Agar
M1469 - HiFluoro <i>Pseudomonas</i> Agar Base
M1826- Coliform Broth w/SLS
M1540 - L.mono Differential Agar Base
M1354 - M-CP Agar Base
M1465 - Rapid HiColiform Agar
M1453 - Rapid HiColiform Broth
M1082 - <i>Salmonella</i> Differential Agar, Modified (Twin Pack)
M1078 - <i>Salmonella</i> Differential Agar (Twin Pack) (RajHans Medium)

HiCrome™ Veg-Peptone Based media

#MCD1651 - HiCrome™ Bacillus HiCynth™ Agar
#MV1297A- HiCrome™ Candida Differential HiVeg™ Agar
MV1456A- HiCrome™ Candida Differential HiVeg™ Agar Base, Modified
#MCD1991I - HiCrome™ Chromogenic Coliform HiCynth™ Agar
MV1300 - HiCrome™ Coliform HiVeg™ Agar w/SLS
MV1293 - HiCrome™ ECC HiVeg™ Agar
MV1294 - HiCrome™ ECC Selective HiVeg™ Agar Base
MV1488 - HiCrome™ ECD HiVeg™ Agar w/ MUG
MV1295 - HiCrome™ <i>E. coli</i> HiVeg™ Agar
MV1575A - HiCrome™ EC 0157:H7 Selective HiVeg™ Agar Base
MV1577- HiCrome™ <i>Enterobacter sakazakii</i> HiVeg™ Agar
MV1641 - HiCrome™ <i>Enterobacter sakazakii</i> HiVeg™ Agar, Modified
MV1376 - HiCrome™ Enterococci HiVeg™ Broth
MV1580- HiCrome™ <i>Enterococcus faecium</i> HiVeg™ Agar Base
#MV1466 - HiCrome™ Improved <i>Salmonella</i> HiVeg™ Agar
MV1573 - HiCrome™ <i>Klebsiella</i> Selective HiVeg™ Agar Base
MV1393 - HiCrome™ MM HiVeg™ Agar
MV1674 - HiCrome™ MeReSa HiVeg™ Agar Base
MV1712 - HiCrome™ Nickels and Leesment HiVeg™ Agar Base
MV1296 - HiCrome™ <i>Salmonella</i> HiVeg™ Agar
#MV1353 - HiCrome™ UTI HiVeg™ Agar
MV1418 - HiCrome™ UTI HiVeg™ Agar, Modified
MV1505 - HiCrome™ UTI Selective HiVeg™ Agar
#MV1682 - HiCrome™ <i>Vibrio</i> HiVeg™ Agar
MV1469 - HiFluoro <i>Pseudomonas</i> HiVeg™ Agar Base
MV1826 - Coliform HiVeg™ Broth w/SLS
#MV1540 - L.mono Differential HiVeg™ Agar Base
MV1354 - M-CP HiVeg™ Agar Base
MV1465 - Rapid HiColiform HiVeg™ Agar
MV1453 - Rapid HiColiform HiVeg™ Broth
MV1082 - <i>Salmonella</i> Differential HiVeg™ Agar, Modified (Twin Pack)
MV1078 - <i>Salmonella</i> Differential HiVeg™ Agar (Twin Pack) (RajHans Medium)

HiCrome™ HiCynth™ Media available



Standard culture media procedures, though a reliable way of detecting micro-organisms lack speed and accuracy of isolation and differentiation so critical in dealing with pathogenic microorganisms.

HiMedia's HiCrome™ range of culture media employing the chromogen technology of visual identification significantly removes the guesswork out of identification and differentiation, thereby obviating the need for subculturing, thus saving time.

The methodology is simple and precisely designed for each bacteria. Organisms are identified through simple enzymatic reactions specific to their species, yielding visually distinct colours. Over the past two decades chromogenic media have been well researched and documented to merit incorporation in standard microbiology lab protocols.

HiMedia have been a part of this global revolution in diagnostic microbiology where we have developed the largest range of chromogenic media, and the research continues.

We also offer the range of these HiCrome™ media in the HiCromeVeg™ and HiCrome™ HiCynth™ versions, wherein the animal based nutrients have been substituted with their veg-based counterparts or chemically defined peptones.

• • • *expect only Quality from us™*



Chromogenic Media Index / Cross References

Equivalent Media of various Brands

Organisms	Page No.	Code	HiMedia	Chromagar	Merck	Oxoid	Remel
<i>E. coli</i> and Total coliforms	6	M1293	HiCrome™ ECC Agar	CHROMagar™ ECC	ECC ChromoSelect Agar	Brilliance™ <i>E. coli</i> / coliform Agar	—
	8	M1294 / M2056	HiCrome™ ECC Selective Agar Base / Modified	—	Chromocult Coliform Agar	Brilliance™ <i>E. coli</i> /coliform Selective Agar	—
	10	M1300 / M1826	HiCrome™ Coliform Agar w/ SLS / Coliform Broth w/SLS	—	Coliform Chromo-Select Agar	—	—
	12	M2011	HiCrome™ Rapid ECC Broth	—	—	—	—
	14	M1951	HiCrome™ M-Coliform Differential Agar Base	—	—	—	—
	16	M2064	HiCrome™ M-Coliconfirm Broth Base	—	—	—	—
	18	M1426	M-E.coli Broth	—	—	—	—
	20	M1569	HiCrome™ M-Lauryl Sulphate Agar	—	—	—	—
	22	M1991I	HiCrome™ Chromogenic Coliform Agar (CCA)	CHROMagar™ CCA	ChromoCult® Coliform Agar	—	—
	24	M1295 / M1295I	HiCrome™ <i>E. coli</i> Agar	—	—	—	—
	26	M1832	HiCrome™ Coliform Agar Modified	—	—	—	—
	64	M1663	HiCrome™ PA Broth	—	—	—	—
	65	M1850	HiColiform™ Broth, Modified	—	—	—	—
	66	M1465 / M1453	Rapid HiColiform Agar / Broth	AquaCHROM™ ECC	Fluorocult LMX Broth Modified (Manafi and Ossmer)	—	—
	28	M1571 / M1713	HiCrome™ M-TEC Agar / HiCrome™ M-TEC Broth	m-TEC Agar	m-TEC ChromoSelect Agar	—	—
<i>E. coli</i>	14	M1295 / M1295I	HiCrome™ <i>E. coli</i> Agar	CHROMagar™ <i>E. coli</i> / CHROMagar™ TBX	ChromoCult® TBX-Agar	Tryptone Bile X-Glucuronide Medium (TBX)	—
	67	M1488	HiCrome™ ECD Agar w/ MUG	—	ECD ChromoSelect Agar w/MUG	—	—
	138	M2073	HiCrome™ EC Broth w/ RUG	—	—	—	—
<i>Escherichia coli</i> O157:H7	30	M1340	HiCrome™ MacConkey Sorbitol Agar Base	CHROMagar™ .0157	MacConkey Sorbitol ChromoSelect Agar	MacConkey Sorbitol Agar (with chromogenic substrate)	Sorbitol MacConkey Agar w/BCIG
	32	M1574A	HiCrome™ EC 0157:H7 Agar	CHROMagar™ .0157	ECO157 ChromoSelect Agar	—	—
	34	M1575A	HiCrome™ EC 0157:H7 Selective Agar Base*	—	ECO157:H7 ChromoSelect Agar Base	—	—
	36	M1598	HiCrome™ Enrichment Broth Base for EC 0157:H7	—	Enrichment ChromoSelect Broth Base for ECO157:H7	—	—
	38	M1862	HiCrome™ Modified ECO157:H7 Selective Agar Base	—	—	—	—
<i>Cronobacter sakazakii</i>	40	M1577 M1641	HiCrome™ Enterobacter sakazakii Agar, Modified (HiCrome™ Cronobacter sakazakii Agar, Modified)	CHROMagar™ <i>E. sakazakii</i>	ChromoCult® Enterobacter sakazakii Agar	Brilliance Enterobacter sakazakii Agar	Chromo <i>E. sakazakii</i> Medium
	42	M2062I	HiCrome™ Cronobacter Isolation Agar (CCI Agar)	—	—	—	—
<i>Salmonella</i> species	44	M1078 / M1082	Salmonella Differential Agar (RajHans Medium / Modified) (Twin pack)	Rambach Agar	Rambach Agar	—	—
	46	M1633/ M1634	HiCrome RajHans Medium/ Modified (Salmonella Agar/ Modified)	Rambach Agar	RajHans Chromo-Select Medium, Modified	Rambach Agar	—
	48	M1296 / M1466	HiCrome™ Salmonella Agar / HiCrome™ Improved Salmonella Agar	CHROMagar™ <i>Salmonella</i>	Salmonella Chromo-Select Agar, Improved	Salmonella Chromogenic Agar Base	Salmonella Chromogenic Agar
	50	M1842	HiCrome™ Selective Salmonella Agar Base	CHROMagar™ <i>Salmonella</i> Plus	—	Brilliance Salmonella Agar Base	—
	52	M1393	HiCrome™ MM Agar	—	—	—	—
	54	M1816	HiCrome™ MM Agar Modified	—	—	—	—
<i>Klebsiella</i> species	56	M1573	HiCrome™ Klebsiella Selective Agar Base	—	Klebsiella Chromo-Select Agar Base	—	—
ESBL/ Carbapenem resistant	58	M1829	HiCrome™ ESBL Agar Base	CHROMagar™ ESBL	ESBL ChromoSelect Agar Base	Brilliance ESBL Agar	—

Organisms	Page No.	Code	HiMedia	Chromagar	Merck	Oxoid	Remel
	60	M1831	HiCrome™ KPC Agar Base	CHROMagar™ KPC	—	Brilliance CRE Agar	—
<i>Vibrio species</i>	62	M1682	HiCrome™ Vibrio Agar	CHROMagar™ Vibrio	Vibrio ChromoSelect Agar	—	—
UTI Infections	64	M1418 / M1505	HiCrome™ UTI Agar, modified / HiCrome™ UTI Selective Agar	—	—	—	—
	66	M1353 / M1353R	HiCrome™ UTI Agar / Modified	CHROMagar™ Orientation	UTI ChromoSelect Agar, Modified	Brilliance UTI Clarity Agar Brilliance UTI Agar	Chromogenic UTI Medium
	68	M1600	HiCrome™ Universal Differential Medium	—	Universal Differential ChromoSelect Medium	—	—
	70	M2010	HiCrome Mueller Hinton Agar	CHROMagar™ Orientation	—	—	—
<i>Enterococcus</i>	78	M1414 / M1376	HiCrome Enterococci Agar / Broth	Aquachrom™ Enterococcus	Chromocult Enterococci Broth / ChromoCult® enterococci Agar	—	—
	80	M1580	HiCrome™ Enterococcus faecium Agar Base	—	Enterococcus faecium ChromoSelect Agar Base	—	—
	82	M1840 / M1966	HiCrome™ Strep B Selective Agar Base / Modified	CHROMagar™ / CHROMagar™ Strep B	—	Brilliance GBS Agar	—
VRE	84	M1830 / M1925	HiCrome™ VRE Agar Base / Modified	CHROMagar™ VRE	VRE ChromoSelect Agar Base	Brilliance VRE Agar	—
<i>Listeria species</i>	86	M1417 / M1417F	HiCrome™ Listeria Agar Base, Modified	—	—	—	—
	88	M1540	L.mono Differential Agar Base	CHROMagar™™ ALOA	ChromoCult® Listeria Selective Agar (ALOA®)	—	—
	90	M1924	HiCrome™ L. mono Rapid Differential Agar Base	CHROMagar™ Listeria	Chromocult Listeria Selective Agar	Brilliance™ Listeria Agar	—
	92	M2009	HiCrome™ L.mono Differential Agar Base	—	—	Brilliance™ Listeria Agar Base	—
<i>Staphylococcus</i>	94	M1468	HiCrome™ Aureus Agar Base	—	Aureus ChromoSelect Agar Base	—	—
	96	M1837	HiCrome™ Staph Agar Base, Modified	CHROMagar™ Staph aureus	—	Brilliance Staph 24	—
	98	M1931	HiCrome™ Staph Selective Agar	—	—	—	—
MRSA/ MRSE	100	M1674	HiCrome™ MeReSa Agar Base	CHROMagar™ MRSA	MeReSa ChromoSelect Agar	—	—
	102	M1953	HiCrome™ MRSA Agar Base, Modified	—	—	Brilliance MRSA Agar	—
	104	M1974	HiCrome™ Rapid MRSA Agar Base	—	—	Brilliance MRSA 2 Agar	—
<i>Bacillus</i>	106	M1651	HiCrome™ Bacillus Agar	CHROMagar™ B.cereus	Bacillus ChromoSelect Agar	Brilliance Bacillus cereus Agar	—
<i>Clostridium species</i>	108	M1354	M-CP Agar Base	—	m-CP Agar Base	—	—
	74	M2026	HiCrome™ Clostridial Agar Base	—	—	—	—
Lactic acid bacteria	110	M1712	HiCrome™ Nickels and Leesment Medium	—	—	—	—
<i>Bifidobacterium species</i>	112	M1960	HiCrome™ Bifidobacterium Agar	—	—	—	—
Acinetobacter	114	M1938	HiCrome™ Acinetobacter Agar Base	CHROMagar™ Acinetobacter	—	—	—
Yeasts and Moulds	118	M1297A / M1456A	HiCrome™ Candida Differential Agar / Base, Modified	CHROMagar™ Candida	—	—	—
	122	M1297AR	HiCrome™ Candida Differential Agar Base	—	—	Brilliance Candida Agar	Chromogenic Candida Agar
	120	M2067	HiCrome™ Mueller Hinton Agar (for Antifungal Testing)	—	—	—	—
	116	M1467	HiCrome™ OGYE Agar Base	—	OGYE ChromoSelect Agar	—	—
	124	M1985	HiCrome™ Malassezia Agar Base (Twin Pack)	CHROMagar™ Malassezia	—	—	—
<i>Pseudomonas</i>	136	M1469	HiFluoro™ Pseudomonas Agar Base	—	—	—	—
<i>Yersinia</i>	72	M2025	HiCrome™ Yersinia Agar Base	—	—	—	—
Lactobacillus	76	M2065	HiCrome™ Lactobacillus Selective Agar Base	—	—	—	—

For identification and differentiation of *E. coli* and Total coliforms

HiCrome™ ECC Agar

Recommended as a differential medium for presumptive identification of *Escherichia coli* and other coliforms in food and environmental samples.

M1293

Composition **

Ingredients	Grams/Litre
Peptone, special	5.00
Yeast extract	3.00
Lactose	2.50
Disodium hydrogen phosphate	3.50
Potassium dihydrogen phosphate	1.50
Sodium chloride	5.00
Chromogenic mixture	20.30
Neutral red	0.03
Agar	15.00

Final pH (at 25°C) 6.8 ± 0.2

** Formula adjusted, standardized to suit performance parameters

Directions

Suspend 55.83 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Mix well and pour into sterile Petri plates.

Principle and Interpretation

Escherichia coli, a member of the family *Enterobacteriaceae* is a part of normal flora of the intestinal tract of humans and a variety of animals. Although most of *E. coli* does not cause gastrointestinal illnesses, certain groups of *E. coli* can cause life-threatening diarrhoea and severe sequelae or disability (1). HiCrome™ ECC Agar is a differential medium recommended for the presumptive identification of *E. coli* and other coliforms in food and environmental samples (2). The medium contains two chromogens. One of the chromogen is cleaved by the enzyme β -glucuronidase produced by *E. coli* to give blue to purple coloured colonies whereas the other chromogen is cleaved by the enzyme galactosidase, produced by majority of coliforms, resulting in the formation of rose-pink coloured colonies (3, 4).

Peptone special, yeast extract provide nitrogenous substances, carbonaceous compounds, long chain amino acids, vitamin B complex and other essential growth nutrients. Lactose is the fermentable carbohydrate, which aids in detecting lactose fermenters with neutral red as an indicator. Disodium hydrogen phosphate and potassium dihydrogen phosphate buffers the medium well. Sodium chloride maintains the osmotic equilibrium. Dry the surface of plate medium

before use.

Dilute the food sample 1: 5 or 1: 10 with 0.1% sterile Peptone Water (M028) and homogenize in a blender or a stomacher. Spread 0.5 ml or 1.0 ml of the homogenate over the agar surface with a sterile glass spreader and incubate the plates at 35-37°C for 18-24 hours. Count the blue/purple colonies and multiply with the dilution factor. The number of *E. coli* is reported per gram of food. The medium should be used only for in-vitro diagnostic purpose. Wear mask while handling the dehydrated product and avoid contact with eyes.

Type of specimen


Food and environmental sample

Specimen Collection and Handling

For food and dairy samples, follow appropriate techniques for sample collection and processing as per guidelines (5, 6).

Warning and Precautions

Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets



M1293 HiCrome™ ECC Agar

HiCromeVeg™ Freedom from BSE / TSE worries
Single Streak Rapid Differentiation Series

HiCrome™ ECC Agar (M1293) is also available as HiCrome™ ECC HiVeg™ Agar (MV1293) wherein all the animal origin nutrients have been replaced by vegetable based nutrients

HiCrome™ ECC Agar

Recommended as a differential medium for presumptive identification of *Escherichia coli* and other coliforms in food and environmental samples.

M1293

Limitations

1. β -glucuronidase is present in 97% of *E. coli* strains, however few *E. coli* may be negative.
2. Certain species of *Shigella* and *Salmonella* are β -glucuronidase positive which may appear as light blue.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

- Appearance of Powder** : Light yellow to pink coloured, homogeneous, free flowing powder.
- Gelling** : Firm, comparable with 1.5% Agar gel.
- Colour and Clarity of prepared medium** : Reddish pink coloured, opaque gel forms in Petri plates.
- Reaction** : Reaction of 5.58% w/v aqueous solution at 25°C. pH: 6.8 ± 0.2
- Cultural Response** : Cultural characteristics observed after an incubation at 35-37°C for 18-24 hours.

Organisms (ATCC)	Inoculum CFU	Growth	Recovery of colony	Colour
<i>Escherichia coli</i> (25922) (00013*)	50-100	luxuriant	>=70%	blue / purple
<i>Klebsiella pneumoniae</i> (13883) (00097*)	50-100	luxuriant	>=70%	rose / pink
<i>Pseudomonas aeruginosa</i> (27853) (00025*)	50-100	good - luxuriant	>=70%	straw
<i>Salmonella</i> Enteritidis (13076) (00030*)	50-100	luxuriant	>=70%	pink

Key : * : corresponding WDCM Numbers

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (7, 8).

References

1. Doyle M. P., (Ed.), 1989, Foodborne Bacterial Pathogens, Marcel Dekker, New York
2. Frampton E.W., Restaino L. and Blaszkowski N., 1988, J. Food Prot., 51:402.
3. Kilian M. and Bülow P., 1976, Acta. Pathol. Microbiol. Scand., Sect. B, 84:245.
4. Kilian M. and Bülow P., 1979, Acta. Pathol. Microbiol. Scand., Sect. B, 87:271.
5. American Public Health Association, Standard Methods for the Examination of Dairy Products, 1978, 14th Ed., Washington D.C.
6. Salfinger Y., and Tortorello M.L. Fifth (Ed.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.
7. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
8. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S. and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1

Ready Prepared Media			
Code	Product Name	Usage	Packing
Category : HiDip Slides			
HD036	HiDip™ HiCrome™ ECC Agar-HiCrome™ Salmonella Agar	for chromogenic screening of <i>E. coli</i> , coliforms and Salmonella on surfaces or food or water	5 tubes / 10 tubes
HD037	HiDip™ HiCrome™ ECC Agar-PCA w/TTC & Neutralizers	for detection of <i>E. coli</i> , coliforms and for total bacterial count with inactivation of disinfectants	5 tubes / 10 tubes
HD038	HiDip™ HiCrome™ ECC Agar-Baird Parker Agar w/Neutralizers	for detection of <i>E. coli</i> , coliforms and <i>S. aureus</i> with inactivation of disinfectants	5 tubes / 10 tubes

For identification and differentiation of *E. coli* and Total coliforms

HiCrome™ ECC Selective Agar Base / Modified

Recommended for detection of *Escherichia coli* and coliforms in water and food samples.

M1294/
M2056

Composition **	M1294	M2056
Ingredients	Grams/Litre	Grams/Litre
Peptone, special	6.00	-
Peptone	-	10.00
Tryptone	3.30	-
Sodium dihydrogen phosphate	0.60	2.20
Disodium hydrogen phosphate	1.00	2.70
Sodium chloride	2.00	5.00
Sodium pyruvate	1.00	1.00
L-Tryptophan	1.00	1.00
Sorbitol	1.00	1.00
Tergitol-7® (Sodium heptadecyl sulphate)	0.15	-
Potassium nitrate	-	1.00
Sodium lauryl sulphate (SLS)	-	0.200
Chromogenic mixture	0.43	0.200
Agar	10.00	15.00
Final pH (at 25°C)	6.8 ± 0.2	7.00 ± 0.2

** Formula adjusted, standardized to suit performance parameters

Directions

Suspend 26.48 grams of M1294 / 39.30 grams of M2056 in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. If desired, selective medium can be prepared by aseptically adding the rehydrated contents of 1 vial of HiCrome™ ECC Selective Supplement (FD190) to M1294. Mix well and pour into sterile Petri plates. Medium may show haziness, but it does not affect the performance of the medium.

Principle and Interpretation

These are selective medium recommended for the simultaneous detection of *Escherichia coli* and total coliforms in water and food samples (1, 3). The chromogenic mixture contains two chromogenic substrates. The enzyme β -D-galactosidase produced by coliforms cleaves one of the chromogen to form salmon to red coloured colonies (4). The enzyme β -D-glucuronidase produced by *E. coli*, cleaves X-glucuronide, the other chromogen (5). *E. coli* gives dark blue to violet coloured colonies due to cleavage of both the chromogens. Addition of L-Tryptophan improves the indole reaction, thereby increasing the detection reliability.

Peptone special, Peptone, Tryptone and sodium pyruvate provide nitrogenous substances, carbonaceous compounds, long chain amino acids, vitamins and other essential growth nutrients for the organisms. Sorbitol is the fermentable carbohydrate. Phosphates buffer the medium. The media formulation helps even the sublethally injured coliforms to recover and grow rapidly. Tergitol in M1294 inhibits gram-

positive as well as some gram-negative bacteria other than coliforms (3). Sodium lauryl sulphate in M2056 inhibits gram-positive organisms. Addition of HiCrome™ ECC Selective Supplement (FD190) in M1294 helps to inhibit the accompanying heterogenous microflora.

The medium is inoculated either by pour plate technique or by spreading the sample on the surface of plated medium. Membrane filter technique can also be used. Other gram negative bacteria forms colourless colonies, except some organisms which are β -glucuronidase positive. β -glucuronidase positive organisms gives light blue to turquoise colonies. Glucuronidase is present in 94-96% of *E. coli* strains and in some *Salmonella*, *Shigella* and *Yersinia* spp (2). To confirm *E. coli*, add a drop of Kovac's reagent on the dark blue to violet colony. Formation of cherry red colour indicates a positive reaction.

Type of specimen

Water samples ; Food samples

Specimen Collection and Handling


For water samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (6).

For food and dairy samples, follow appropriate techniques for sample collection and processing as per guidelines (7, 8).

After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

Read the label before opening the container. Wear protective gloves/ protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets



M1294 – HiCrome™ ECC Selective Agar Base

HiCromeVeg™ Freedom from BSE / TSE worries
Single Streak Rapid Differentiation Series

HiCrome™ ECC Selective Agar Base (M1294) is also available as HiCrome ECC Selective HiVeg™ Agar Base (MV1294) wherein all the animal origin nutrients have been replaced by vegetable based nutrients.

For identification and differentiation of *E. coli* and Total coliforms

HiCrome™ ECC Selective Agar Base / Modified

Recommended for detection of *Escherichia coli* and coliforms in water and food samples.

M1294/
M2056

Limitations

1. β -glucuronidase is present in 97% of *E. coli* strains, however few *E. coli* may be negative.
2. Certain species of *Shigella* and *Salmonella* are β -glucuronidase positive which may appear as light blue.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

Appearance of powder : Light yellow to pink (M1294) / Cream to yellow (M2056) coloured, homogeneous, free flowing powder.

Gelling : Firm, comparable with 1.0% Agar gel (M1294) Firm, comparable with 1.5% Agar gel (M2056)

Colour and Clarity of prepared medium : Light pink coloured (M1294) / Light yellow coloured (M2056), clear to slightly opalescent gel forms in Petri plates

Reaction : Reaction of 2.65% w/v aqueous solution of M1294 at 25°C. pH : 6.8 \pm 0.2.
Reaction of 3.93% w/v aqueous solution of M2056 at 25°C. pH : 7.00 \pm 0.2.

Cultural Response : Cultural characteristics observed after an incubation at 35-37°C for 18-24 hours.

Organisms (ATCC)	Inoculum (CFU)	M1294 Growth w/FD190	M2056 Growth	Recovery	Colour of colony (M1294)	Colour of colony (M2056)
<i>Escherichia coli</i> (25922) (00013*)	50-100	good - luxuriant	good - luxuriant	\geq 50%	dark blue to violet ●	dark blue
<i>Escherichia coli</i> O157:H7 (NCTC 12900)	50-100	luxuriant	-	\geq 50%	salmon to red ●	-
# <i>Klebsiella aerogenes</i> (13048) (00175*)	50-100	luxuriant	luxuriant	\geq 50%	salmon to red	Pink
<i>Klebsiella pneumoniae</i> ATCC 13883 (00097*)	50-100	luxuriant	luxuriant	\geq 50%	-	Pink
<i>Citrobacter freundii</i> (8090)	50-100	luxuriant	luxuriant	\geq 50%	salmon to red (big)	Pink

<i>Salmonella</i> Enteritidis (13076) (00030*)	50-100	good	-	40-50%	colourless	-
<i>Shigella flexneri</i> (29508)	50-100	good	-	40-50%	light blue to turquoise	-
<i>Enterococcus faecalis</i> (29212) (00087*)	\geq 10 ³	inhibited	inhibited	0%	-	-

Key : ● : positive reaction, confirmation of red colour around the colony by addition of Kovac's reagent (R008)

* : corresponding WDCM Numbers

: Formerly known as *Enterobacter aerogenes*

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (9, 10).

References

1. Frampton E. W., Restaino L. and Blaszkowski N., 1988, J. Food Prof., 51:402.
2. Hartman, P.A., 1989 Turano, A. (Ed.), Brixia Academic Press, Brescia, Italy, pp. 290-308.
3. Kilian M. and Bulow P., 1976, Acta. Pathol. Microbiol. Scand Sect. B, 84:245.
4. LeMinor L. and Hamida F., 1962, Ann. Inst. Pasteur 102:267.
5. Manafi M. and Kneifel W., 1989, Zentralbl. Hyg., 189:225.
6. Baird R.B., Eaton A.D., and Rice E.W., (Eds.), 2015, Standard Methods for the Examination of Water and Wastewater, 23rd ed., APHA, Washington, D.C
7. American Public Health Association, Standard Methods for the Examination of Dairy Products, 1978, 14th Ed., Washington D.C.
8. Salfinger Y., and Tortorello M.L. Fifth (Ed.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.
9. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
10. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

Ready Prepared Media

Code	Product Name	Usage	Packing
Category : HiTouch™ FlexiPlates			
FL022	HiTouch™ ECC Count Flexi Plate™	for enumeration (count) of <i>Escherichia coli</i> and coliforms.	50 plts
Category : DriFilter Membrane Nutrient Pad			
MF028	ECC Selective Medium (without Membrane Filter)	for detection and enumeration of total coliforms and <i>E. coli</i> based on chromogenic differentiation.	50 plts

For identification and differentiation of *E. coli* and Total coliforms

HiCrome™ Coliform Agar w/ SLS/ Coliform Broth w/ SLS

Recommended for the simultaneous detection of *Escherichia coli* and total coliforms in water and food samples.

M1300/
M1826

Composition **	M1300	M1826
Ingredients	Grams/Litre	Grams/Litre
Peptone, special	3.00	3.00
Sodium chloride	5.00	5.00
Dipotassium hydrogen phosphate	3.00	3.00
Potassium dihydrogen phosphate	1.70	1.70
Sodium pyruvate	1.00	1.00
L-Tryptophan	1.00	1.00
Sodium lauryl sulphate (SLS)	0.10	0.10
Chromogenic mixture	0.20	-
Chromogenic substrate	-	0.30
Agar	12.00	-
Final pH (at 25°C)	6.8 ± 0.2	6.8 ± 0.2

** Formula adjusted, standardized to suit performance parameters

Directions

Suspend 27 grams of M1300 or 15.10 grams of M1826 in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. For M1826 dispense into tubes or flasks or as desired. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Add 5mg/ l novobiocin before autoclaving the medium, when a high number of gram-positive accompanying bacteria are expected. Cool to 45-50°C. Mix well and pour into sterile Petri plates for M1300.

Principle and Interpretation

HiCrome™ Coliform Agar / Broth w/ SLS is a selective medium recommended for the simultaneous detection of *Escherichia coli* and total coliforms in water and food samples (5). Peptone special and Sodium pyruvate provides nitrogenous, carbonaceous, compounds, long chain amino acid and other essential growth nutrients. The phosphates buffer the medium well. The medium composition helps even the sublethally injured coliforms to grow rapidly. Sodium lauryl sulphate inhibits gram-positive organisms.

The chromogenic mixture contains two chromogenic substrates. The enzyme β -galactosidase produced by coliforms cleaves one chromogen, resulting in the salmon to red colouration of coliform colonies. The enzyme β -glucuronidase produced by *E. coli*, cleaves X-glucuronide. *E. coli* forms dark blue to violet coloured colonies due to cleavage of both the chromogens (1, 3, 4). The addition of L-Tryptophan improves the indole reaction, thereby increasing detection reliability in combination with the two chromogens. In M1826 chromogenic substrate in the medium help to detect β -glucuronidase positive *E. coli*. Other gram negative bacteria forms colourless colonies, except some organisms which are β -glucuronidase positive. β -glucuronidase positive organisms gives light blue to turquoise colonies. GUD is present in 94–96% of *E. coli* strains and in some *Salmonella*, *Shigella* and *Yersinia* spp (2). To confirm *E. coli*, add a drop of Kovac's reagent (R008) on the dark-blue

to violet colony/Blue colour broth. Formation of cherry-red colour indicates positive reaction.

Type of specimen

Water samples ; Food samples.

Specimen Collection and Handling


For water samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (6).

For food samples, follow appropriate techniques for sample collection and processing as per guidelines (7).


After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

Read the label before opening the container. Wear protective gloves/



M1300 HiCrome™ Coliform Agar w/SLS



M1826 Coliform Broth w/SLS

- Control
- Citrobacter freundii* ATCC 8090
- Klebsiella pneumoniae* ATCC 13883
- Salmonella* Enteritidis ATCC 13076
- Shigella flexneri* ATCC 12022
- Escherichia coli* ATCC 25922 (00013*)
- Enterococcus faecalis* ATCC 29212 (00087*)
- Staphylococcus aureus* subsp aureus ATCC 25923 (00034*)

HiCromeVeg™ Freedom from BSE / TSE worries
Single Streak Rapid Differentiation Series

HiCrome™ Coliform Agar w/ SLS (M1300) is also available as HiCrome™ Coliform HiVeg™ Agar w/ SLS (MV1300) wherein all the animal origin nutrients have been replaced by vegetable based nutrients.

HiCrome™ Coliform Agar w/ SLS/ Coliform Broth w/ SLS

Recommended for the simultaneous detection of *Escherichia coli* and total coliforms in water and food samples.

M1300/
M1826

protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets

Limitations

1. β -glucuronidase is present in 97% of *E. coli* strains, however few *E. coli* may be negative.
2. Certain species of *Shigella* and *Salmonella* are β -glucuronidase positive which may appear as light blue.
3. Further biochemical and serological test must be carried out for confirmation of suspected colonies.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

Appearance of powder : Light yellow to beige coloured (M1300)/Cream to yellow coloured (M1826), homogeneous, free flowing powder.

Gelling : Firm, comparable with 1.2% Agar gel (M1300).

Colour and Clarity of prepared medium : Colourless, clear to very slightly opalescent gel forms in Petri plates (M1300).
Cream, clear to slightly opalescent solution, may have slight precipitate (M1826).

Reaction : Reaction of 2.70% w/v (M1300) / 1.51% w/v (M1826) aqueous solution of M1300 at 25°C. pH: 6.8 ± 0.2.

Cultural Response : Cultural characteristics observed after an incubation at 35-37°C for 18-24 hours (48 hours if necessary).

Organisms (ATCC)	Inoculum (CFU)	Growth	Recovery (M1300)	Colour of colony (M1300)	Colour of medium (M1826)
<i>Escherichia coli</i> (25922) (00013*)	50-100	good - luxuriant	>=50%	dark blue / violet*	blue
<i>Enterobacter cloacae</i> (23355) (00082*)	50-100	good - luxuriant	>=50%	salmon to red	-
<i>Citrobacter freundii</i> (8090)	50-100	good - luxuriant	>=50%	salmon to red	colourless
<i>Escherichia coli</i> (35218)	50-100	-	>=50%	-	blue
<i>Klebsiella pneumoniae</i> (13883) (00097*)	50-100	good - luxuriant	>=50%	light pink	colourless
<i>Salmonella</i> Enteritidis (13076) (00030*)	50-100	good	40-50%	colourless	colourless

<i>Shigella flexneri</i> (12022) (00126*)	50-100	good	40-50%	colourless	colourless
<i>Enterococcus faecalis</i> (29212) (00087*)	>=10 ³	inhibited	0%	-	-
<i>Staphylococcus aureus</i> subsp <i>aureus</i> (25923) (00034*)	>=10 ³	-	0%	-	-
<i>Staphylococcus aureus</i> subsp <i>aureus</i> (6538) (00032*)	>=10 ³	-	0%	-	-

Key : * : positive reaction, confirmation of red colour around the colony by addition of Kovac's reagent (R008) - : negative reaction.

* : corresponding WDCM Numbers

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (8, 9).

References

1. Frampton E.W., Restaino L. and Blaszkowski N., (1988), J. Food Prot., 51:402.
2. Hartman, P.A., 1989 Turano, A. (Ed.), Brixia Academic Press, Brescia, Italy, pp. 290-308.
3. Kilian M. and Bülow P., (1976), Acta. Pathol. Microbiol. Scand., Sect. B, 84:245.
4. LeMinor L. and Hamida F., (1962), Ann. Inst. Pasteur (Paris), 102:267.
5. Manafi M. and Kneifel W., (1989), Zentralbl. Hyg., 189:225.
6. Baird R.B., Eaton A.D., and Rice E.W., (Eds.), 2015, Standard Methods for the Examination of Water and Wastewater, 23rd ed., APHA, Washington, D.C
7. Salfinger Y., and Tortorello M.L. Fifth (Ed.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.
8. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
9. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S. and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

Ready Prepared Media

Code	Product Name	Usage	Packing
Category : DriFilter Membrane Nutrient Pad			
MF026E	#Coliform Medium w/ SLS (Economy Pack)(without Membrane Filter)	for detection and enumeration of total coliforms and <i>E. coli</i> based on chromogenic differentiation	50 pcs
MF026F	Coliform Medium w/ SLS w/ Sterile Membrane Filter	for detection and enumeration of total coliforms and <i>E. coli</i> based on chromogenic differentiation	50 plts

HiCrome™ Rapid ECC Broth

Recommended for rapid detection of *Escherichia coli* and other *Enterobacteriaceae* from water samples.

M2011

Composition **

Ingredients	Grams/Litre
Peptone special	24.00
Sodium chloride	5.00
Disodium hydrogen phosphate	1.00
Sodium thiosulphate	5.00
Ferric citrate	1.00
Lactose	5.00
Phenol red	0.018
Selective mix	1.50
Chromogenic substrate	3.83

Final pH (at 25°C) 7.4 ± 0.2

** Formula adjusted, standardized to suit performance parameters

Directions

Suspend 46.35 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. DO NOT AUTOCLAVE. Mix well and dispense into sterile tubes or flasks as desired.

Principle and Interpretation

HiCrome™ Rapid ECC Broth is designed for detection and confirmation of *Escherichia coli* and other coliforms from water samples. The major microbial water contaminants are coliforms include *Escherichia coli*, *Klebsiella pneumoniae*, *Salmonella*, *Citrobacter*, *Vibrio*, and *Pseudomonas* (1). This test was designed for the rapid detection and differentiation of these organisms.

Peptone special provides nitrogen and carbon source, long chain amino acids, vitamins and other essential growth nutrients. Phosphates buffer the medium. Lactose is the fermentable carbohydrate and phenol red is the indicator. Lactose fermenting organisms gives yellow colour to the medium while lactose non-fermentors gives pink to red colour. The chromogenic substrate is used to detect the presence of β-D-glucuronidase produced by *E.coli* thus imparting blue colour to the medium. However since *E.coli* also ferments lactose, the presence of *E.coli* is indicated by bluish green to green colour. The detection of H₂S production is enhanced by the presence of specific H₂S detectors. The medium turns black in case of H₂S producers such as *Salmonella*, *Citrobacter* etc are present. The phosphate salts provide buffering action for rapid growth of coliforms. Sodium chloride helps to maintain the osmotic balance. Selective mix present in the medium suppresses the growth of gram positive microorganisms. Recovery of these pathogens is faster and reliable.

Type of specimen

Water and waste water samples

Specimen Collection and Handling

For water samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (1, 2, 3).

After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets

Limitations


1. β-glucuronidase is present in 97% of *E.coli* strains, however few *E.coli* may be negative.
2. β-glucuronidase negative *E.coli* will impart yellow colour to the medium.
3. Further confirmation must be carried out for *E.coli* by indole test.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

- Appearance of Powder :** Light yellow to pink homogeneous free flowing powder
- Colour and Clarity :** Red coloured clear solution in tubes
- Reaction :** Reaction of 4.63% w/v aqueous solution at 25°C. pH : 7.4±0.2
- Cultural Response :** Cultural characteristics observed after an incubation at 35-37°C for 12-18 hours.



1. Control
2. *Escherichia coli* ATCC 25922 (00013*)
3. *Klebsiella pneumoniae* ATCC (13883) (00097*)
4. *Citrobacter freundii* ATCC 8090
5. *Salmonella* Typhimurium ATCC 14028 (00031*)

HiCrome™ Rapid ECC Broth

Recommended for rapid detection of *Escherichia coli* and other *Enterobacteriaceae* from water samples.

M2011

Organism (ATCC)	Inoculum (CFU)	Growth	Colour change in medium
<i>Escherichia coli</i> ATCC 25922 (00013*)	50-100	luxuriant	green
<i>Klebsiella pneumoniae</i> ATCC (13883) (00097*)	50-100	luxuriant	yellow
<i>Citrobacter freundii</i> ATCC 8090	50-100	luxuriant	black
<i>Salmonella</i> Typhimurium ATCC 14028 (00031*)	50-100	luxuriant	black
<i>Enterococcus faecalis</i> ATCC 29212 (00087*)	>=10 ³	inhibited	
<i>Staphylococcus aureus</i> subsp <i>aureus</i> (25923) (00034*)	>=10 ³	inhibited	

Key : * = corresponding WDCM Numbers

Storage and Shelf-life

Store dehydrated powder and prepared medium at 2-8°C in tightly closed container. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation

due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label. Product performance is best if used within stated expiry period.

Disposal


User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (4, 5).

References

1. Methods for Examination of Waters and Associated Materials, Environment Agency, 1998, Standing Committee of Analysts.
2. Greenberg A. E., Clesceri L. S. and Eaton A. D., (Eds.), 2005, Standard Methods for the Examination of Water and Wastewater, 21st ed., APHA, Washington, D.C.
3. Baird R.B., Eaton A.D., and Rice E.W., (Eds.), 2015, Standard Methods for the Examination of Water and Wastewater, 23rd ed., APHA, Washington, D.C.
4. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
5. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.



Rapid - Reliable



HiCrome™ M-Coliform Differential Agar Base

Recommended as a selective and differential agar for the detection of coliform bacteria using membrane filtration technique.

M1951

Composition **

Ingredients	Grams/Litre
Peptone	5.00
Tryptone	10.00
Yeast extract	3.00
Lactose	12.50
Sodium deoxycholate	0.15
Aniline Blue	0.10
Chromogenic substrate	0.50
Agar	15.00

Final pH (at 25°C) 7.2 ± 0.2

** Formula adjusted, standardized to suit performance parameters

Directions

Suspend 46.25 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. DO NOT AUTOCLAVE. Cool to 45- 50°C. Aseptically add the rehydrated contents of one vial of Monensin Selective supplement (FD309). Mix well and pour into sterile Petri plates.

Principle and Interpretation

HiCrome™ M-Coliform Differential Agar Base is based on coliform enumeration medium, M-FC Agar (1). This medium was modified for detection and enumeration of total coliforms by addition of Monensin supplement to improve the recovery of injured coliforms (2).

Peptone, Tryptone and Yeast extract provides nitrogenous compounds, carbonaceous compounds, long chain amino acids and other growth nutrients and vitamins. Lactose is the fermentable carbohydrate. Monensin and sodium deoxycholate acts as selective agents, inhibiting Gram-positive bacteria. Aniline blue forms the indicator system of the medium. The chromogenic substrate is utilized by *E.coli* which detects the presence of β -glucuronidase. The medium helps injured coliforms to grow in the presence of selective agents.

Type of specimen

Water samples

Specimen Collection and Handling

For water samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (3).

After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions


Read the label before opening the container. Wear protective gloves/ protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets

Limitations

1. β -glucuronidase is present in 97% of *E.coli* strains, however few *E.coli* may be negative.
2. Since the medium is highly selective, some strains may show poor growth due to nutritional variations.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.



M1951 – HiCrome™ M-Coliform Differential Agar Base

HiCrome™ M-Coliform Differential Agar Base

Recommended as a selective and differential agar for the detection of coliform bacteria using membrane filtration technique.

M1951

Quality Control

- Appearance of Powder** : Light yellow to greyish yellow homogeneous free flowing powder
- Gelling** : Firm, comparable with 1.5% Agar gel.
- Colour and Clarity of prepared medium** : Light yellow, clear to slightly opalescent gel forms in Petri plates
- Reaction** : Reaction of 4.63% w/v aqueous solution at 25°C. pH : 7.2 ± 0.2
- Cultural Response** : Cultural characteristics observed after an incubation at 35-37°C for 18-24 hours

Organisms (ATCC)	Inoculum (CFU)	Growth	Recovery	Colour of colony (On membrane filter)
<i>Escherichia coli</i> (25922) (00013*)	50-100	good-luxuriant	≥50%	blue
<i>Proteus vulgaris</i> (13315)	50-100	good-luxuriant	≥50%	tan
<i>Bacillus spizizenii sub spizizenii</i> (6633) (00003*)	≥10 ³	inhibited	0%	

Key : * : corresponding WDCM Numbers

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (4, 5).

References

1. Brodsky, M. H., P. Entis, A. N. Sharpe, and G. A. Jarvis. 1982. Enumeration of indicator organisms in foods using the automated hydrophobic membrane filter technique. *J. Food Prod.* 45:292-296.
2. Entis, P., and P. Boleszczuk. 1990. Direct enumeration of coliforms and *Escherichia coli* by hydrophobic grid membrane filter in 24 hours using MUG. *J. Food Prot.* 53:948-952.
3. Baird R.B., Eaton A.D., and Rice E.W., (Eds.), 2015, Standard Methods for the Examination of Water and Wastewater, 23rd ed., APHA, Washington, D.C
4. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
5. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1

For the detection of *E. coli* and total coliforms

HiCrome™ M-Coliconfirm Broth Base

Recommended for detection of *E. coli* and other total coliforms in water samples by membrane filtration.

M2064

Composition **

Ingredients	Grams/Litre
Tryptone	8.00
Yeast extract	0.50
Lactose	0.60
Sodium chloride	3.00
Dipotassium hydrogen phosphate	1.75
Potassium dihydrogen phosphate	1.25
Sodium pyruvate	1.00
Octyphenol ethoxylate	0.50
Magnesium sulphate	0.30
Sodium azide	0.02
L-Methionine	0.10
Methylene blue	0.016
Cyclohexylammonium salt	0.20
Chromogenic mixture	0.20

Final pH (at 25°C) 7.0±0.2

** Formula adjusted, standardized to suit performance parameters

Directions

Suspend 17.43 grams in 1000 ml distilled water. Heat if necessary to dissolve the medium completely. DO NOT AUTOCLAVE. Cool to 45-50°C. Aseptically add the rehydrated contents of ECC Selective Supplement (FD344) and 7ml of TTC Solution, 1% (FD057). Mix well and aseptically add desired quantity (2 to 5 ml) of broth on sterile absorbent cotton pad or sterile filter paper for saturation. The nutrient pad should be used within 24 hours of saturation.

Principle and Interpretation

This is a selective medium recommended for the simultaneous detection of *Escherichia coli* and total coliforms in water (1). The water sample is filtered through membranes and then placed on pad saturated with medium and incubated at 35 ±5°C for 24 hours in sealed Petri plates. Tryptone provides nitrogenous and carbonaceous compounds, long chain amino acids, vitamins and other essential nutrients. Yeast extract serves as a source of vitamins. Lactose is the fermentable carbohydrate. The phosphates in the medium buffers the medium. Sodium chloride maintains the osmotic balance. The enzyme beta-glucuronidase produced by *E. coli* utilizes the chromogenic substrate to produce blue-

purple coloured colonies. Coliforms other than *Escherichia coli* turn red as they reduce TTC (2,3,5-triphenyl tetrazolium chloride). Thus, the resulting colour distinction allows simple interpretation of test without further confirmation. Methylene blue and ECC selective supplement containing imparts selectivity to the medium. Non-coliforms usually give white coloured colonies.

Type of specimen

Water samples

Specimen Collection and Handling

For water samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (1).

After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

Read the label before opening the container. Wear protective gloves/ protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets



M2064 HiCrome™ M-Coliconfirm Broth Base

For the detection of *E. coli* and total coliforms

HiCrome™ M-Coliconfirm Broth Base

Recommended for detection of *E. coli* and other total coliforms in water samples by membrane filtration..

M2064

Limitations

1. β -glucuronidase is present in 97% of *E. coli* strains, however few *E. coli* may be negative.
2. Since the medium is highly selective, some strains may show poor growth due to nutritional variations.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

- Appearance of Powder** : Cream to yellow homogeneous free flowing powder.
- Colour and Clarity** : Cream, clear to slightly opalescent solution, may have slight precipitate.
- Reaction** : Reaction of 1.74% w/v aqueous solution at 25°C. pH : 7.0±0.2
- Cultural Response** : Cultural characteristics observed after an incubation at 34.5-35.5°C for 24 hours.

Organism (ATCC)	Inoculum (CFU)	Growth	Colour of colony on membrane
<i>Citrobacter freundii</i> ATCC 8090	50-100	luxuriant	red
<i>Escherichia coli</i> ATCC 25922 (00013*)	50-100	luxuriant	blue
<i>Escherichia coli</i> ATCC 35218	50-100	luxuriant	blue
<i>Enterococcus faecalis</i> ATCC 29212 (00087*)	>=10 ³	inhibited	-
<i>Klebsiella pneumoniae</i> ATCC 13883 (00097*)	50-100	luxuriant	red
<i>Staphylococcus aureus</i> ATCC 25923 (00034*)	>=10 ³	inhibited	-
<i>Staphylococcus aureus</i> ATCC 6538 (00032*)	>=10 ³	inhibited	-

Key : (*) Corresponding WDCM numbers

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.


Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (2,3).

References

1. Baird R.B., Eaton A.D., and Rice E.W., (Eds.), 2015, Standard Methods for the Examination of Water and Wastewater, 23rd ed., APHA, Washington, D.C.
2. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
3. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.



M2064 HiCrome™ M-Coliconfirm Broth Base

For identification and differentiation of *E. coli* and Total coliforms

M-E. coli Broth

Recommended for the detection, differentiation and enumeration of *Escherichia coli* and coliforms in water samples by membrane filtration technique.

M1426

Composition **

Ingredients	Grams/Litre
Tryptone	20.00
Bile salts mixture	1.50
Chromogenic mixture	0.175

Final pH (at 25°C) 7.2 ± 0.2

** Formula adjusted, standardized to suit performance parameters

Directions

Suspend 21.67 grams in 1000 ml distilled water. Heat if necessary to dissolve the medium completely. DO NOT AUTOCLAVE. Cool to 45-50°C. Aseptically add desired quantity (2 to 5 ml) of broth on sterile absorbent cotton pad or sterile filter paper for saturation. The medium should be used within 24 hours of rehydration.

Principle and Interpretation

M-E.coli Broth is used for detection and differentiation of *Escherichia coli* and coliforms in water samples using membrane filter technique. It is based on Tryptone Bile Agar used for detection of *Escherichia coli* in foods (1) where recovery of *Escherichia coli* is faster, more reliable and accurate.

The water sample is filtered through membranes and then placed on pad saturated with M-E.coli Broth and incubated at 37°C in sealed Petri plates. Glucuronidase test is used increasingly for detection of *E. coli* in water and food microbiology as *E. coli* is an important indicator of fecal contamination in samples from the food processing and water purification plants. Other *Escherichia* spp. do not produce this enzyme (3). The medium contains chromogenic mixture which helps to detect glucuronidase activity of *Escherichia coli* (2). This specific enzyme differentiates *Escherichia coli* from other coliforms. *Escherichia coli* cells split the chromogenic mixture with the help of the enzyme glucuronidase to give blue to green colouration to the colonies. Coliforms other than *Escherichia coli* turn red as they reduce TTC (2,3,5-triphenyl tetrazolium chloride). Thus, the resulting colour distinction allows simple interpretation of test without further confirmation.

Tryptone provides nitrogen and carbon source, long chain amino acids, vitamins and other essential growth nutrients to the organisms. Bile salt mixture inhibit gram-positive organisms.

Type of specimen

Water samples


Specimen Collection and Handling

For water samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (4).

After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

Read the label before opening the container. Wear protective gloves/ protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets



M1426 M-E. coli Broth

M-E. coli Broth

M1426

Recommended for the detection, differentiation and enumeration of *Escherichia coli* and coliforms in water samples by membrane filtration technique.

Limitations

1. β -glucuronidase is present in 97% of *E. coli* strains, however few *E. coli* may be negative.
2. Since the medium is highly selective, some strains may show poor growth due to nutritional variations.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

- Appearance of Powder** : Light yellow to beige coloured, homogeneous, free flowing powder.
- Colour and Clarity of prepared medium** : Light yellow coloured, clear solution without any precipitate
- Reaction** : Reaction of 2.17% w/v aqueous solution at 25°C. pH : 7.2 ± 0.2.
- Cultural Response** : Cultural characteristics observed after an incubation at 35-37°C for 18-24 hours.

Organisms (ATCC)	Inoculum (CFU)	Growth	Colour of colony on membrane filter
<i>Escherichia coli</i> (25922) (00013*)	50-100	luxuriant	blue
# <i>Klebsiella aerogenes</i> (13048) (00175*)	50-100	luxuriant	red
<i>Staphylococcus aureus</i> subsp <i>aureus</i> (25923) (00034*)	≥10 ³	inhibited	—

Key : * : corresponding WDCM Numbers

: Formerly known as *Enterobacter aerogenes*

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (5, 6).

References

1. Anderson J. M. and Baird Parker A.C., (1975), J. Appl. Bact., 39:111.
2. Hansen W. and Yourassawsky E., (1984), J. Clin. Microbiol. 20:1177.
3. Rice, E.W., Allen, M.J., Brenner, D.J., Edberg, S.C., 1991. Appl. Environ. Microbiol. 57, 592–593.
4. Baird R.B., Eaton A.D., and Rice E.W., (Eds.), 2015, Standard Methods for the Examination of Water and Wastewater, 23rd ed., APHA, Washington, D.C
5. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
6. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

Ready Prepared Media			
Code	Product Name	Usage	Packing
Category : DriFilter Membrane Nutrient Pad			
MF027	M-E.coli Medium (without Membrane Filter)	for detection and enumeration of total coliforms and <i>E. coli</i> based on chromogenic differentiation	50 plts
MF027E	M-E.coli Medium (Economy Pack) (without Membrane Filter)	for detection and enumeration of total coliforms and <i>E. coli</i> based on chromogenic differentiation	50 pcs
MF027F	M-E. coli Medium w/ Sterile Membrane Filter	for detection and enumeration of total coliforms and <i>E. coli</i> based on chromogenic differentiation	50 plts

For identification and differentiation of *E. coli* and Total coliforms

HiCrome™ M-Lauryl Sulphate Agar

Recommended for the differentiation and enumeration of *Escherichia coli* and other coliforms by a single membrane filtration technique

M1569

Composition **

Ingredients	Grams/Litre
Peptone	40.00
Yeast extract	6.00
Lactose	30.00
Phenol red	0.20
Sodium lauryl sulphate (SLS)	1.00
Sodium pyruvate	0.50
Chromogen	0.20
Agar	10.00

Final pH (at 25°C) 7.4 ± 0.2

** Formula adjusted, standardized to suit performance parameters

Directions

Suspend 87.9 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45- 50°C. Mix well and pour into sterile Petri plates.

Principle and Interpretation

HiCrome™ M-Lauryl Sulphate Agar is a modification of the Lauryl Tryptose Broth, formulated by Mallman and Darby, (1). This chromogenic medium is recommended for the presumptive identification and differentiation of *Escherichia coli* and other coliforms by a single membrane filtration technique (2, 3). The incorporation of chromogen X-glucuronide and the dye phenol red favours the differentiation of *E. coli* and other coliforms on the basis of colour.

Peptone provide nitrogenous and carbonaceous compounds, long chain amino acids and other essential growth nutrients to the organisms. Yeast extract serves as a source of vitamins especially group B vitamins. Lactose acts as a source of fermentable sugar while sodium lauryl sulphate inhibits gram positive organisms. The enzyme β-D-glucuronidase produced by *E. coli*, cleaves X-glucuronide, imparting a green colour to the colonies. Lactose fermentation is detected by phenol red indicator. Other lactose fermentors not possessing β-D-Glucuronidase enzyme will show yellow colonies, whereas lactose non-fermentors will exhibit pink coloured colonies.

Type of specimen

Water samples

Specimen Collection and Handling


For water samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (4). After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

Read the label before opening the container. Wear protective gloves/ protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets

Limitations

1. β-glucuronidase is present in 97% of *E. coli* strains, however few *E. coli* may be negative.
2. Some species may show poor growth due to nutritional variations.



M1569 – HiCrome™ M-Lauryl Sulphate Agar

For identification and differentiation of *E. coli* and Total coliforms

HiCrome™ M-Lauryl Sulphate Agar

Recommended for the differentiation and enumeration of *Escherichia coli* and other coliforms by a single membrane filtration technique

M1569

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

- Appearance of Powder** : Light yellow to pink coloured, homogeneous, free flowing powder.
- Gelling** : Firm, comparable with 1.0% Agar gel.
- Colour and Clarity of prepared medium** : Red coloured, clear to slightly opalescent gel forms in Petri plates.
- Reaction** : Reaction of 8.8% w/v aqueous solution at 25°C. pH:7.4 ± 0.2 .
- Cultural Response** : Cultural characteristics observed after an incubation at 35-37°C for 18-24 hours.

Organisms (ATCC)	Inoculum (CFU)	Growth	Colour of colony
<i>Escherichia coli</i> (25922) (00013*)	50-100	luxuriant	green
<i>Klebsiella pneumoniae</i> (13883) (00097*)	50-100	good	yellow, mucoid
<i>Staphylococcus aureus</i> subsp. aureus (25923) (00034*)	≥10 ³	inhibited	-
<i>Salmonella</i> Enteritidis (13076) (00030*)	50-100	good	pink

Key : * : corresponding WDCM Numbers

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.


Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (5, 6).

References

- Mallman and Darby, 1941, Am. J. Public Health, 31:127.
- Methods for Examination of Waters and Associated Materials, Environment Agency, 1998, Standing Committee of Analysts.
- Sartory D.P. and Howard L, 1992, Lett Appl. Microbiol. 15:273-276.
- Baird R.B., Eaton A.D., and Rice E.W., (Eds.), 2015, Standard Methods for the Examination of Water and Wastewater, 23rd ed., APHA, Washington, D.C
- Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
- Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1



M1569 – HiCrome™ M-Lauryl Sulphate Agar

Ready Prepared Media			
Code	Product Name	Usage	Packing
Category : 55 mm Scored Polystyrene Plates			
SP1569	HiCrome™ M-Lauryl Sulphate Agar Plate	for differentiation and enumeration of <i>Escherichia coli</i> and other coliforms by single membrane filtration	100 plts

For the detection of *E. coli* and coliforms

HiCrome™ Chromogenic Coliform Agar (CCA)

Recommended for detection of *Escherichia coli* and coliforms in water samples. The composition and performance criteria of this medium are as per the specifications laid down in ISO 9308-1:2014.

M1991I

Composition **	
Ingredients	Grams/Litre
Tryptone#	1.000
Yeast extract	2.000
Sodium chloride	5.000
Sodium dihydrogen phosphate, dihydrate	2.200
Disodium hydrogen phosphate	2.700
Sodium pyruvate	1.000
Sorbitol	1.000
Tryptophan	1.000
Tergitol-7	0.150
6-chloro-3-indoxyl-β-D-galactopyranoside	0.200
5-bromo-4-chloro-3-indoxyl-β-D-glucuronic acid cyclohexamine ammonium salt, monohydrate	0.100
IPTG (Isopropyl-β-D-thiogalactopyranoside)	0.100
Agar	15.000

Final pH (at 25°C) 6.8 ± 0.2

** Formula adjusted, standardized to suit performance parameters

Enzymatic digest of casein

Directions

Suspend 30.92 grams (the equivalent weight of dehydrated medium per litre) in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. DO NOT AUTOCLAVE. DO NOT OVERHEAT. Cool to 45-50°C. Mix well and pour into sterile Petri plates.

Principle and Interpretation

HiCrome™ Chromogenic Coliform Agar is a selective medium recommended for the simultaneous detection of *Escherichia coli* and total coliforms in water samples (1). The medium contains three chromogenic substrates. The enzyme β-D-galactosidase produced by coliforms cleaves 6-chloro-3-indoxyl-β-D-galactopyranoside to form pink to red coloured colonies (3). The enzyme β-D-glucuronidase produced by *E. coli*, cleaves 5-bromo-4-chloro-3-indoxyl-β-D-glucuronic acid (2). Colonies of *E. coli* give dark blue to violet coloured colonies due to cleavage of both the chromogens. The presence of the third chromogen IPTG enhances the colour reaction. Addition of L-Tryptophan improves the indole reaction thereby increasing the detection reliability. Tryptone, yeast extract, sodium pyruvate and sorbitol provide nitrogenous substances, fermentable carbohydrate and other essential growth nutrients for the organisms. Phosphates buffer the medium. The media formulation helps even sublethally injured coliforms to recover

and grow rapidly. Tergitol-7 inhibits gram-positive as well as some gram-negative bacteria other than coliforms (3).

The medium is inoculated either by pour plate technique or by spreading the sample on the surface of plated medium. Membrane filter technique can also be used. To confirm *E. coli*, add a drop of Kovacs reagent on the dark blue to violet colony. Formation of cherry red colour indicates a positive reaction.

Type of specimen

Water samples


Specimen Collection and Handling

For water samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (4).

After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets



M1991I- HiCrome™ Chromogenic Coliform Agar (CCA)

For the detection of *E. coli* and coliforms

HiCrome™ Chromogenic Coliform Agar (CCA)

Recommended for detection of *Escherichia coli* and coliforms in water samples. The composition and performance criteria of this medium are as per the specifications laid down in ISO 9308-1:2014.

M1991I

Limitations

- β-glucuronidase is present in 97% of *E. coli* strains, however few *E. coli* may be negative hence *E. coli* species may show pink to red colonies.
- Certain species of *Shigella* and *Salmonella* are β-glucuronidase positive, hence they appear light blue to turquoise colonies.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

- Appearance of Powder** : Cream to yellow homogeneous free flowing powder
- Gelling** : Firm, comparable with 1.5% Agar gel.
- Colour and Clarity of prepared medium** : Light yellow coloured opalescent gel forms in Petri plates
- Reaction** : Reaction of 3.09% w/v aqueous solution at 25°C. pH : 6.8 ± 0.2
- Cultural Response** : Cultural characteristics observed after an incubation at 34-38°C for 24 hours.

Organism (ATCC)	Inoculum (CFU)	Growth	Recovery	Colour of Colony
<i>Citrobacter freundii</i> (43864) (00006*)	50-100	luxuriant	≥70%	pink to red
# <i>Klebsiella aerogenes</i> (13048) (00175*)	50-100	luxuriant	≥70%	pink to red
<i>Escherichia coli</i> (25922) (00013*)	50-100	luxuriant	≥70%	dark blue to violet
<i>Enterococcus faecalis</i> (29212) (00087*)	≥10 ³	inhibited	0%	-
<i>Pseudomonas aeruginosa</i> (27853) (00025*)	50-100	luxuriant	≥70%	colourless

Key: * : corresponding WDCM Numbers

: Formerly known as *Enterobacter aerogenes*

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (5, 6).

References

- International Organization for Standardization. Water quality: Enumeration of *E. coli* and coliform bacteria. Part I Membrane filtration methods for bacteria with low bacterial background flora. ISO 9308-1:2014.
- Kilian M. and Bülow P., 1976, Acta. Pathol. Microbiol. Scand Sect. B, 84:245.
- Manafi M. and Kneifel W., 1989, Zentralbl. Hyg., 189:225.
- Baird R.B., Eaton A.D., and Rice E.W., (Eds.), 2015, Standard Methods for the Examination of Water and Wastewater, 23rd ed., APHA, Washington, D.C.
- Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
- Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

Ready Prepared Media

Code	Product Name	Usage	Packing
Category : 90 mm Ready prepared Plates			
MP1991I	HiCrome™ Chromogenic Coliform Agar Plate (CCA Plate)	for detection of <i>Escherichia coli</i> and coliforms in water samples. The composition and performance criteria of this medium are as per the specifications laid down in ISO 9308-1:2014.	20 plts 50 plts
Category : 55 mm Scored Polystyrene Plates			
SP1991I	HiCrome™ Chromogenic Coliform Agar Plate	for detection of <i>Escherichia coli</i> and coliforms in water samples. The composition and performance criteria of this medium are as per the specifications laid down in ISO 9308-1:2014.	100 plts
Category : Drifilter™ Membrane Nutrient Pad			
MF034	Chromogenic Coliform Medium (without membrane filter)	for detection of <i>E. coli</i> and coliforms in water samples	20 plts 50 plts
MF034F	Chromogenic Coliform Medium w/sterile membrane filter	for detection of <i>E. coli</i> and coliforms in water samples	20 plts 50 plts

HiCrome™ E. coli Agar

Recommended for the detection and enumeration of *Escherichia coli* in foods and water without further confirmation on membrane filter or by indole reagent.

M1295/
M1295I

Composition **

	M1295	M1295I
Ingredients	Grams/Litre	Grams/Litre
Tryptone	14.00	20.00
Peptone, special	5.00	—
Bile salts mixture	1.50	1.50
Disodium hydrogen phosphate	1.00	—
Sodium dihydrogen phosphate	0.60	—
Sodium chloride	2.40	—
X-Glucuronide	0.075	0.075
Agar	12.00	15.00

Final pH (at 25°C) 7.2 ± 0.2

** Formula adjusted, standardized to suit performance parameters

Directions

Suspend 36.57 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Mix well and pour into sterile Petri plates.

Principle and Interpretation

HiCrome™ *E. coli* Agar is based on Tryptone Bile Agar to detect *Escherichia coli* in foods (1), where recovery of *E. coli* is faster, more reliable and accurate. Most of the *E. coli* strains can be differentiated from other coliforms by the presence of enzyme glucuronidase, which is highly specific for *E. coli* (2). Glucuronidase test is used increasingly for detection of *E. coli* in water and food microbiology as *E. coli* is an important indicator of fecal contamination in samples from the food processing and water purification plants. Other *Escherichia* spp. do not produce this enzyme (4). The chromogenic agent X-glucuronide used in this medium helps to detect glucuronidase activity of *E. coli*. *E. coli* cells absorb X-glucuronide and the intracellular glucuronidase enzyme splits the bond between the chromophore and the glucuronide. The released chromophore gives bluish green colouration to the *E. coli* colonies. Formulation of M1295I is in accordance with ISO (3).

Tryptone and peptone special provide the nitrogenous compounds, carbon, amino acids, vitamins and other essential growth nutrients to the organisms. Bile salts mixture inhibits gram-positive organisms. Sodium chloride and phosphates maintain osmotic balance and buffering action respectively.

The surface of the plated medium is dried before use. Dilute food samples 1:5 or 1:10 with 0.1% (w/v) sterile Peptone Water (M028) and homogenize in a blender or a stomacher. Pipette 0.5 ml or 1.0 ml of the homogenized food sample on to the plate and spread with sterile glass spreader. Incubate the plates at 30°C for 4 hours and then at 44°C for 18 hours.

Type of specimen


Water samples ; Food samples

Specimen Collection and Handling

For water samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (5).

For food samples, follow appropriate techniques for sample collection and processing as per guidelines (6).

After use, contaminated materials must be sterilized by autoclaving before discarding.



M1295 – HiCrome™ *E. coli* Agar

HiCromeVeg™ Freedom from BSE / TSE worries
Single Streak Rapid Differentiation Series

HiCrome™ *E. coli* Agar (M1295) is also available as HiCrome™ *E. coli* HiVeg™ Agar (MV1295) wherein all the animal origin nutrients have been replaced by vegetable based nutrients.

HiCrome™ E. coli Agar

Recommended for the detection and enumeration of *Escherichia coli* in foods and water without further confirmation on membrane filter or by indole reagent.

M1295/
M1295I

Warning and Precautions

Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets

Limitations

1. β -glucuronidase is present in 97% of *E. coli* strains, however few *E. coli* may be negative hence *E. coli* species may show pink to red colonies.
2. Certain species of *Shigella* and *Salmonella* are β -glucuronidase positive, hence they appear light blue to turquoise colonies.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

- Appearance of powder** : Cream to yellow coloured, homogeneous, free flowing powder.
- Gelling** : Firm, comparable with 1.2% Agar gel of M1295 or 1.5% Agar gel of M1295I .
- Colour and Clarity of prepared medium** : Light yellow coloured, clear to slightly opalescent gel forms in Petri plates.
- Reaction** : Reaction of 3.66% w/v aqueous solution at 25°C. pH:7.2 \pm 0.2.
- Cultural Response** : Cultural characteristics observed after an incubation at 44°C for 18-24 hours.

Organisms (ATCC)	Inoculum (CFU)	Growth	Recovery	Colour of colony
<i>Escherichia coli</i> (25922) (00013*)	50-100	luxuriant	\geq 50%	bluish green
<i>Salmonella</i> Enteritidis (13076) (00030*)	50-100	luxuriant	\geq 50%	colourless

<i>Staphylococcus aureus</i> subsp aureus (25923) (00034*)	$\geq 10^3$	inhibited	0%	-
--	-------------	-----------	----	---

Key : * : corresponding WDCM Numbers

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (7, 8).

References

1. Anderson J. M. and Baird-Parker A. C., 1975, J. Appl. Bacteriol., 39:111.
2. Hansen W. and Yourassawsky E., 1984, J. Clin. Microbiol., 20:1177.
3. International Standard ISO 16649:2001. Microbiology of food and animal feeding stuff - horizontal method for the enumeration of presumptive *Escherichia coli*; Part 2: Colony count technique at 44°C using 5-bromo-4-chloro-3-indolyl- β -D-glucornic acid.
4. Rice, E.W., Allen, M.J., Brenner, D.J., Edberg, S.C., 1991. Appl. Environ. Microbiol. 57, 592–593.
5. Baird R.B., Eaton A.D., and Rice E.W., (Eds.), 2015, Standard Methods for the Examination of Water and Wastewater, 23rd ed., APHA, Washington, D.C.
6. Salfinger Y., and Tortorello M.L. Fifth (Ed.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.
7. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
8. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S. and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

Ready Prepared Media			
Code	Product Name	Usage	Packing
Category : 90 mm Ready prepared Plates			
MP1295	HiCrome™ E.coli Agar Plate	for the detection and enumeration of <i>Escherichia coli</i> in foods without further confirmation on membrane filtration or by indole reagent	20 plts 50 plts
Category : HiTouch™ FlexiPlates			
FL002	HiTouch™ E.coli/Coliform Count Flexi Plate™	for enumeration (count) of all coliforms along with differential count of <i>E. coli</i>	50 plts

For the detection of *E. coli* and thermotolerant coliforms

HiCrome™ Coliform Agar Modified

Recommended for the simultaneous detection of *Escherichia coli* and thermotolerant coliforms in water, milk, dairy products and other food samples.

M1832

Composition **

Ingredients	Grams/Litre
Peptone, special	8.000
Sodium chloride	1.000
Yeast extract	3.000
Potassium dihydrogen phosphate	0.200
Dipotassium phosphate	0.600
Bile salts	0.800
Magnesium sulphate	0.200
Chromogenic mixture	0.200
Agar	10.000

Final pH (at 25°C) 7.2 ± 0.2

** Formula adjusted, standardized to suit performance parameters

Directions

Suspend 24 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Mix well and pour into sterile Petri plates.

Principle and Interpretation

HiCrome™ Coliform Agar Modified is a selective medium recommended for the simultaneous detection of *E. coli* and thermotolerant coliforms in water and food samples (4). Peptone special and yeast extract provide the necessary nitrogen compound, carbon, vitamins and also some trace ingredients required for the growth of organisms. The phosphates buffer the medium well. Magnesium sulphate helps colour development. Bile salts inhibits gram-positive organisms. Sodium chloride maintains osmotic balance.

The chromogenic mixture contains two chromogenic substrates, which enables the detection of two specific enzymes, β -galactosidase and β -glucuronidase. β -galactosidase produced by coliforms cleaves one chromogen, resulting in the pink colouration of coliform colonies. The enzyme β -glucuronidase produced by *E. coli*, cleaves X-glucuronide. *E. coli* forms dark blue to violet coloured colonies due to cleavage of both the chromogens (1, 2, 3). *E. coli* strains that are β -glucuronidase negative (serotype O157:H7) produce pink coloured colonies. Other gram negative bacteria able to grow at (44±0.5)°C produce white or colourless colonies.

Transfer 1 ml of product to analyse and its tenfold dilutions to sterile Petri plates. Pour 12 ml of medium, mix well and allow to solidify. Overlay with 4 ml of medium, allow to solidify and incubate at 43-45°C for 18-24 hours.

Type of specimen

Water, Food and Dairy samples

Specimen Collection and Handling

For water samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (5).

For food and dairy samples, follow appropriate techniques for sample collection and processing as per guidelines (6, 7, 8).


After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets

Limitations

1. β -glucuronidase is present in 97% of *E. coli* strains, however few *E. coli* may be negative.
2. Some species may show poor growth due to nutritional requirements.



M1832 – HiCrome™ Coliform Agar Modified

For the detection of *E. coli* and thermotolerant coliforms

HiCrome™ Coliform Agar Modified

Recommended for the simultaneous detection of *Escherichia coli* and thermotolerant coliforms in water, milk, dairy products and other food samples.

M1832

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

- Appearance of Powder** : Light yellow to beige homogeneous free flowing powder
- Gelling** : Firm, comparable with 1.0% Agar gel.
- Colour and Clarity of prepared medium** : Light yellow to slightly opalescent gel forms in Petri plates
- Reaction** : Reaction of 2.4% w/v aqueous solution at 25°C. pH : 7.2 ± 0.2
- Cultural Response** : Cultural characteristics observed after an incubation at 43-45°C for 24 hours (48 hours if necessary).

Organism (ATCC)	Inoculum (CFU)	Growth	Recovery	Colour of Colony
<i>Escherichia coli</i> (10536)	50-100	good - luxuriant	≥50%	dark blue/violet
<i>Escherichia coli</i> (25922) (00013*)	50-100	good - luxuriant	≥50%	dark blue/violet
<i>Enterobacter cloacae</i> (23355) (00082*)	50-100	good - luxuriant	≥50%	pink
<i>Klebsiella pneumoniae</i> (13883) (00097*)	50-100	good - luxuriant	≥50%	light pink
<i>Enterococcus faecalis</i> (29212) (00087*)	≥10 ³	inhibited	0%	-
<i>Staphylococcus aureus</i> subsp <i>aureus</i> (25923) (00034*)	≥10 ³	inhibited	0%	-

Key : * : corresponding WDCM Numbers

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (9, 10).

References

1. Frampton E. W., Restaino L. and Blaszkowski N., 1988, J. Food Prot., 51:402.
2. Kilian M. and Bülow P., 1976, Acta. Pathol. Microbiol. Scand., Sect. B, 84:245.
3. LeMinor L. and Hamida F., 1962, Ann. Inst. Pasteur (Paris), 102:267.
4. Manafi M. and Kneifel W., 1989, Zentralbl. Hyg., 189:225.
5. Baird R.B., Eaton A.D., and Rice E.W., (Eds.), 2015, Standard Methods for the Examination of Water and Wastewater, 23rd ed., APHA, Washington, D.C
6. American Public Health Association, Standard Methods for the Examination of Dairy Products, 1978, 14th Ed., Washington D.C.
7. Salfinger Y., and Tortorello M.L. Fifth (Ed.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.
8. Wehr H. M. and Frank J. H., 2004, Standard Methods for the Microbiological Examination of Dairy Products, 17th Ed APHA Inc., Washington, D.C.
9. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
10. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S. and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1

Ready Prepared Media			
Code	Product Name	Usage	Packing
Category : 90 mm Ready prepared Plates			
MP1832	HiCrome™ Coliform Agar Plate, Modified	it is a selective agar recommended for the simultaneous detection of <i>Escherichia coli</i> and total coliforms in water and food samples.	20 plts 50 plts

For Identification of thermotolerant *E. coli*

HiCrome™ M-TEC Agar / HiCrome M-TEC Broth

Recommended by the U.S. Environmental Protection Agency (USEPA) for differentiation and enumeration of thermotolerant *Escherichia coli* in water by the membrane filtration technique.

M1571/
M1713

Composition **

	M1571	M1713
Ingredients	Grams/Litre	Grams/Litre
Proteose peptone	5.00	5.00
Yeast extract	3.00	3.00
Lactose	10.00	10.00
Sodium chloride	7.50	7.50
Dipotassium phosphate	3.30	3.30
Potassium dihydrogen phosphate	1.00	1.00
Sodium lauryl sulphate (SLS)	0.20	0.20
Sodium deoxycholate	0.10	0.10
Chromogen	0.50	0.50
Agar	15.00	—

Final pH (at 25°C) 7.3 ± 0.2

** Formula adjusted, standardized to suit performance parameters

Directions

Suspend 45.6 grams of M1571 or 30.6 grams of M1713 in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C Mix well & pour into sterile Petri plates (M1571). In M1713 aseptically add desired quantity (2-5 ml broth) on sterile absorbent pad for saturation in a sterile Petri plate. The medium should be used within 24 hours after rehydration.

Principle and Interpretation

HiCrome™ M-TEC Agar/Broth are the chromogenic media used for detection and enumeration of thermotolerant *Escherichia coli* (TEC) in water by membrane filtration (2). HiCrome™ M-TEC Broth is a modification of the M-TEC Agar developed by Dufour (1). The modified medium contains the chromogen, 5-bromo-6-chloro-3-indolyl-β-D-glucuronide that is cleaved by enzyme β-D-glucuronidase to yield glucuronic acid, produced by *E. coli* strains. This imparts a purple-magenta colour to the colonies of *E. coli* only.

Proteose peptone and yeast extract provides nitrogenous carbonaceous compounds, amino acids and long chain peptides for the growth of microorganisms. Lactose is the fermentable carbohydrate. Sodium chloride maintains osmotic equilibrium. Potassium dihydrogen phosphate and dipotassium phosphate provide strong buffering system to control the pH in the presence of fermentative action. Sodium lauryl sulphate and sodium deoxycholate make the medium more selective by inhibiting gram-positive bacteria. Saturate a sterile cotton absorbent pad with about 2ml of HiCrome™ M-TEC Broth (M1713). Membrane filter through which water sample has been passed is aseptically placed on

the saturated absorbent cotton pad face upwards. This absorbent pad is then incubated at 44.5 ± 0.2°C for 22-24 hours. Following incubation *E. coli* will form purple to magenta coloured colonies on the membrane filters.

Type of specimen

Water samples


Specimen Collection and Handling

For water samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (3).


After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets



M1571 - HiCrome™ M-TEC Agar



M1713 - HiCrome™ M-TEC Broth

For Identification of thermotolerant *E. coli*

HiCrome™ M-TEC Agar / HiCrome M-TEC Broth

Recommended by the U.S. Environmental Protection Agency (USEPA) for differentiation and enumeration of thermotolerant *Escherichia coli* in water by the membrane filtration technique.

M1571/
M1713

Limitations

1. β -glucuronidase is present in 97% of *E.coli* strains, however few *E.coli* may be negative.
2. Overgrowth of non-coliform organisms may interfere with the total coliform organisms.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

- Appearance of Powder** : Cream to yellow homogeneous free flowing powder
- Gelling** : Firm, comparable with 1.5% Agar gel of M1571
- Colour and Clarity of prepared medium** : Light amber coloured, clear to slightly opalescent gel forms in Petri plates (M1571) / clear solution in tubes (M1713).
- Reaction** : Reaction of 4.56% w/v aqueous solution of M1571 and 3.06% w/v of M1713 at 25°C. pH : 7.3 ± 0.2 .
- Cultural Response** : Cultural characteristics observed after an incubation at $44.5 \pm 0.2^\circ\text{C}$ for 22-24 hours.

Organisms (ATCC)	Inoculum (CFU)	Growth	Colour of Colony
<i>Escherichia coli</i> (25922) (00013*)	50-100	good to luxuriant	purple / magenta
<i>Proteus mirabilis</i> (25933)	50-100	good	colourless-light brown
<i>Klebsiella pneumoniae</i> (13883) (00097*)	50-100	good	colourless-tan
<i>Enterococcus faecalis</i> (29212) (00087*)	$\geq 10^3$	inhibited	—

Key : * : corresponding WDCM Numbers

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (4, 5).

References

1. Dufour, Strickland and Cabelli, 1981, Appl. Environ. Microbiol. 41: 1152.
2. U.S. Environmental Protection Agency, 2002, Method 1603; Publication EPA-821-R-02-023.
3. Baird R.B., Eaton A.D., and Rice E.W., (Eds.), 2015, Standard Methods for the Examination of Water and Wastewater, 23rd ed., APHA, Washington, D.C
4. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
5. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1



Escherichia coli

M1571 – HiCrome™ M-TEC Agar

HiCrome™ MacConkey Sorbitol Agar Base

Recommended for the selective isolation of *Escherichia coli* O157:H7 from clinical food and animal feeding stuff.

M1340

Composition **

Ingredients	Grams/Litre
Tryptone	17.00
Proteose peptone	3.00
Sorbitol	10.00
Bile salts mixture	1.50
Sodium chloride	5.00
Crystal violet	0.001
Neutral red	0.03
B.C. Indicator	0.10
Agar	13.50

Final pH (at 25°C) 7.1 ± 0.2

** Formula adjusted, standardized to suit performance parameters

Directions

Suspend 25.06 grams in 495 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Mix well and pour into sterile Petri plates. If desired rehydrated contents of 1 vial of Tellurite-Cefixime Supplement (FD147) may be added aseptically to 495 ml sterile molten, cooled (45-50°C) medium before pouring into sterile Petri plates.

Principle and Interpretation

Sorbitol MacConkey Agar is based on the formulation described by Rappaport and Henigh (4). The medium contains sorbitol instead of lactose and it is recommended for the detection of enteropathogenic strains of *Escherichia coli* O157:H7 that ferments lactose but does not ferment sorbitol (2) and hence produce colourless colonies. Sorbitol fermenting strains of *Escherichia coli* produce pink-red colonies. The red colour is due to production of acid from sorbitol, absorption of neutral red and a subsequent colour change of the dye when pH of the medium falls below 6.8. *Escherichia coli* O157:H7 has been recognised as a cause of hemorrhagic colitis (2). March and Ratnam (3) reported that the detection of *Escherichia coli* O157:H7 had a sensitivity of 100% and specificity of 85% on Sorbitol MacConkey Agar and they recommended this medium as reliable means of screening *Escherichia coli* O157:H7. B.C. indicator is added to detect the presence of the enzyme β -D-glucuronidase which is specific for *Escherichia coli*. (1). Strains of *Escherichia coli* fermenting sorbitol and possessing β -D-glucuronidase appear as blue - purple coloured colonies on the medium. Enteropathogenic strains of *Escherichia coli* O157:H7 do not possess β -D-glucuronidase activity (5) and thus produce colourless colonies.

Tryptone and proteose peptone provide carbonaceous, nitrogenous and other essential growth nutrients. Most of the gram-positive organisms are inhibited by crystal violet and bile salts. Sodium chloride maintains the osmotic equilibrium.

Addition of Tellurite-Cefixime Supplement makes the medium selective (6). Potassium Tellurite selects the serogroups and inhibits *Aeromonas* species and *Providencia* species. Cefixime inhibits *Proteus* species. *Pseudomonas* if present produces colourless colonies on this medium. For confirmation oxidase test may be performed with suspected colonies and results should be noted within 5-10 seconds.

Type of specimen

Clinical, Food and animal feeding stuff, Dairy samples.


Specimen Collection and Handling

For clinical samples follow appropriate techniques for handling specimens as per established guidelines (10, 11).

For food and dairy samples, follow appropriate techniques for sample collection and processing as per guidelines (7, 8, 9).

Warning and Precautions

In Vitro Diagnostic use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets.



M1340 HiCrome™ MacConkey Sorbitol Agar Base

For identification of *E. coli* 0157:H7

HiCrome™ MacConkey Sorbitol Agar Base

Recommended for the selective isolation of *Escherichia coli* 0157:H7 from clinical food and animal feeding stuff.

M1340

Limitations

1. *Pseudomonas* if present produces colourless colonies on this medium. For confirmation oxidase test may be performed with suspected colonies and results should be noted within 5-10 seconds.
2. Some species may show poor growth due to nutritional requirements.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

Appearance of powder : Light yellow to pink coloured, homogeneous, free flowing powder.

Gelling : Firm, comparable with 1.35% Agar gel.

Colour and Clarity of prepared medium : Purplish red coloured, clear to slightly opalescent gel forms in Petri plates.

Reaction : Reaction of 5.01% w/v aqueous solution at 25°C. pH : 7.1 ± 0.2.

Cultural Response : Cultural characteristics observed after an incubation at 35-37°C for 18-24 hours (48 hours if necessary).

Organisms (ATCC)	Inoculum (CFU)	Growth	Recovery	#Colour of colony	Oxidase
<i>Escherichia coli</i> 0157:H7 (NCTC 12900)	50-100	good - luxuriant	>50%	colourless	-
<i>Escherichia coli</i> (25922) (00013*)	50-100	good	40-50%	blue-green	-
<i>Pseudomonas aeruginosa</i> (27853) (00025*)	50-100	fair-good	30-40%	colourless	+
<i>Klebsiella pneumoniae</i> (13883) (00097*)	50-100	good	40-50%	pink-red	-

Key : # = Colour of the colony without addition of Tellurite-Cefixime Supplement (FD147)

+ = positive reaction deep-purple blue colour develops within 10 seconds

- = negative reaction

* = corresponding WDCM Numbers

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (10, 11).

References

1. Hansen W. and Yourassawsky E., (1984), J. Clin. Microbiol., 20:1177.
2. Karmali M.A., Petric M., Lim C., et al, (1985), J. Infect. Dis., 151-775.
3. March S.B. and Ratnam S., (1986) : J. Clin. Microbiol. 23, 869-872.
4. Rappaport F. and Henigh E., (1952), J. Clin. Path., 5:361.
5. Thompson et al. (1990). J. Clin. Microbiol. 29, 2165-2168.
6. Zadik P.M., Chapman P.A. and Siddons C.A., (1993), J. Med. Microbiol., 39, 155-158.
7. American Public Health Association, Standard Methods for the Examination of Dairy Products, 1978, 14th Ed., Washington D.C.
8. Salfinger Y., and Tortorello M.L. Fifth (Ed.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.
9. Wehr H. M. and Frank J. H., 2004, Standard Methods for the Microbiological Examination of Dairy Products, 17th Ed APHA Inc., Washington, D.C.
10. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
11. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S. and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

HiCrome™ EC 0157: H7 Agar, Modified

Recommended for isolation and differentiation of *Escherichia coli* O157:H7 from food and environmental samples.

M1574A

Composition **

Ingredients	Grams/Litre
Tryptone	8.000
Sorbitol	7.00
Bile salts mixture	1.50
Sodium Lauryl Sulphate	0.10
Chromogenic mixture	0.25
Agar	12.00

Final pH (at 25°C) 6.8 ± 0.2

** Formula adjusted, standardized to suit performance parameters

Directions

Suspend 28.85 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Mix well and pour into sterile Petri plates. This medium can be made more selective by aseptically adding 0.25 ml of rehydrated contents of one vial of 1% Potassium Tellurite Solution (FD052) to 1000 ml molten and cooled medium (45-50°C).

Principle and Interpretation

Escherichia coli O157:H7 belongs to the Enterohemorrhagic *Escherichia coli* (EHEC) group and it predominates as a food borne pathogen. *E. coli* O157: H7 was first recognized as a human pathogen in 1982 when two outbreaks of hemorrhagic colitis were associated with consumption of undercooked ground beef that has been contaminated with this organism (1).

HiCrome™ EC 0157:H7 Agar, Modified is a chromogenic medium recommended for the isolation and differentiation of *E. coli* O157:H7 from food and environmental samples. HiCrome™ ECO157:H7 Agar, Modified is based on the formulation described by Rappaport and Henigh (2). The medium contains sorbitol and a proprietary chromogenic mixture instead of lactose and indicator dyes respectively, as is conventionally used. The chromogenic substrate is specifically and selectively cleaved by *E. coli* O157: H7 resulting in a dark purple to magenta coloured moiety. *E. coli* gives bluish green coloured colonies.

Tryptone and yeast extract provides carbonaceous, nitrogenous and growth nutrients. Bile salts mixture and SLS inhibits gram-positive organisms. Potassium tellurite selects the serogroups and inhibits *Aeromonas* species and *Providencia* species.

Type of specimen

Food and Environmental samples.

Specimen Collection and Handling


For food and dairy samples, follow appropriate techniques for sample collection and processing as per guidelines (3, 4, 5).

Warning and Precautions

Read the label before opening the container. Wear protective gloves/ protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets

Limitations

1. Some species may show poor growth due to varying nutritional requirements.
2. Further biochemical test must be carried out for confirmation.



M1574A HiCrome™ EC 0157: H7 Agar, Modified

For identification of *E. coli* O157:H7

HiCrome™ EC 0157: H7 Agar, Modified

Recommended for isolation and differentiation of *Escherichia coli* O157:H7 from food and environmental samples.

M1574A

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

- Appearance of Powder** : Cream to yellow coloured, homogeneous, free flowing powder.
- Gelling** : Firm, comparable with 1.2% Agar gel.
- Colour and Clarity of prepared medium** : Light amber coloured, clear to slightly opalescent gel forms in Petri plates.
- Reaction** : Reaction of 2.88% w/v aqueous solution at 25°C. pH : 6.8 ± 0.2.
- Cultural Response** : Cultural characteristics observed after an incubation at 35-37°C for 18-24 hours.

Organisms (ATCC)	Inoculum (CFU)	Growth	Recovery	Colour of colony
<i>Escherichia coli</i> O157:H7 (NCTC 12900)	50-100	luxuriant	≥50%	dark purple-magenta
<i>Escherichia coli</i> (25922) (00013*)	50-100	luxuriant	≥50%	bluish green
<i>Klebsiella pneumoniae</i> (13883) (00097*)	50-100	luxuriant	≥50%	colourless mauve (mucoid)
<i>Pseudomonas aeruginosa</i> (27853) (00025*)	50-100	luxuriant	≥50%	colourless
<i>Bacillus subtilis</i> sub spizizenii (6633) (00003*)	≥10 ³	inhibited	0%	-
<i>Staphylococcus aureus</i> subsp aureus (25923) (00034*)	≥10 ³	inhibited	0%	-

Key : * : corresponding WDCM Numbers

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition. Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (6, 7).

References

- Downes F. P. and Ito K., (Ed.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 4th Ed., American Public Health Association, Washington, D.C.
- Rappaport F. and Henigh E., 1952, J. Clin. Pathol., 5:361.
- American Public Health Association, Standard Methods for the Examination of Dairy Products, 1978, 14th Ed., Washington D.C.
- Salfinger Y., and Tortorello M.L. Fifth (Ed.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.
- Wehr H. M. and Frank J. H., 2004, Standard Methods for the Microbiological Examination of Dairy Products, 17th Ed APHA Inc., Washington, D.C.
- Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
- Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

Ready Prepared Media			
Code	Product Name	Usage	Packing
Category : Dirfilter Membrane Nutrient Pad			
MF019	EC O157 : H7 Filter Membrane Medium (without Membrane Filter)	for detection and enumeration of enterohaemorrhagic <i>E. coli</i> based on chromogenic differentiation.	20 plts 50 plts
MF019E	EC O157 : H7 Filter Membrane Medium (Economy pack) (without Membrane Filter)	for detection and enumeration of enterohaemorrhagic <i>E. coli</i> based on chromogenic differentiation.	20 plts 50 plts
MF019F	EC O157 : H7 Filter Membrane Medium w/ Sterile Membrane Filter	for detection and enumeration of enterohaemorrhagic <i>E. coli</i> based on chromogenic differentiation.	20 plts 50 plts

For identification of *E. coli* O157:H7

HiCrome™ EC 0157:H7 Selective Agar Base, Modified

Recommended for selective isolation and easy detection of *Escherichia coli* O157:H7 from food samples.

M1575A

Composition **

Ingredients	Grams/Litre
Tryptone	5.00
Yeast extract	3.00
Sorbitol	7.00
Bile salts mixture	1.50
Sodium lauryl sulphate (SLS)	0.10
Chromogenic mixture	0.25
Agar	15.00

Final pH (at 25°C) 6.8 ± 0.2

** Formula adjusted, standardized to suit performance parameters

Directions

Suspend 31.85 grams in 990 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Add rehydrated contents of 1 vial of HiCrome EC 0157:H7 Selective Supplement (FD187) aseptically. Mix well and pour into sterile Petri plates.

Principle and Interpretation

Enterohaemorrhagic *E. coli* strains are also termed as verocytotoxin-producing *E. coli* (VTEC/ EHEC). Although many different serotypes of *Escherichia coli* are known to produce verocytotoxin (2) those of *Escherichia coli* O157:H7 and O157:H are so far the common types causing human infections. O157 VTEC strains have several unusual biochemical characters that are exploited in methods for their laboratory identification. They belong to the minority of *E. coli* that are β-glucuronidase negative and do not ferment sorbitol or rhamnose within 24 hours. These can be isolated from faecal specimens by plating on media containing D-sorbitol instead of lactose.

HiCrome™ EC 0157:H7 Selective Agar Base, Modified is based on the formulation described by Rappaport and Henigh (1). The medium contains sorbitol and a proprietary chromogenic mixture instead of lactose and indicator dyes respectively. The chromogenic substrate is specifically and selectively cleaved by *Escherichia coli* O157:H7 resulting in a dark purple to magenta coloured moiety. *E. coli* forms bluish green coloured colonies.


Tryptone and yeast extract provides carbonaceous, nitrogenous and growth nutrients. Addition of HiCrome™ EC 0157:H7 Selective Supplement (FD187) makes the medium selective (3). Potassium tellurite selectively inhibits *Aeromonas* and *Providencia* species. Novobiocin inhibits gram-positive bacteria. Sodium lauryl sulphate helps to inhibit the accompanying gram-positive flora.

Type of specimen

Food and dairy samples.

Specimen Collection and Handling

For food and dairy samples, follow appropriate techniques for sample collection and processing as per guidelines (4, 5, 6).



M1575A – HiCrome™ EC 0157:H7 Selective Agar Base, Modified

HiCromeVeg™ Freedom from BSE / TSE worries
Single Streak Rapid Differentiation Series

HiCrome™ EC 0157:H7 Selective Agar Base (M1575A) is also available as HiCrome™ EC 0157:H7 Selective HiVeg Agar Base (MV1575A) wherein all the animal origin nutrients have been replaced by vegetable based nutrients.

HiCrome™ EC 0157:H7 Selective Agar Base, Modified

Recommended for selective isolation and easy detection of *Escherichia coli* O157:H7 from food samples.

M1575A

Warning and Precautions

Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets

Limitations

1. Some species may show poor growth due to varying nutritional requirements.
2. Further biochemical test must be carried out for confirmation.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

- Appearance of Powder** : Cream to yellow coloured, homogeneous, free flowing powder.
- Gelling** : Firm, comparable with 1.5% Agar gel.
- Colour and Clarity of prepared medium** : Light amber coloured, clear to slightly opalescent gel forms in Petri plates.
- Reaction** : Reaction of 3.18% w/v aqueous solution at 25°C. pH : 6.8 ± 0.2.
- Cultural Response** : Cultural characteristics observed with added HiCrome™ EC 0157:H7 Selective Supplement (FD187). after an incubation at 35-37°C for 18-24 hours

Organisms (ATCC)	Inoculum (CFU)	Growth	Recovery	Colour of colony
<i>Escherichia coli</i> O157:H7 (NCTC 12900)	50-100	luxuriant	≥50%	dark purple-magenta
<i>Escherichia coli</i> (25922) (00013*)	50-100	none to poor	≤10%	bluish green
<i>Pseudomonas aeruginosa</i> (27853) (00025*)	50-100	fair to good	30-40%	colourless
<i>Klebsiella pneumoniae</i> (13883) (00097*)	>10 ³	fair to good	30-40%	colourless - mauve (mucoid)
<i>Bacillus subtilis</i> subsp. spizizenii (6633) (00003*)	≥10 ³	inhibited	0%	-
<i>Staphylococcus aureus</i> subsp. aureus (25923) (00034*)	≥10 ³	inhibited	0%	-

Key : * : corresponding WDCM Numbers

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (7, 8).

References

1. Rappaport F. and Henigh E., 1952, J. Clin. Pathol., 5:361.
2. Smith and Scotland, 1988, J. Med. Microbiol., 26:77-85.
3. Zadik P. M., Chapman P. A. and Siddons C. A., 1993, J. Med. Microbiol., 39, 155-158.
4. American Public Health Association, Standard Methods for the Examination of Dairy Products, 1978, 14th Ed., Washington D.C.
5. Salfinger Y., and Tortorello M.L. Fifth (Ed.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.
6. Wehr H. M. and Frank J. H., 2004, Standard Methods for the Microbiological Examination of Dairy Products, 17th Ed APHA Inc., Washington, D.C.
7. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
8. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

HiCrome™ Enrichment Broth Base for EC 0157:H7

Recommended for isolation and selective differentiation of *Escherichia coli* O157:H7 from food and environmental samples by chromogenic method.

M1598

Composition **

Ingredients	Grams/Litre
Tryptone	10.00
Sorbitol	10.00
Bile salts mixture	1.50
Chromogenic mixture	1.30

Final pH (at 25°C) 7.1 ± 0.2

** Formula adjusted, standardized to suit performance parameters

Directions

Suspend 11.4 grams in 500 ml distilled water. Heat if necessary to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. For selective isolation of *E. coli* O157:H7, aseptically add the rehydrated contents of 1vial of HiCrome™ ECO157:H7 Selective Supplement I (FD230). Mix well and dispense into sterile test tubes or flasks as desired.

Principle and Interpretation

March and Ratnam (1) reported the inability of *Escherichia coli* O157:H7 to ferment sorbitol while developing Sorbitol MacConkey medium. Subsequently Thomson et al (2) observed the absence of β -glucuronidase activity in *E. coli* O157:H7 from a variety of samples by direct culture.

The medium contains Tryptone that provides nitrogenous, carbonaceous compounds and other essential growth nutrients. Sorbitol is the fermentable carbohydrate, bile salt mixture inhibits most of the gram-positive organisms. Addition of tellurite (FD230) makes the medium more specific and selective. The bluish colour development of *E. coli* and *Klebsiella* in the medium is due to the enzymes β -D-galactosidase and β -D-glucuronidase respectively that cleaves the chromogenic substrates present in chromogenic mixture. However *E. coli* O157:H7 gives a purple colour to the medium due to the absence of β -glucuronidase and its inability to ferment sorbitol.

Type of specimen

Food and Environmental samples

Specimen Collection and Handling

For food and dairy samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (3, 4, 5). After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets

Limitations

1. Certain species of *Shigella* and *Salmonella* are β -glucuronidase positive which may appear as light blue.
2. Further biochemical test must be carried out for confirmation.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.


Quality Control

Appearance of Powder : Cream to yellow coloured, homogeneous, free flowing powder.

Colour and Clarity of prepared medium : Light yellow coloured, clear solution without any precipitate.

Reaction : Reaction of 2.28% w/v aqueous solution at 25°C. pH : 7.1 ± 0.2.

Cultural Response : Cultural characteristics observed with added HiCrome™ EC O157:H7 Selective Supplement I (FD230) after an incubation at 35-37°C for 18-24 hours.



M1598 HiCrome™ Enrichment Broth Base for EC 0157:H7

1. Control
2. *E. coli* O157:H7 (NCTC 12900)
3. *Escherichia coli* (ATCC 25922 (00013*))
4. *Cronobacter sakazakii* (ATCC 12868)
5. *Klebsiella pneumoniae* (ATCC 13883)

For identification of *E. coli* O157:H7

HiCrome™ Enrichment Broth Base for EC 0157:H7

Recommended for isolation and selective differentiation of *Escherichia coli* O157:H7 from food and environmental samples by chromogenic method.

M1598

Organisms (ATCC)	Inoculum (CFU)	Growth	Colour of Medium	Growth«	Colour of Medium«
<i>Escherichia coli</i> O157:H7 (NCTC 12900)	50-100	good-luxuriant	purple#	good-luxuriant	purple#
<i>Escherichia coli</i> (25922) (00013*)	50-100	good-luxuriant	blue#	inhibited	-
** <i>Cronobacter sakazakii</i> (12868)	50-100	good-luxuriant	white#	none-poor	colourless#
<i>Klebsiella pneumoniae</i> (13883) (00097*)	50-100	good-luxuriant	bluish green	good	bluish green#
<i>Salmonella</i> Enteritidis (13076) (00030*)	50-100	good-luxuriant	colourless#	good	colourless#
<i>Shigella flexneri</i> (12022) (00126*)	50-100	good-luxuriant	colourless	inhibited	-
<i>Enterococcus faecalis</i> (29212) (00087*)	50-100	good	-	inhibited	-
<i>Staphylococcus aureus</i> subsp <i>aureus</i> (25923) (00034*)	≥10 ³	inhibited	-	inhibited	-

KEY : « : after addition of HiCrome™ ECO157:H7 Selective Supplement I (FD230)

: may show slight precipitation of growth

** : Formerly known as *Enterobacter sakazakii*

* = corresponding WDCM Numbers

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (6, 7).

References

1. March S. B. and Ratnam S., (1986), J. Clin. Microbiol. 23, 869 - 872.
2. Thompson et al. (1990), J. Clin. Microbiol. 29, 2165 - 2168.
3. American Public Health Association, Standard Methods for the Examination of Dairy Products, 1978, 14th Ed., Washington D.C.
4. Salfinger Y., and Tortorello M.L. Fifth (Ed.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.
5. Wehr H. M. and Frank J. H., 2004, Standard Methods for the Microbiological Examination of Dairy Products, 17th Ed., APHA Inc., Washington, D.C.
6. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
7. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S. and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

HiCrome™ M-Modified ECO157:H7 Selective Agar Base

Recommended for presumptive enumeration of *Escherichia coli* O157:H7 by membrane filtration technique.

M1862

Composition **	
Ingredients	Grams/Litre
Peptone	5.000
Yeast extract	3.000
Sodium chloride	5.000
Lysine	10.000
Sorbitol	20.000
Dextrose (Glucose)	2.500
Magnesium sulphate	1.500
Sodium deoxycholate	0.150
Sodium glucuronate	0.500
Phenol red	0.120
Chromogenic mixture	0.050
Agar	15.000

Final pH (at 25°C) 7.2 ± 0.2

** Formula adjusted, standardized to suit performance parameters

Directions

Suspend 62.82 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. DO NOT AUTOCLAVE. Cool to 45-50°C and aseptically add rehydrated contents of one vial of HiCrome™ ECO157: H7 Selective Supplement, Modified (FD295). Mix well and pour in to sterile Petri plates.

Principle and Interpretation

Escherichia coli O157:H7 belongs to the Enterohemorrhagic *Escherichia coli* (EHEC) group and it predominates as a food borne pathogen. *E. coli* O157:H7 was first recognized as a human pathogen in 1982 when two outbreaks of hemorrhagic colitis were associated with consumption of undercooked ground beef that has been contaminated with this organism (3) that results from the action of a shiga-like toxin (SLT) (1, 7). This medium is recommended for isolation of enteropathogenic *Escherichia coli* O157:H7 in meats, poultry, dairy foods, infant formula, liquid eggs, mayonnaise and apple cider (4, 5). The medium is based on three differential biochemical reactions - lysine decarboxylase (positive for typical EHEC O157 strains), sorbitol fermentation and beta-glucuronidase (2). This medium is also used for the enumeration of β-glucuronidase-positive *E. coli* from foods (6).

Peptone and yeast extract provides nitrogenous and carbonaceous compounds, long chain amino acids and other essential growth nutrients. Sodium chloride maintains the osmotic environment of the medium. Bacteria which were able to grow on this medium will ferment dextrose first. Once dextrose has been depleted, sorbitol positive bacteria will begin to ferment sorbitol, producing a drop in pH of medium, which produces yellow colour to the colony due to phenol red which is a pH indicator. Glucuronidase positive *E. coli* will break down

X-Gluc, resulting in the production of an insoluble blue precipitate in the colony. This will combine with the colour of the pH indicator dye to produce a green colony in case of sorbitol positive or lysine negative bacteria. This medium also contains lysine, lysine positive organisms decarboxylates lysine which produces an increase in pH of medium, hence produces pink coloured colonies. Selectivity is achieved through the use of monensin (FD295) which inhibits gram positive bacteria and incubation at 44 - 44.5°C inhibits gram negative bacteria. Most of the other organisms are unable to grow and if any develop yellow colonies.

Type of specimen

Food samples ; Dairy samples


Specimen Collection and Handling

For food samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (3, 8, 9).

After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

Read the label before opening the container. Wear protective gloves/ protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets



M1862 – HiCrome™ M-Modified ECO157:H7 Selective Agar Base

HiCrome™ M-Modified ECO157:H7 Selective Agar Base

Recommended for presumptive enumeration of *Escherichia coli* O157:H7 by membrane filtration technique.

M1862

Limitations

1. β -glucuronidase is present in 97% of *E. coli* strains, however few *E. coli* may be negative.
2. Certain species of *Shigella* and *Salmonella* are β -glucuronidase positive which may appear as light blue.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

- Appearance of Powder** : Light yellow to pink homogeneous free flowing powder
- Gelling** : Firm, comparable with 1.5 % Agar gel.
- Colour and Clarity of prepared medium** : Red coloured, clear to slightly opalescent gel forms in Petri plates
- Reaction** : Reaction of 6.28% w/v aqueous solution at 25°C. pH : 7.2 \pm 0.2
- Cultural Response** : Cultural characteristics observed with added HiCrome™ ECO157:H7 Selective Supplement, Modified (FD295), after an incubation at 44 - 44.5°C for 18-24 hours.

Organism (ATCC)	Inoculum (CFU)	Growth	Recovery	Colour of Colony (On membrane filter)
<i>Escherichia coli</i> (25922) (00013*)	50-100	luxuriant	>50%	green
<i>Escherichia coli</i> O157:H7 (NCTC 12900)	50-100	luxuriant	\geq 50%	pink
<i>Klebsiella pneumoniae</i> (13883) (00097*)	50-100	fair	20-30%	yellow
<i>Staphylococcus aureus</i> subsp aureus (25923) (00034*)	$\geq 10^3$	inhibited	0%	-
<i>Enterococcus faecalis</i> (29212) (00087*)	$\geq 10^3$	inhibited	0%	-

Key : * = corresponding WDCM Numbers

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (10, 11).

References

1. Centre for Diseases Control, 1991, Morbid. Mortal, Weekly Rep 40:265.
2. Corry J.E.L, Curtis G.D.W., Baird R.M., Culture Media for Food Microbiology, Progress in Industrial Microbiology, Volume 37.
3. Salfinger Y, and Tortorello M.L. Fifth (Ed.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 4th Ed., APHA, Washington, D.C.
4. Entis, P., and I. Lerner. 1997. 24-hour presumptive enumeration of *Escherichia coli* O157:H7 in food using the ISO-GRID method with SD-39 agar. J. Food Prot. 60:883-890.
5. Entis, P. 1998. Direct 24-hour presumptive enumeration of *Escherichia coli* O157:H7 in food using the hydrophobic grid membrane filter, followed by serological confirmation : collaborative study. J. AOAC Int. 81:403-418.
6. Entis, P., and I. Lerner. 1998. Enumeration of β -glucuronidase positive *E. coli* in foods by using the ISO-GRID method with SD-39 agar. J. Food Prot. 61:913-916.
7. March S. B. and Ratnam S., 1986, J. Clin. Microbiol., 23:869.
8. American Public Health Association, Standard Methods for the Examination of Dairy Products, 1978, 14th Ed., Washington D.C.
9. Wehr H. M. and Frank J. H., 2004, Standard Methods for the Microbiological Examination of Dairy Products, 17th Ed., APHA Inc., Washington, D.C.
10. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
11. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

HiCrome™ *Enterobacter sakazakii* Agar / Modified

Recommended for the isolation and identification of *Enterobacter sakazakii* from food, milk and dairy products (*Enterobacter sakazakii* now referred as *Cronobacter sakazakii*)

M1577/
M1641

Composition **	M1577	M1641
	Grams/Litre	Grams/Litre
Ingredients		
Tryptone	15.00	7.00
Soya peptone	5.00	-
Yeast extract	-	3.00
Sodium chloride	5.00	5.00
Sodium deoxycholate	0.50	0.60
Sodium thiosulphate	1.00	-
Chromogenic mixture	10.17	-
Chromogenic substrate	-	0.15
Crystal violet	-	0.002
Agar	15.00	15.00
Final pH (at 25°C)	7.3 ± 0.2	7.0 ± 0.2

** Formula adjusted, standardized to suit performance parameters

Directions

Suspend 51.67 grams of M1577 or 30.75 grams of M1641 in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Mix well and pour into sterile Petri plates.

Principle and Interpretation

Enterobacter species are widely distributed in nature occurring in fresh water, soil, sewage, plants, vegetables, animal and human faeces.**Cronobacter sakazakii* has been closely associated with neonatal meningitis and sepsis (3). The chromogenic substrate in HiCrome™ *Enterobacter sakazakii* Agar (M1577) is cleaved specifically (2) by the glucosidase enzyme possessed by *Enterobacter* species resulting in formation of blue-green colonies. Other organisms, which do not cleave this substrate, produce yellow coloured colonies. Incorporation of the chromogenic mixture in the media renders an intense blue colour to **C.sakazakii* colonies and light blue green colour to other *Enterobacter* species. HiCrome™ *Enterobacter sakazakii* Agar, Modified is recommended by ISO Committee for the isolation and identification of **C.sakazakii* (1). The chromogenic substrate is cleaved specifically (2) by **C.sakazakii* resulting in the formation of blue green colonies. Other organisms, which do not cleave this substrate, produce colorless to slightly violet coloured colonies.

Tryptone, soya peptone and yeast extract provide the essential growth nutrients along with nitrogenous and carbonaceous compounds, long chain amino acids and vitamins. Sodium chloride helps in maintaining the osmotic equilibrium of the medium. Sodium deoxycholate and crystal violet (in M1641) inhibits the accompanying gram-positive flora.

*: Formerly known as *Enterobacter sakazakii*

HiCromeVeg™ Freedom from BSE / TSE worries
Single Streak Rapid Differentiation Series

HiCrome™ Ent. sakazakii Agar (M1577) is also available as HiCrome™ Ent. sakazakii HiVeg™ Agar (MV1577) wherein all the animal origin nutrients have been replaced by vegetable based nutrients.

Type of specimen

Food, milk & dairy products

Specimen Collection and Handling


For food and dairy samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (4, 5, 6). After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions


Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets.

Limitations

1. Slight variation in colour may be observed depending on enzyme production by organism and substrate utilization from the medium
2. Some species may show poor growth due to nutritional variations.
3. Further biochemical tests must be carried out for confirmation.



M1577 HiCrome™ *Enterobacter sakazakii* Agar,



M1641 HiCrome™ *Enterobacter sakazakii* Agar, Modified

For identification of *Cronobacter (Enterobacter) sakazakii*

HiCrome™ *Enterobacter sakazakii* Agar / Modified

Recommended for the isolation and identification of *Enterobacter sakazakii* from food, milk and dairy products (*Enterobacter sakazakii* now referred as *Cronobacter sakazakii*)

M1577/
M1641

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

- Appearance of Powder** : Light yellow to pink coloured, homogeneous, free flowing powder.
- Gelling** : Firm, comparable with 1.5% Agar gel.
- Colour and Clarity of prepared medium** : Purple coloured (M1577) or light purple coloured (M1641), clear to slightly opalescent gel forms in Petri plates.
- Reaction** : Reaction of 5.16% w/v aqueous solution of M1577 at 25°C. pH : 7.3 ± 0.2 .
Reaction of 3.07% w/v aqueous solution of M1641 at 25°C. pH : 7.0 ± 0.2.
- Cultural Response** : M1577 Cultural characteristics observed after an incubation at 35-37°C for 18-24 hours.

Organisms (ATCC)	Inoculum (CFU)	Growth	Recovery	Colour of colony
<i>Escherichia coli</i> (25922) (00013*)	50-100	good-luxuriant	≥50%	yellow
▲ <i>Klebsiella aerogenes</i> (13048) (00175*)	50-100	good-luxuriant	≥50%	green
# <i>Cronobacter sakazakii</i> (12868)	50-100	good-luxuriant	≥50%	blue
<i>Klebsiella pneumoniae</i> (13883) (00097*)	50-100	good-luxuriant	≥50%	green (mucoid)
<i>Staphylococcus aureus</i> sub-sp <i>aureus</i> (25923) (00034*)	≥10 ³	inhibited	0%	—
<i>Enterococcus faecalis</i> (29212) (00087*)	≥10 ³	inhibited	0%	—

(M1641) : Cultural characteristics observed after an incubation at 44±1°C for 18-24 hours.

Organisms (ATCC)	Inoculum (CFU)	Growth	Recovery	Colour of colony
<i>Escherichia coli</i> (25922) (00013*)	50-100	good-luxuriant	≥50%	colourless with blue centre

* <i>Cronobacter sakazakii</i> (12868)	50-100	good-luxuriant	≥50%	blue - green
<i>Staphylococcus aureus</i> sub-sp <i>aureus</i> (25923) (00034*)	≥10 ³	inhibited	0%	—
<i>Enterococcus faecalis</i> (29212) (00087*)	≥10 ³	inhibited	0%	—

Key : * = corresponding WDCM Numbers
▲ : Formerly known as *Enterobacter aerogenes*
: Formerly known as *Enterobacter sakazakii*

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (7, 8).

References

1. International Organization for Standardization Draft ISO/ TS 22964, 2006 (E).
2. Isenberg, (Ed.), 1992, Clinical Microbiology Procedures Handbook, Vol. 1, American Society for Microbiology, Washington, D. C.
3. Muytjens H. L., Zanen H. C., Sonderkamp H. J. et al, J. Clin Microbiol 18:115-120,1983.
4. American Public Health Association, Standard Methods for the Examination of Dairy Products, 1978, 14th Ed., Washington D.C.
5. Salfinger Y., and Tortorello M.L. Fifth (Ed.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.
6. Wehr H. M. and Frank J. H., 2004, Standard Methods for the Microbiological Examination of Dairy Products, 17th Ed., APHA Inc., Washington, D.C.
7. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
8. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

Ready Prepared Media

Code	Product Name	Usage	Packing
Category : HiTouch™ FlexiPlate™			
FL036	HiTouch™- HiCrome Ent. Sakazakii Agar Flexi Plate™	for enumeration (count) and differentiation of <i>Enterobacter sakazakii</i> (<i>Cronobacter sakazakii</i>).	50 plts

For the isolation and identification of *Cronobacter sakazakii*

HiCrome™ Cronobacter Isolation Agar (CCI Agar)

Recommended for the isolation and identification of *Cronobacter sakazakii* from food products. The composition and performance of this media are as per specifications laid down in in ISO /TS 22964: 2017

M2062I

Composition **

Ingredients	Grams/Litre
Tryptone#	7.00
Yeast extract	3.00
Sodium chloride	5.00
Sodium deoxycholate	0.25
5-Bromo-4-chloro-3-indolyl α -D-glucopyranoside	1.50
Ammonium iron(III) citrate	1.00
Sodium thiosulfate	1.00
Agar	15.00

Final pH (at 25°C) 7.3 \pm 0.2

** Formula adjusted, standardized to suit performance parameters

- Equivalent to Tryptic digest of casein

Directions

Suspend 32.4 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Mix well and pour into sterile Petri plates.

Principle and Interpretation

Enterobacter species are widely distributed in nature occurring in fresh water, soil, sewage, plants, vegetables, animal and human faeces. *Cronobacter sakazakii* has been closely associated with neonatal meningitis and sepsis (1). HiCrome™ Cronobacter isolation Agar is recommended by ISO Committee for the isolation and identification of *C.sakazakii* from food samples (2).

The chromogenic substrate (5-Bromo-4-chloro-3-indolyl α -D-glucopyranoside) is cleaved specifically (3) by *C.sakazakii* resulting in the formation of blue green colonies. Other organisms, which do not cleave this substrate, produce colourless colonies.

Tryptone and yeast extract provides nitrogenous and carbonaceous compounds, long chain amino acids, vitamins and other essential growth nutrients. Sodium chloride helps in maintaining the osmotic equilibrium of the medium. Sodium deoxycholate inhibits the accompanying gram-positive flora.

Key: #: Formerly known as *Enterobacter sakazakii*.

Type of specimen

Food samples

Specimen Collection and Handling

For food and dairy samples, follow appropriate techniques for sample collection and processing as per guidelines (2, 3).


After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets

Limitations

1. Slight variation in colour may be observed depending on enzyme production by organism and substrate utilization from the medium
2. Some species may show poor growth due to nutritional variations.
3. Further biochemical tests must be carried out for confirmation.



M2062I HiCrome™ Cronobacter Isolation Agar (CCI Agar)

HiCrome™ Cronobacter Isolation Agar (CCI Agar)

Recommended for the isolation and identification of *Cronobacter sakazakii* from food products. The composition and performance of this media are as per specifications laid down in in ISO /TS 22964: 2017

M2062I

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

- Appearance of Powder** : Cream to yellow to pink homogeneous free flowing powder.
- Gelling** : Firm, comparable with 1.5% Agar gel
- Colour and Clarity** : Yellow coloured, clear to slightly opalescent gel forms in Petri plates
- Reaction** : Reaction of 3.24% w/v aqueous solution at 25°C. pH : 7.3±0.2
- Cultural Response** : Cultural characteristics observed after an incubation at 41.5±1°C for 24±2 hours.

Organism (ATCC)	Inoculum (CFU)	Growth	Recovery	Colour of colony
# <i>Cronobacter sakazakii</i> (29544) (00214*)	50-100	good-luxuriant	>=50%	blue-green
<i>Cronobacter muytjensii</i> (51329) (00213*)	50-100	good-luxuriant	>=50%	blue-green
<i>Enterobacter cloacae</i> (13047) (00083*)	50-100	good-luxuriant	>=50%	colourless without green or blue green colour
<i>Staphylococcus aureus subsp aureus</i> (25923) (00034*)	>=10 ³	inhibited	0%	
<i>Staphylococcus aureus subsp aureus</i> (6538) (00032*)	>=10 ³	inhibited	0%	

Key: * Corresponding WDCM numbers
#: Formerly known as *Enterobacter sakazakii*

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (4, 5).

References

- Muytjens H. L., Zanen H. C., Sonderkamp H. J. et al, J. Clin Microbiol 18:115-120, 1983.
- International Organization for Standardization. Microbiology of the food chain- Horizontal method for the detection of *Cronobacter* spp. Draft ISO/ TS 22964, 20176 (E).
- Salinger Y., and Tortorello M.L. Fifth (Ed.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.
- Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
- Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

For identification and differentiation of *Salmonella* species

Salmonella Differential Agar / Modified (Twin pack) (RajHans Medium)

Recommended for identification and differentiation of *Salmonella* species from members of *Enterobacteriaceae*, especially *Proteus* species.

M1078/
M1082

Composition **	M1078	M1082
Ingredients	Grams/Litre	Grams/Litre
Part A :		
Peptone, special	8.00	8.00
Yeast extract	2.00	3.00
Sodium deoxycholate	1.00	1.00
Sodium chloride	—	5.00
B.C. indicator	2.00	2.00
Agar	12.00	12.00
Part B :		
Propylene glycol	10.00	10.00

Final pH (at 25°C) 7.3 ± 0.2

** Formula adjusted, standardized to suit performance parameters

Directions

Suspend 10 grams of fluid Part B in 1000 ml distilled water. Add 25 grams of Part A (M1078) or 31 grams of Part A (M1082). Mix well and heat to boiling to dissolve the medium completely. DO NOT AUTOCLAVE. Cool to 45- 50°C. Mix well before and pour into sterile Petri plates.

Principle and Interpretation

Salmonella Differential Agar media are slight modification of original formulation of Rambach (3) used for differentiation of *Salmonella* species from *Proteus* species and other enteric bacteria. Production of acid from propylene glycol is a novel characteristic of *Salmonella* species and is utilized in these media. Many of the media such as SS Agar, XLD Agar recommended for the identification and differentiation of *Salmonella* species (1) are based on lactose fermentation and hydrogen sulphide production.

Peptone special and yeast extract provides nitrogenous and carbonaceous compounds, long chain amino acid, vitamins and other essential growth nutrients. Sodium deoxycholate inhibits gram-positive organisms rendering the medium selective for enteric microorganisms. The BC indicator turns pink in presence of acid produced from propylene glycol. Lactose fermenting ability is determined by using an indicator, which can detect the presence of enzyme β -galactosidase. Lactose fermenting (β -galactosidase producing) bacteria yield blue violet coloured colony (2). Salmonellae produce acid from propylene glycol and on combining with the BC indicator gives typical pink

red colonies. Other enteric gram-negative bacteria form colourless colonies. *Salmonella* Typhimurium and *Salmonella* Enteritidis produce pink to red colonies. Specimen should be enriched in an appropriate selective enrichment broth. This enriched culture is then inoculated on Salmonella Differential Agar/ Salmonella Differential Agar, Modified and incubated at 35-37°C for 24-48 hours.

Type of specimen

Clinical : faeces, urine; Water samples and Food samples

Specimen Collection and Handling

For clinical samples follow appropriate techniques for handling specimens as per established guidelines (4, 5).


For food samples, follow appropriate techniques for sample collection and processing as per guidelines (6).

For water samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (7).

After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

In Vitro diagnostic use only. Read the label before opening the container. Wear protective gloves/protective clothing/ eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets.



M1082 – *Salmonella* Differential Agar Modified (Twin pack) (RajHans Medium)

HiCromeVeg™ Freedom from BSE / TSE worries
Single Streak Rapid Differentiation Series

Salmonella Differential Agar / Modified (Twin pack) (RajHans Medium)(M1078/M1082) is also available as *Salmonella* Differential HiVeg™ Agar / Modified (Twin pack) (RajHans Medium) (MV1078/MV1082) wherein all the animal origin nutrients have been replaced by vegetable based nutrients.

For identification and differentiation of *Salmonella* species

Salmonella Differential Agar / Modified (Twin pack) (RajHans Medium)

Recommended for identification and differentiation of *Salmonella* species from members of *Enterobacteriaceae*, especially *Proteus* species.

M1078/
M1082

Limitations

1. The medium is selective for *Salmonella* and may not support the growth of other microorganisms.
2. Most of the *Salmonella* strains shows pink-red colonies except few which may show colourless colonies.
3. Due to nutritional variations, some strains may show poor growth.
4. Final confirmation of suspected colonies must be carried out by serological and biochemical tests.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

Appearance of Powder : Part A : Light yellow to light pink coloured, homogeneous, free flowing powder.
Part B : Colourless, viscous, solution.

Gelling : Firm, comparable with 1.2% Agar gel.

Colour and Clarity of prepared medium : Light orange coloured, clear to slightly opalescent gel forms in Petri plates.

Reaction : Reaction of 2.5% w/v Part A of M1078 or 3.1% w/v of Part A of M1082 aqueous solution at 25°C. pH : 7.3 ± 0.2.

Cultural Response : Cultural characteristics observed after an incubation at 35-37°C for 24 - 48 hours.

Organisms (ATCC)	Inoculum (CFU)	Growth	Recovery	Colour of colony
<i>Salmonella</i> Enteritidis (13076) (00030*)	50-100	luxuriant	≥50%	pink-red
<i>Salmonella</i> Typhimurium (14028) (00031*)	50-100	luxuriant	≥50%	pink-red
<i>Salmonella</i> Typhi (6539)	50-100	luxuriant	≥50%	colourless
<i>Escherichia coli</i> (25922) (00013*)	50-100	luxuriant	≥50%	blue-green
<i>Klebsiella pneumoniae</i> (13883) (00097*)	50-100	luxuriant	≥50%	blue-violet
<i>Proteus mirabilis</i> (25933)	50-100	luxuriant	≥50%	colourless
<i>Shigella flexneri</i> (12022) (00126*)	50-100	luxuriant	≥50%	colourless
<i>Staphylococcus aureus</i> subsp <i>aureus</i> (25923) (00034*)	≥10 ³	inhibited	0%	-

Key : * = corresponding WDCM Numbers

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.


Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (4, 5).

References

1. Rambach A., 1990, Appl Environ. Microbiol., 56:301.
2. Eaton A. D., Clesceri L. S., Rice E. W. and Greenberg A W., (Eds.), 2005, Standard Methods for the Examination of Water and Wastewater, 21st Ed., APHA, Washington, D.C.
3. Greenwald R., Henderson R.W. and Yappaw S., 1991, J. Clin. Microbiol. 29:2354.
4. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
5. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S. and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.
6. Salfinger Y., and Tortorello M.L. Fifth (Ed.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.
7. Baird R.B., Eaton A.D., and Rice E.W., (Eds.), 2015, Standard Methods for the Examination of Water and Wastewater, 23rd ed., APHA, Washington, D.C.



M1078 – *Salmonella* Differential Agar (Twin pack) (RajHans Medium)

For identification and differentiation of *Salmonella* species

HiCrome™ RajHans Medium/Modified (Salmonella Agar/Modified)

Recommended For identification and differentiation of *Salmonella* species from among the members of *Enterobacteriaceae*, especially *Proteus* species.

M1633/
M1634

Composition **

	M1633	M1634
Ingredients	Grams/Litre	Grams/Litre
Tryptone	8.00	8.00
Yeast extract	5.00	5.00
Peptone	4.00	4.00
Sodium chloride	5.00	5.00
Sodium deoxycholate	1.00	1.00
Neutral red	0.02	0.02
Lactose	3.00	3.00
Chromogenic mixture	7.30	4.32
Agar	13.50	12.00

Final pH (at 25°C) 7.3 ± 0.2

** Formula adjusted, standardized to suit performance parameters

Directions

Suspend 46.82 grams of M1633 and 42.34 grams of M1634 in 1000 ml distilled water. Mix well and heat to boiling to dissolve the medium completely. DO NOT AUTOCLAVE. Cool to 45 - 50°C. Mix well and pour into sterile Petri plates.

Principle and Interpretation

HiCrome™ RajHans Medium/Modified is a modification of the original formulation of Rambach (2), used for differentiation of *Salmonella* species from *Proteus* species and other enteric bacteria. The original formulation is based on the novel characteristic of *Salmonella* species to produce acid from propylene glycol, which is detected by indicators present in the medium. These media are unique, because it is not based on acid production by propylene glycol. These media like many other media such as SS Agar, XLD Agar, recommended for the identification and differentiation of *Salmonella* species are based on lactose fermentation (1).

Tryptone, peptone and yeast extract supports the luxuriant growth of bacteria by providing carbonaceous and nitrogenous compounds, long chain amino acids, vitamin B complex and other essential nutrients. Sodium deoxycholate inhibits gram-positive organisms rendering the medium selective for enteric microorganisms. The chromogenic mixture incorporated in the medium yields pink to red colonies of *Salmonella*. Lactose fermenting organisms form light purple to blue violet colonies. Other enteric gram-negative bacteria form colourless colonies.

Type of specimen

Clinical: faeces, urine; Water samples and Food samples

Specimen Collection and Handling

For clinical samples follow appropriate techniques for handling specimens as per established guidelines (3, 4).

For food samples, follow appropriate techniques for sample collection and processing as per guidelines (5).

For water samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (1).


After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

In Vitro diagnostic use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets

Limitations

1. The medium is selective for *Salmonella* may not support the growth of other microorganisms.
2. Most of the *Salmonella* strains shows pink-red colonies except few which may show colourless colonies.
3. Due to nutritional variations, some strains may show poor growth.
4. Final confirmation of suspected colonies must be carried out by serological and biochemical tests.



M1633 – HiCrome™ RajHans Medium (*Salmonella* Agar)

For identification and differentiation of *Salmonella* species

HiCrome™ RajHans Medium/Modified (Salmonella Agar/Modified)

Recommended For identification and differentiation of *Salmonella* species from among the members of *Enterobacteriaceae*, especially *Proteus* species.

M1633/
M1634

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

Appearance of Powder : Light yellow to beige coloured, homogeneous, free flowing powder.

Gelling : Firm, comparable with 1.35% Agar gel of M1633 and 1.2% Agar gel of M1634.

Colour and Clarity of prepared medium : Light orange coloured, clear to slightly opalescent gel forms in Petri plates.

Reaction : Reaction of 4.68% w/v of M1633 and 4.23% w/v of M1634 aqueous solution at 25°C. pH:7.3 ± 0.2.

Cultural Response : Cultural characteristics observed after an incubation at 35-37°C for 24 - 48 hours.

Organisms (ATCC)	Inoculum (CFU)	Growth	Recovery	Colour of colony
<i>Escherichia coli</i> (25922) (00013*)	50-100	luxuriant	≥50%	light purple
<i>Klebsiella pneumoniae</i> (13883) (00097*)	50-100	luxuriant	≥50%	blue-violet
<i>Proteus mirabilis</i> (25933)	50-100	luxuriant	≥50%	colourless
<i>Salmonella</i> Typhi (6539)	50-100	luxuriant	≥50%	colourless
<i>Salmonella</i> Typhimurium (14028) (00031*)	50-100	luxuriant	≥50%	pink-red
<i>Salmonella</i> Enteritidis (13076) (00030*)	50-100	luxuriant	≥50%	pink-red
<i>Shigella flexneri</i> (12022) (00126*)	50-100	luxuriant	≥50%	colourless

<i>Staphylococcus aureus</i> subsp aureus (25923) (00034*)	≥10 ³	inhibited	0%	-
--	------------------	-----------	----	---

Key : * = corresponding WDCM Numbers

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (3, 4).

References

- Baird R.B., Eaton A.D., and Rice E.W., (Eds.), 2015, Standard Methods for the Examination of Water and Wastewater, 23rd ed., APHA, Washington, D.C.
- Rambach A., 1990, Environment. Microbiol, 56:301.
- Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
- Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.
- Salfinger Y., and Tortorello M.L. Fifth (Ed.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.



Klebsiella pneumoniae (13883) (00097*)



Salmonella Enteritidis (13076) (00030*)

For identification and differentiation of *Salmonella* species

HiCrome™ Salmonella Agar /HiCrome™ Improved Salmonella Agar

Recommended for the simultaneous detection of *Salmonella* and *Escherichia coli* from food, water and clinical samples.

M1296/
M1466

Composition **

	M1296	M1466
Ingredients	Grams/Litre	Grams/Litre
Peptone	6.00	—
Peptone special	—	8.00
Yeast extract	2.50	2.00
Bile salts mixture	1.00	—
Sodium deoxycholate	—	1.00
Chromogenic mixture	5.40	3.25
Agar	13.00	12.00

Final pH (at 25°C) 7.7 ± 0.2 7.3 ± 0.2

** Formula adjusted, standardized to suit performance parameters

Directions

Suspend 27.9 grams of M1296 or 26.25 grams of M1466 in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. DO NOT AUTOCLAVE. Cool to 45-50°C. Mix well and pour into sterile Petri plates.

Principle and Interpretation

Salmonella species have been isolated from humans and almost all animals throughout the world. They cause many types of infections from mild, self-limiting gastroenteritis to life threatening typhoid fever. *Salmonella* Typhi and *Salmonella* Paratyphi A & B cause gastroenteritis, bacteremia and enteric fever, *Salmonella* Choleraesuis causes gastroenteritis and enteric fever, especially in children. *Salmonella* Typhimurium is the most frequently isolated serotype of *Salmonella* (2). HiCrome™ *Salmonella* Agar medium is a modification of the original formulation of Rambach (3) and is used for the differentiation of *Salmonella* species from other enteric bacteria. Rambach formulation differentiates *Salmonella* based on propylene glycol utilization and presence of a chromogenic indicator. However, HiCrome™ *Salmonella* Agar medium uses only a chromogenic mixture for identification and differentiation of *Salmonella* species.

Peptone, peptone special and yeast extract provides nitrogenous and carbonaceous compounds, long chain amino acids, vitamins and other essential growth nutrients.

Escherichia coli and *Salmonella* are easily distinguishable due to their colony characteristics. *Salmonella* forms light purple coloured colonies with a purple halo on (M1296) and pink to red colonies on (M1466). *E. coli* exhibits a characteristic blue colour, due to presence of the enzyme β -glucuronidase. Other organisms form colourless colonies. The characteristic light purple and blue colour is due to the chromogenic mixture (1). Bile salts mixture or sodium deoxycholate inhibits gram-positive organisms.

Type of specimen

Clinical: faeces, urine; Water samples and Food samples

Specimen Collection and Handling

For clinical samples follow appropriate techniques for handling specimens as per established guidelines (4,5).

For food samples, follow appropriate techniques for sample collection and processing as per guidelines (6).

For water samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (7).


After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions


In Vitro Diagnostic use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets

Limitations

1. The medium is selective for *Salmonella* may not support the growth of other microorganisms.



M1296-HiCrome™ *Salmonella* Agar



M1466-HiCrome™ Improved *Salmonella* Agar

HiCromeVeg™ Freedom from BSE / TSE worries
Single Streak Rapid Differentiation Series

HiCrome™ *Salmonella* Agar / HiCrome™ Improved *Salmonella* Agar (M1296/M1466) is also available as HiCrome™ *Salmonella* HiVeg™ Agar / HiCrome™ Improved *Salmonella* HiVeg™ Agar (MV1296/MV1466) HiCrome™ Improved *Salmonella* HiCynth™ Agar (MCD1466) wherein all the animal origin nutrients have been replaced by vegetable based nutrients / or chemical defined nutrients.

For identification and differentiation of *Salmonella* species

HiCrome™ Salmonella Agar /HiCrome™ Improved Salmonella Agar

Recommended for the simultaneous detection of *Salmonella* and *Escherichia coli* from food, water and clinical samples.

M1296/
M1466

- Most of the *Salmonella* strains shows purple (M1296) or pink-red (M1466) colonies except few which may show colourless colonies.
- Due to nutritional variations, some strains may show poor growth.
- Final confirmation of suspected colonies must be carried out by serological and biochemical tests.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

- Appearance of powder** : Cream to yellow coloured, homogeneous, free flowing powder.
- Gelling** : Firm, comparable with 1.3% Agar gel of M1296 or 1.2% Agar gel of M1466.
- Colour and Clarity of prepared medium** : Light amber coloured (M1296) or reddish pink coloured (M1466), slightly opalescent gel forms in Petri plates.
- Reaction** : Reaction of 2.79% w/v aqueous solution of M1296 at 25°C. pH:7.7 ± 0.2.
Reaction of 2.62% w/v aqueous solution of M1466 at 25°C. pH : 7.3 ± 0.2.
- Cultural Response** : Cultural characteristics observed after an incubation at 35-37°C for 24-48 hours.

Organisms (ATCC)	Inoculum (CFU)	Growth	Recovery	Colour of colony (M1296)	Colour of colony (M1466)
<i>Escherichia coli</i> (25922) (00013*)	50-100	luxuriant	>50%	blue	blue to purple
<i>Salmonella</i> Enteritidis (13076) (00030*)	50-100	luxuriant	>50%	light purple with halo	pink to red
<i>Salmonella</i> Typhi (6539)	50-100	luxuriant	>50%	light purple with halo	light pink
<i>Salmonella</i> Typhimurium (14028) (00031*)	50-100	luxuriant	>50%	light purple with halo	pink to red
<i>Proteus vulgaris</i> (13315)	50-100	luxuriant	40-50%	colourless	light brown

<i>Staphylococcus aureus</i> subsp <i>aureus</i> (25923) (00034*)	>10 ³	inhibited	0%	-	-
<i>Bacillus subtilis</i> sub <i>spizizenii</i> (6633) (00003*)	>10 ³	inhibited	0%	-	-

Key : * = corresponding WDCM Numbers

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (4, 5).

References

- Greenwald R., Henderson R. W. and Yappan S., 1991, J. Clin. Microbiol., 29:2354.
- Murray P. R., Baron J. H., Pfaller M. A., Jorgensen J. H. and Tenover F. C., (Eds.), 2003, Manual of Clinical Microbiology, 8th Ed., American Society for Microbiology, Washington, D.C.
- Rambach A., 1990, Appl. Environ. Microbiol., 56:301.
- Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
- Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.
- Salinger Y., and Tortorello M.L. Fifth (Ed.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.
- Baird R.B., Eaton A.D., and Rice E.W., (Eds.), 2015, Standard Methods for the Examination of Water and Wastewater, 23rd ed., APHA, Washington, D.C.

Ready Prepared Media

Code	Product Name	Usage	Packing
Category : HiDip™ Slides			
HD036	HiDip™ Hicrome™ ECC Agar- Hicrome™ Salmonella Agar	for chromogenic screening of <i>E.coli</i> , coliforms and <i>Salmonella</i> on surfaces or food or water	5 tubes 10 tubes

HiCrome™ Selective Salmonella Agar Base

Recommended for the selective isolation of *Salmonella* species from food samples

M1842

Composition **

Ingredients	Grams/Litre
HI powder#	12.000
Yeast hydrolysate	5.000
Tryptose	5.000
Sodium cholate	3.000
Sodium taurocholate	5.000
Sodium deoxycholate	1.000
Chromogenic mixture	8.000
Agar	15.000

Final pH (at 25°C) 7.3 ± 0.2

** Formula adjusted, standardized to suit performance parameters

#Equivalent to Heart Infusion powder

Directions

Suspend 54.00 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. DO NOT AUTOCLAVE. Cool to 45-50°C. Aseptically add the rehydrated contents of one vial of HiCrome™ Selective *Salmonella* Agar Supplement (FD274). Mix well and pour into sterile Petri plates.

Principle and Interpretation

Salmonella species have been isolated from humans and almost all animals throughout the world. They cause many types of infections from mild, self-limiting gastroenteritis to life threatening typhoid fever. *Salmonella* Typhi and *Salmonella* Paratyphi A & B cause gastroenteritis, bacteremia and enteric fever, *Salmonella* Choleraesuis causes gastroenteritis and enteric fever, especially in children. *Salmonella* Typhimurium is the most frequently isolated serotype of *Salmonella*. *Salmonella* species are the major cause of food poisoning (1) Various chromogenic media are available for the differentiation of *Salmonella* species. The original media formulated by Rambach (2) differentiates *Salmonella* based on propylene glycol utilization *Enterobacter* and presence of a chromogenic indicator. However HiCrome™ Selective Salmonella Agar Base uses chromogenic mixture for identification and differentiation of *Salmonella* species. Sodium cholate, Sodium taurocholate and Sodium deoxycholate in the medium helps to restrict the growth of other organisms. Besides the selective supplement added to the medium inhibits competing microorganisms.

HI powder, yeast hydrolysate and tryptose in the medium provides nitrogenous, carbonaceous compounds, long chain amino acids, vitamins and other essential growth nutrients. Due to the presence of chromogenic mix in the medium *Salmonella* are easily distinguishable and forms purple coloured colonies while some *Enterobacteriaceae* like *Klebsiella* and *Enterobacter* forms blue to dark blue coloured colonies. Conventional method employs the H₂S production property for *Salmonella* detection which is also exhibited by other non *Salmonella*

species such as *Citrobacter*, *Proteus*, etc. Hence further biochemical confirmation is required for further identification.

This medium is specially employed for food samples where the sample is initially enriched in Salmonella Selective Enrichment Broth (M1843) and then isolated on HiCrome™ Selective Salmonella Agar Base. *Salmonella* species give purple coloured colonies due to the enzyme specificity.

Type of specimen

Food samples


Specimen Collection and Handling

For food samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (3).

After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets



M1842 – HiCrome™ Selective *Salmonella* Agar Base

HiCrome™ Selective Salmonella Agar Base

Recommended for the selective isolation of *Salmonella* species from food samples

M1842

Limitations

1. Being highly selective, some strains may show poor growth.
2. Most of the *Salmonella* strains shows purple colonies except few.
3. Final confirmation of suspected colonies must be carried out by serological and biochemical tests.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

- Appearance of Powder** : Light yellow to beige homogeneous free flowing powder
- Gelling** : Firm, comparable with 1.5 % Agar gel.
- Colour and Clarity of prepared medium** : Whitish cream coloured, opalescent gel forms in Petri plates
- Reaction** : Reaction of 5.4% w/v aqueous solution at 25°C. pH : 7.3 ± 0.2.
- Cultural Response** : Cultural characteristics observed with added HiCrome™ Selective Salmonella Agar Supplement (FD274), after an incubation at 35-37°C for 22-24 hours.

Organism (ATCC)	Inoculum (CFU)	Growth	Recovery	Colour of colony
<i>Salmonella</i> Typhimurium (14028) (00031*)	50 -100	good-luxuriant	≥50%	purple
<i>Salmonella</i> Enteritidis (13076) (00030*)	50 -100	good-luxuriant	≥50%	≥50%
<i>Klebsiella pneumoniae</i> (13883) (00097*)	50 -100	good	40-50%	blue
<i>Enterococcus faecalis</i> (29212) (00087*)	>10 ³	inhibited	0 %	-
<i>Staphylococcus aureus</i> subsp aureus (25923) (00034*)	>10 ³	inhibited	0 %	-

Key : * = corresponding WDCM Numbers

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle

in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (4, 5).

References

1. Murray P. R., Baron J. H., Pfaller M. A., Tenover J. C. and Tenover F. C., (Ed.), 2003, Manual of Clinical Microbiology, 8th Ed., American Society for Microbiology, Washington, D.C.
2. Rambach A., 1990, Appl. Environ. Microbiol., 56:301.
3. Salfinger Y., and Tortorello M.L. Fifth (Ed.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.
4. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
5. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S. and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

For identification and differentiation of *Salmonella* species

HiCrome™ MM Agar

Recommended For identification and differentiation of *Salmonella* and non-*Salmonella* like *Citrobacter* from food, water and clinical samples.

M1393

Composition **

Ingredients	Grams/Litre
Peptone	10.00
HM peptone B#	2.00
D-Cellobiose	3.00
Lactose	10.00
D-Mannitol	1.20
D-Trehalose	1.33
Chromogenic mixture	6.60
Agar	15.00

Final pH (at 25°C) 7.6 ± 0.2

** Formula adjusted, standardized to suit performance parameters

#Equivalent to Beef extract

Directions

Suspend 49.13 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Mix well and pour into sterile Petri plates.

Principle and Interpretation

HiCrome™ MM Agar was formulated by Miller and Mallison (1) for specific isolation and detection of *Salmonellae*. This medium is superior to XLT4 Agar in supporting growth of *Salmonella* due to the presence of appropriate proportion of four sugars. Most differential and selective media are formulated with one or more sugars and pH indicators respectively. The utilization of sugars by organisms results in pH-changes. This is used as a means of distinguishing *Salmonella* from competing bacteria on the basis of colony colour. *Salmonella* usually are unable to ferment the sugars (2) that support growth of competing bacteria. Thus other bacteria tend to overgrow *Salmonellae*, masking their presence. The inclusion of sugars like mannitol, cellobiose and trehalose stimulate the better initial growth of *Salmonella* cells. However, the low concentrations of these sugars do not interfere with the utilization of protein and H₂S production. Presence of lactose suppresses H₂S production by non-*Salmonellae* like *Citrobacter freundii*. A chromogenic mixture, present in this medium helps to differentiate between lactose fermenters and nonfermenters. Lactose fermenters

give bluish green coloured colonies, which would have been impossible to differentiate with an indicator based on pH change. Inclusion of tergitol 4 (included in chromogenic mixture) in the medium suppresses the presence of *Proteus* and *Providencia* colonies. Peptone and HM peptone B provide nitrogenous and carbonaceous compounds, long chain amino acids, vitamins and other essential growth nutrients.

Type of specimen


Clinical: faeces, urine; Water samples and Food samples

Specimen Collection and Handling

For clinical samples follow appropriate techniques for handling specimens as per established guidelines (4, 5). For food samples, follow appropriate techniques for sample collection and processing as per guidelines (6). For water samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (7). After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

In Vitro diagnostic use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets.



M1393 HiCrome™ MM Agar

HiCromeVeg™ Freedom from BSE / TSE worries
Single Streak Rapid Differentiation Series

HiCrome™ MM Agar (M1393) is also available as HiCrome™ MM HiVeg™ Agar (MV1393) wherein all the animal origin nutrients have been replaced by vegetable based nutrients.

HiCrome™ MM Agar

Recommended For identification and differentiation of *Salmonella* and non-*Salmonella* like *Citrobacter* from food, water and clinical samples.

M1393

Limitations

1. Due to nutritional variations, some strains may show poor growth.
2. Though most of the *Salmonella* produce H₂S, certain non H₂S producing *Salmonella* species may appear as colourless colonies.
3. Due to nutritional variations, some strains may show poor growth.
4. Final confirmation of suspected colonies must be carried out by serological and biochemical tests.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

Appearance of Powder : Cream to yellow coloured, homogeneous, free flowing powder.

Gelling : Firm, comparable with 1.5% Agar gel.

Colour and Clarity of prepared medium : Light amber coloured, slightly opalescent gel forms in Petri plates.

Reaction : Reaction of 4.91% w/v aqueous solution at 25°C. pH:7.6 ± 0.2.

Cultural Response : Cultural characteristics observed after an incubation at 35-37°C for 18-24 hours.

Organisms (ATCC)	Inoculum (CFU)	Growth	Recovery	Colour of colony
<i>Escherichia coli</i> (25922) (00013*)	50-100	luxuriant	≥50%	light blue
<i>Salmonella</i> Enteritidis (13076) (00030*)	50-100	luxuriant	≥50%	black centered
<i>Salmonella</i> Typhimurium (14028) (00031*)	50-100	luxuriant	≥50%	black centered
<i>Citrobacter freundii</i> (8090)	50-100	good-luxuriant	≥50%	colourless#
<i>Pseudomonas aeruginosa</i> (27853) (00025*)	50-100	good-luxuriant	≥50%	colourless
<i>Enterococcus faecalis</i> (29212) (00087*)	≥10 ³	inhibited	0%	-

Key : # = may show bluish green colour on prolonged incubation

* = corresponding WDCM Numbers

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label. Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (4, 5).

References

1. Miller R.G. and Mallison E.T., 2000, J. Food Protection, 63(10), 1443-46.
2. Miller R.G., Tate C.R., Mallinson E.T. and Scherrer J.A., 1991, Pault Sa 70:2429-32.
3. Greenwald R., Henderson R.W. and Yappaw S., 1991, J. Clin. Microbiol. 29:2354.
4. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
5. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S. and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.
6. Salfinger Y., and Tortorello M.L. Fifth (Ed.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.
7. Baird R.B., Eaton A.D., and Rice E.W., (Eds.), 2015, Standard Methods for the Examination of Water and Wastewater, 23rd ed., APHA, Washington, D.C.

For the detection of *Salmonella* species

HiCrome™ MM Agar, Modified

Recommended For identification and differentiation of *Salmonella* and non-*Salmonella* like *Citrobacter* from food, water and clinical samples.

M1816

Composition **	
Ingredients	Grams/Litre
Proteose peptone	6.00
Yeast extract	10.00
L-Lysine hydrochloride	5.00
D-Cellobiose	10.00
Lactose	10.00
Sucrose	10.00
D-Xylose	3.75
Ferric ammonium citrate	0.80
Sodium thiosulphate	6.80
Chromogenic mixture	0.20
Bromothymol blue	0.10
Agar	18.00

Final pH (at 25°C) 7.6 ± 0.2

** Formula adjusted, standardized to suit performance parameters

Directions

Suspend 80.65 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. DO NOT AUTOCLAVE OR OVERHEAT. Cool to 45-50°C. Mix well and pour into sterile Petri plates.

Principle and Interpretation

HiCrome™ MM Agar was formulated by Miller and Mallison (1) for specific isolation and detection of *Salmonellae*. This medium is superior to XLT4 Agar in supporting growth of *Salmonella* due to the presence of appropriate proportion of four sugars. HiCrome™ MM Agar, Modified is a slight modification of HiCrome™ MM Agar and designed to differentiate *Enterobacteriaceae* especially *Salmonella* from *Proteus* and *Citrobacter* group. The utilization of sugars by organisms results in pH-changes. This is used as a means of distinguishing *Salmonella* from competing bacteria on the basis of colony colour.

Salmonella are gram negative rods in the family *Enterobacteriaceae* present in the stomach and intestinal tissues of human & animals and are found in their wastes. *Salmonella* usually are unable to ferment the sugars (2) that support growth of competing bacteria. Thus other bacteria tend to overgrow *Salmonellae*, masking their presence. Proteose peptone is a source of carbon, nitrogen and other essential amino acid and growth factor. Yeast extract provides vitamin especially Group B vitamins required for growth. To add to the differentiating ability of the formulation, an H₂S indicator system, consisting of sodium thiosulphate and ferric ammonium citrate, is included for the visualization of hydrogen sulphide produced, resulting in the formation of colonies with black centers. Bromothymol blue act as a pH indicator. The inclusion of sugars like lactose, sucrose, xylose and cellobiose

provides source of fermentable carbohydrate which stimulate the better initial growth of *Salmonella* cells. Presence of lactose suppresses H₂S production by non-*Salmonellae* like *Citrobacter freundii*. A chromogenic mixture, present in this medium helps to differentiate between lactose fermenters and nonfermenters. Lactose fermenters give bluish green coloured colonies, which would have been impossible to differentiate with an indicator based on pH change.

Type of specimen

Clinical samples : Faeces, water and food samples

Specimen Collection and Handling

For clinical samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (3, 4).


For food samples, follow appropriate techniques for sample collection and processing as per guidelines (5).

For water samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (6).

After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

In Vitro diagnostic use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets



M1816 HiCrome™ MM Agar Modified

For the detection of *Salmonella* species

HiCrome™ MM Agar, Modified

Recommended For identification and differentiation of *Salmonella* and non-*Salmonella* like *Citrobacter* from food, water and clinical samples.

M1816

Limitations

1. Due to nutritional variations, some strains may show poor growth.
2. Though most of the *Salmonella* produce H₂S, certain non H₂S producing *Salmonella* species may appear as colourless colonies.
3. Certain *Salmonella* species which are lactose fermenters may show as bluish green coloured colonies
4. Further confirmation may be carried out on suspected colonies.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

- Appearance of Powder** : Cream to yellow homogeneous free flowing powder
- Gelling** : Firm, comparable with 1.8% Agar gel
- Colour and Clarity of prepared medium** : Bluish green coloured, clear to slightly opalescent gel forms in Petri plates
- Reaction** : Reaction of 8.07 % w/v aqueous solution at 25°C. pH : 7.6 ± 0.2.
- Cultural Response** : Cultural characteristics observed after an incubation at 35-37°C for 18 - 24 hours.

Organism (ATCC)	Inoculum (CFU)	Growth	Recovery	Colour of Colony
<i>Citrobacter freundii</i> (8090)	50-100	good - luxuriant	≥50%	Yellow# coloured
<i>Escherichia coli</i> (25922) (00013*)	50-100	luxuriant	≥50%	Bluish green
<i>Salmonella</i> Typhimurium (14028) (00031*)	50-100	luxuriant	≥50%	Black centered
<i>Salmonella</i> Enteritidis (13076) (00030*)	50-100	luxuriant	≥50%	Black centered with yellow zone
<i>Salmonella</i> Typhi (6539)	50-100	good - luxuriant	≥50%	Black centered
<i>Proteus mirabilis</i> (25933)	50-100	good - luxuriant	≥50%	Gray coloured
<i>Klebsiella pneumoniae</i> (13883) (00097*)	50-100	luxuriant	≥50%	Yellowish green, mucoid

key # : may show bluish green colour on prolonged incubation.

* = corresponding WDCM Numbers

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (3, 4).

References

1. Miller R.G. and Mallison E.T., 2000, J. Food Protection, 63(10), 1443-46.
2. Miller R.G., Tate C.R., Mallinson E.T. and Scherrer J.A., 1991, Pault Sa 70:2429-32.
3. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
4. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S. and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.
5. Salfinger Y., and Tortorello M.L. Fifth (Ed.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.
6. Baird R.B., Eaton A.D., and Rice E.W., (Eds.), 2015, Standard Methods for the Examination of Water and Wastewater, 23rd ed., APHA, Washington, D.C.

For identification of *Klebsiella pneumoniae*

HiCrome™ *Klebsiella* Selective Agar Base

Recommended for the isolation and detection of *Klebsiella* species from water and other sources. This medium can also be used in membrane filtration procedure.

M1573

Composition **

Ingredients	Grams/Litre
Peptone, special	12.00
Yeast extract	7.00
Sodium chloride	5.00
Bile salts mixture	1.50
Sodium lauryl sulphate (SLS)	0.10
Chromogenic mixture	0.20
Agar	15.00

Final pH (at 25°C) 7.1 ± 0.2

** Formula adjusted, standardized to suit performance parameters

Directions

Suspend 20.4 grams in 500 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C and aseptically add rehydrated contents of one vial of *Klebsiella* Selective Supplement (FD225). Mix well and pour into sterile Petri plates.

Principle and Interpretation

HiCrome™ *Klebsiella* Selective Agar Base is recommended for isolation and enumeration of *Klebsiella* species based on chromogenic differentiation. *Klebsiella pneumoniae* strains are widely distributed in the environment and contribute to biochemical and geochemical process (1).

K. pneumoniae causes severe often fatal pneumonia. It also proves to be the source of lung infections that generally occur in patients with debilitating conditions such as alcoholism, diabetes mellitus, and chronic obstructive pulmonary disease (2). The chromogenic substrate incorporated in the media is cleaved specifically by *Klebsiella* species. *K. pneumoniae*, the causative agent of pneumonia, produces a purple-magenta coloured colony thereby aiding in the easy detection of the organisms. Most of the frequently encountered gram-negative faecal contaminants are inhibited on this media using a selective supplement.

Peptone special and yeast extract provide nitrogenous and carbonaceous compounds, long chain amino acids, vitamins and other essential nutrients required for the growth of the organism. Sodium chloride maintains the osmotic equilibrium of the medium. Bile salts mixture and sodium lauryl sulphate (SLS) inhibits most of the accompanying flora. Addition of the selective supplement further increases the selectivity of the medium.


Type of specimen

Water samples

Specimen Collection and Handling

For water samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (3).

After use, contaminated materials must be sterilized by autoclaving before discarding.



M1573 HiCrome™ *Klebsiella* Selective Agar Base

HiCromeVeg™ Freedom from BSE / TSE worries
Single Streak Rapid Differentiation Series

HiCrome™ *Klebsiella* Selective Agar Base (M1573) is also available as HiCrome™ *Klebsiella* Selective HiVeg Agar Base (MV1573) wherein all the animal origin nutrients have been replaced by vegetable based nutrients.

For identification of *Klebsiella pneumoniae*

HiCrome™ *Klebsiella* Selective Agar Base

Recommended for the isolation and detection of *Klebsiella* species from water and other sources. This medium can also be used in membrane filtration procedure.

M1573

Warning and Precautions

Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets

Limitations

Some organisms may show poor growth due to nutritional variation.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

- Appearance of Powder** : Cream to yellow coloured, homogeneous, free flowing powder.
- Gelling** : Firm, comparable with 1.5% Agar gel.
- Colour and Clarity of prepared medium** : Light amber coloured, clear to slightly opalescent gel forms in Petri plates.
- Reaction** : Reaction of 4.08% w/v aqueous solution at 25°C. pH:7.1 ± 0.2.
- Cultural Response** : Cultural characteristics observed with added *Klebsiella* Selective Supplement (FD225) after an incubation at 35-37°C for 18-24 hours.

Organisms (ATCC)	Inoculum (CFU)	Growth	Recovery	Colour of Colony
<i>Klebsiella pneumoniae</i> (13883) (00097*)	50-100	luxuriant	≥50%	purple-magenta (mucoid)
# <i>Klebsiella aerogenes</i> (13048) (00175*)	≥10 ³	inhibited	0%	-
<i>Escherichia coli</i> (25922) (00013*)	≥10 ³	inhibited	0%	-
<i>Serratia marcescens</i> (8100)	≥10 ³	inhibited	0%	-
<i>Salmonella</i> Typhi (6539)	≥10 ³	inhibited	0%	-

Key : * = corresponding WDCM Numbers

: Formerly known as *Enterobacter aerogenes*

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (4, 5).

References

- Krieg, N. R., and J. G. Holt, (Eds.), 1984, Bergey's Manual of Systematic Bacteriology, Vol. 1, p. 408 - 516. The Williams and Wilkins Co., Baltimore, Md.
- Wyngaarden J. B., Smith L. H., (Eds.), Cecil Text book of Medicine, 16th Ed, pp 1430 -1432, Philadelphia, W. B. Saunders, 1982.
- Baird R.B., Eaton A.D., and Rice E.W., (Eds.), 2015, Standard Methods for the Examination of Water and Wastewater, 23rd ed., APHA, Washington, D.C.
- Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
- Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

For the selective isolation of ESBL producing *Enterobacteriaceae*

HiCrome™ ESBL Agar Base

Recommended for selective isolation of Extended-Spectrum β -lactamase-producing *Enterobacteriaceae*.

M1829

Composition **

Ingredients	Grams/Litre
Peptone mix	12.000
Chromogenic mixture	4.000
Sodium chloride	5.000
Buffer mix	4.000
Agar	15.000

Final pH (at 25°C) 6.8 ± 0.2

** Formula adjusted, standardized to suit performance parameters

Directions

Suspend 40 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C and add rehydrated contents of two vials of HiCrome™ ESBL Agar Supplement (FD278). Mix well and pour into sterile Petri plates.

Principle and Interpretation

Extended-spectrum β -lactamase (ESBL)-producing organisms are an increasing challenge for healthcare practitioners fighting healthcare-associated infections (HAIs). *Escherichia coli*, *Klebsiella pneumoniae* and *Klebsiella oxytoca* are the most common ESBL-producing pathogens (1). ESBL-producing organisms are generally resistant to many classes of antibiotics, including aminoglycosides and fluoroquinolones; ESBL-producing organisms are able to attack newer cepheims and monobactams as well as narrow-spectrum cephalosporins and anti gram-negative penicillins (1). They are associated with increased mortality and are difficult to detect and treat. The widespread use of extended-spectrum, third-generation cephalosporins, introduced in the 1980s to treat antibiotic-resistant bacteria, is believed to be a major contributor to the emergence of ESBL-producing organisms.

HiCrome™ ESBL Agar Base is chromogenic screening medium for the selective isolation of ESBL producing organisms. It contains peptone mix which serves as the carbon and nitrogen sources, long chain amino acids, vitamins and other growth nutrients. Chromogenic mixture is used to differentiate the ESBL producing organisms on the basis of colour. HiCrome™ ESBL Agar Supplement (FD278) helps in inhibition of other contaminating organisms. ESBL producing *E.coli* grow as pink to purple colonies. ESBL producing members of the KESC group produce bluish green colonies; *Proteus*, *Morganella* and *Providencia* do not utilize any chromogen resulting in colourless to light brown colonies. This medium can be inoculated with liquid suspension equivalent to 0.5 McFarland turbidity, prepared from rectal screening swabs, faecal samples or from isolated colony. Isolated colonies should not be directly plated on to this medium, because the high level inoculum may cause false positive results. Further confirmation using biochemical identification tests is recommended.

Type of specimen

Clinical samples

Specimen Collection and Handling

For clinical samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (2, 3).


After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

In Vitro diagnostic use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets

Limitations

1. Some species may show poor growth due to nutritional variations.
2. Slight colour variation may be observed depending upon strains.
3. Isolated colonies should not be directly plated on to this medium, because the high level inoculum may cause false positive results. Further confirmation using biochemical identification tests is recommended.



M1829 – HiCrome™ ESBL Agar Base

For the selective isolation of ESBL producing *Enterobacteriaceae*

HiCrome™ ESBL Agar Base

Recommended for selective isolation of Extended-Spectrum β -lactamase-producing *Enterobacteriaceae*.

M1829

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

- Appearance of Powder** : Cream to yellow homogeneous free flowing powder
- Gelling** : Firm, comparable with 1.5% Agar gel
- Colour and Clarity of prepared medium** : Yellow coloured opalescent gel forms in Petri plates.
- Reaction** : Reaction of 4.0% w/v aqueous solution at 25°C. pH : 6.8 \pm 0.2.
- Cultural Response** : Cultural characteristics observed with added HiCrome™ ESBL Agar Supplement (FD278) after an incubation at 35-37°C for 18-24 hours.

Organism (ATCC)	Inoculum (CFU)	Growth	Recovery	Colour of Colony
<i>Escherichia coli</i> (NCTC 13351)	50-100	luxuriant	\geq 50%	pink to purple
<i>Klebsiella pneumoniae</i> (700603)	50-100	luxuriant	\geq 50%	bluish green
<i>Enterobacter cloacae</i> (23355) (00082*)	$\geq 10^3$	inhibited	0%	—
<i>Citrobacter freundii</i> (8090)	$\geq 10^3$	inhibited	0%	—
<i>Candida albicans</i> (10231) (00054*)	$\geq 10^3$	inhibited	0%	—

Key : * = corresponding WDCM Numbers

Storage and Shelf-life


Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label. Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (2, 3).

References

1. Journal of Clinical Microbiology, February 2007, Page 501-505, Vol. 45, No. 2 pp 1430 -1432, Philadelphia, W. B. Saunders, 1982.
2. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
3. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.



M1829 – HiCrome™ ESBL Agar Base

For the detection of Carbapenem Resistant *Enterobacteriaceae*

HiCrome™ KPC Agar Base

Recommended for the detection of gram negative bacteria with a reduced susceptibility to a carbapenem agents.

M1831

Composition **

Ingredients	Grams/Litre
Peptone	15.00
Chromogenic mixture	3.00
Agar	15.00

Final pH (at 25°C) 7.0 ± 0.2

** Formula adjusted, standardized to suit performance parameters

Directions

Suspend 16.50 grams in 500 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C and aseptically add rehydrated contents of 1 vial of HiCrome™ KPC Agar Supplement (FD279). Mix well and pour into sterile Petri plates.

Principle and Interpretation

HiCrome™ KPC Agar Base is a chromogenic medium designed for the detection and differentiation of KPC producing gram negative bacterial species without selective pre-enrichment. Carbapenems are the last line of defense against invasive or serious infections and are used to treat these life threatening infections that are caused by gram negative, drug resistant pathogens (2). Production of carbapenemase enzyme results in resistance to penicillins, cephalosporins (i.e. cefepime, ceftriaxone), carbapenems (i.e. meropenem, ertapenem) and aztreonam there by making these pathogens multi drug resistant.

Most carbapenemase producing bacteria are included in the family *Enterobacteriaceae* and are thus termed as carbapenem resistant *Enterobacteriaceae* (CRE). Besides the *Enterobacteriaceae* family, rare strains of *Pseudomonas aeruginosa* and *Acinetobacter baumannii* have also found to produce carbapenemase (1, 2, 3).

Peptone provides nitrogenous and carbonaceous compounds long chain amino acids and other essential growth nutrients. This medium can be made selective by supplementation with antibiotics for detecting microorganisms associated with hospital borne infections. Selective supplement have been added to inhibit the growth of yeast, gram positive organisms and gram negative organisms that do not produce carbapenemase.

This medium is intended to be used as a screening medium. Isolates should be tested further for carbapenem susceptibility following CLSI guidelines. Indole test may be perform for the confirmation of carbapenem resistant *E. coli* because some rare strains of *C. freundii* may produce small pink to magenta coloured colonies similar to *E. coli*. Carbapenem resistant strains of *Klebsiella*, *Enterobacter* and *Serratia* species produce bluish green colonies. *Acinetobacter* and *Salmonella* species produce smooth, colourless colonies. *Pseudomonas* species

produce colourless to light yellowish green, translucent colonies with wrinkled edges. Further biochemical tests may be needed for complete identification.

Type of specimen

Clinical samples

Specimen Collection and Handling

For clinical samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (1).


After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

In Vitro diagnostic use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets

Limitations

1. Some species may show poor growth due to nutritional variations and resistance to antibiotics.
2. Slight colour variation may be observed depending upon strains.



Klebsiella pneumoniae
ATCC BAA 1705

M1831 – HiCrome™ KPC Agar Base

For the detection of Carbapenem Resistant *Enterobacteriaceae*

HiCrome™ KPC Agar Base

Recommended for the detection of gram negative bacteria with a reduced susceptibility to a carbapenem agents.

M1831

3. Indole test may be performed for the confirmation of carbapenem resistant *E. coli* because some rare strains of *C. freundii* may produce small pink to magenta coloured colonies similar to *E. coli*.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

- Appearance of Powder** : Cream to yellow homogeneous free flowing powder.
- Gelling** : Firm, comparable with 1.5% Agar gel
- Colour and Clarity of prepared medium** : Light amber coloured, clear to slightly opalescent gel forms in Petri plates.
- Reaction** : Reaction of 3.3% w/v aqueous solution at 25°C. pH : 7.0 ± 0.2.
- Cultural Response** : Cultural characteristics observed with added HiCrome™ KPC Agar Supplement (FD279) after an incubation at 35-37°C for 18-24 hours.

Organism (ATCC)	Inoculum (CFU)	Growth	Recovery	Colour of Colony
<i>Enterococcus faecalis</i> (29212) (00087*)	≥10 ³	inhibited	0%	—
<i>Klebsiella pneumoniae</i> (BAA 1705)	50-100	luxuriant	≥50%	bluish green
<i>Klebsiella pneumoniae</i> (13883) (00097*)	≥10 ³	inhibited	0%	—
<i>Candida albicans</i> (10231) (00054*)	≥10 ³	inhibited	0%	—
<i>Staphylococcus aureus sub-sp aureus</i> (25923) (00034*)	≥10 ³	inhibited	0%	—

Key : * = corresponding WDCM Numbers

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition. Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (4, 5).

References

- Hindiyeth, M., et. al. 2008, J. Clin. Microbiol.; Vol. 46, p.2879 -2883
- Pillai D.R. et.al. 2009. Emerg. Infect. Dis; Vol. 15, P.827-829
- Samra, Z., 2008, J. Clin. Microbiol; Vol. 146, P.3110-3111.
- Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
- Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S. and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

Ready Prepared Media			
Code	Product Name	Usage	Packing
Category : 90 mm Sterile Ready Prepared Plate			
MP1831	HiCrome™ KPC Agar Plate	for detection of Gram-negative bacteria with a reduced susceptibility to carbapenem agents	20 plts 50 plts

HiCrome™ Vibrio Agar

Recommended for the isolation and selective chromogenic differentiation of *Vibrio* species from seafood.

M1682

Composition **

Ingredients	Grams/Litre
Peptone	10.00
Sodium chloride	25.00
Sodium thiosulphate	5.00
Sodium citrate	6.00
Sodium cholate	1.00
Chromogenic mixture	5.50
Agar	15.00

Final pH (at 25°C) 8.5 ± 0.2

** Formula adjusted, standardized to suit performance parameters

Directions

Suspend 67.5 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. DO NOT AUTOCLAVE. Cool to 45-50°C. Mix well and pour into sterile Petri plates.

Principle and Interpretation

Vibrio's have played a significant role in human history. Outbreaks of cholera, caused by *Vibrio cholerae*, can be traced back in time to early recorded descriptions of enteric infections. The *Vibrios* have also received the attention of marine microbiologists who observed that the readily cultured bacterial population in near-shore waters and those associated with fish and shell fish were predominantly *Vibrio* species (4). *Vibrio* species are mainly responsible for causing cholera and food poisoning in humans. *Vibrio cholerae* causes cholera due to the intake of contaminated food such as raw oysters. *Vibrio parahaemolyticus* is a major cause of food borne infections, causing food poisoning (1). Since *Vibrio* species naturally occur in sea-water, worth special mention is their need for sodium chloride, although some species can grow with minimum sodium chloride concentration (4). The widely used media for *Vibrio* isolation are TCBS Agar and Alkaline Peptone Water (2). However accompanying sucrose-fermenting bacteria pose a problem in the identification of *Vibrio* species on TCBS Agar. On HiCrome™ Vibrio Agar, the colour development by *Vibrio* species is not affected by the presence of colonies of other bacteria. This is because, the amount of colour developed depends on the reaction of the bacterial beta-galactosidase with the substrate contained in the media (3).

Peptone provides carbonaceous and nitrogenous compounds, long chain amino acids and essential nutrients to the organisms. High concentration of sodium chloride in addition to maintaining the osmotic equilibrium also has an inhibitory action on the accompanying microflora. Sodium thiosulphate, sodium citrate and sodium cholate are used in the formulation because they can inhibit the growth of gram positive and some gram negative bacteria, but not members of *Enterobacteriaceae*. The proprietary chromogenic mixture incorporated in the medium helps in the chromogenic differentiation of *Vibrio cholerae* and *Vibrio parahaemolyticus*. The high (alkaline) pH of the medium helps in selective isolation of *Vibrio* species.


Type of specimen

Food samples

Specimen Collection and Handling

For food samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (1).

After use, contaminated materials must be sterilized by autoclaving before discarding.



M1682 HiCrome™ Vibrio Agar

HiCromeVeg™ Freedom from BSE / TSE worries
Single Streak Rapid Differentiation Series

HiCrome™ Vibrio Agar (M1682) is also available as HiCrome™ Vibrio HiVeg™ Agar (MV1682) & HiCrome™ Vibrio HiCynth™ Agar (MCD1682) wherein all the animal origin nutrients have been replaced by vegetable based nutrients and chemically defined peptones respectively.

HiCrome™ Vibrio Agar

Recommended for the isolation and selective chromogenic differentiation of *Vibrio* species from seafood.

M1682

Warning and Precautions

Read the label before opening the container. Wear protective gloves/ protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets

Limitations

1. Being highly selective, some species may show poor growth due to nutritional variations.
2. Slight colour variation may be observed depending upon strains.
3. Further biochemical tests must be carried out for confirmation.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

- Appearance of Powder** : Light yellow to light tan coloured, homogeneous, free flowing powder.
- Gelling** : Firm, comparable with 1.5% Agar gel.
- Colour and Clarity of prepared medium** : Light yellow coloured, clear to slightly opalescent gel forms in Petri plates.
- Reaction** : Reaction of 6.75% w/v aqueous solution at 25°C. pH:8.5 ± 0.2.
- Cultural Response** : Cultural characteristics observed after an incubation at 35-37°C for 18-24 hours.

Organisms (ATCC)	Inoculum (CFU)	Growth	Recovery	Colour of colony
<i>Vibrio cholerae</i> (15748)	50-100	good-luxuriant	≥50%	purple
<i>Vibrio parahaemolyticus</i> (17802) (00037*)	50-100	good-luxuriant	≥50%	bluish-green
<i>Enterococcus faecalis</i> (29212) (00087*)	≥10 ³	inhibited	0%	
<i>Staphylococcus aureus</i> subsp aureus (25923) (00034*)	≥10 ³	inhibited	0%	
<i>Escherichia coli</i> (25922) (00013*)	≥10 ³	inhibited	0%	

Key * = corresponding WDCM Numbers

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (6, 7).

References

1. Alcamo, E.I., 2001. Fundamentals of Microbiology, 6th ed, Jones and Bartlett Publishers, Inc. pg 254, 244.
2. Clesceri, Greenberg and Eaton (ed.), 1998. Standard Method for the examination of Water and Waste water, 20th ed. American Public Health Association, Washington, D. C.
3. Kudo, H. Y et al, 2001. Improved Method for Detection of *Vibrio parahaemolyticus* in Seafood. ASM. Vol 67, No. 12, pg 5819-5823.
4. Thompson et al (ed.). 2006. The Biology of Vibrios, ASM Press, chapter1, pg 3.
5. Salfinger Y., and Tortorello M.L. Fifth (Ed.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.
6. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
7. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

For identification and differentiation of UTI causing organisms

HiCrome™ UTI Agar, Modified / HiCrome™ UTI Selective Agar

A chromogenic differential medium for identification, differentiation and confirmation of enteric bacteria from specimens such as urine which may contain large number of *Proteus* species as well potentially pathogenic gram - postive organisms.

M1418/
M1505

Composition **

	M1418	M1505
Ingredients	Grams/Litre	Grams/Litre
Peptone	18.00	18.00
Tryptone	4.00	4.00
HM peptone B#	6.00	6.00
Chromogenic mixture	12.44	12.44
Bile salts	-	1.50
Agar	15.00	15.00

Final pH (at 25°C) 7.2 ± 0.2

** Formula adjusted, standardized to suit performance parameters

#Equivalent to Beef extract

Directions

Suspend 56.94 grams of M1505 or 55.44 grams of M1418 in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Mix well and pour into sterile Petri plates.

Principle and Interpretation

HiCrome™ UTI Agar, Modified/ HiCrome™ UTI Selective Agar is formulated on the basis of work carried out by Pezzlo (4), Wilkie et al (6), Friedman et al (1), Murray et al (3), Soriano and Ponte (5) and Merlino et al (2). These media are the modifications of HiCrome™ UTI Agar (M1353), which can be used in place of MacConkey Agar for isolation, and confirmation of various microorganisms. It facilitates and expedites the identification of some gram-negative bacteria and some gram-positive bacteria on the basis of different contrasted colony colours produced by reactions of genus or species specific enzymes with two chromogenic substrates.

Enzymes produced by *Enterococcus* species, *Escherichia coli* and coliforms cleave the chromogenic substrates incorporated in the medium. Presence of rich source of phenylalanine and tryptophan from peptone and tryptone provides an indication of tryptophan deaminase activity, revealed with TDA Reagent (R036) indicating the presence of *Proteus* species, *Morganella* species and *Providencia* species, which appear brown. One chromogenic substrate is cleaved by β -glucosidase possessed by enterococci resulting in formation of blue colonies.


E. coli produce purple to magenta colonies due to the enzyme β -D-galactosidase which cleaves the other chromogenic substrate. Further confirmation of *E. coli* can be done by performing indole test using DMACA Reagent (R035). Also, some strains of *Enterobacter cloacae* lacking β -glucosidase show pink colonies indistinguishable from *E. coli*. The DMACA reagent for indole test (should be performed on filter paper) distinguishes between *E. coli* and *Enterobacter* and TDA reagent between *Proteus mirabilis* and other species. Coliforms produce purple coloured colonies due to cleavage of both the chromogenic substrates. Peptone, HM peptone B and Tryptone provides nitrogenous and carbonaceous compounds, long chain amino acids and other essential growth nutrients. HiCrome™ UTI Selective Agar is made selective by the addition of bile salts, which selectively inhibits gram-positive bacteria.

Type of specimen

Clinical samples

Specimen Collection and Handling

For clinical samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (7, 8). After use, contaminated materials must be sterilized by autoclaving before discarding.



M1418 HiCrome™ UTI Agar, Modified

HiCromeVeg™ Freedom from BSE / TSE worries
Single Streak Rapid Differentiation Series

HiCrome™ UTI Agar, Modified / HiCrome™ UTI Selective Agar (M1418/M1505) is also available as HiCrome™ UTI HiVeg™ Agar, Modified / HiCrome™ UTI Selective HiVeg™ Agar (MV1418/MV1505) wherein all the animal origin nutrients have been replaced by vegetable based nutrients.

For identification and differentiation of UTI causing organisms

HiCrome™ UTI Agar, Modified / HiCrome™ UTI Selective Agar

A chromogenic differential medium for identification, differentiation and confirmation of enteric bacteria from specimens such as urine which may contain large number of *Proteus* species as well potentially pathogenic gram - positive organisms.

M1418/
M1505

Warning and Precautions

In Vitro diagnostic use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets

Limitations

1. Some strains of *Enterobacter cloacae* lacking β -glucosidase show pink colonies indistinguishable from *E. coli*.
2. TDA reagent between *Proteus mirabilis* and other species.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

Appearance of Powder : Cream to yellow coloured, homogeneous, free flowing powder.

Gelling : Firm, comparable with 1.5% Agar gel.

Colour and Clarity of prepared medium : Light amber coloured, clear to slightly opalescent gel forms in Petri plates.

Reaction : Reaction of 5.54% w/v of M1418 or 5.69% w/v of M1505 aqueous solution at 25°C. pH:7.2 ± 0.2.

Cultural Response : Cultural characteristics observed after an incubation at 35-37°C for 18-24 hours.

Organisms (ATCC)	Inoculum (CFU)	Growth «	Recovery «	Growth ««	Recovery ««	Colour of colony	TDA #	DMACA ##
<i>Escherichia coli</i> (25922) (00013*)	50-100	luxuriant	≥70%	luxuriant	≥50%	purple-magenta	-	+
<i>Proteus mirabilis</i> (12453)	50-100	luxuriant	≥70%	luxuriant	≥50%	light brown	+	-
<i>Klebsiella pneumoniae</i> (13883) (00097*)	50-100	luxuriant	≥70%	luxuriant	≥50%	blue to purple, mucoid	-	-
<i>Pseudomonas aeruginosa</i> (27853) (00025*)	50-100	luxuriant	≥70%	luxuriant	≥50%	colourless greenish pigment may be observed	-	-
<i>Enterococcus faecalis</i> (29212) (00087*)	50-100	luxuriant	≥70%	fair	20-30%	blue - blue green (small)	-	-
<i>Staphylococcus aureus</i> subsp <i>aureus</i> (25923) (00034*)	50-100	luxuriant	≥70%	inhibited	0%	golden yellow*	-	-

Key : TDA + : Tryptophan deaminase present , TDA - : Tryptophan deaminase absent

DMACA + : Indole positive, DMACA - : Indole negative,

« : on HiCrome UTI Agar, Modified «« : on HiCrome UTI Selective Agar

: Add 1-2 drops of TDA reagent directly on suspected colony. Brown colouration-positive.

: Transfer suspected colony on filter paper, dipped in DMACA reagent. Bluish purple colouration- positive.

* = corresponding WDCM Numbers

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (8, 9).

References

1. Friedman M. P. et al, 1991, J. Clin. Microbiol., 29:2385-2389.
2. Merlino et al, 1995, Abstr. Austr. Microbiol., 16(4):17-3.
3. Murray P. R., Tenover P. C. and Tenover F. C., 1992, J. Clin. Microbiol., 30:1600-1601.
4. Pezzlo M., 1998, Clin. Microbiol. Rev., 1:268-280.
5. Soriano F. and Ponte C., 1992, J. Clin. Microbiol., 30:3033-3034.
6. Wilkie M. E., Almond M. K. and Marsh F. P., 1992, British Medical Journal, 305:1137-1141.
7. Lapage S., Shelton J. and Mitchell T., 1970, Methods in Microbiology, Norris J. and Ribbons D., (Eds.), Vol. 3A, Academic Press, London.
8. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S. and Tenover, F.C. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.
9. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.

For identification and differentiation of UTI causing organisms

HiCrome™ UTI Agar

A differential medium recommended for presumptive identification and confirmation of microorganisms mainly causing urinary tract infections. Can also be used for testing water, food, environmental and other clinical samples.

M1353 /
M1353R

Composition **

	M1353	M1353R
Ingredients	Grams/Litre	Grams/Litre
Peptone, special	15.00	-
Peptone	-	15.00
Chromogenic mixture	2.45	26.80
Agar	15.00	15.00

Final pH (at 25°C) 6.8 ± 0.2

** Formula adjusted, standardized to suit performance parameters

Directions

Suspend 32.45 grams of M1353 or 56.8 grams of M1353R in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Mix well and pour into sterile Petri plates.

Principle and Interpretation

Urinary tract infections are bacterial infections affecting parts of urinary tract. The common symptoms of urinary tract infection are urgency and frequency of micturition, with associated discomfort or pain. The common condition is cystitis, due to infection of the bladder with a uropathogenic bacterium, which most frequently is *Escherichia coli*, but sometimes *Staphylococcus saprophyticus* or especially in hospital-acquired infections, *Klebsiella* species, *Proteus mirabilis*, other coliforms, *Pseudomonas aeruginosa* or *Enterococcus faecalis* (1). HiCrome™ UTI Agar is formulated on basis of work carried out by Pezzlo (5) Wilkie et al (7), Friedman et al (2), Murray et al (4), Soriano and Ponte (6) and Merlino et al (3). These media are recommended for the detection of urinary tract pathogens where HiCrome™ UTI Agar has broader application as a general nutrient agar for isolation of various microorganisms. It facilitates and expedites the identification of some gram-negative bacteria and some gram-positive bacteria on the basis of different contrasted colony colours produced by reactions of genus or species specific enzymes with two chromogenic substrates. The chromogenic substrates are specifically cleaved by enzymes produced by Enterococcus species, *E. coli* and coliforms. Presence of amino acids like phenylalanine and tryptophan from peptones helps for detection of tryptophan deaminase activity, indicating the presence of *Proteus* species, *Morganella* species and *Providencia* species.

One of the chromogenic substrate is cleaved by β -glucosidase possessed by enterococci resulting in formation of blue colonies. *E. coli* produce pink-purple colonies due to the enzyme β -D-galactosidase that cleaves the other chromogenic substrate. Further confirmation of *E. coli* can be done by performing the indole test. Coliforms produce purple coloured colonies due to cleavage of both the chromogenic substrate. Colonies of *Proteus*, *Morganella* and *Providencia* species appear brown because of tryptophan deaminase activity. Peptone or

peptone special provides nitrogenous, carbonaceous compounds long chain amino acids and other essential growth nutrients.

This medium can be made selective by supplementation with antibiotics for detecting microorganisms associated with hospital borne infections.

Type of specimen

Clinical, Food & Water samples

Specimen Collection and Handling

For clinical samples follow appropriate techniques for handling specimens as per established guidelines (8, 9).


For food and dairy samples, follow appropriate techniques for sample collection and processing as per guidelines (10).

For water samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards.(11).

After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

In Vitro diagnostic use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets



M1353R HiCrome™ UTI Agar

HiCromeVeg™ Freedom from BSE / TSE worries
Single Streak Rapid Differentiation Series

HiCrome™ UTI Agar (M1353) is also available as HiCrome™ UTI HiVeg™ Agar (MV1353) & HiCrome™ UTI HiCynth™ Agar (MCD1353) wherein all the animal origin nutrients have been replaced by vegetable based nutrients & Chemically defined peptones respectively.

For identification and differentiation of UTI causing organisms

HiCrome™ UTI Agar

A differential medium recommended for presumptive identification and confirmation of microorganisms mainly causing urinary tract infections. Can also be used for testing water, food, environmental and other clinical samples.

M1353 /
M1353R

Limitations

1. Some species may show poor growth due to nutritional variations.
2. Slight colour variation may be observed depending upon strains.
3. Further confirmation of *E. coli* can be done by performing the indole test.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

- Appearance of Powder** : Cream to yellow (M1353) or white to cream (M1353R) coloured, homogeneous, free flowing powder.
- Gelling** : Firm, comparable with 1.5% Agar gel.
- Colour and Clarity** : Light amber coloured, clear to slightly opalescent of prepared medium gel forms in Petri plates of M1353 or white coloured opaque gel forms with precipitate in Petri plates of M1353R.
- Reaction** : Reaction of 3.24% w/v of M1353 or 5.68% w/v of M1353R aqueous solution at 25°C.
pH : 6.8 ± 0.2.
- Cultural Response** : Cultural characteristics observed after an incubation at 35-37°C for 18-24 hours.

Organisms (ATCC)	Inoculum	Growth	Recovery	Colour of colony
<i>Escherichia coli</i> (25922) (00013*)	50-100	luxuriant	≥70%	pink-purple
<i>Proteus mirabilis</i> (12453)	50-100	luxuriant	≥70%	light brown
<i>Klebsiella pneumoniae</i> (13883) (00097*)	50-100	luxuriant	≥70%	blue to purple, mucoid
<i>Pseudomonas aeruginosa</i> (27853) (00025*)	50-100	luxuriant	≥70%	colourless, greenish pigment may be observed
<i>Staphylococcus aureus subsp aureus</i> (25923) (00034*)	50-100	luxuriant	≥70%	golden yellow
<i>Enterococcus faecalis</i> (29212) (00087*)	50-100	luxuriant	≥70%	blue, small

Key : * = corresponding WDCM Numbers

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (8, 9).

References

1. Collee J. G., Fraser A. G., Marmion B. P., Simmons A., (Eds.), Mackie and McCartney, Practical Medical Microbiology, 1996, 14th Edition, Churchill Livingstone
2. Friedman M. P. et al, 1991, J. Clin. Microbiol., 29:2385-2389.
3. Merlino et al, 1995, Abstr. Austr. Microbiol. 16(4):17-3.
4. Murray P., Traynor P. Hopson D., 1992, J. Clin. Microbiol., 30:1600-1601.
5. Pezzlo M., 1998, Clin. Microbiol. Rev., 1:268-280
6. Soriano F., Ponte C., 1992, J. Clin. Microbiol., 30:3033-3034.
7. Wilkie M. E., Almond M. K., Marsh F. P., 1992, British Medical Journal 305:1137-1141.
8. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
9. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.
10. Salfinger Y., and Tortorello M.L. Fifth (Ed.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.
11. Baird R.B., Eaton A.D., and Rice E.W., (Eds.), 2015, Standard Methods for the Examination of Water and Wastewater, 23rd ed., APHA, Washington, D.C.

Ready Prepared Media			
Code	Product Name	Usage	Packing
Category : 90 mm Ready prepared Plates			
MP1353	HiCrome™ UTI Agar Plate	for presumptive identification of microorganisms mainly causing urinary tract infections	20 plts 50 plts
Category : Ready Prepared Solid Media in Glass Bottles			
SM1353	HiCrome™ UTI Agar	for presumptive identification & confirmation of microorganisms mainly causing urinary tract infections & other clinical samples.	5X100 ml
Category : HiTouch™ Flexi Plate™			
FL031	HiTouch™ HexaCrome Flexiplate™	for differentiation of six pathogenic organisms - <i>E.coli</i> , <i>Klebsiella</i> , <i>Enterococcus</i> , <i>Proteus</i> , <i>S. aureus</i> & <i>Pseudomonas</i>	50 plts

For identification and differentiation of UTI causing organisms

HiCrome™ Universal Differential Medium

Recommended for presumptive identification of microorganisms from clinical and non-clinical specimens.

M1600

Composition **

Ingredients	Grams/Litre
Peptone	15.00
Tryptone	4.00
Chromogenic mixture	2.50
Agar	13.50

Final pH (at 25°C) 7.2 ± 0.2

** Formula adjusted, standardized to suit performance parameters

Directions

Suspend 35.00 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Mix well and pour into sterile Petri plates.

Principle and Interpretation

HiCrome™ Universal Differential Medium is a modification of the medium formulated on basis of the work carried out by Pezzlo (4), Wilkie et al (6), Friedman et al (1), Murray et al (3), Soriano and Ponte (5) and Merlino et al (2). HiCrome™ Universal Differential Medium is recommended for the presumptive identification of microorganisms from clinical and non-clinical specimens where the medium has broader application as a general nutrient agar for isolation of various microorganisms. This medium helps in the identification of some gram-positive bacteria and gram-negative bacteria on the basis of different colony colours exhibited by them. These colours are formed due to the reactions of genus or species specific enzymes with the two chromogenic substrates incorporated in the medium. Enterococcus species, *Escherichia coli* and coliforms produce enzymes which specifically cleave these chromogenic substrates to give characteristically distinctive colony colours. Peptone in the medium serve as sources of amino acids like phenylalanine and tryptophan which aids in indicating tryptophan deaminase activity, thereby facilitating the identification of *Proteus* species, *Morganella* species and *Providencia* species. One of the chromogenic substrate is cleaved by β -glucosidase enzyme possessed by Enterococci resulting in the formation of bluish green colonies. *Escherichia coli* possesses the enzyme β -galactosidase which specifically cleaves the other chromogenic substrate resulting in the formation of purple coloured colonies. *Escherichia coli* can be differentiated and confirmed from other similar coloured colonies, by performing the indole test.

Coliforms cleave both the chromogenic substrates forming blue to purple coloured colonies. Colonies of *Proteus*, *Morganella* and *Providencia* species appear brown due to tryptophan deaminase activity. Peptone and Tryptone provide nitrogenous and carbonaceous compounds, essential growth nutrients and also serve as a source of amino acids.

Type of specimen

Clinical samples , Food samples.

Specimen Collection and Handling

For clinical samples follow appropriate techniques for handling specimens as per established guidelines (7, 8).


For food samples, follow appropriate techniques for sample collection and processing as per guidelines (9).

For water samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards.(10)

After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

In Vitro diagnostic use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets



M1600 HiCrome™ Universal Differential Medium

For identification and differentiation of UTI causing organisms

HiCrome™ Universal Differential Medium

Recommended for presumptive identification of microorganisms from clinical and non-clinical specimens.

M1600

Limitations

1. Some species may show poor growth due to nutritional variations.
2. Slight colour variation may be observed depending upon strains.
3. *Escherichia coli* can be differentiated and confirmed from other similar coloured colonies, by performing the indole test.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

- Appearance of Powder** : Cream to yellow homogeneous free flowing powder
- Gelling** : Firm, comparable with 1.35% Agar gel.
- Colour and Clarity** : Light amber coloured, clear to slightly of prepared medium opalescent gel forms in Petri plates
- Reaction** : Reaction of 3.5% w/v aqueous solution at 25°C. pH : 7.2 ± 0.2.
- Cultural Response** : Cultural characteristics observed after an incubation at 35-37°C for 18-24 hours.

Organism (ATCC)	Inoculum (CFU)	Growth	Recovery	Colour of colony
<i>Enterococcus faecalis</i> (29212) (00087*)	50-100	luxuriant	≥70%	blue, small
<i>Escherichia coli</i> (25922) (00013*)	50-100	luxuriant	≥70%	purple
<i>Klebsiella pneumoniae</i> (13883) (00097*)	50-100	luxuriant	≥70%	blue - green, mucoid
<i>Pseudomonas aeruginosa</i> (27853) (00025*)	50-100	luxuriant	≥70%	colourless (greenish pigment may be observed)
<i>Proteus mirabilis</i> (12453)	50-100	luxuriant	≥70%	light brown
<i>Staphylococcus aureus</i> subsp aureus (25923) (00034*)	50-100	luxuriant	≥70%	golden yellow

<i>Salmonella</i> Typhi (6539)	50-100	luxuriant	≥70%	colourless
<i>Salmonella</i> Typhimurium (14028) (00031*)	50-100	luxuriant	≥70%	colourless

Key : * = corresponding WDCM Numbers

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (7, 8).

References

1. Friedman M.P. et al (1991), Journal of Clinical Microbiology, 29:2385-2389.
2. Merlino et al (1995) Abstr. Austr. Microbiol. 16(4):17-3.
3. Murray P., Traynor P. Hopson D., (1992), Journal of Clinical Microbiology 30:1600-1601.
4. Pezzlo M (1998), Clinical Microbiology Reviews 1:268-280.
5. Soriano F., Ponte C., (1992), Journal of Clinical Microbiology 30:3033-3034.
6. Wilkie M.E., Almond M.K., Marsh F.P. (1992), British Medical Journal 305:1137-1141.
7. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
8. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.
9. Salfinger Y., and Tortorello M.L. Fifth (Ed.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.
10. Baird R.B., Eaton A.D., and Rice E.W., (Eds.), 2015, Standard Methods for the Examination of Water and Wastewater, 23rd ed., APHA, Washington, D.C.

Ready Prepared Media			
Code	Product Name	Usage	Packing
Category : 90 MM Ready Prepared Petri Plates			
MP1600	HiCrome™ Universal Agar Plate	for presumptive identification of microorganisms from clinical & non clinical specimens	20 plts
Category : HiTouch™ FlexiPlates			
FL042	HiTouch™ HiCrome™ Universal Agar Flexi Plate™	for presumptive identification of microorganisms from clinical & non clinical specimens	50 plts
Category : HiDip™ slides			
HD041	HiDip™ HiCrome™ Universal Agar-PCA	for differential & presumptive identification of microorganisms and total bacteria count.	50 plts

For antimicrobial susceptibility

HiCrome™ Mueller Hinton Agar

Recommended for differentiation of organisms based on chromogenic differentiation and determination of susceptibility of microorganisms to antimicrobial agents.

M2010

Composition **

Ingredients	Grams/Litre
Acicase#	20.00
Chromogenic mixture	1.50
Agar	17.00

Final pH (at 25°C) 7.3±0.1

** Formula adjusted, standardized to suit performance parameters

Casein acid hydrolysate

Directions

Suspend 38.50 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Mix well and pour into sterile Petri plates.

Principle and Interpretation

The Mueller Hinton formulation was originally developed for the cultivation of pathogenic *Neisseria* species (1). Other media were subsequently developed that replaced the use of Mueller Hinton Agar for the cultivation of pathogenic *Neisseria* species, but it became widely used in the determination of sulfonamide resistance of gonococci and other organisms. Mueller Hinton Agar is now used as a test medium for antimicrobial susceptibility testing (2).

Acicase provide nitrogenous compounds, carbon, sulphur and other essential nutrients. These ingredients are selected for low thymine and thymidine content as determined by MIC values for *Enterococcus faecalis* with sulfamethoxazole trimethoprim (SXT). Chromogenic mixture incorporated helps in colour differentiation. One of the chromogenic substrate is cleaved by β -glucosidase possessed by Enterococci resulting in formation of blue colonies. *E.coli* produce pink to purple colonies due to the enzyme β -D-galactosidase that cleaves the other chromogenic substrate. *Staphylococcus aureus* produces colourless colonies. *Pseudomonas aeruginosa* produces greenish pigmentation. *Klebsiella* and *Enterobacter* species produces metallic blue colored colonies. Colonies of *Proteus*, *Morganella* and *Providencia* species appear brown. This medium can be employed in screening urinary tract pathogens wherein organisms can be differentiated based on colour and simultaneously the antibiotic sensitivity can be determined.

Type of specimen

Clinical samples


Specimen Collection and Handling

For clinical samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (3,4).

After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

In Vitro diagnostic use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets



M2010 HiCrome™ Mueller Hinton Agar

Mixture of *Klebsiella*, *S. faecalis* and *E. coli*

HiCrome™ Mueller Hinton Agar

Recommended for differentiation of organisms based on chromogenic differentiation and determination of susceptibility of microorganisms to antimicrobial agents.

M2010

Limitations

1. Inoculum density may effect the zone size. Heavy inoculum may result in smaller zones or too less inoculum may result in bigger zones.
2. As antimicrobial susceptibility is carried with antibiotic disc, proper storage of the disc is desired which may effect the potency of the disc.
3. Under certain circumstances, the in vitro results of antibiotic susceptibility may not show the same in vivo.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

- Appearance of Powder :** Cream to yellow homogeneous free flowing powder
- Gelling :** Firm, comparable with 1.7% agar gel.
- Colour and Clarity :** Light amber coloured clear to slightly opalescent gel forms in Petri plates
- Reaction :** Reaction of 3.85% w/v aqueous solution at 25°C. pH : 7.3±0.1
- Cultural Response :** Cultural characteristics observed after an incubation at 35 - 37°C for 18 - 24 hours.

Organisms (ATCC)	Inoculum (CFU)	Growth	Recovery	Colour of colony
<i>Escherichia coli</i> (25922) (00013*)	50-100	luxuriant	≥70%	pink-purple
<i>Pseudomonas aeruginosa</i> (27853) (00025*)	50-100	luxuriant	≥70%	greenish pigment may be observed
<i>Staphylococcus aureus</i> (25923) (00034*)	50-100	luxuriant	≥70%	colourless-golden yellow
<i>Enterococcus faecalis</i> (29212) (00087*)	50-100	luxuriant	≥70%	blue
<i>Klebsiella pneumoniae</i> (13883) (00097*)	50-100	luxuriant	≥70%	metallic blue

Key : * = corresponding WDCM Numbers

Storage and Shelf-life


Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label. Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (3,4).

References

1. Mueller J. H. and Hinton J., 1941, Proc. Soc. Exp. Biol. Med., 48:330.
2. National Committee for Clinical Laboratory Standards, 2000, Approved Standard: M7-A5. Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria that grow aerobically, 5th Ed., NCCLS, Wayne, Pa.
3. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
4. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S. and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.



M2010 HiCrome™ Mueller Hinton Agar

Mixture of *Klebsiella*, *S. faecalis* and *E. coli*

For detection and isolation of pathogenic *Yersinia enterocolitica*

HiCrome™ *Yersinia* Agar Base

Recommended for isolation of pathogenic *Yersinia enterocolitica* from clinical specimens and food samples.

M2025

Composition **

Ingredients	Grams/Litre
Peptone mix	24.24
Selective mix	7.74
Chromogenic mixture	10.45
Growth factor	3.00
Agar	12.50

Final pH (at 25°C) 7.4±0.2

** Formula adjusted, standardized to suit performance parameters

Directions

Suspend 57.93 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45 - 50°C and aseptically add reconstituted contents of 1 vial of *Yersinia* Selective Supplement (FD034). Mix well before pouring into sterile Petri plates.

Principle and Interpretation

Yersinia enterocolitica is widely distributed in lakes and reservoirs ranging from intestinal tracts of numerous mammals. Environmental isolates are avirulent, however, isolates recovered from porcine sources contain human pathogens. Epizootic outbreaks of diarrhea, lymphadenopathy, pneumonia and spontaneous abortions occur in various animals due to these pathogenic serogroups. *Y. enterocolitica* is biochemically more active at room temperature than at 37°C. *Yersinia* Selective Agar Base with added *Yersinia* Selective Supplement is used to isolate *Y. enterocolitica* from clinical and food samples (3). *Yersinia* Selective Agar Base is recommended for selective isolation of *Yersinia* (1, 2) with modification of chromogenic identification.

Peptone mix and growth factor provides nitrogen and carbon source, long chain amino acids, vitamins and other essential growth nutrients. The medium is selective due to the presence of selective mix, which inhibit gram-positive and a number of gram-negative bacteria. Addition of antibiotic supplement makes it highly selective for *Yersinia*, thus imparting additional selectivity. One of the chromogen is split by *Yersinia* species and results in purple coloured colonies. Other organisms are either inhibited or results in colourless colonies.

For the isolation of *Y. enterocolitica* by direct plating and pour plating, inoculate the specimen directly onto the medium. Incubate at 22-32°C for 24-48 hours or suspend the sample (food, faeces, etc.) in sterile Phosphate Buffer Saline and incubate for upto 21 days (4) at 4°C.

Type of specimen

Clinical, food samples

Specimen Collection and Handling


For Clinical samples, follow appropriate techniques for handling specimens as per established guidelines (5, 6).

For food sample follow appropriate techniques for sample collection, processing as per guidelines and local standards (7).

After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

In Vitro diagnostic use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets



M2025 HiCrome™ *Yersinia* Agar Base

For detection and isolation of pathogenic *Yersinia enterocolitica*

HiCrome™ Yersinia Agar Base

Recommended for isolation of pathogenic *Yersinia enterocolitica* from clinical specimens and food samples.

M2025

Limitations

1. Some species may show poor growth due to nutritional variations.
2. Slight colour variation may be observed depending upon strains.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature..

Quality Control

- Appearance of Powder** : Light yellow to greenish yellow homogeneous free flowing powder.
- Gelling** : Firm, comparable with 1.25% Agar gel.
- Colour and Clarity** : Reddish purple coloured clear to slightly opalescent gel forms in Petri plates.
- Reaction** : Reaction of 5.8% w/v aqueous solution at 25°C. pH : 7.4±0.2
- Cultural Response** : Cultural characteristics observed with added *Yersinia* Selective Supplement (FD034) after an incubation at 22-32°C for 24-48 hours.

Organism (ATCC)	Inoculum (CFU)	Growth	Colour of colony
<i>Escherichia coli</i> O157:H7 (NCTC 12900)	>=10 ³	inhibited	
<i>Salmonella</i> Typhimurium ATCC 14028 (00031*)	>=10 ³	inhibited	
<i>Listeria monocytogenes</i> ATCC 19112	>=10 ³	inhibited	
<i>Campylobacter jejuni</i> ATCC 29428	>=10 ³	inhibited	

<i>Yersinia enterocolitica</i> ATCC 27729	50-100	good - luxuriant	Purple
<i>Escherichia coli</i> ATCC 25922 (00013*)	>=10 ³	inhibited	
<i>Enterococcus faecalis</i> ATCC 29212 (00087*)	>=10 ³	inhibited	

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (5, 6).

References

1. Schiemann D. A., 1979, Can. J. Microbiol., 25: 1298.
2. Schiemann D. A., 1980, Can. J. Microbiol., 26: 1232.
3. International Organization for Standardization (ISO), 1994, Draft ISO/DIS 10273.
4. Weissfeild and Sonnenwirth, 1982, J. Clin. Microbiol. 15 :508.
5. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
6. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.
7. Salfinger Y., and Tortorello M.L. Fifth (Ed.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.

For selective isolation and presumptive identification of *Clostridium* species

HiCrome™ Clostridial Agar Base

Recommended for selective isolation and presumptive identification of *Clostridium* species

M2026

Composition **

Ingredients	Grams/Litre
Tryptone	15.00
Yeast extract	10.00
Dextrose (Glucose)	1.00
Sodium chloride	5.00
Sodium thioglycollate	0.50
Chromogenic mixture	3.31
Agar	13.00

Final pH (at 25°C) 7.1 ± 0.2

** Formula adjusted, standardized to suit performance parameters

Directions

Suspend 47.81 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Aseptically add rehydrated contents of one vial of Perfringens Supplement II (FD012). Mix well and pour into sterile Petri plates.

Principle and Interpretation

One of the major species of anaerobic bacteria to cause disease in humans is *Clostridium*. *Clostridium* species cause tetanus and gas gangrene that ultimately leads to tissue damage. Another *Clostridium* species produces the lethal botulinum toxin, the causative agent of botulism (1). Clostridial Agar formulated by Vera is recommended for the selective isolation of pathogenic Clostridia from mixed flora (2). HiCrome is the modification for chromogenic differentiation.

Tryptone provide the essential nutrients, mainly the nitrogen compounds. Yeast extract serves as source of vitamins especially of the B group. Dextrose acts as fermentable carbohydrate source. Sodium thioglycollate is the reducing agents that help to create low oxidation-reduction potential enabling the growth of Clostridia. Also the media is well supplemented to support luxuriant growth of *Clostridium* species. The selective supplements inhibits other enteric bacteria.

The ideal method of inoculation of Clostridial Agar is direct inoculation of sterile, cooled medium with the specimen (in tubes). Alternatively agar plates of the medium can also be inoculated by streaking.

Type of specimen

Clinical samples

Specimen Collection and Handling

For clinical samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (3,4). After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

In Vitro diagnostic use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets

Limitations

1. Some species may show poor growth due to nutritional variations.
2. Slight colour variation may be observed depending upon strains.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature..



M2026 HiCrome™ Clostridial Agar Base

HiCrome™ Clostridial Agar Base

Recommended for selective isolation and presumptive identification of *Clostridium* species

M2026

Quality Control

- Appearance of Powder** : Cream to beige homogeneous free flowing powder.
- Gelling** : Firm, comparable with 1.3% Agar gel
- Colour and Clarity** : Yellow coloured, clear to slightly opalescent gel forms in Petri plates
- Reaction** : Reaction of 4.78% w/v aqueous solution at 25°C. pH : 7.1±0.2
- Cultural Response** : Cultural characteristics observed after an incubation at 35-37°C for 24-48hours(under anaerobic condition).

Organism (ATCC)	Inoculum (CFU)	Growth	Recovery	Colour of colony
<i>Clostridium perfringens</i> ATCC 13124 (00007*)	50-100	luxuriant	>=50%	Pale yellowish green
<i>Clostridium sporogenes</i> ATCC 11437	50-100	luxuriant	>=50%	Pale green-bluish green
<i>Clostridium sporogenes</i> ATCC 19404 (00008*)	50-100	luxuriant	>=50%	Pale green-bluish green
<i>Escherichia coli</i> ATCC 25922 (00013*)	>=10 ³	inhibited	0%	
<i>Staphylococcus aureus</i> ATCC 25923 (00034*)	>=10 ³	inhibited	0%	

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (3, 4).

References

1. Alcamo E. I., 2001, Fundamentals of Microbiology, 6th Ed., Jones and Bartlett Publishers.
2. Vera, 1962, Presented Pa. Soc. Med. Tech., York, Pa.
3. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
4. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

For isolation and differentiation between various species of *Lactobacillus*

HiCrome™ Lactobacillus Selective Agar Base

Recommended for isolation and differentiation between various species of *Lactobacillus* from a mixed culture by chromogenic method

M2065

Composition **

Ingredients	Grams/Litre
Peptone	10.00
HM Extract #	1.00
M-Protein powder ##	5.00
D-Mannitol	10.00
Sodium chloride	10.00
Chromogenic mixture	3.20
Phenol red	0.025
Agar	15.00

Final pH (at 25°C) 7.1 ± 0.2

** Formula adjusted, standardized to suit performance parameters

Equivalent to Meat Extract

Equivalent to Milk Protein

Directions

Suspend 54.22 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Aseptically add rehydrated contents of 1 vial of Ciprofloxacin Supplement (FD345). Mix well and pour into sterile Petri plates.

Principle and Interpretation

Lactobacillus is a genus of Gram-positive, facultative anaerobic or microaerophilic, rod-shaped, non-spore-forming bacteria. They are a major part of the lactic acid bacteria group. As more LABs have been developed and sold in mixed forms as probiotics, it is necessary to develop a method for counting each LAB in a mixture(1)

The medium contains peptone and HM extract, which provide nitrogenous and carbonaceous compounds, long chain amino acids, vitamins and other essential nutrients. Mannitol serves as the fermentable carbohydrate, fermentation of which can be detected by phenol red. M-protein aids in detecting casein hydrolysis activity. The chromogenic mixture present in the medium is cleaved by the enzyme β -glucosidase resulting in greenish blue to blue coloured colonies.

For selective isolation of *Lactobacillus*, Ciprofloxacin supplement is added (FD345) which inhibits the accompanying bacteria.

Type of specimen

Dairy samples: milk and milk products

Specimen Collection and Handling

For dairy samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (2, 3). After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions


Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets

Limitations

1. Some species may show poor growth due to nutritional variations.
2. Slight colour variation may be observed depending upon the utilization of the substrate by the organism.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature..



M2065 HiCrome™ Lactobacillus Selective Agar Base

HiCrome™ Lactobacillus Selective Agar Base

Recommended for isolation and differentiation between various species of *Lactobacillus* from a mixed culture by chromogenic method

M2065

Quality Control

- Appearance of Powder** : Light yellow to pink homogeneous free flowing powder.
- Gelling** : Firm, comparable with 1.5% Agar gel
- Colour and Clarity** : Red coloured, clear to slightly opalescent gel forms in Petri plates
- Reaction** : Reaction of 5.42% w/v aqueous solution at 25°C. pH : 7.1±0.2
- Cultural Response** : Cultural characteristics observed with addition of Ciprofloxacin supplement (FD345) after an incubation at 25-30°C for 24-48 hours (with 5% CO₂).

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (4, 5).

References

1. De Man, J.C., Rogosa, M. and Sharpe, E.M. (1960) A medium for the cultivation of lactobacilli. *J Appl Bacteriol* 23, 30-35.
2. American Public Health Association, Standard Methods for the Examination of Dairy Products, 1978, 14th Ed., Washington, D.C.
3. Wehr H. M. and Frank J. H., 2004, Standard Methods for the Microbiological Examination of Dairy Products, 17th Ed., APHA Inc., Washington, D.C.
4. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
5. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S. and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

Organism (ATCC)	Inoculum (CFU)	Growth	Recovery	Colour of colony
<i>Lactobacillus acidophilus</i> ATCC 4356 (WDCM 00098)	50-100	good - luxuriant	>=50%	Pale pink - pink
<i>Lactobacillus rhamnosus</i> ATCC 9595	50-100	good	>=50%	Light green
<i>Lactobacillus fermentum</i> ATCC 9338	50-100	good - luxuriant	>=50%	Yellow
<i>Lactobacillus plantarum</i> ATCC 8014 (WDCM 00016)	50-100	good - luxuriant	>=50%	Light green-green colonies w/ hazy background
<i>Lactococcus lactis subsp. lactis</i> ATCC 19435	50-100	good - luxuriant	>=50%	Light green-green colonies w/ hazy background
<i>Bacillus spizizenii subsp. spizizenii</i> ATCC 6633 (WDCM 00003)	>=10 ³	inhibited	0%	
<i>Bacillus cereus</i> ATCC 10876	>=10 ³	inhibited	0%	
<i>Staphylococcus aureus subsp. aureus</i> ATCC 6538 (WDCM 00032)	>=10 ³	inhibited	0%	

For Identification of *Enterococcus* species

HiCrome™ Enterococci Agar / Broth

Recommended for the identification and differentiation of Enterococci from water samples.

M1414 /
M1376

Composition **

	M1414	M1376
Ingredients	Grams/Litre	Grams/Litre
Peptone, special	10.00	10.00
Sodium chloride	5.00	5.00
Sodium azide	0.30	0.30
Chromogenic substrate	0.06	0.040
Polysorbate 80 (Tween 80)	2.00	2.00
Disodium hydrogen phosphate	1.25	1.25
Agar	15.00	-

Final pH (at 25°C) 7.5 ± 0.2

** Formula adjusted, standardized to suit performance parameters

Directions

Suspend 33.61 gm of M1414 and 18.59 grams (single strength) or 37.18 grams (double strength) of M1376 in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Mix well and pour/dispense into sterile Petri plates / tubes.

Principle and Interpretation

HiCrome™ Enterococci media are formulated on the basis of the work carried out by Althous et al (1), Amoras (2), Litsky et al (3), and Manafi and Sommer (4) and Snyder and Lichstein (5). These media are recommended for the rapid detection of Enterococci from water samples. The presence of *Enterococcus* group, which is a subgroup of the faecal streptococci, serves as a valuable bacterial indicator for determining the extent of faecal contamination (1, 6) and it is more specific than the detection of coliforms, which may originate from non-faecal sources. The enzyme β -glucosidase produced by Enterococci cleaves the chromogenic substrate, resulting in a bluish green colour. The medium contains peptone special, which provides nitrogenous and carbonaceous compounds, long chain amino acids, vitamins and other essential nutrients. Sodium chloride maintains the osmotic balance of the medium. Sodium azide inhibits the accompanying microflora, especially gram-negative organisms. Polysorbate 80 (Tween 80) acts as a source of fatty acids.

Type of specimen

Water samples

Specimen Collection and Handling

For water samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (7). After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions


Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets

Limitations

1. Some species may show poor growth due to nutritional variations.
2. Slight colour variation may be observed depending upon the utilization of the substrate by the organism.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.



M1414 – HiCrome™ Enterococci Agar

HiCromeVeg™ Freedom from BSE / TSE worries
Single Streak Rapid Differentiation Series

HiCrome™ Enterococci Broth (M1376) is also available as HiCrome™ Enterococci HiVeg™ Broth (MV1376) wherein all the animal origin nutrients have been replaced by vegetable based nutrients

HiCrome™ Enterococci Agar / Broth

Recommended for the identification and differentiation of Enterococci from water samples.

M1414 /
M1376

Quality Control

- Appearance of Powder** : Cream to yellow coloured, homogeneous, free flowing powder.
- Gelling** : Firm, comparable with 1.5% Agar gel of M1414.
- Colour and Clarity of prepared medium** : Light amber coloured, clear to slightly opalescent gel forms in Petri plates / clear solution in tubes.
- Reaction** : Reaction of 3.36% w/v of M1414 and 1.86% w/v of M1376 aqueous solution at 25°C. pH:7.5 ± 0.2.
- Cultural Response** : Cultural characteristics observed after an incubation at 35-37°C for 18-24 hours.

Organism (ATCC)	Inoculum (CFU)	Growth	Recovery on M1414	Colour of colony (M1414)/	Colour of Medium (M1376)
<i>Enterococcus faecalis</i> (29212) (00087*)	50-100	good	40-50%	blue green	light blue-green
<i>Staphylococcus aureus subsp aureus</i> (25923) (00034*)	50-100	good	40-50%	colourless	light yellow
<i>Escherichia coli</i> (25922) (00013*)	50-100	none to poor	≤10%	-	light yellow
<i>Pseudomonas aeruginosa</i> (27853) (00025*)	50-100	none to poor	≤10%	-	light yellow

Key : * = corresponding WDCM Numbers

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (8, 9).

References

1. Althous, H., Dott, W., Havemeister, G, Muller, H.E, a. Sacre',C., 1982, Zbl. Bakt. Hyg. I. Abt. Orig. A. 252:154-165.
2. Amoras I, 1995, Poster Präsentation congress of Spanish Society of Microbiology, Madrid.
3. Litsky, W., Mallmann, W.L., a Fifield, C.W. 1953, Amer. J. Pbl. Hlth. 43:873-879.
4. Manafi M., a. Sommer R, 1993, Wat. Sci. Tech. 27:271-274.
5. Snyder M.L., and Lichstein, H.C. 1940, J. Infect. Dis. 67. 113-115.
6. Standard Methods for the Examination of Water and Wastewater, 20th Edition, Edited by L.S.Clesceri, A.E., Greenberg and A.D. Eaton, Published by APHA, AWWA and WEF (1998).
7. Baird R.B., Eaton A.D., and Rice E.W., (Eds.), 2015, Standard Methods for the Examination of Water and Wastewater, 23rd ed., APHA, Washington, D.C.
8. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
9. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

HiCrome™ *Enterococcus faecium* Agar Base

Recommended for the identification of *Enterococcus faecium* from water, faeces and sewage samples.

M1580

Composition **

Ingredients	Grams/Litre
Peptone, special	23.00
Corn starch	1.00
Sodium chloride	5.00
Chromogenic substrate	0.10
Arabinose	10.00
Phenol red	0.10
Agar	15.00

Final pH (at 25°C) 7.8 ± 0.2

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 27.1 grams in 500 ml distilled water. Heat to boiling to dissolve the medium completely. DO NOT AUTOCLAVE. Cool to 45-50°C and aseptically add the rehydrated contents of 1 vial of *Enterococcus faecium* Selective Supplement (FD226). Mix well and pour into sterile Petri plates.

Principle and Interpretation

HiCrome™ *Enterococcus faecium* Agar Base is recommended for the chromogenic detection of *Enterococcus faecium* from urine, faeces, soil, food, water, plants and animals. *E. faecium* is commonly found in the gastrointestinal tracts of humans (1). The resistance exhibited by *Enterococcus* species to various antimicrobials has led them to being a major cause of human infections including nosocomial infections (2). *E. faecalis* causes 80-90% of infection while *E. faecium* causes the majority of the remainder (3). The use of selective media for the isolation of Enterococci has been previously reviewed, including those containing chromogenic substrates (4) and media containing cephalixin-aztreonam supplements. Enterococcus species possess the enzyme glucosidase, which specifically cleaves the chromogenic substrate to produce blue coloured colonies. *E. faecium* ferment arabinose; and cleaves the chromogenic substrate present in the media to produce green coloured colonies along with yellow colouration to the medium. *E. faecalis* does not ferment arabinose and therefore retains the blue colour.

Peptone special serves as a source of carbon, nitrogen and essential growth nutrients. Corn starch neutralizes the toxic metabolites while sodium chloride maintains the osmotic equilibrium. Phenol red serves as a pH indicator with arabinose being the fermentable carbohydrate.

Type of specimen

Clinical samples - Urine, faeces, Food samples, Water samples.

Specimen Collection and Handling

For clinical samples follow appropriate techniques for handling specimens as per established guidelines (5, 6).

For food samples, follow appropriate techniques for sample collection and processing as per guidelines (7).

For water samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards.(8)


After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

In Vitro diagnostic use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets

Limitations

1. Some species may show poor growth due to nutritional variations.
2. Slight colour variation may be observed depending upon the utilization of the substrate by the organism.



M1580 HiCrome™ *Enterococcus faecium* Agar Base (Mixture)

HiCromeVeg™ Freedom from BSE / TSE worries
Single Streak Rapid Differentiation Series

HiCrome™ *Enterococcus faecium* Agar Base (M1580) is also available as HiCrome™ *Enterococcus faecium* HiVeg™ Agar Base (MV1580) wherein all the animal origin nutrients have been replaced by vegetable based nutrients.

For Identification of *Enterococcus* species

HiCrome™ *Enterococcus faecium* Agar Base

Recommended for the identification of *Enterococcus faecium* from water, faeces and sewage samples.

M1580

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

- Appearance of Powder** : Light yellow to pinkish beige coloured, homogeneous, free flowing powder.
- Gelling** : Firm, comparable with 1.5% Agar gel.
- Colour and Clarity of prepared medium** : Red coloured, clear to slightly opalescent gel forms in Petri plates.
- Reaction** : Reaction of 5.42% w/v aqueous solution at 25°C. pH:7.8 ± 0.2 .
- Cultural Response** : Cultural characteristics observed (with added *Enterococcus faecium* Selective Supplement (FD226) after an incubation at 35-37°C for 24-48 hours.

Organisms (ATCC)	Inoculum (CFU)	Growth	Recovery	Colour of colony
<i>Enterococcus faecium</i> (19434) (00010*)	50-100	luxuriant	≥50%	green
<i>Enterococcus faecalis</i> (29212) (00087*)	50-100	luxuriant	≥50%	blue
<i>Enterococcus hirae</i> (10541)	50-100	luxuriant	≥50%	blue
<i>Escherichia coli</i> (25922) (00013*)	≥10 ³	inhibited	0%	-
<i>Pseudomonas aeruginosa</i> (27853) (00025*)	≥10 ³	inhibited	0%	-
<i>Staphylococcus aureus</i> subsp aureus (25923) (00034*)	≥10 ³	inhibited	0%	-

Key : * = corresponding WDCM Numbers

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.


Product performance is best if used within stated expiry period.

Disposal


User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (5, 6).

References

1. Chenoweth C., Schaberg D., The Epidemiology of Enterococci, Eur. J. Clin. Microbiol. Infect. Dis., 9:80-89, 1990.
2. Moellering R. C., 1992, Clin. Infect. Dis. 14: 1173.
3. Skinner F. A. and Quesnel L. B., (Ed.), 1978, Streptococci. Academic Press, Inc. (London) Ltd., London, United Kingdom, p. 245 261
4. Willinger B. and Manafi M., 1995, Lett. Appl. Microbiol., 20: 300-302.
5. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
6. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S. and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.
7. Salfinger Y., and Tortorello M.L. Fifth (Ed.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.
8. Baird R.B., Eaton A.D., and Rice E.W., (Eds.), 2015, Standard Methods for the Examination of Water and Wastewater, 23rd ed., APHA, Washington, D.C.



M1580 HiCrome™ *Enterococcus faecium* Agar Base



M1580 HiCrome™ *Enterococcus faecium* Agar Base

For the detection of Group B Streptococci

HiCrome™ Strep B Selective Agar Base / Modified

Recommended for selective isolation of Group B streptococci from clinical samples.

M1840/
M1966

Quality Control

- Appearance of Powder** : Cream to yellow (M1840), light yellow to pink (M1966) homogeneous free flowing powder
- Gelling** : Firm, comparable with 1.5% Agar gel
- Colour and Clarity of prepared medium** : Yellow coloured (M1840) opaque gel forms in Petri plates or red coloured (M1966) clear to slightly opalescent gel forms in Petri plates.
- Reaction** : Reaction of 3.77% w/v (M1840) and 3.68% w/v (M1966) aqueous solution at 25°C.
pH : 7.3 ± 0.2.
- Cultural Response** : Cultural characteristics observed with added HiCrome™ Strep B Selective Supplement (FD273), after an incubation at 35-37°C for 18 - 24 hours.

Organism (ATCC)	Inoculum (CFU)	Growth	Recovery	Colour of Colony (M1840)	Colour of Colony (M1966)
<i>Streptococcus agalactiae</i> (13813)	50-100	luxuriant	≥50%	blue	purple
<i>Escherichia coli</i> (25922) (00013*)	≥10 ³	inhibited	0%	—	—
<i>Neisseria meningitidis</i> (13090)	≥10 ³	inhibited	0%	—	—
<i>Staphylococcus aureus</i> subsp <i>aureus</i> (25923) (00034*)	≥10 ³	inhibited	0%	—	—
<i>Enterococcus faecalis</i> ATCC 29212 (00087*)	50-100	luxuriant	≥50%	—	bluish green
<i>Enterococcus faecium</i> ATCC 19434 (00010*)	50-100	luxuriant	≥50%	—	green w/ yellow background

Key: * = corresponding WDCM Numbers

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.


Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (5, 6).

References

1. Anthony BF, Okada DM, Hobel CJ. Epidemiology of group B Streptococcus: longitudinal observations during pregnancy. J.Infect Dis 1978; 137:524-30.
2. Murray P.R., Baron J.H., Manual of Clinical Microbiology Murray P. R., Baron J. H. Pfaller M. A., Jorgensen J. H. and Tenover F. C., (Eds.), 2003, Manual of Clinical Microbiology, 8th Ed., American Society for Microbiology, Washington, D.C.
3. NHS Processing swabs for Group B Streptococcal carriage Issue no.2.1,2006.
4. Prevention of perinatal group B Streptococcal disease: a public health perspective. Centres for Disease control and Prevention. MMWR Recomm Rep 1996; 51:1-22.
5. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
6. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S. and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.



M1966 HiCrome™ Strep B Selective Agar Base / Modified
Streptococcus agalactiae ATCC 13813

For the detection of Vancomycin Resistant Enterococci

HiCrome™ VRE Agar Base / Modified

Recommended for identification of Vancomycin Resistant Enterococci from clinical specimens

M1830/
M1925

Limitations

1. Some Intermediate species may show poor growth due to nutritional variations and tolerance to vancomycin.
2. Slight colour variation may be observed depending upon the utilization of the substrate by the organism.
3. Further confirmation has to be carried using sensitivity testing.
4. For M1830, interspecies differentiation between *Enterococcus faecalis* and *Enterococcus faecium* cannot be confirmed.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

- Appearance of Powder** : Cream to yellow (M1830) or Light yellow to pink (M1925) homogeneous free flowing powder
- Gelling** : Firm, comparable with 1.5% Agar gel.
- Colour and Clarity of prepared medium** : Off white coloured (M1830) or red coloured (M1925) opaque gel forms in Petri plates.
- Reaction** : Reaction of 5.1% w/v (M1830) or 5.37% w/v (M1925) aqueous solution at 25°C.
pH : 6.5 ± 0.2 (M1830) pH : 7.8 ± 0.2 (M1925)
- Cultural Response** : Cultural characteristics observed with added HiCrome™ VRE Agar Supplement (FD277) after an incubation at 35-37°C for 24-48 hours.

Organism (ATCC)	Inoculum (CFU)	Growth	Recovery	Colour of colony (M1830)	Colour of colony (M1925)
<i>Enterococcus faecalis</i> (VRE) (51299) (00085*) (00152*)	50-100	luxuriant	≥50%	bluish green	blue
<i>Enterococcus faecium</i> (VRE) (700221)	50-100	luxuriant	≥50%	—	green w/ yellow background
<i>Enterococcus faecalis</i> (29212) (00087*)	≥10 ³	inhibited	0%	—	—

<i>Staphylococcus aureus</i> subsp <i>aureus</i> (25923) (00034*)	≥10 ³	inhibited	0%	—	—
---	------------------	-----------	----	---	---

Key : * = corresponding WDCM Numbers

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.


Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (3, 4).

References

1. Mara D., Horan NJ: The Handbook of water, wastewater and microbiology, Amsterdam, The Netherlands, Academic Press; 2003.
2. Mascini EM, Bonten MJ: Vancomycin- resistant enterococci: consequences for therapy and infection control. Clin Microbiol Infect.2005,11 (Suppl.4) :43-56
3. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
4. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.



M1830 HiCrome™ VRE Agar Base / Modified

For Identification of *Listeria* species

HiCrome™ Listeria Agar Base / Modified

Recommended as a selective and differential agar medium recommended for rapid and direct identification of *Listeria* species.

M1417F/
M1417

Limitations

1. Some species may show poor growth due to nutritional variations.
2. Slight colour variation may be observed depending upon the utilization of the substrate by the organism.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

Appearance of powder : Light yellow to pink coloured, homogeneous, free flowing powder.

Gelling : Firm, comparable with 1.3% Agar gel.

Colour and Clarity of prepared medium : Red coloured, clear to slightly opalescent gel forms in Petri plates.

Reaction : Reaction of 6.72% w/v aqueous solution of M1417 at 25°C. pH : 7.3 ± 0.2.
Reaction of 7.32% w/v aqueous solution of M1417F at 25°C. pH : 7.3 ± 0.1.

Cultural Response : Cultural characteristics observed with added HiCrome™ Listeria Selective Supplement (FD181) after incubation at 35-37°C for 24-48 hours.

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (5, 6).


References

1. FDA U.S., Bacteriological Analytical Manual 8 ed. Gaithersburg, MD, AOAC international, 1998.
2. Mengaud J., Braun-Breton C. and Cossart P., (1991), Molecular Microbiology, (2): 367-372.
3. Notermans S.H. and Dufrenne J., (1991), Applied and Environmental Microbiology, 57(09):2666-70.
4. Salfinger Y., and Tortorello M.L. Fifth (Ed.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.
5. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
6. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

Organisms (ATCC)	Inoculum (CFU)	Growth	Recovery	Colour of colony	(M1417) Rhamnose fermentation	(M1417F) Xylose fermentation
<i>Listeria monocytogenes</i> (19118)	50-100	luxuriant	≥50%	bluish-green	+	-
<i>Listeria ivanovii</i> (19119) (00018*)	50-100	luxuriant	≥50%	bluish-green	-	+
<i>Listeria innocua</i> (33090) (00017*)	50-100	luxuriant	≥50%	bluish-green	+	-
<i>Escherichia coli</i> (25922) (00013*)	≥10 ³	inhibited	0%			-
<i>Bacillus spizizenii</i> sub <i>spizizenii</i> (6633) (00003*)	≥10 ³	inhibited	0%			-
<i>Pseudomonas aeruginosa</i> (27853) (00025*)	≥10 ³	inhibited	0%			-
<i>Candida albicans</i> (10231) (00054*)	≥10 ³	inhibited	0%			-

Key : + = positive reaction, - = negative reaction.

* = corresponding WDCM Numbers



M1417 HiCrome™ Listeria Agar Base

For Identification of *Listeria* species

L. Mono Differential Agar Base

Recommended for the selective and differential isolation of *Listeria monocytogenes* based on PIPLC activity.

M1540

Composition **

Ingredients	Grams/Litre
HM peptone#	18.00
Tryptone	6.00
Yeast extract	10.00
Sodium pyruvate	2.00
Glucose (Dextrose)	2.00
Magnesium glycerophosphate	1.00
Magnesium sulphate	0.50
Sodium chloride	5.00
Lithium chloride	10.00
Disodium hydrogen phosphate	2.50
Chromogenic substrate	0.05
Agar	15.00

Final pH (at 25°C) 7.2 ± 0.2

** Formula adjusted, standardized to suit performance parameters

Equivalent to Meat peptone

Directions

Suspend 36.02 grams in 460 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Aseptically add sterile rehydrated contents of 1 vial of L. mono Selective Supplement I (FD212) and L. mono Selective Supplement II (FD213). For enrichment add sterile content of L. mono Enrichment Supplement I (FD214). Mix well and pour into sterile Petri plates.

Principle and Interpretation

L. mono Differential Agar Base is based on the formulation of Ottoviani and Agosti (2, 3) for the selective and differential isolation of *Listeria monocytogenes* from food and animal feeds which is adopted by ISO Committee (1).

HM peptone, Tryptone, yeast extract and sodium pyruvate provide essential growth nutrients and nitrogenous, carbonaceous compounds, long chain amino acids and vitamin B complex. Glucose is the fermentable carbohydrate. Sodium chloride maintains osmotic equilibrium. Phosphate buffers the medium. Lithium chloride and added selective supplements (FD212 and FD213) inhibit accompanying microflora and allow the growth of *Listeria* species. *Listeria* species hydrolyse the chromogenic substrate which produces greenish-blue

coloured colonies. Differentiation of *Listeria monocytogenes* from other *Listeria* species is based on phosphatidylinositol-specific phospholipase C (PIPLC) activity. Phospholipase C enzyme hydrolyses the purified substrate (FD214) added to the medium resulting in an opaque halo around *Listeria monocytogenes* colonies.

Type of specimen

Food samples

Specimen Collection and Handling

For Food samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (4).


After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

Read the label before opening the container. Wear protective gloves/ protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets

Limitations

1. Due to variable nutritional requirements, some strains may show poor growth on this medium.
2. Further biochemical tests must be carried out to differentiate between *L. monocytogenes* and *L. ivannovi*.



M1540 L. Mono Differential Agar Base

HiCromeVeg™ Freedom from BSE / TSE worries
Single Streak Rapid Differentiation Series

L. Mono Differential Agar Base(M1540) is also available as L. Mono Differential HiVeg™ Agar Base (MV1540) & L. Mono Differential HiCynth™ Agar Base (MCD1540) wherein all the animal origin nutrients have been replaced by vegetable based nutrients and chemically defined peptones respectively.

L. Mono Differential Agar Base

Recommended for the selective and differential isolation of *Listeria monocytogenes* based on PIPLC activity.

M1540

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

- Appearance of Powder** : Cream to yellow coloured, homogenous free flowing powder.
- Gelling** : Firm, comparable with 1.5% Agar gel.
- Colour and Clarity of prepared medium** : Light amber coloured, opalescent gel forms in Petri plates.
- Reaction** : Reaction of 7.2% w/v aqueous solution at 25°C.
- pH** : 7.2 ± 0.2
- Cultural Response** : Cultural characteristics observed with added L. mono Selective Supplement I (FD212), L. mono Selective Supplement II (FD213) and L. mono Enrichment Supplement I (FD214) after an incubation at 35-37°C for 24-48 hours.

Organisms (ATCC)	Inoculum (CFU)	Growth	Recovery	Colour of the Colony	PIPLC activity #
<i>Candida albicans</i> (10231) (00054*)	≥10 ³	inhibited	0%	-	-
<i>Enterococcus faecalis</i> (29212) (00087*)	≥10 ³	inhibited	0%	-	-
<i>Escherichia coli</i> (25922) (00013*)	≥10 ³	inhibited	0%	-	-
<i>Listeria innocua</i> (33090) (00017*)	50-100	luxuriant	≥50%	greeish-blue	-
<i>Listeria grayi</i> (19120)	50-100	luxuriant	≥50%	greeish-blue	-
<i>Listeria ivanovii</i> (19119) (00018*)	50-100	luxuriant	≥50%	greeish-blue	+
<i>Listeria monocytogenes</i> (19112)	50-100	luxuriant	≥50%	greeish-blue	+
<i>Listeria seeligeri</i> (35967)	50-100	luxuriant	≥50%	greeish-blue	-
<i>Listeria welshimeri</i> (43549)	50-100	luxuriant	≥50%	greeish-blue	-
<i>Pseudomonas aeruginosa</i> (27853) (00025*)	≥10 ³	inhibited	0%	-	-

Key : * : Corresponds to WDCM number

PIPLC activity # : opaque halo around the colony exhibiting phosphatidylinositol - specific phospholipase activity.

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (5, 6).

References

1. Draft Amendment ISO 11290-2:1996/DAM 1.
2. Ottaviani F., Ottaviani M., and Agosti M. (1997 a), *Industrie Alimentari* 36, 1-3.
3. Ottaviani F., Ottaviani M., and Agosti M. (1997 b), *Quimper Froid Symposium Proceedings* p. 6, A.D.R.I.A. Quimper, France, 16-18 June 1997.
4. Salfinger Y., and Tortorello M.L. Fifth (Ed.), 2001, *Compendium of Methods for the Microbiological Examination of Foods*, 5th Ed., American Public Health Association, Washington, D.C.
5. Isenberg, H.D. *Clinical Microbiology Procedures Handbook*. 2nd Edition.
6. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock, D.W. (2015) *Manual of Clinical Microbiology*, 11th Edition. Vol. 1



M1540 L. Mono Differential Agar Base

HiCrome™ L. mono Rapid Differential Agar Base

Recommended for the rapid identification and differentiation of *Listeria monocytogenes* from other *Listeria* species based on rhamnose fermentation and PIPLC activity.

M1924

Composition **

Ingredients	Grams/Litre
Peptone special	23.00
Tryptone	10.00
Soya peptone	2.00
Sodium chloride	4.00
Lithium chloride	5.00
Chromogenic mixture	1.16
Rhamnose	10.00
Phenol red	0.12
Agar	15.00

Final pH (at 25°C) 7.4 ± 0.2

** Formula adjusted, standardized to suit performance parameters

Directions

Suspend 35.14 grams in 470 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Aseptically add sterile contents of 1 vial of L. mono Enrichment Supplement I (FD214) and sterile rehydrated contents of 1 vial of HiCrome™ Listeria Selective Supplement (FD181). Mix well and pour into sterile Petri plates.

Principle and Interpretation

Listeria monocytogenes is a gram-positive foodborne human pathogen responsible for serious infections in pregnant women that may ultimately result in abortion, stillbirth, birth of a child with neonatal listeriosis and meningitis or primary bacteremia in adults and juveniles. The pathogenicity of *Listeria ivanovii* for humans is uncertain (5). Since *L. monocytogenes* and *L. innocua* have similar biochemical properties, they cannot be differentiated on traditional media (PALCAM, Oxford). This medium is based on the specific chromogenic detection of β-glucosidase activity, rhamnose fermentation and PIPLC activity. *Listeria* species hydrolyse the purified chromogenic substrate in the medium giving blue coloured colonies. Since β-glucosidase activity is specific for *Listeria* species, other organisms cannot utilize the chromogenic substrate and therefore give white colonies. Differentiation between *Listeria* species is based on the property of rhamnose fermentation and PIPLC activity. The colonies of *L. monocytogenes* appear bluish green with a yellow halo (rhamnose positive) while the colonies of *L. ivanovii* appear bluish green without a yellow halo (Rhamnose negative) (1, 2). The differentiation of *L. mono* and *L. innocua* is based on PIPLC (phosphatidylinositol-specific phospholipase C) activity. Phospholipase C enzyme hydrolyses the purified substrate (FD214) added to the medium resulting in an opaque halo around *Listeria monocytogenes* colonies. *L. ivanovii* also demonstrates PIPLC activity however since it does not ferment

rhamnose it can be easily distinguished from *L. monocytogenes* (3, 4). Peptone special, tryptone and soya peptone provide nitrogenous compounds, carbon, long chain amino acids vitamin B complex and other essential growth nutrients. Rhamnose is the fermentable carbohydrate with phenol red as an indicator. Sodium chloride maintains the osmotic equilibrium. The added lithium chloride and HiCrome™ Listeria Selective Supplement (FD181) inhibit growth of most gram- positive bacteria, gram-negative bacteria, yeasts and moulds. Phospholipase C enzyme hydrolyses the purified substrate (FD214) added to the medium resulting in an opaque halo around *Listeria monocytogenes* colonies demonstrating PIPLC activity.

Type of specimen

Food samples


Specimen Collection and Handling

For Food samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (6).

After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

Read the label before opening the container. Wear protective gloves/ protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets



M1924 – HiCrome™ L. mono Rapid Differential Agar Base

HiCrome™ L. mono Rapid Differential Agar Base

Recommended for the rapid identification and differentiation of *Listeria monocytogenes* from other *Listeria* species based on rhamnose fermentation and PIPLC activity.

M1924

Limitations

1. Due to variable nutritional requirements, some strains may show poor growth on this medium.
2. Slight colour variation may be observed depending upon the utilization of the substrate by the organism.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

- Appearance of Powder** : Light yellow to pink homogeneous free flowing powder
- Gelling** : Firm, comparable with 1.5% Agar gel.
- Colour and Clarity of prepared medium** : Red coloured, opalescent gel forms in Petri plates
- Reaction** : Reaction of 7.03% w/v aqueous solution at 25°C. pH : 7.4±0.2.
- Cultural Response** : Cultural characteristics observed w/added HiCrome™ *Listeria* Selective Supplement (FD181) and L.mono Enrichment supplement I (FD214), after an incubation at 35-37°C for 24-48 hours.

Organism (ATCC)	Inoculum (CFU)	Growth	Recovery	Colour of colony	Rhamnose fermentation	PIPLC Activity
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> (6633) (00003*)	≥10 ³	inhibited	0%			
<i>Candida albicans</i> (10231) (00054*)	≥10 ³	inhibited	0%			
<i>Escherichia coli</i> (25922) (00013*)	≥10 ³	inhibited	0%			
<i>Listeria innocua</i> (33090) (00017*)	50-100	luxuriant	≥50%	bluish green (yellow background)	positive reaction,	negative reaction
<i>Listeria ivanovii</i> (19119) (00018*)	50-100	luxuriant	≥50%	bluish green	negative reaction	**positive
<i>Listeria monocytogenes</i> (19118)	50-100	luxuriant	≥50%	bluish green (yellow background)	positive reaction,	**positive,
<i>Pseudomonas aeruginosa</i> (27853) (00025*)	≥10 ³	inhibited	0%			

Key : * : Corresponds to WDCM number

** : opaque halo around the colony exhibiting phosphatidyl inositol specific phospholipase activity.

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (7, 8).

References

1. Mengaud J., Braun-Breton C. and Cossart P., (1991), *Molecular Microbiology*, 5(2): 367-372.P
2. Notermans S.H. and Dufrenne J., (1991), *Applied and Environmental Microbiology*, 57(09): 2666-70.
3. Ottaviani F., Ottaviani M., and Agosti M. (1997 a), *Industrie Alimentari* 36, 1-3.
4. Ottaviani F., Ottaviani M., and Agosti M. (1997 b), *Quimper Froid Symposium Proceedings* p. 6, A.D.R.I.A. Quimper, France, 16-18 June 1997.
5. Schlech WF, Lavigne PM, Bortolussi RA, et al. (January 1983). "Epidemic listeriosis-evidence for transmission by food". *N. Engl. J. Med.* 308(4): 203-6. doi:10.1056.
6. Salfinger Y., and Tortorello M.L. Fifth (Ed.), 2001, *Compendium of Methods for the Microbiological Examination of Foods*, 5th Ed., American Public Health Association, Washington, D.C.
7. Isenberg, H.D. *Clinical Microbiology Procedures Handbook*. 2nd Edition.
8. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock, D.W. (2015) *Manual of Clinical Microbiology*, 11th Edition. Vol. 1

For identification of *Listeria* species

HiCrome™ L. mono Differential Agar Base

Recommended for the selective and differential isolation, enumeration and identification of *Listeria monocytogenes* and *Listeria* species based on PCPLC activity.

M2009

Composition **

Ingredients	Grams/Litre
Peptone	15.00
Tryptone	6.00
Yeast extract	10.00
Sodium pyruvate	2.00
Maltose	4.00
Magnesium glycerophosphate	1.00
Magnesium sulphate	0.50
Sodium chloride	5.00
Lithium chloride	5.00
Disodium hydrogen phosphate	2.50
Chromogenic substrate	2.20
Agar	14.00

Final pH (at 25°C) 7.2±0.2

** Formula adjusted, standardized to suit performance parameters

Directions

Suspend 33.60 grams in 480 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Aseptically add sterile contents of 1 vial of Lecithin solution (FD332) and sterile rehydrated contents of Modified L.mono Selective Supplement (FD333). Mix well and pour into sterile Petri plates.

Principle and Interpretation

Listeria monocytogenes is a gram-positive foodborne human pathogen responsible for serious infections in pregnant women that may ultimately result in abortion, stillbirth, birth of a child with neonatal listeriosis and meningitis or primary bacteremia in adults and juveniles (2). The pathogenicity of *Listeria ivanovii* for humans is uncertain. Since *L.monocytogenes* and *L.innocua* have similar biochemical properties, they cannot be differentiated on traditional media (PALCAM, Oxford). HiCrome *L.mono* Differential Agar Base is based on for the selective and differential isolation of *Listeria* species on the basis of utilization of chromogenic substrate and lecithinase activity [Phosphotidylcholine hospholipase C (PCPLC)] (3). PI-PLC and PC-PLC, the major virulence factors, are only produced by pathogenic *L. monocytogenes* and *Listeria ivanovii* (1)

Peptone, tryptone, yeast extract and sodium pyruvate provide nitrogenous and carbonaceous compounds, long chain amino acids, vitamins and essential growth nutrients . Maltose is the

fermentable carbohydrate. Sodium chloride maintains osmotic equilibrium. Phosphate buffers the medium. Lithium chloride and added selective supplement (FD333) inhibit accompanying microflora and allow the growth of *Listeria* species. *Listeria* species hydrolyse the chromogenic substrate and produces green coloured colonies. Lecithin solution (FD332) helps in detecting PCPLC activity. Differentiation of *Listeria* species is based on phosphatidylcholine phospholipase C (PCPLC) activity. *L. monocytogenes* and *L.ivanovii* exhibits PCPLC activity which is seen as opaque halo around the colony.

Type of specimen


Food samples

Specimen Collection and Handling

For Food samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (4). After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

Read the label before opening the container. Wear protective gloves/ protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets



M2009 - HiCrome L. mono Differential Agar Base

HiCrome™ L. mono Differential Agar Base

Recommended for the selective and differential isolation, enumeration and identification of *Listeria monocytogenes* and *Listeria species* based on PCPLC activity.

M2009

Limitations

1. Due to variable nutritional requirements, some strains may show poor growth on this medium.
2. Slight colour variation may be observed depending upon the utilization of the substrate by the organism.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

- Appearance of Powder** : Cream to yellow homogeneous free flowing powder
- Gelling** : Firm, comparable with 1.4% Agar gel
- Colour and Clarity of prepared medium** : Light amber coloured, opalescent gel forms in Petri plates.
- Reaction** : Reaction of 6.72% w/v aqueous solution at 25°C. pH : 7.2±0.2
- pH** : 7.00-7.40
- Cultural Response** : Cultural characteristics observed with added sterile Modified L.mono Selective Supplement (FD333) and Lecithin solution (FD332) after an incubation at 35 - 37°C for 24 - 48 hours.

Organism (ATCC)	Inoculum (CFU)	Growth	Recovery	Colour of colony	PIPLC Activity#
<i>Enterococcus faecalis</i> (29212) (00087*)	≥10 ³	inhibited	0%	-	-
<i>Listeria innocua</i> (33090) (00017*)	50-100	luxuriant	≥50%	greeish-blue	negative
<i>Listeria ivanovii</i> (19119) (00018*)	50-100	luxuriant	≥50%	greeish-blue	positive#
<i>Listeria monocytogenes</i>	50-100	luxuriant	≥50%	greeish-blue	positive#

Key: # : opaque halo around the colony exhibiting phophatidylcholine phospholipase activity

* : corresponds to WDCM numbers

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (5, 6).

References

1. Mengaud J, Braun-Breton C, Cossart P 1991. Identification of phosphatidylinositol-specific phospholipase C activity in *Listeria monocytogenes* : a novel type of virulence factor. *Mol. Microbiol.* 5:367-372. doi:10.1111/j.1365-2958.1991
2. Painter J, Slutsker L. 2007. Listeriosis in humans, p 85-109. In Ryser ET, Marth EH (ed), *Listeria, listeriosis, and food safety*. Marcel Dekker, New York, NY.
3. Sang-Hyun Park, Pahn-Shick Chang, Sangryeol Ryu and Dong-Hyun Kang. Development of a Novel Selective and Differential Medium for the Isolation of *Listeria monocytogenes*. *Applied and Environmental Microbiology* 2014.
4. Salfinger Y., and Tortorello M.L. Fifth (Ed.), 2001, *Compendium of Methods for the Microbiological Examination of Foods*, 5th Ed., American Public Health Association, Washington, D.C.
5. Isenberg, H.D. *Clinical Microbiology Procedures Handbook*. 2nd Edition.
6. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) *Manual of Clinical Microbiology*, 11th Edition. Vol. 1

HiCrome™ Aureus Agar Base

Recommended for isolation and identification of *Staphylococci* from environmental samples.

M1468

Composition **

Ingredients	Grams/Litre
Tryptone	12.00
Gelatin peptone	3.00
HM peptone B#	6.00
Yeast extract	5.00
Sodium pyruvate	10.00
Lithium chloride	5.00
Chromogenic mixture	2.10
Agar	20.00

Final pH (at 25°C) 7.0 ± 0.2

** Formula adjusted, standardized to suit performance parameters

#Equivalent to Beef extract

Directions

Suspend 63.1 grams in 950 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C and aseptically add 50 ml concentrated Egg Yolk Tellurite Emulsion (FD046). Mix well and pour into sterile Petri plates.

Principle and Interpretation

HiCrome™ Aureus Agar Base is recommended for isolation and enumeration of coagulase positive *Staphylococcus aureus* from environment samples. Coagulase positive *S. aureus* gives brown black colonies with clear zone around the colony whereas *S. epidermidis* gives slightly brownish colonies. Other organisms give either colourless colonies or bluish coloured colonies due to the presence of chromogen. *Listeria monocytogenes* colonies are bluish in colour whereas *Bacillus*, *E. coli* and *Micrococcus* give colourless colonies.

Tryptone, gelatin peptone, HM peptone B and yeast extract provide nitrogenous, carbonaceous compounds, long chain amino acids, vitamin B complex and other essential growth nutrients. Sodium pyruvate protects injured cells, helps recovery and enhances growth of *Staphylococcus*. Lithium chloride and potassium tellurite inhibit most of the contaminating microflora except *Staphylococcus aureus* (1). Due to addition of egg yolk, proteolytic bacteria produce a clear zone around colony (1).

Type of specimen

Clinical samples, Environmental samples

Specimen Collection and Handling

For clinical samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (2, 3). After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions


In Vitro diagnostic use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets

Limitations

1. Due to variable nutritional requirements, some strains may show poor growth on this medium.
2. Slight colour variation may be observed depending upon the utilization of the substrate by the organism.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.



M1468 HiCrome™ Aureus Agar Base

HiCrome™ Aureus Agar Base

Recommended for isolation and identification of *Staphylococci* from environmental samples.

M1468

Quality Control

- Appearance of Powder** : Cream to yellow homogeneous free flowing powder
- Gelling** : Firm, comparable with 2.0 % Agar gel.
- Colour and Clarity of prepared medium** : Yellow coloured opaque gel forms in Petri plates.
- Reaction** : Reaction of 6.31% w/v aqueous solution at 25°C. pH:7.0±0.2
- Cultural Response** : Cultural characteristics observed with added Egg Yolk Tellurite Emulsion (FD046) after an incubation at 35-37°C for 24-48 hours.

Organism (ATCC)	Inoculum (CFU)	Growth	Recovery	Colour of colony	Lecithinase activity
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> (6633) (00003*)	50-100	none to poor	≤10%	colourless	Negative reaction
<i>Escherichia coli</i> (25922) (00013*)	50-100	none to poor	≤10%	colourless	Negative reaction
<i>Listeria monocytogenes</i> (19112)	50-100	fair - good	30-40%	bluish	Negative reaction
<i>Micrococcus luteus</i> (10240)	50-100	none to poor	≤10%	colourless	Negative reaction
<i>Staphylococcus aureus</i> subsp <i>aureus</i> (25923) (00034*)	50-100	good-luxuriant	≥50%	brown-black halo or clear zone around the colony	Positive reaction
<i>Staphylococcus epidermidis</i> (12228) (00036*)	50-100	none to poor	≤10%	yellow-slight brownish	Negative reaction

Key : * : Corresponds to WDCM number

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (2, 3).

References

1. Baird Parker, Ac (1962) J appl. Bact., 25:12.
2. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
3. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1

For Identification of *Staphylococcus aureus*.

HiCrome™ Staph Agar Base, Modified

Recommended as a selective medium for the isolation and enumeration of *Staphylococcus aureus*.

M1837

Composition **

Ingredients	Grams/Litre
Peptone special	23.000
Sodium pyruvate	4.000
Sodium chloride	40.000
Lithium chloride	5.000
Chromogenic mixture	5.300
Agar	15.000

Final pH (at 25°C) 7.2 ± 0.2

** Formula adjusted, standardized to suit performance parameters

Directions

Suspend 46.15 grams in 500 ml distilled water. Heat to boiling to dissolve the medium completely. DO NOT AUTOCLAVE. Cool to 45-50°C and aseptically add rehydrated contents of 1 vial of Polymyxin B Selective supplement (FD003). Mix well and pour into sterile Petri plates.

Principle and Interpretation

Staphylococci are widespread in nature, though they are mainly found living on the skin, skin glands and mucous membranes of mammals and birds. Humans and animals are the primary source of this organism. Because of its widespread nature it is easily transferred to food and a cause of food poisoning if not handled properly.(1)

The coagulase positive species *S. aureus* is well documented as a human opportunistic pathogen. *Staphylococcus* species are a major cause of food poisoning and produces a wide variety of enterotoxins, thus causing various types of disease symptoms. The ability to clot plasma continues to be the most widely used and accepted criterion for the identification of pathogenic staphylococci associated with acute infections (2).

This medium is a selective chromogenic medium recommended for the isolation and enumeration of coagulase positive staphylococci in foods within 24 hours. This medium has an advantage over the traditional media which requires 48 hours. Peptone special in the medium supplies the essential nitrogenous, carbonaceous compounds long chain aminoacids, vitamins and other essential growth nutrients required for the growth. The chromogenic mixture incorporated in the medium is specifically cleaved by *Staphylococcus aureus* to give bluish green coloured colonies which are clearly visible against the opaque background. Sodium pyruvate enhances the growth of *Staphylococcus* species. Sodium chloride in the medium helps to maintain the osmotic equilibrium of the medium. High concentration of sodium chloride also helps in inhibiting the accompanying microflora. Lithium chloride inhibits most of the contaminating microflora. Addition of Polymyxin B Sulphate (FD003) helps to restrict growth of gram-negative bacteria such as *Escherichia coli* and *Pseudomonas aeruginosa*.

Type of specimen

Clinical and food samples

Specimen Collection and Handling

For Clinical samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (3, 4).

For food and dairy samples, follow appropriate techniques for sample collection and processing as per guidelines (5)

After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions


In Vitro diagnostic use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets

Limitations

1. Due to variable nutritional requirements, some strains may show poor growth on this medium.
2. Slight colour variation may be observed depending upon the utilization of the substrate by the organism.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.



M1837 – HiCrome™ Staph Agar Base, Modified

For identification of *Staphylococcus aureus*

HiCrome™ Staph Selective Agar

Recommended as a selective medium for the isolation and enumeration of *Staphylococcus aureus*.

M1931

Composition **

Ingredients	Grams/Litre
Peptone special	25.000
Sodium chloride	50.000
Chromogenic mixture	3.200
Selective mixture	2.800
D-Mannitol	10.000
Phenol red	0.025
Agar	12.000

Final pH (at 25°C) 7.4 ± 0.2

** Formula adjusted, standardized to suit performance parameters

Directions

Suspend 103.03 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. DO NOT AUTOCLAVE. Cool to 45-50°C. Mix well and pour into sterile Petri plates.

Principle and Interpretation

Staphylococci are widespread in nature, though they are mainly found living on the skin, skin glands and mucous membranes of mammals and birds. Humans and animals are the primary source of this organism. Because of its widespread nature it is easily transferred to food and a cause of food poisoning if not handled properly.(1)

The coagulase positive species *S. aureus* is well documented as a human opportunistic pathogen. *Staphylococcus species* are a major cause of food poisoning and produces a wide variety of enterotoxins, thus causing various types of disease symptoms. The ability to clot plasma continues to be the most widely used and accepted criterion for the identification of pathogenic staphylococci associated with acute infections (2).

This medium is a selective chromogenic medium recommended for the isolation and enumeration of coagulase positive staphylococci in foods within 24 hours. Peptone special in the medium supplies the essential nitrogenous, carbonaceous required for the growth. Phenol red is pH indicator. Mannitol in the medium is fermented by *Staphylococcus aureus* and the chromogenic mixture incorporated in the medium is specifically cleaved by *Staphylococcus aureus* to give greenish coloured colonies which are easily distinguishable. *Staphylococcus epidermidis* does not ferment mannitol hence blue coloured colonies are observed. Sodium chloride in the medium helps to maintain the osmotic equilibrium of the medium. High concentration of sodium chloride also helps in inhibiting the accompanying microflora.

Type of specimen

Clinical and food samples

Specimen Collection and Handling

For Clinical samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (3, 4).

For food and dairy samples, follow appropriate techniques for sample collection and processing as per guidelines (5)

After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions


In Vitro diagnostic use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets

Limitations

1. Due to variable nutritional requirements, some strains may show poor growth on this medium.
2. Slight colour variation may be observed depending upon the utilization of the substrate by the organism.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.



M1931 HiCrome™ Staph Selective Agar

For Identification of Methicillin Resistant *Staphylococcus aureus* (MRSA)

HiCrome™ MeReSa Agar Base

For the isolation and selective identification of Methicillin Resistant *Staphylococcus aureus* (MRSA) from clinical isolates.

M1674

Composition **

Ingredients	Grams/Litre
Tryptone	13.00
Yeast extract	2.50
HM peptone B#	2.50
Sodium pyruvate	5.00
Sodium chloride	40.00
Chromogenic mixture	5.30
Agar	15.00

Final pH (at 25°C) 7.0 ± 0.2

** Formula adjusted, standardized to suit performance parameters

#Equivalent to Beef extract

Directions

Suspend 41.65 grams in 500 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Aseptically add sterile rehydrated contents of 1 vial of MeReSa Selective Supplement (FD229) and Cefoxitin supplement (FD259) for selectivity. Mix well and pour into sterile Petri plates.

Principle and Interpretation

Staphylococcus aureus is an invasive pathogen that can cause disease in almost any tissue or organ in the human body, primarily in compromised individuals (3). Staphylococcal infections were earlier treated using Penicillin. But over the year's resistance to this drug developed. Methicillin was the next drug of choice. While methicillin is very effective in treating most *Staphylococcus* infections some strains have developed resistance to methicillin and can no longer be killed by this antibiotic. These resistant bacteria are called Methicillin Resistant *Staphylococcus aureus* (MRSA) (4). Patients with breaks in their skin due to wounds, indwelling catheters or burns are those with certain risk of developing MRSA infection (2). Spread of MRSA infections can be controlled to a

great extent by maintaining personal hygiene after interaction with an MRSA infected person (4).

CLSI recommends the usage of cefoxitin instead of oxacillin for determination of resistance against Methicillin for *S. aureus* (1). To increase the sensitivity for the detection of heterogeneously resistant MRSA strains, cefoxitin is used which selectively inhibits the susceptible strains.

Tryptone, HM peptone B and yeast extract provide the essential nutrients along with carbonaceous, nitrogenous and Vitamin B complex nutrients. The proprietary chromogenic mixture incorporated in the medium is specifically cleaved by *Staphylococcus aureus* to form bluish green coloured colonies. Sodium pyruvate enhances the growth of *Staphylococcus* species. Sodium chloride in the medium helps to maintain the osmotic equilibrium of the medium. High concentration of sodium chloride also helps in inhibiting the accompanying microflora. The medium is made selective for MRSA by the addition of MeReSa Selective Supplement (FD229) and Cefoxitin Supplement (FD259).

Type of specimen

Clinical samples

Specimen Collection and Handling

For clinical samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (5, 6). After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

In Vitro diagnostic use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets



M1674 – HiCrome™ MeReSa Agar Base

HiCromeVeg™
Single Streak Rapid Differentiation Series
Freedom from BSE / TSE worries

HiCrome™ MeReSa Agar Base (M1674) is also available as HiCrome™ MeReSa HiVeg™ Agar Base (MV1674) wherein all the animal origin nutrients have been replaced by vegetable based nutrients.

HiCrome™ MeReSa Agar Base

For the isolation and selective identification of Methicillin Resistant *Staphylococcus aureus* (MRSA) from clinical isolates.

M1674

Limitations

1. Some intermediate strains may show poor growth due to nutritional variations and resistance to methicillin/cefoxitin.
2. Slight colour variation may be observed depending upon the utilization of the substrate by the organism.
3. Further confirmation must be carried out by sensitivity testing.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

- Appearance of Powder** : Cream to yellow coloured, homogeneous, free flowing powder.
- Gelling** : Firm, comparable with 1.5% Agar gel.
- Colour and Clarity of prepared medium** : Light yellow coloured, opaque gel forms in Petri plates.
- Reaction** : Reaction of 8.33% w/v aqueous solution at 25°C. pH:7.0 ± 0.2.
- Cultural Response** : Cultural characteristics observed with added MeReSa Selective Supplement (FD229) and Cefoxitin Supplement (FD259) after an incubation at 30-35°C for 18-48 hours.

Organisms (ATCC)	Inoculum CFU	Growth w/FD229 & FD259	Recovery w/FD229 & FD259	Colour of colony
<i>Staphylococcus aureus subsp aureus</i> (25923) (00034*)	≥10 ³	inhibited	0%	-
<i>Staphylococcus aureus</i> (MRSA) (43300)	50-100	luxuriant	≥50%	bluish - green
<i>Staphylococcus epidermidis</i> (12228) (00036*)	≥10 ³	inhibited	0%	-
<i>Escherichia coli</i> (25922) (00013*)	≥10 ³	inhibited	0%	-
<i>Enterococcus faecalis</i> (29212) (00087*)	≥10 ³	inhibited	0%	-

<i>Staphylococcus aureus subsp aureus</i> (6538) (00032*)	≥10 ³	inhibited	0%	-
<i>Staphylococcus xylosus</i> (29971)	≥10 ³	inhibited	0%	-

Key : * : Corresponds to WDCM number

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (5, 6).

References

1. CLSI 2010. Performance Standard for antimicrobial disk susceptibility testing, Twentieth Informational Supplem.
2. Dr. Alan Johnson, methicillin resistant *Staphylococcus aureus* (MRSA) infection. The Support group for MSRA sufferers and Dependents, Aug 1st, 2005.ent.
3. DWorkin M et. al 2006. The Prokaryotes (a Handbook on the Biology of Bacteria) 3rd ed, Vol. 2, page 345.
4. Methicillin Resistant *Staphylococcus aureus* Copyright ©1997-2005 Canadian Centre for Occupational Health and Safety, Sept 19th, 2005.
5. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
6. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1

Ready Prepared Media

Code	Product Name	Usage	Packing
Category : 90 mm Ready Prepared Petri Plates			
MP1674	HiCrome™ MeReSa Agar Plate	for isolation and selective identification of Methicillin Resistant <i>Staphylococcus aureus</i> (MRSA) from clinical specimens.	50 plts

For identification of MRSA and MRSE

HiCrome™ MRSA Agar Base, Modified

Recommended for the differentiation and identification of MRSA and MRSE *Staphylococcus* species.

M1953

Composition **

Ingredients	Grams/Litre
Peptone	23.000
Sodium chloride	10.000
Sodium pyruvate	5.000
Chromogenic substrate	0.770
Inhibitor mixture	7.000
Agar	15.00

Final pH (at 25°C) 7.2±0.2

** Formula adjusted, standardized to suit performance parameters

Directions

Suspend 30.38 grams in 500 ml distilled water. Heat to boiling to dissolve the medium completely. DO NOT AUTOCLAVE. Cool to 45-50°C. Aseptically add sterile rehydrated contents of 1 vial of MeReSa Selective Supplement (FD229) and Cefoxitin supplement (FD259), both in combination for more selectivity as desired. Mix well and pour into sterile Petri plates.

Principle and Interpretation

MRSA is a resistant variation of the common bacterium *Staphylococcus aureus* and MRSE is a resistant variation of the common bacterium *Staphylococcus epidermidis*. *Staphylococcus aureus* is an invasive pathogen that can cause disease in almost any tissue or organ in the human body, primarily in compromised individuals (1). Staphylococcal infections were earlier treated using Penicillin. But over the years resistance to this drug developed. Methicillin was the next drug of choice. While methicillin is very effective in treating most *Staphylococcus* infections some strains have developed resistance to methicillin and can no longer be killed by this antibiotic. These resistant bacteria are called Methicillin Resistant *Staphylococcus aureus* (MRSA) (2). Patients with breaks in their skin due to wounds, indwelling catheters or burns are those with certain risk of developing MRSA infection (3). Spread of MRSA infections can be controlled to a great extent by maintaining personal hygiene after interaction with an MRSA infected person (2).

Peptone provide the essential nutrients along with carbonaceous, nitrogenous and Vitamin B complex nutrients. The chromogenic mixture incorporated in the medium is specifically cleaved by *Staphylococcus aureus* to give green coloured colonies whereas Methicillin Resistant *Staphylococcus epidermidis* gives blue coloured colonies. This medium helps in identification and differentiation of MRSA and MRSE. Sodium pyruvate enhances the growth of *Staphylococcus* species. Sodium chloride in the medium helps to maintain the osmotic equilibrium of the medium. High concentration of sodium chloride also helps in inhibiting

the accompanying microflora. Inhibitor mixture imparts selectivity to the medium. Cefoxitin is recommended to use for selective isolation of MRSA. The medium is made selective for MRSA by the addition of MeReSa Selective Supplement (FD229) and Cefoxitin Supplement (FD259) in combination.

Type of specimen

Clinical samples

Specimen Collection and Handling

For clinical samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (4, 5). After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

In Vitro diagnostic use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets.



M1953 HiCrome™ MRSA Agar Base, Modified (Mixture)

For identification of MRSA and MRSE

HiCrome™ MRSA Agar Base, Modified

Recommended for the differentiation and identification of MRSA and MRSE *Staphylococcus* species.

M1953

Limitations

1. Some intermediate strains may show poor growth due to nutritional variations and resistance to methicillin/cefoxitin.
2. Slight colour variation may be observed depending upon the utilization of the substrate by the organism.
3. Further confirmation must be carried out by sensitivity testing.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

- Appearance of Powder** : Cream to beige homogeneous free flowing powder
- Gelling** : Firm, comparable with 1.5% Agar gel
- Colour and Clarity** : Light purple coloured, clear to slightly
- of prepared medium** opalescent gel forms in Petri plates
- Reaction** : Reaction of 6.08% w/v aqueous solution at 25°C. pH : 7.2±0.2.
- Cultural Response** : Cultural characteristics observed with added MeReSa Selective Supplement (FD229) and Cefoxitin Supplement (FD259) after an incubation at 30-35°C for 18-48 hours.

Organism (ATCC)	Inoculum (CFU)	Growth	Recovery	Colour of colony
<i>Escherichia coli</i> (25922) (00013*)	≥10 ³	inhibited	0%	-
<i>Enterococcus faecalis</i> (29212) (00087*)	≥10 ³	inhibited	0%	-
<i>Staphylococcus aureus</i> , MRSA (43300)	50-100	luxuriant	≥50%	green

<i>Staphylococcus epidermidis</i> , MRSE	50-100	luxuriant	≥50%	blue
<i>Staphylococcus xylosus</i> (29971)	≥10 ³	inhibited	0%	-

Key : * : Corresponds to WDCM number

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.


Product performance is best if used within stated expiry period.

Disposal


User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (4, 5).

References

1. DWorkin M et. al 2006. The Prokaryotes (a Handbook on the Biology of Bacteria) 3rd ed, Vol. 2, page 345.
2. Methicillin Resistant *Staphylococcus aureus* Copyright © 1997-2005 Canadian Centre for Occupational Health and Safety, Sept 19th, 2005.
3. Dr. Alan Johnson, methicillin resistant *Staphylococcus aureus* (MRSA) infection. The Support group for MSRA sufferers and Dependents, Aug 1st, 2005.
4. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
5. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1



M1953 HiCrome™ MRSA Agar Base, Modified



M1953 HiCrome™ MRSA Agar Base, Modified

For Identification of MRSA and MRSE

HiCrome™ Rapid MRSA Agar Base

It is recommended for rapid isolation and identification of Methicillin Resistant *Staphylococcus aureus* (MRSA).

M1974

Composition **

Ingredients	Grams/Litre
Special peptone	20.00
Casitose#	20.00
Sodium chloride	8.50
Carbohydrate	14.00
Phenol red	0.025
Chromogenic mix	6.50
Amino-Vitamin mix	1.20
Agar	15.00

Final pH (at 25°C) 7.4 ± 0.2

** Formula adjusted, standardized to suit performance parameters

Equivalent to Casein peptone

Directions

Suspend 85.23 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Cool to 45-50°C. Aseptically add sterile rehydrated contents of 1 vial of MRSA Selective Supplement (FD319). Mix well and pour into sterile Petri plates. DO NOT AUTOCLAVE.

Principle and Interpretation

MRSA is a resistant variation of the common bacterium *Staphylococcus aureus*. It is an invasive pathogen that can cause disease in almost any tissue or organ in the human body, primarily in compromised individuals (1). Staphylococcal infections were earlier treated using Penicillin. But over the years resistance to this drug developed. Methicillin was the next drug of choice. While methicillin is very effective in treating most *Staphylococcus* infections some strains have developed resistance to methicillin and can no longer be killed by this antibiotic. These resistant bacteria are called Methicillin Resistant *Staphylococcus aureus* (MRSA) (2). Patients with breaks in their skin due to wound, indwelling catheters or burns are those with certain risk of developing MRSA infection (3).

Special peptone, Casitose and amino-vitamin mix provides essential nutrients for growth. Carbohydrate is the source of carbon and energy. Phenol red is the pH indicator. The chromogenic mixture incorporated in the medium is specifically cleaved by *Staphylococcus aureus* (MRSA) to give greenish yellow coloured colonies. Sodium chloride in the medium helps to maintain the osmotic equilibrium of the medium. High concentration of sodium chloride also helps in inhibiting the accompanying microflora. Agar acts as solidifying agent.

Type of specimen


Clinical samples

Specimen Collection and Handling

For Clinical samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (4, 5). After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions


In Vitro diagnostic use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets



M1974 HiCrome™ Rapid MRSA Agar Base

HiCrome™ Rapid MRSA Agar Base

It is recommended for rapid isolation and identification of Methicillin Resistant *Staphylococcus aureus* (MRSA).



Limitations

1. Some intermediate strains may show poor growth due to nutritional variations and resistance to methicillin/cefoxitin.
2. Slight colour variation may be observed depending upon the utilization of the substrate by the organism.
3. Further confirmation must be carried out by sensitivity testing.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

- Appearance of Powder** : Cream to beige homogeneous free flowing powder
- Gelling** : Firm, comparable with 1.5% Agar gel.
- Colour and Clarity of prepared medium** : Red coloured, clear to slightly opalescent gel forms in Petri plates.
- Reaction** : Reaction of 8.52% w/v aqueous solution at 25°C. pH: 7.4 ± 0.2.
- Culture Response** : Cultural characteristics observed with added MRSA Selective Supplement (FD319) after an incubation at 30-35°C for 18-24 hours

Organism (ATCC)	Inoculum (CFU)	Growth	Recovery	Colour of colony
<i>Staphylococcus aureus</i> , MRSA ATCC 43300	50-100	luxuriant	≥50%	greenish yellow (Note: Green colour may develop after 48 hours)
<i>Staphylococcus epidermidis</i> , MRSE	50-100	luxuriant	≥50%	blue
<i>Staphylococcus aureus subsp aureus</i> ATCC 25923 (00034*)	≥10 ³	inhibited	0%	-

<i>Staphylococcus aureus subsp aureus</i> ATCC 6538 (00032*)	≥10 ³	inhibited	0%	-
<i>Escherichia coli</i> ATCC 25922 (00013*)	≥10 ³	inhibited	0%	-
<i>Candida albicans</i> ATCC 10231 (00054*)	≥10 ³	inhibited	0%	-

Key : * : Corresponds to WDCM number

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (4, 5).

References

1. DWorkin M et. al 2006. The Prokaryotes (a Handbook on the Biology of Bacteria) 3rd ed, Vol. 2, page 345.
2. Methicillin Resistant *Staphylococcus aureus* Copyright © 1997-2005 Canadian Centre for Occupational Health and Safety, Sept 19th, 2005.
3. Dr. Alan Johnson, methicillin resistant *Staphylococcus aureus* (MRSA) infection. The Support group for MSRA sufferers and Dependents, Aug 1st, 2005.
4. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
5. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1

Ready Prepared Media

Code	Product Name	Usage	Packing
Category : 90 mm Ready Prepared Petri Plates			
MP1974	HiCrome™ Rapid MRSA Agar Plate	for rapid isolation and identification of Methicillin Resistant <i>Staphylococcus aureus</i> (MRSA) from clinical specimens.	50 plts

For identification of *Bacillus* species

HiCrome™ Bacillus Agar

Recommended for isolation and differentiation between various species of *Bacillus* from a mixed culture by chromogenic method.

M1651

Composition **

Ingredients	Grams/Litre
Peptone	10.00
HM extract#	1.00
D-Mannitol	10.00
Sodium chloride	10.00
Chromogenic mixture	3.20
Phenol red	0.025
Agar	15.00

Final pH (at 25°C) 7.1 ± 0.2

** Formula adjusted, standardized to suit performance parameters

#Equivalent to Meat extract

Directions

Suspend 49.22 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C and aseptically add rehydrated contents of 1 vial of Bacillus Selective Supplement (FD324) if desired. Mix well and pour into sterile Petri plates.

Principle and Interpretation

Majority of *Bacillus* species apparently have little or no pathogenic potential and are rarely associated with disease in humans or lower animals. The principal exception to this are *Bacillus anthracis*, the agent of anthrax, and *Bacillus cereus*. However a number of other species, particularly those of the *B. subtilis* group, have been implicated in food poisoning and other human and animal infections (4). *Bacillus cereus* causes food poisoning due to consumption of contaminated rice (2, 1, 5) other starchy foods such as potato, pasta and cheese have also been implicated, eye infections and a wide range of other clinical conditions like abscess formation, meningitis, septicemia and wound infection.

HiCrome™ *Bacillus* Agar is based on the formulation of MYP Agar formulated by Mossel et al (2) used for enumeration of *Bacillus cereus* and *Bacillus thuringiensis* when present in large number in certain foodstuffs. The medium contains peptone and HM extract, which provide nitrogenous and carbonaceous compounds, long chain amino acids, vitamins and other growth nutrients. Mannitol serves as the fermentable carbohydrate, fermentation of which can be detected by phenol red. Mannitol fermenting organisms like *B. megaterium* yield yellow coloured colonies. The chromogenic mixture present in the medium is cleaved by the enzyme b-galactosidase found in *B. cereus* resulting in the formation of blue colonies. *B. thuringiensis* also grows as blue/green colonies on this medium as *B. cereus* and *B. thuringiensis* are biochemically identical. If selective isolation of *B. cereus* or *B. thuringiensis* is required aseptically add Bacillus Selective Supplement (FD324).

Type of specimen

Clinical, food samples

Specimen Collection and Handling


For Clinical samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (6, 7).

For Food samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (8)

After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

In Vitro diagnostic use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets



M1651 HiCrome™ Bacillus Agar

HiCrome™ Bacillus Agar

Recommended for isolation and differentiation between various species of *Bacillus* from a mixed culture by chromogenic method.

M1651

Limitations

1. Due to variable nutritional requirements, some strains may show poor growth on this medium.
2. Slight colour variation may be observed depending upon the utilization of the substrate by the organism.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

- Appearance of Powder** : Light yellow to pink coloured, homogeneous, free flowing powder.
- Gelling** : Firm, comparable with 1.5% Agar gel.
- Colour and Clarity of prepared medium** : Red coloured, clear to slightly opalescent gel forms in Petri plates.
- Reaction** : Reaction of 4.92% w/v aqueous solution at 25°C. pH : 7.1 ± 0.2.
- Cultural Response** : Cultural characteristics observed after an incubation at 30°C for 24-48 hours.

Organisms (ATCC)	Inoculum (CFU)	Growth **	Recovery **	Growth ***	Recovery ***	Colour of colony
<i>Bacillus subtilis</i> (6633) (00003*)	50-100	fair	20-30%	inhibited	0%	yellowish green to green
<i>Bacillus cereus</i> (10876)	50-100	good-luxuriant	≥50%	good-luxuriant	≥50%	light blue,# large, flat with blue centre
<i>Bacillus thuringiensis</i> (10792)	50-100	good-luxuriant	≥50%	good-luxuriant	≥50%	light blue, large, flat with irregular margins
<i>Bacillus megaterium</i> (14581)	50-100	good-luxuriant	≥50%	inhibited	0%	yellow, mucoid colonies
<i>Bacillus coagulans</i> (7050) (00002*)	50-100	good-luxuriant	≥50%	inhibited	0%	pink, small, raised colonies
<i>Bacillus pumilis</i> (14884)	50-100	good-luxuriant	≥50%	poor	10-20%	light green to green colonies

<i>Staphylococcus aureus</i> subsp <i>aureus</i> (25923) (00034*)	50-100	luxuriant	≥50%	inhibited	0%	yellow colonies
<i>Enterococcus faecalis</i> (29212) (00087*)	50-100	luxuriant	≥50%	inhibited	0%	light green to green colonies

Key : * : Corresponds to WDCM number
 ** : Growth without addition of FD324 *** : Growth with addition of FD324
 # : Colony surrounded by pink halo

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (6, 7).

References

1. Bouza E., Grant S., Jordan C. et al, 1979, Arch. Ophthalmol., 97:488.
2. Mortimer P. R. and McCann G., 1974, Lancet, 1043.
3. Mossel D. A. A., Koopman M. J. and Jongerium E., 1967, Appl. Microbiol., 15:650.
4. Murray P. R., Baron J. H., Pfaller M. A., Jorgensen J. H. and Tenover F. C., (Eds.), 2003, Manual of Clinical Microbiology, 8th Ed. American Society for Microbiology, Washington, D.C.
5. Wohlgenuth K., Kirkbride C. A., Bicknell E. J. and Ellis R. P., 1972 Am. Vet. Met. Ass. 161:1691.
6. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1
7. Salfinger Y., and Tortorello M.L. Fifth (Ed.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.

Ready Prepared Media			
Code	Product Name	Usage	Packing
Category : 90 mm Ready Prepared Petri Plates			
MP1651	HiCrome™ Bacillus Agar Plate	for isolation and differentiation between various species of <i>Bacillus</i> from a mixed cultures.	50 plts

For Identification of *Clostridium perfringens*

M-CP Agar Base

Recommended by the Directive of the Council of the European Union 98/83/EC for isolation and enumeration of *Clostridium perfringens* from water sample using membrane filtration technique.

M1354

Composition **	
Ingredients	Grams/Litre
Tryptose	30.00
Yeast extract	20.00
Sucrose	5.00
L-Cysteine hydrochloride	1.00
Magnesium sulphate, 7H ₂ O	0.10
Bromo cresol purple	0.04
Ferric chloride, 6H ₂ O	0.09
Indoxyl-β-D-glucoside	0.06
Agar	15.00

Final pH (at 25°C) 7.6 ± 0.2

** Formula adjusted, standardized to suit performance parameters

Directions


Suspend 35.60 grams (the equivalent weight of dehydrated powder per 485 ml) in 485 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Aseptically add the rehydrated contents of 1 vial of M-CP Selective Supplement I (FD153) and 1 vial of M-CP Selective Supplement II (FD154) or rehydrated contents of 1 vial of M-CP Selective Supplement II Modified (FD154A). Mix well and pour into sterile Petri plates.

Principle and Interpretation

Clostridial species are one of the major causes of food poisoning/ gastro-intestinal illnesses. They are gram-positive, spore-forming rods that occur naturally in the soil (2). Among the family are: *Clostridium botulinum* which produces one of the most potent toxins in existence; *Clostridium tetani*, causative agent of tetanus; and *Clostridium perfringens* commonly found in wound infections and diarrhoea cases. The use of toxins to damage the host is a method deployed by many bacterial pathogens. The major virulence factor of *C. perfringens* is the CPE enterotoxin, which is secreted upon invasion of the host gut, and contributes to food poisoning and other gastrointestinal illnesses (2). Several solid media have been devised for quantitation of *C. perfringens*. The selectivity of the media is achieved by incorporation of one or more antibiotics that inhibit certain anaerobes or facultative anaerobes. M-CP Agar Base is prepared as per the formula of Armon

and Payment (1). It is also recommended by the Directive of the Council of the European Union 98/83/EC (3) for isolation and enumeration of *Clostridium perfringens* from water sample using membrane filtration technique.

Tryptose, yeast extract provide nitrogenous and carbonaceous compounds, long chain amino acids, vitamin B complex and other essential growth nutrients compounds while sucrose is the fermentable carbohydrate. Bromocresol purple serves as a pH indicator. Indoxyl-β-D-glucoside is a chromogenic substrate for β-D-glucosidase or cellobiose and phenolphthalein diphosphate for the detection of acid phosphatase. The addition of D-cycloserine and polymyxin B (FD153) makes the medium inhibitory to accompanying non-clostridial microflora and thus allows analysis of both clostridial vegetative cells and spores. Further selectivity is provided by incubation under anaerobic conditions. Yellow (cellobiose-negative) colonies becoming old rose to pink-red upon exposure to ammonia fumes for 30 seconds are considered to be presumptive *C. perfringens*. Colour differentiation on M-CP Agar Base is sometimes difficult, so typical colonies (yellow turning into pink) as well as atypical colonies (green or those that remain yellow upon exposure to ammonia fumes) are picked for confirmation. Presumptive *C. perfringens* can be confirmed by sulphite reduction, gram-positive, sporulating rods, non-motile, reduction of nitrate, gelatine liquefaction, lactose fermentation and other biochemical tests (4).




M1354 M-CP Agar Base

HiCromeVeg™ Freedom from BSE / TSE worries
Single Streak Rapid Differentiation Series

M-CP Agar Base (M1354) is also available as M-CP HiVeg™ Agar Base (MV1354) wherein all the animal origin nutrients have been replaced by vegetable based nutrients.

M-CP Agar Base

Recommended by the Directive of the Council of the European Union 98/83/EC for isolation and enumeration of *Clostridium perfringens* from water sample using membrane filtration technique.



M1354

Type of specimen

Food samples

Specimen Collection and Handling

For Food samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (5, 6). After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets

Limitations

1. Colour differentiation on M-CP Agar Base is sometimes difficult, so typical colonies (yellow turning into pink) as well as atypical colonies (green or those that remain yellow upon exposure to ammonia fumes) are picked for confirmation.
2. Due to variable nutritional requirements, some strains may show poor growth on this medium.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

- Appearance of Powder** : Light yellow to light green coloured, homogeneous, free flowing powder.
- Gelling** : Firm, comparable with 1.5% Agar gel.
- Colour and Clarity of prepared medium** : Purple coloured, clear to slightly opalescent gel forms in Petri plates.
- Reaction** : Reaction of 7.12% w/v aqueous solution at 25°C. pH: 7.6 ± 0.2.
- Culture Response** : Cultural characteristics observed after an incubation at 44°C for 24-48 hours with added contents of 1 vial of M-CP Selective Supplement I (FD153) and 1 vial of M-CP Selective Supplement II (FD154) or rehydrated contents of 1 vial of M-CP Selective Supplement II Modified (FD154A) under anaerobic conditions.

Organisms (ATCC)	Inoculum (CFU)	Growth	Colour of Colonies
<i>Clostridium perfringens</i> (12924)	50-100	good	yellow**
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> (25923) (00034*)	≥10 ³	inhibited	-
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> (6633) (00003*)	≥10 ³	inhibited	-
<i>Salmonella</i> Typhi (6539)	≥10 ³	inhibited	-

Key : * : Corresponds to WDCM number

** : colonies becomes old rose to light pink-red upon exposure to ammonia fumes for 30 seconds.

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition. Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (6, 7).

References

1. Armon R. and Payment P., 1988, Can. J. Microbiol., 34:78-79.
2. Czeizulin J. R., Hanna P. C., McClane B. A., 1993, Infect. Immun., 61: 3429-3439.
3. Directive of the Council of the European Union 98/83/EC Sartory D. P., Field M., Curbishley S. M., Pritchard A. M., 1998, Lett. Appl. Microbiol., 27:323-327.
4. Sartory D.P., Field M, Curbishley S.M., Pritchard A.M., 1998, Lett. Appl. Microbiol., 27:323-327.
5. Salfinger Y., and Tortorello M.L. Fifth (Ed.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.
6. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
7. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S. and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1

For the enumeration of citrate-fermenting lactic acid bacteria

HiCrome™ Nickels and Leesment Medium

Recommended for the enumeration of citrate-fermenting lactic acid bacteria from milk, milk products and mesophilic starter cultures.

M1712

Composition **

Ingredients	Grams/Litre
Tryptone	18.00
Yeast extract	4.50
Gelatine	2.25
Glucose (Dextrose)	4.50
Lactose	4.50
Sodium chloride	3.60
Trisodium citrate dihydrate	1.80
Calcium lactate pentahydrate	8.00
Tricalcium dicitrate tetrahydrate	6.65
Carboxymethyl cellulose (CMC)	0.40
Chromogenic substrate (X-gal)	0.20
Agar	15.00

Final pH (at 25°C) 6.65 ± 0.05

** Formula adjusted, standardized to suit performance parameters

Directions

Suspend 66.0 grams (the equivalent weight of dehydrated medium per litre) in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. If desired, add rehydrated contents of 2 vials of HiCrome™ Nickels and Leesment Selective Supplement (FD245). Mix well and pour into sterile Petri plates.

Principle and Interpretation

Lactic acid bacteria are widespread in nature and are best known for their activities in major food such as dairy, meat and vegetable products (1). Testing for lactic acid bacteria in dairy products may be useful for various reasons like evaluating lactic starter cultures; determining the cause of acid defects in milk products, controlling the quality of cured cheese, cultured milks and uncultured products containing added cultures (2). HiCrome™ Nickels and Leesment Medium is a modification of Modified Nickels and Leesment Medium formulated as per APHA (1) and is used for the enumeration of citrate-fermenting lactic acid bacteria using colony count technique at 25°C. Tryptone and yeast extract serve as carbon and nitrogen sources, long chain amino acids, vitamin B complex and other essential growth nutrients. Lactose and Glucose (Dextrose) are the carbohydrate source in the medium.

X-gal differentiates between *Lactococcus lactis* subsp. *lactis* and *Leuconostoc* species. *Lactococcus lactis* subsp. *lactis biovar diacetylactis* colonies are white with a clear zone. *Lactococcus lactis* subsp. *lactis* and *Lactococcus lactis* subsp. *cremoris* colonies are white without a clear zone. *Leuconostoc* species are blue, with or without a clear zone. HiCrome™ Nickels and Leesment Medium with the addition of HiCrome™ Nickels and Leesment Supplement (FD245) can be used for enumeration of *Leuconostoc* (1). Vancomycin acts as a supplement for the selective isolation of *Leuconostoc* from a mix flora of lactic acid bacteria. Sodium chloride maintains osmotic equilibrium and various salts provide essential ions.

Type of specimen

Dairy : milk and milk product samples


Specimen Collection and Handling

For dairy samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (3, 4).

After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

Read the label before opening the container. Wear protective gloves/ protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture.



M1712 HiCrome™ Nickels and Leesment Medium

HiCromeVeg™ Freedom from BSE / TSE worries
Single Streak Rapid Differentiation Series
HiCrome™ Nickels and Leesment Medium (M1712) is also available as HiCrome™ Nickels and Leesment HiVeg™ Medium (MV1712) wherein all the animal origin nutrients have been replaced by vegetable based nutrients.

For the enumeration of citrate-fermenting lactic acid bacteria

HiCrome™ Nickels and Leesment Medium

Recommended for the enumeration of citrate-fermenting lactic acid bacteria from milk, milk products and mesophilic starter cultures.

M1712

Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets

Limitations

1. Due to variable nutritional requirements, some strains may show poor growth on this medium.
2. Slight colour variation may be observed depending upon the utilization of the substrate by the organism.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

- Appearance of Powder** : Cream to light yellow coloured homogeneous, free flowing powder.
- Gelling** : Firm, comparable with 1.5% Agar gel.
- Colour and Clarity of prepared medium** : White coloured opaque gel forms in Petri plates containing white precipitate.
- Reaction** : Reaction of 6.6% w/v aqueous solution at 25°C. pH : 6.65 ± 0.05.
- Cultural Response** : Cultural characteristics observed after an incubation at 25-30°C for 48-72 hours.

Organisms (ATCC)	Growth	Growth**	Colour of colony
<i>Lactococcus lactis</i> subsp <i>lactis</i> biovar <i>diacetylactis</i>	good-luxuriant	inhibited	white with a clear zone
<i>Lactococcus lactis</i> subsp <i>lactis</i> (19435) (00016*)	good-luxuriant	inhibited	white without a clear zone
<i>Lactococcus lactis</i> subsp <i>cremoris</i> (19257)	good-luxuriant	inhibited	white without a clear zone
<i>Leuconostoc mesenteroides</i> (9135) (00108*)	good-luxuriant	good-luxuriant	blue without a clear zone

Key : * : Corresponds to WDCM number

** = with the addition of HiCrome Nickels and Leesment Selective Supplement (FD245)

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (4, 5).

References

1. Salfinger Y., and Tortorello M.L. Fifth (Ed.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.
2. Marshall R.T., 1992, Standard Methods for the Examination of Dairy products, 16th Ed, American Public Health Association, Washington D.C.
3. American Public Health Association, Standard Methods for the Examination of Dairy Products, 1978, 14th Ed., Washington D.C
4. Wehr H. M. and Frank J. H., 2004, Standard Methods for the Microbiological Examination of Dairy Products, 17th Ed., APHA Inc., Washington, D.C.
5. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
6. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1

HiCrome™ Bifidobacterium Agar

Recommended for the differentiation of *Bifidobacterium* and *LactoBacillus* species.

M1960

Composition **

Ingredients	Grams/Litre
Peptone special	23.00
Sodium chloride	5.00
M Protein powder#	5.00
Chromogenic mixture	10.48
Agar	16.00

Final pH (at 25°C) 7.2 ± 0.2

** Formula adjusted, standardized to suit performance parameters
#Equivalent to Milk Protein

Directions

Suspend 59.48 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Mix well and pour into sterile Petri plates.

Principle and Interpretation

The genus *Bifidobacterium* is the third most numerous bacterial populations found in the human intestine after *Bacteroides* and *Eubacterium*. It is an anaerobic bacteria that makes up the gut microbial flora. It resides in the colon and have health benefits for their hosts. Bifidobacteria are also associated with lower incidences of allergies (1, 2). Bifidobacterium Agar is used for the cultivation and maintenance of *Bifidobacterium* species (3).

Peptone special provides nitrogenous and carbonaceous compounds, long chain amino acids, vitamins and other essential growth nutrients. Sodium chloride maintains osmotic balance. M Protein powder aids in detecting casein hydrolysis activity which is exhibited by *Bifidobacterium breve*. A halo zone is observed around the colony in case of casein hydrolysis. The indicator system in the chromogenic mixture helps in distinguishing between *Lactobacillus* and *Bifidobacterium* species. *Lactobacillus* species usually produce green colonies with opaque background. *Bifidobacterium infantis* produces dark blue to bluish green colonies. Agar serves as a solidifying agent.

Type of specimen

Clinical, Dairy : Milk & Milk products samples

Specimen Collection and Handling

For dairy samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (4, 5).

For Clinical samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (6, 7).

After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions


In Vitro diagnostic use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets.

Limitations

1. Due to variable nutritional requirements, some strains may show poor growth on this medium.
2. Slight colour variation may be observed depending upon the utilization of the substrate by the organism.
3. *Bifidobacterium* species are strict anaerobes, hence condition must be appropriately maintained

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.



M1960 - HiCrome™ Bifidobacterium Agar

HiCrome™ Bifidobacterium Agar

Recommended for the differentiation of *Bifidobacterium* and *LactoBacillus* species.

M1960

Quality Control

- Appearance of Powder** : Cream to yellow homogeneous free flowing powder
- Gelling** : Firm, comparable with 1.6% Agar gel
- Colour and Clarity of prepared medium** : Reddish orange coloured clear to slightly opalescent gel forms in Petri plates
- Reaction** : Reaction of 5.95% w/v aqueous solution at 25°C. pH :7.2±0.2
- Cultural Response** : Cultural characteristics observed after an incubation at 35-37°C for 48 hours in an anaerobic conditions.

Organisms (ATCC)	Inoculum (CFU)	Growth	Recovery	Colour of the colony
<i>Bifidobacterium infantis</i> (25962)	50-100	good-luxuriant	≥50%	Dark blue - bluish green
<i>Bifidobacterium breve</i> (15698)	50-100	good-luxuriant	≥50%	Red-pink with halo zone
<i>Lactobacillus plantarum</i> (8014)	50-100	good-luxuriant	≥50%	Green colonies w/ hazy background
<i>Lactobacillus fermentum</i> (9338)	50-100	good-luxuriant	≥50%	Pink without halo zone

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal


User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (6, 7).

References

1. Björkstén B., Sepp E., Julge K., Voor T., and Mikelsaar M., 2001, J. Allergy Clin. Microbiol., Volume 108, Issue 4, 516-520.
2. Guarner F., and Malagelada J. R., 2003, The Lancet, Vol. 361, Issue 9356, 8 February 2003, 512-519.
3. Atlas R. M. 2004, 3rd Edi. Handbook of Microbiological Media, Parks, L. C. (Ed.), CRC Press, Boca Raton.
4. American Public Health Association, Standard Methods for the Examination of Dairy Products, 1978, 14th Ed., Washington, D.C
5. Wehr H. M. and Frank J. H., 2004, Standard Methods for the Microbiological Examination of Dairy Products, 17th Ed., APHA Inc., Washington, D.C.
6. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
7. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1



M1960 - HiCrome™ Bifidobacterium Agar



M1960 - HiCrome™ Bifidobacterium Agar

For Identification of Acinetobacter

HiCrome™ Acinetobacter Agar Base

Recommended for selective isolation of *Acinetobacter* species from environmental and clinical samples.

M1938

Composition **

Ingredients	Grams/Litre
Peptone special	9.00
Sodium chloride	5.00
Selective mix	0.50
Chromogenic mixture	1.35
Agar	15.00

Final pH (at 25°C) 7.0 ± 0.2

** Formula adjusted, standardized to suit performance parameters

Directions

Suspend 30.85 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. DO NOT AUTOCLAVE. Cool to 45-50°C and add the rehydrated contents of two vials of MDR Selective Supplement (FD271) or two vials of Leeds Acinetobacter Selective Supplement (FD335). Mix well and pour into sterile Petri plates.

Principle and Interpretation

Acinetobacter species are gram negative ubiquitous bacteria that have been isolated from patients with nosocomial infection, environment, soil, and water. Acinetobacter is mostly found in every type of infections (3). There is an alarming situation as Acinetobacter baumannii is found to be resistant to most commonly used antibiotics which includes beta-lactams and aminoglycosides (2,3). Immunocompromised patients requiring mechanical respirations are at more risk of infection by Acinetobacter species.(1)

Peptone special provides nitrogenous and carbonaceous compounds, long chain amino acids and vitamins to the organisms. Sodium chloride maintains the osmotic balance. Selective mix inhibits gram positive organisms. The chromogenic mixture in the medium allows the differentiation of Acinetobacter species from other organisms.

Type of specimen

Clinical Samples, environment

Specimen Collection and Handling

For clinical samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (4, 5). After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

In Vitro diagnostic use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets

Limitations

1. Due to variable nutritional requirements, some strains may show poor growth on this medium.
2. Slight colour variation may be observed depending upon the utilization of the substrate by the organism.

Performance and Evaluation


Performance of the medium is expected when used as per the direction on the label within the recommended temperature.



M1938 HiCrome™ Acinetobacter Agar Base

HiCrome™ Acinetobacter Agar Base

Recommended for selective isolation of *Acinetobacter* species from environmental and clinical samples.



M1938

Quality Control

- Appearance of Powder** : Light yellow to yellow homogeneous free flowing powder.
- Gelling** : Firm, comparable with 1.5% Agar gel.
- Colour and Clarity of prepared medium** : Yellow coloured, clear to slightly opalescent gel forms in Petri plates.
- Reaction** : Reaction of 3.09% w/v aqueous solution at 25°C. pH:7.0 ± 0.2.
- Cultural Response** : Cultural characteristics observed after an incubation at 35-37°C for 24-48 hours.

Organism (ATCC)	Inoculum (CFU)	Growth	Recovery with FD271	Colour of colony
<i>Acinetobacter baumannii</i> (BAA-1605)	50-100	luxuriant	≥50%	Light purple with halo
<i>Acinetobacter baumannii</i> (BAA-747)	≥10 ³	inhibited	0%	-
<i>Acinetobacter baumannii</i> (19606)	≥10 ³	inhibited	0%	-
<i>Acinetobacter Iwofii</i> (15309)	≥10 ³	inhibited	0%	-
<i>Acinetobacter haemolyticus</i> (19002)	≥10 ³	inhibited	0%	-
<i>Escherichia coli</i> (25922) (00013*)	≥10 ³	inhibited	0%	-
<i>Enterococcus faecalis</i> (29212) (00087*)	≥10 ³	inhibited	0%	-

Key : * : Corresponds to WDCM number

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (4, 5).

References

- Bergogne-Berezin, E., m. L. Joly-Guillou, and J.F. Vieu. 1987. Epidemiology of nosocomial infections due to *Acinetobacter calcoaceticus*. J. Hosp. Infect. 10:105-113
- Montefour, K., et.al.2008. *Acinetobacter baumannii* : An EmergingMultidrug Resistant pathogen in critical care Nurse; 28:15-25
- Valentine, S.C., et.al. 2008 Phenotypic and molecular characterization of *Acinetobacter baumannii*. Clinical isolates from nosocomial outbreaks in Los Angeles Country, California. J.Clin. Microbiology.; 46:2499-2507
- Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
- Jorgensen,J.H., Pfaller , M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1

Ready Prepared Media

Code	Product Name	Usage	Packing
Category : 90 mm Ready Prepared Petri Plates			
MP1938	HiCrome™ MDR Acinetobacter Agar Plate	for selective isolation of <i>Acinetobacter</i> species from enviromental and clinical samples.	50 plts

HiCrome™ OGYE Agar Base

Recommended for isolation and enumeration of yeasts and moulds from milk and milk products by chromogenic method.

M1467

Composition **

Ingredients	Grams/Litre
Yeast extract	4.00
Dextrose (Glucose)	20.00
Chromogenic mixture	1.10
Agar	12.00

Final pH (at 25°C) 7.0 ± 0.2

** Formula adjusted, standardized to suit performance parameters

Directions

Suspend 18.55 grams in 500 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C and aseptically add reconstituted contents of one vial of Oxytetra Selective Supplement (FD032). Mix well and pour into sterile Petri plates.

Principle and Interpretation

OGYE Agar Media were originally formulated by Mossel et al (1, 2) for the isolation and enumeration of yeasts and moulds from foodstuffs. Mossel et al (3) further added Oxytetracycline as a selective agent and found that the use of Oxytetracycline in a medium with a neutral pH gives increased counts of yeasts and moulds as compared to media having a low pH to suppress bacterial growth. HiCrome™ OGYE Agar is a selective and differential medium, which facilitates rapid isolation of yeasts and moulds from milk and milk products.

Yeast extract provides essential growth nutrients. Dextrose acts as carbon and energy source. Although the low pH helps to reduce the bacterial flora, Oxytetracycline makes the medium, more selective by inhibiting the growth of lactobacilli encountered in milk and milk-products at low pH. Incorporation of chromogenic compounds into the growth medium helps in identification of yeasts and moulds isolates directly on primary isolation. **Aspergillus brasiliensis* appear as light blue coloured colonies with black spores due to presence of chromogenic mixture, *C.albicans* shows green coloured colonies and *Saccharomyces cerevisiae* forms colourless colonies.

Type of specimen

Dairy : Milk & Milk products samples

Specimen Collection and Handling

For dairy samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (4, 5).


After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

Read the label before opening the container. Wear protective gloves/ protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets

Limitations

1. Due to variable nutritional requirements, some strains may show poor growth on this medium.
2. Slight colour variation may be observed depending upon the utilization of the substrate by the organism.



M1467 – HiCrome™ OGYE Agar Base
*Formerly known as *Aspergillus niger*

HiCrome™ OGYE Agar Base

Recommended for isolation and enumeration of yeasts and moulds from milk and milk products by chromogenic method.

M1467

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

- Appearance of Powder** : Cream to yellow coloured homogeneous, free flowing powder.
- Gelling** : Firm, comparable with 1.2% Agar gel.
- Colour and Clarity of prepared medium** : Light amber coloured, clear to slightly opalescent gel forms in Petri plates.
- Reaction** : Reaction of 3.71% w/v aqueous solution at 25°C. pH:7.0 ± 0.2.
- Cultural Response** : Cultural characteristics observed with added Oxytetra Selective Supplement (FD032) after an incubation at 25-30°C for 2-3 days.

Organisms (ATCC)	Inoculum (CFU)	Growth	Recovery	Colour of colony
# <i>Aspergillus brasiliensis</i> (16404) (00053*)	50-100	luxuriant	—	light blue with black spores
<i>Candida albicans</i> (10231) (00054*)	50-100	luxuriant	≥50%	green
<i>Escherichia coli</i> (25922) (00013*)	≥10 ³	inhibited	0%	—
<i>Saccharomyces cerevisiae</i> (9763) (00058*)	50-100	luxuriant	≥50%	colourless

Key : * : Corresponds to WDCM number
: Formerly known as *Aspergillus niger*

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (6, 7).

References

- Mossel D.A.A. et al, 1970, J. Appl. Bact., 33:454.
- Mossel D.A.A., Harrewijn G.A. and Elzebrock J.M., 1973, UNICEF.
- Mossel D.A.A., Visser M. and Mengerink W.H.J., 1962, Lab. Prac. 11:109.
- American Public Health Association, Standard Methods for the Examination of Dairy Products, 1978, 14th Ed., Washington, D.C
- Wehr H. M. and Frank J. H., 2004, Standard Methods for the Microbiological Examination of Dairy Products, 17th Ed., APHA Inc., Washington, D.C.
- Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
- Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S. and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1

For Identification and Differentiation of *Candida* species

HiCrome™ *Candida* Differential Agar / Base, Modified

Recommended for rapid isolation and identification of *Candida* species from mixed cultures.

M1297A /
M1456A

Composition **	M1297A	M1456A
Ingredients	Grams/Litre	Grams/Litre
Peptone	—	5.00
Peptone special	15.00	—
Yeast extract	4.00	3.00
Malt extract	—	3.00
Dipotassium hydrogen phosphate	1.00	—
Glucose (Dextrose)	—	10.00
Chromogenic mixture	7.22	3.00
Chloramphenicol	0.50	0.05
Agar	15.00	18.00
Final pH (at 25°C)	6.3 ± 0.2	7.2 ± 0.2

** Formula adjusted, standardized to suit performance parameters

Directions

Suspend 42.72 grams of M1297A in 1000 ml distilled water and 21.02 grams of M1456A in 500 ml distilled water. Heat to boiling to dissolve the medium completely. DO NOT AUTOCLAVE. Cool to 45-50°C and aseptically add rehydrated contents of 1 vial of HiCrome™ *Candida* Selective Supplement (FD192) to M1456A. Mix well and pour into sterile Petri plates.

Principle and Interpretation

Perry and Miller (1) reported that *Candida albicans* produces an enzyme β -N-acetyl-galactosaminidase and according to Rousselle et al (2) incorporation of chromogenic or fluorogenic hexosaminidase substrates into the growth medium helps in identification of *C. albicans* isolates directly on primary isolation. HiCrome™ *Candida* Differential Agar is a selective and differential medium, which facilitates rapid isolation of yeasts from mixed cultures and allows differentiation of *Candida* species namely *C. albicans*, *C. krusei*, *C. tropicalis* and *C. glabrata* on the basis of colouration and colony morphology. On this medium, results are obtained within 48 hours and it is useful for the rapid and presumptive identification of common yeasts in Mycology and Clinical Microbiology Laboratory.

Peptone, peptone special, malt extract and yeast extract provides nitrogenous, carbonaceous compounds and other essential growth nutrients. Phosphate buffers the medium well. Chloramphenicol suppresses the accompanying bacterial flora. *C. albicans* appear as light green coloured smooth colonies, *C. tropicalis* appear as blue to purple coloured raised colonies. *C. glabrata* colonies appear as cream to white smooth colonies, while *C. krusei* appear as purple fuzzy

colonies. *C. glabrata*, *C. kefyr*, *C. parapsilosis* colonies appear as cream to white, beige/yellow due to natural pigmentation and some alkaline phosphatase activity. The use of HiCrome Selective Supplement (FD192) in M1456A imparts additional selectivity to the medium.

Type of specimen

Food samples, Clinical samples

Specimen Collection and Handling

For clinical samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (3, 4).

Food samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (5).


After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

In Vitro diagnostic use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets

Limitations

1. Due to variable nutritional requirements, some strains may show poor growth on this medium.
2. Slight colour variation may be observed depending on the presence of enzyme in the organism and substrate utilization provided in the medium.



M1297A HiCrome™ *Candida* Differential Agar

HiCromeVeg™ Freedom from BSE / TSE worries
Single Streak Rapid Differentiation Series

HiCrome™ *Candida* Differential Agar / Base, Modified (M1297A)/(M1456A) is also available as HiCrome™ *Candida* Differential HiVeg™ Agar / Base, Modified (MV1297A)/(MV1456A) / HiCrome™ *Candida* Differential HiCynth™ Agar Base (MCD1297A) wherein all the animal origin nutrients have been replaced by vegetable based nutrients & Chemically defined peptones respectively.

HiCrome™ *Candida* Differential Agar / Base, Modified

Recommended for rapid isolation and identification of *Candida* species from mixed cultures.

M1297A
M1456A

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

- Appearance of Powder** : Cream to beige coloured, homogeneous, free flowing powder.
- Gelling** : Firm, comparable with 1.5% Agar gel of M1297A or 1.8% Agar gel of M1456A.
- Colour and Clarity of prepared medium** : Light amber coloured, clear to slightly opalescent gel forms in Petri plates.
- Reaction** : Reaction of 4.27% w/v aqueous solution of M1297A at 25°C. pH:6.3 ± 0.2
Reaction of 4.20% w/v aqueous solution of M1456A at 25°C. pH:7.2 ± 0.2
- Cultural Response** : Cultural characteristics observed after an incubation at 30°C for 40-48 hours on addition of HiCrome™ *Candida* Selective Supplement (FD192) in M1456A.

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (3, 4).

References

- Perry J. L. and Miller G. R., 1987, J. Clin. Microbiol., 25:2424-2425.
- Rousselle P., Freydiere A., Couillerot P., de Montclos H. and Gille Y., 1994, J.Clin. Microbiol., 32:3034-3036.
- Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
- Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1
- Salfinger Y., and Tortorello M.L. Fifth (Ed.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.

Organisms (ATCC)	Inoculum (CFU)	Growth	Recovery	Colour of the colony
<i>Candida albicans</i> (10231) (00054*)	50-100	good-luxuriant	≥50%	light green
<i>Candida tropicalis</i> (750)	50-100	good-luxuriant	≥50%	blue to purple
<i>Candida krusei</i> (24408)	50-100	good-luxuriant	≥50%	purple, fuzzy
<i>Candida glabrata</i> (15126)	50-100	good-luxuriant	≥50%	cream to white
<i>Candida parapsilosis</i> (22019)	50-100	good-luxuriant	≥50%	cream to white (may have mauve center)
<i>Escherichia coli</i> (25922) (00013*)	≥10 ³	inhibited	0%	—
<i>Staphylococcus aureus subsp aureus</i> (25923) (00034*)	≥10 ³	inhibited	0%	—

Key : * : Corresponds to WDCM number

Ready Prepared Media			
Code	Product Name	Usage	Packing
Category : 90 mm Ready Prepared Petri Plates			
MP1297A	HiCrome™ <i>Candida</i> Differential Agar Plate	for rapid isolation and identification of <i>Candida</i> species from mixed cultures.	50 plts
Category : Ready Prepared Solid Media in Glass Bottles			
SM1297A	HiCrome™ <i>Candida</i> Differential Agar	for rapid isolation and identification of <i>Candida</i> species from mixed cultures.	5X100ml

For Identification and Differentiation of *Candida* species

HiCrome™ Mueller Hinton Agar (For Antifungal testing)

Recommended for the chromogenic differentiation of yeasts from clinical samples and determination of susceptibility to antifungal agents.

M2067

Composition **

Ingredients	Grams/Litre
Acicase #	14.00
Dextrose (Glucose)	20.00
Chromogenic mixture	1.80
Agar	17.00
Final pH (at 25°C)	7.3 ± 0.1

** Formula adjusted, standardized to suit performance parameters

#Equivalent to Casein acid hydrolysate

Directions

Suspend 52.8 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. DO NOT AUTOCLAVE. Cool to 45-50°C. Mix well and pour into sterile Petri plates. The performance of this batch has been tested and standardized as per the current CLSI (formerly NCCLS) document M44-A2 in Method for Antifungal Disk Diffusion susceptibility Testing of yeasts.

Principle and Interpretation

The Mueller Hinton formulation was originally developed as a simple, transparent agar medium for the cultivation of pathogenic species (1). Mueller Hinton Agar, Modified (as per CLSI for antifungal) is recommended for the diffusion of antifungal agents impregnated on paper disc through an agar gel as described in CLSI Approved Standard (2). When supplemented with glucose to a final concentration of 2%, it provides for suitable fungal growth. Kirby-Bauer et al recommended Mueller Hinton Agar for performing antibiotic susceptibility tests using a single disc of high concentration (4). WHO Committee on Standardization of Susceptibility Testing has accepted Mueller Hinton Agar for determining the susceptibility of microorganisms because of its reproducibility (3). Mueller Hinton Agar with 5% sheep blood and Mueller Hinton Agar with Haemoglobin have been recommended for antimicrobial susceptibility testing of *Streptococcus pneumoniae* and *Haemophilus influenzae*. Similarly Mueller Hinton Agar, Modified (as per CLSI for antifungal) is recommended for antifungal susceptibility testing of discs. Perry and Miller (1) reported that *Candida albicans* produces an enzyme β -N-acetyl-galactosaminidase and according to Rousselle et al (2) incorporation of chromogenic or fluorogenic hexosaminidase substrates into the growth medium helps in identification of *C. albicans* isolates directly on primary isolation.

Acicase provide nitrogenous and carbonaceous compounds, long chain amino acids, sulphur and other essential nutrients. Starch acts as a protective colloid against toxic substances present in the medium. Starch hydrolysis yields dextrose, which serves as a source of energy.

Glucose serves as an energy source for fungal cultures. *C. albicans* appear as light green coloured smooth colonies, *C. tropicalis* appear as blue to metallic blue coloured raised colonies. *C. glabrata* colonies appear as cream to white smooth colonies, while *C. krusei* appear as purple fuzzy colonies.

Type of specimen

Clinical samples : Pure cultures isolated from urine , stool, blood etc.


Specimen Collection and Handling

For clinical samples follow appropriate techniques for handling specimens as per established guidelines (5, 6).

After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

In Vitro diagnostic use only. Read the label before opening the container. Wear protective gloves/protective clothing/ eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets.



M2067 HiCrome™ Mueller Hinton Agar (For Antifungal testing) Mixture

HiCromeVeg™ Freedom from BSE / TSE worries
Single Streak Rapid Differentiation Series

HiCrome™ *Candida* Differential Agar / Base, Modified (M1297A)/(M1456A) is also available as HiCrome™ *Candida* Differential HiVeg™ Agar / Base, Modified (MV1297A)/(MV1456A) / HiCrome™ *Candida* Differential HiCynth™ Agar Base (MCD1297A) wherein all the animal origin nutrients have been replaced by vegetable based nutrients & Chemically defined peptones respectively.

HiCrome™ Mueller Hinton Agar (For Antifungal testing)

Recommended for the chromogenic differentiation of yeasts from clinical samples and determination of susceptibility to antifungal agents.

M2067

Limitations

1. This medium is recommended for susceptibility testing of pure cultures only.
2. Inoculum density may effect the zone size. Heavy inoculum may result in smaller zones or too less inoculum may result in bigger zones.
3. Fastidious organisms may not grow on this medium.
4. Certain species of *Candida* may show variation in colour intensity depending on the presence of enzyme.
5. As antimicrobial susceptibility is carried with antibiotic disc, proper storage of the disc is desired which may effect the potency of the disc. Under certain circumstances, the in vitro results of antibiotic susceptibility may not show the same in vivo.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

- Appearance of Powder** : Cream to yellow homogeneous free flowing powder.
- Gelling** : Firm, comparable with 1.5% Agar gel
- Colour and Clarity of prepared medium** : Yellow to amber coloured opalescent gel forms in Petri plates
- Reaction** : Reaction of 5.28% w/v aqueous solution at 25°C. pH : 7.3±0.2
- Cultural Response** : A luxuriant growth of test organisms was observed in 48 hours at 33-37°C along with inhibition zones with respective antibiotic concentrations.

Organisms (ATCC)	Growth	Colour of the colony	Sensitivity testing	
			Fluconazole FLC (25 mcg)	Voriconazole VRC (1 mcg)
<i>Candida albicans</i> (90028)	luxuriant	light green	28-39 mm	31-42mm
<i>Candida parapsilosis</i> 22019	luxuriant	cream	22-33 mm	28-37 mm
<i>Candida tropicalis</i> (750)	luxuriant	blue to purple	26-37mm	-
<i>Candida krusei</i> (6258)	luxuriant	purple	-	16-25mm

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (5, 6).

References

1. Mueller J. H. and Hinton J., 1941, Proc. Soc. Exp. Biol. Med., 48:330.
2. Method for Antifungal Disk Diffusion Susceptibility Testing of yeasts; Approved Guideline Second Edition M44-A2 Vol.24 No.17.
3. Present Status and Future Work, WHO Sponsored collaborative study, Chicago, Oct. 1967.
4. Bauer A. W., Kirby W. M., Sherris J. L. and Turck M., 1966, Am. J. Clin. Pathol., 45:493.
5. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
6. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol

HiCrome™ Candida Differential Agar Base

Recommended as a selective and differential medium for rapid isolation and identification of *Candida* species from mixed cultures.

M1297AR

Composition **

Ingredients	Grams/Litre
Peptone	4.000
Chromogenic mixture	13.600
Agar	13.600

Final pH (at 25°C) 6.0±0.2

** Formula adjusted, standardized to suit performance parameters

Directions

Suspend 15.6 grams in 500 ml distilled water. Add the rehydrated contents of one vial of HiCrome™ Candida Differential Selective Supplement (FD283R). Heat to boiling with frequent agitation to dissolve the medium completely. DO NOT AUTOCLAVE. Cool to 45-50°C. Mix well and pour into sterile Petri plates.

Principle and Interpretation

Perry and Miller (1) reported that *Candida albicans* produces an enzyme β -N-acetyl-galactosaminidase and according to Rousselle et al (2) incorporation of chromogenic or fluorogenic hexosaminidase substrates into the growth medium helps in identification of *C. albicans* isolates directly on primary isolation. HiCrome™ Candida Differential Agar Base incorporates two chromogenes X-NAG which detects the activity of hexosaminidase and BCIP which detects phosphatase activity.

HiCrome™ Candida Differential Agar Base is a selective and differential medium, which facilitates rapid isolation of yeasts from mixed cultures and allows differentiation of *Candida* species namely *C. albicans*, *C. krusei*, *C. tropicalis* and *C. glabrata* on the basis of colouration and colony morphology. On this medium results are obtained within 48 hours and it is useful for the rapid and presumptive identification of common yeasts in Mycology and Clinical Microbiology Laboratory.

Peptone provides nitrogenous, carbonaceous compounds and other essential growth nutrients. Chloramphenicol from the supplement suppresses the accompanying bacterial flora. *C. albicans* appear as light green coloured smooth colonies, *C. tropicalis* appear as blue to metallic blue coloured raised colonies. *C. glabrata*, *C. kefyr*, *C. parapsilosis* colonies appear as cream to white, beige/yellow due to natural pigmentation and some alkaline phosphatase activity, while *C. krusei* appear as pink-purple, fuzzy, dry colonies.

Type of specimen

Clinical, Food samples

Specimen Collection and Handling

For clinical samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (3, 4).

Food samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (5).


After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

In Vitro diagnostic use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets

Limitations

1. Due to variable nutritional requirements, some strains may show poor growth on this medium.
2. Slight colour variation may be observed depending on the presence of enzyme in the organism and substrate utilization provided in the medium.



M1297AR – HiCrome™ Candida Differential Agar

HiCrome™ *Candida* Differential Agar Base

Recommended as a selective and differential medium for rapid isolation and identification of *Candida* species from mixed cultures.

M1297AR

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

- Appearance of Powder** : Cream to beige homogeneous free flowing powder
- Gelling** : Firm, comparable with 1.36% Agar gel
- Colour and Clarity of prepared medium** : Light amber coloured, opaque gel forms in Petri plates
- Reaction** : Reaction of 3.12% w/v aqueous solution at 25°C. pH : 6.0±0.2
- Cultural Response** : Cultural characteristics observed with added HiCrome™ *Candida* Differential Selective Supplement (FD283R) after an incubation at 20-25°C for 40-48 hours.

Organism (ATCC)	Inoculum (CFU)	Growth	Recovery	Colour of Colony
<i>Candida albicans</i> (10231) (00054*)	50-100	good-luxuriant	≥50%	light green
<i>Candida krusei</i> (24408)	50-100	good-luxuriant	≥50%	Purple, fuzzy
<i>Candida tropicalis</i> (750)	50-100	good-luxuriant	≥50%	Blue to purple
<i>Candida kefyr</i> (66028)	50-100	good-luxuriant	≥50%	Cream to white
<i>Candida parapsilosis</i> (22019)	50-100	good-luxuriant	≥50%	Cream to white
<i>Candida glabrata</i> (15126)	50-100	good-luxuriant	≥50%	Cream to white
<i>Escherichia coli</i> (25922) (00013*)	≥10 ³	inhibited	0%	
<i>Escherichia coli</i> (8739) (00012*)	≥10 ³	inhibited	0%	

Key : * : Corresponds to WDCM number

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (3, 4).

References

- Perry J. L. and Miller G. R., 1987, J. Clin. Microbiol., 25: 2424 -2425.
- Rousselle P., Freydiere A., Couillerot P., de Montclos H. and GilleY., 1994, J. Clin. Microbiol. 32:3034-3036.
- Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
- Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1
- Salfinger Y., and Tortorello M.L. Fifth (Ed.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.



M1297AR – HiCrome™ *Candida* Differential Agar

For Identification of Malassezia

HiCrome Malassezia Agar (Twin Pack)

Recommended for isolation, cultivation and identification of *Malassezia furfur*

M1985

Composition **

Ingredients	Grams/Litre
Part A -	
Peptone special	30.00
Chromogenic mixture	1.40
Agar	15.00
Part B -	
Tween 40 (Polysorbate 40)	10.00
Glycerol mono-oleate	5.00

Final pH (at 25°C) 5.80±0.2

** Formula adjusted, standardized to suit performance parameters

Directions

Suspend 15ml of fluid Part B in 1000 ml distilled/purified water. Add 46.4 grams of Part A. Mix well and heat to boiling to dissolve the medium completely. DO NOT AUTOCLAVE. Cool to 45-50°C. Mix well and pour into sterile Petri plates .

Principle and Interpretation

Malassezia is a genus of fungi, naturally found on the skin surfaces of many animals, including humans. Media based on malt extract is appreciated by many microbiologists due to their richness and nutrient balance especially for the cultivation of fastidious microorganisms. With acidic pH, they are used for the isolation, cultivation and maintenance of yeast and moulds.

M. furfur is a lipophilic yeast, therefore in vitro growth must be stimulated by natural oils or other fatty substances. Peptone special provides nitrogenous and carbonaceous compounds, long chain amino acids, vitamins and other essential growth nutrients. Low pH favours fungal growth and inhibits contaminating bacteria from test samples (1). Tween 40, Glycerol monooleate enhances the growth of *Malassezia* species as it is a lipophilic yeast. Some pathogenic fungi may produce infective spores which are easily dispersed in air, so examination should be carried out in safety cabinet. For heavily contaminated samples, the plate must be supplemented with inhibitory agents for inhibiting bacterial growth with lower pH.

Type of specimen

Clinical - skin samples

Specimen Collection and Handling


For clinical samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (2, 3). After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

In Vitro diagnostic use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets

Limitations


1. Due to variable nutritional requirements, some strains may show poor growth on this medium.
2. Slight colour variation may be observed depending upon the utilization of the substrate by the organism.
3. For heavily contaminated samples, the plate must be supplemented with inhibitory agents for inhibiting bacterial growth with lower pH.



M1985 HiCrome Malassezia Agar

HiCrome Malassezia Agar (Twin Pack)

Recommended for isolation, cultivation and identification of *Malassezia furfur*



M1985

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

Appearance of Powder : Part A : Cream to yellow homogeneous free flowing powder Part B: Colourless to pale yellow viscous solution

Gelling : Firm, comparable with 1.5% Agar gel.

Colour and Clarity of prepared medium : Yellow coloured, opalescent gel with scum forms in Petri plates

Reaction : Reaction of 4.64% w/v aqueous solution of Part A and 1.5% v/v of Part B at 25°C. pH : 5.80±0.2

Cultural Response : Cultural characteristics observed after an incubation at 35-37°C for 48-72 hours.

Organism (ATCC)	Inoculum (CFU)	Growth	Recovery	Colour of colony
<i>Malassezia furfur</i> (14521)	50-100	good-luxuriant	≥50%	mauve, small
<i>Candida albicans</i> (10231) (00054*)	50-100	good-luxuriant	≥50%	pale green to green
<i>Candida glabrata</i> (15126)	50-100	good-luxuriant	≥50%	colourless
<i>Candida krusei</i> (24408)	50-100	good-luxuriant	≥50%	purple
<i>Candida tropicalis</i> (750)	50-100	good-luxuriant	≥50%	metallic blue

Key : * : Corresponds to WDCM number

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal


User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (2, 3).

References

1. Murray PR, Baren EJ, Jorgensen JH, Pfaller MA, Tenover FC, Tenover FC (editors) 2003, Manual of clinical Microbiology, 8th ed., ASM, Washington, D.C.
2. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
3. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1

2015
iNNOVATION
IN MICROBIOLOGY MEDIA

M1991I - HiCrome™ Chromogenic Coliform Agar




Rapid & Reliable Method of Detection

HiCrome™ CHROMOGENIC COLIFORM AGAR (M1991I)

- » Recommended by ISO 9308 for water testing
- » Recommended for simultaneous detection of *E.coli* and coliforms in water testing
- » Distinct colours allows easy differentiation
 - *E.coli* : dark blue to violet
 - Other Enterobacteriaceae : pink-pink-red or colourless
 - Gram positive bacteria : inhibited

Also available in Chemically Defined Media : MCD1991I


M1974 - HiCrome™ Rapid MRSA Agar



HiCrome™ RAPID MRSA AGAR BASE (M1974)

- » Selective Detection of Methicillin Resistant Staphylococcus
- » Inhibits Methicillin Sensitive Staphylococcus and other gram negative bacteria
- » Differentiates between Methicillin Resistant Staphylococcus aureus (greenish yellow) and Methicillin Resistant Staphylococcus epidermidis (blue)

M1985 - HiCrome™ Malessezia Agar



HiCrome™ MALESSEZIA AGAR (M1985)

- » Selective media for the detection of Malessezia
- » Can be isolated from clinical and veterinary samples



HiCrome™
Fluorogenic
Media

HiCrome™ PA Broth

Recommended for the detection of presence and absence of coliform bacteria in water.

M1663

Composition **

Ingredients	Grams/Litre
Tryptone	20.000
Lactose	5.000
Bile salts mixture	1.500
Dipotassium hydrogen phosphate	3.000
Potassium dihydrogen phosphate	1.500
Sodium chloride	5.000
2-Nitrophenyl β -D-galactopyranoside (ONPG)	1.250
4-methylumbelliferyl β -D-glucuronide (MUG)	0.100

Final pH (at 25°C) 7.0±0.2

** Formula adjusted, standardized to suit performance parameters

Directions

Suspend 37.35 grams in 1000 ml distilled water. Heat if necessary to dissolve the medium completely. DO NOT AUTOCLAVE. Cool to 45-50°C. Dispense into sterile test tubes or as desired.

Principle and Interpretation

Examination of water for the presence of marker groups such as coliforms is one of the most common tests in food microbiology laboratory, partly because of the relative ease and speed with which these tests can be accomplished. Where it is claimed that water has been processed for safety, the finding of such organism demonstrates a failure of the process (1) HiCrome™ PA Broth is a modification of the medium originally devised by Hajna and Perry (4) and is used for the detection of presence and absence of coliform bacteria in water.

The fluorogenic compound 4-Methylumbelliferyl β -D-glucuronide (MUG) is incorporated in the medium for the fluorogenic detection of *Escherichia coli*, the main indicator organism for the faecal contamination of water. The enzyme β -glucuronidase possessed by *Escherichia coli* hydrolyses MUG to yield a fluorescent end product 4-Methylumbelliferone; which can be detected when the medium is observed for fluorescence under UV light (3,7) MUG also detects anaerogenic strains which may not be detected in the conventional procedure (3). ONPG test is used to determine the presence or absence of β -galactosidase in organisms (6) and is also important in differentiating *Enterobacteriaceae* which are commonly classified according to their ability to ferment lactose. ONPG is similar in structure to lactose. The presence of two enzymes, permease and β -D-galactosidase are required to demonstrate lactose fermentation. True lactose non fermenters do not possess either of these enzymes. Late lactose fermenting organisms do not have permease but do possess β -galactosidase. If β -galactosidase is present, the colourless ONPG is split into galactose and o-nitrophenol, a yellow compound (5).

Tryptone provides nitrogenous, carbonaceous compounds and other essential nutrients. Lactose is the fermentable carbohydrate, sodium chloride maintains osmotic equilibrium. The medium has a strong buffering system to control the pH in the presence of fermentative action. Bile salts mixture inhibit gram-positive bacteria especially *Bacillus* species and faecal Streptococci. Mostly β -glucuronidase activity occurs within 4 hours but some weakly β -glucuronidase positive strains require overnight incubation (2).

Type of specimen

Water, waste water, sewage samples and Clinical samples

Specimen Collection and Handling


For water samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (2).

For clinical samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (5, 8).

After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

In Vitro diagnostic use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets



M1663 — HiCrome™ PA Broth

- | | | |
|---------------------------------|----------------------------------|----------------------------------|
| 1. Control | 2. <i>Escherichia coli</i> | 3. <i>Enterobacter aerogenes</i> |
| 4. <i>Klebsiella pneumoniae</i> | 5. <i>Salmonella Typhimurium</i> | 6. <i>Proteus mirabilis</i> |
| 7. <i>Staphylococcus aureus</i> | 8. <i>Enterococcus faecalis</i> | |

HiCrome™ PA Broth

Recommended for the detection of presence and absence of coliform bacteria in water.

M1663

Limitations

1. Due to variable nutritional requirements, some strains may show poor growth in this medium.
2. Slight colour variation may be observed depending upon the utilization of the substrate by the organism.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

- Appearance of Powder** : Cream to yellow homogeneous free flowing powder
- Colour and Clarity of prepared medium** : Light amber coloured, clear solution without any precipitate
- Reaction** : Reaction of 3.7% w/v aqueous solution at 25°C. pH : 7.0±0.2
- Cultural Response** : Cultural characteristics observed after an incubation at 35-37°C for 18-24 hours.

Organism (ATCC)	Inoculum (CFU)	Growth	ONPG	Fluorescence at 366nm
<i>Escherichia coli</i> (25922) (00013*)	50-100	luxuriant	positive reaction, yellow colour	positive, throughout the tube
# <i>Klebsiella aerogenes</i> (13048) (00075*)	50-100	luxuriant	positive reaction, yellow colour	negative
<i>Klebsiella pneumoniae</i> (13883) (00097*)	50-100	luxuriant	positive reaction, yellow colour	negative
<i>Proteus mirabilis</i> (25933)	50-100	luxuriant	negative reaction, yellow colour	negative
<i>Salmonella</i> Typhimurium (14028) (00031*)	50-100	luxuriant	negative reaction, no yellow colour or colourless	negative
<i>Staphylococcus aureus subsp aureus</i> (25923) (00034*)	≥10 ³	inhibited	negative reaction, no yellow colour or colourless	negative
<i>Enterococcus faecalis</i> (29212) (00087*)	≥10 ³	inhibited	negative reaction, no yellow colour	negative

Key : * : Corresponds to WDCM number

: Formerly known as *Enterobacter aerogenes*

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (5, 8).

References

1. Corry J. E. L., Curtis G. D. W., and Baird R. M., Culture Media For Food Microbiology, Vol. 34, Progress in industrial Microbiology, 1995, Elsevier, Amsterdam.
2. Baird R.B., Eaton A.D., and Rice E.W., (Eds.), 2015, Standard Methods for the Examination of Water and Wastewater, 23rd ed., APHA, Washington, D.C.
3. Feng P. C. S and Hartman P. A. S., 1982, Appl. Environ. Microbiol., 43 :132.
4. Hajna A. A. and Perry C. A., 1943, Am. J. Public Health, 33:550.
5. Isenberg H. D., (Eds.), Clinical Microbiology Procedures Handbook, Vol. 1, Washington D.C. American Society for Microbiology; 1992, p.l. 19.20-1.19.22.
6. MacFaddin J. F, 2000, Biochemical Tests for Identification of Medical Bacteria. 3rd Ed. Philadelphia, Lippincott Williams and Wilkins, p. 160-9.
7. Robinson B. J., 1984, Appl. Environ. Microbiol., 48:285.
8. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1

For Identification *E. coli* and Total coliforms

HiColiform™ Broth, Modified

Recommended for the detection and confirmation of *Escherichia coli* and total coliforms from water samples, using a combination of chromogenic and fluorogenic substrates.

M1850

Composition **

Ingredients	Grams/Litre
Peptone	5.000
Sodium chloride	5.000
Potassium sulfate	1.000
Dipotassium hydrogen phosphate	4.000
Potassium dihydrogen phosphate	1.000
Sodium lauryl sulphate (SLS)	0.100
Sodium puruvate	1.000
Chromogenic substrate	0.100
Fluorogenic substrate	0.100
IPTG	0.100

Final pH (at 25°C) 6.8±0.2

** Formula adjusted, standardized to suit performance parameters

Directions

Suspend 17.4 grams in 1000 ml distilled water. Heat if necessary to dissolve the medium completely. Mix well and dispense in tubes or flasks or as desired. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes.

Principle and Interpretation

HiColiform™ Broth, Modified was designed for detection and confirmation of *Escherichia coli* and total coliforms from water samples using a combination of chromogenic and fluorogenic substrates. *Escherichia coli* can be distinguished from other coliforms by its unique ability to fluoresce in the presence of fluorogenic substrate (1, 2). The fluorogenic substrate is split by enzyme β -glucuronidase especially present in *Escherichia coli*. The reaction is indicated by the development of a blue fluorescence under UV light. The presence of total coliforms is indicated by blue-green colourations due to the cleavage of the chromogenic substrate. IPTG amplifies enzyme synthesis and increases the activity of β -galactosidase.

Peptone provides essential growth nutrients and is useful for the simultaneous detection of indole production. The phosphate salts provide buffering action for rapid growth of coliforms. Sodium chloride helps to maintain the osmotic balance. Sodium lauryl sulphate makes the medium selective by inhibiting accompanying microflora, especially the gram-positive organisms.

Type of specimen

Water samples

Specimen Collection and Handling

For water samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (3).


After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets

Limitations

1. β -glucuronidase especially present in 97% of *Escherichia coli*, however few *E. coli* may be negative.
2. Some species may show poor growth due to nutritional variations.



1 – HiColiform™ Broth, Modified
1. *Escherichia coli* 2. *Enterobacter aerogenes* 3. Control

HiColiform™ Broth, Modified

Recommended for the detection and confirmation of *Escherichia coli* and total coliforms from water samples, using a combination of chromogenic and fluorogenic substrates.

M1850

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

- Appearance of Powder** : Cream to yellow homogeneous free flowing powder
- Colour and Clarity of prepared medium** : Light yellow coloured, clear to slightly opalescent solution in tubes
- Reaction** : Reaction of 1.74% w/v aqueous solution at 25°C. pH : 6.8±0.2
- Cultural Response** : Cultural characteristics observed after an incubation at 35-37°C for 18-24 hours.

Organism (ATCC)	Inoculum	Growth	Colour of medium	Fluorescence (under uv)
# <i>Klebsiella aerogenes</i> (13048) (00075*)	50-100	luxuriant	blue-green	negative reaction
<i>Escherichia coli</i> (25922) (00013*)	50-100	luxuriant	blue-green	positive reaction

Key : * : Corresponds to WDCM number

: Formerly known as *Enterobacter aerogenes*

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition. Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (3, 4).

References

1. Feng P.C.S. and Hartman P.A., 1982, J. Appl. Environmental Microbiol. 43. 1320-1323.
2. Harsen W., and Yourassowsky, 1984, J. Clin. Microbiol. 20. 1177-1179.
3. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
4. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S. and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.
5. Baird R.B., Eaton A.D., and Rice E.W., (Eds.), 2015, Standard Methods for the Examination of Water and Wastewater, 23rd ed., APHA, Washington, D.C.

For Fluorogenic Identification and Differentiation of *E. coli* and Coliforms

Rapid HiColiform™ Agar / Broth

For detection and confirmation of *Escherichia coli* and total coliforms on the basis of enzyme substrate reaction from water samples, using a combination of chromogenic and fluorogenic substrates.

M1465 /
M1453

Composition **	M1465	M1453
Ingredients	Grams/Litre	Grams/Litre
Peptone, special	5.00	5.00
Sodium chloride	5.00	5.00
Sorbitol	1.00	1.00
Dipotassium hydrogen phosphate	2.70	2.70
Potassium dihydrogen phosphate	2.00	2.00
Sodium lauryl sulphate (SLS)	0.10	0.10
Chromogenic substrate	0.08	0.08
Fluorogenic substrate	0.05	0.05
IPTG (Isopropyl-β-D-thiogalactopyranoside)	0.10	0.10
Agar	15.00	—

Final pH (at 25°C) 6.8 ± 0.2

** Formula adjusted, standardized to suit performance parameters

Directions

Suspend 31.03 grams of M1465 and 16.03 grams of M1453 in 1000 ml distilled water. For double strength broth use 32.06 grams of M1453 in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Mix well and dispense as desired (for M1453) or pour into sterile Petri plates (for M1465).

Principle and Interpretation

The Rapid HiColiform™ Agar is modification of LMX Broth described by Manafi and Kneifel (2). These media are useful for the detection and confirmation of *Escherichia coli* and total coliforms in water samples on the basis of chromogenic and fluorogenic substrates (1-6).

Peptone special, which is rich in tryptophan provides essential growth nutrients and is useful for the simultaneous detection of indole production. Sorbitol provides the carbon source. The phosphate salts provide buffering action for rapid growth of coliforms. Sodium lauryl sulphate makes the medium selective by inhibiting accompanying microflora, especially the gram-positive organisms. The fluorogenic

substrate is split by enzyme β-D-glucuronidase, which is specifically found in *E. coli*. The reaction is indicated by a blue fluorescence under UV light. The presence of total coliforms is indicated by a blue-green colour of the broth due to cleavage of chromogenic substrate. IPTG, a highly stable synthetic analog of lactose induces synthesis of β-D-glucuronidase. In agar medium, 2-3 drops of Kovac's reagent is added over the suspected colonies. Change in the colour of colony to red confirms *E. coli*. Broth medium is overlaid with Kovac's reagent and formation of red ring confirms *E. coli*. If fluorescence is negative after 24 hours of incubation, continue incubation for another 24 hours without performing the indole test.

Type of specimen

Water samples


Specimen Collection and Handling

For water samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (7).

After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets



Other coliforms

E. coli under uv

M1465 Rapid HiColiform™ Agar

HiCromeVeg™ Freedom from BSE / TSE worries
Single Streak Rapid Differentiation Series

Rapid HiColiform™ Agar / Broth (M1465 / M1453) is also available as Rapid HiColiform™ HiVeg™ Agar / Broth (MV1465 / MV1453) wherein all the animal origin nutrients have been replaced by vegetable based nutrients.

Rapid HiColiform™ Agar / Broth

For detection and confirmation of *Escherichia coli* and total coliforms on the basis of enzyme substrate reaction from water samples, using a combination of chromogenic and fluorogenic substrates.

M1465 /
M1453

Limitations

1. β -glucuronidase especially present in 97% of *Escherichia coli*, however few *E. coli* may be negative.
2. Slight colour variation may be observed depending upon the utilization of the substrate by the organism.
3. Further confirmation of *E.coli* must be confirmed by addition of Kovacs reagent to the suspected colony after visualization of fluorescence

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

- Appearance of powder** : Cream to yellow coloured, homogeneous, free flowing powder.
- Gelling** : Firm, comparable with 1.5% Agar gel of M1465
- Colour and Clarity of prepared medium** : Light yellow coloured, clear to slightly opalescent gel forms in Petri plates (M1465)/ clear solution having slight precipitate in tubes (M1453).
- Reaction** : Reaction of 3.1% w/v of M1465 or 1.6 % w/v of M1453 aqueous solution at 25°C. pH:6.8 ± 0.2
- Cultural Response** : M1453 : Cultural characteristics observed after an incubation at 35-37°C for 18-24 hours.

Organisms (ATCC)	Inoculum (CFU)	Growth	Colour change in medium	Fluorescence (Under UV light)	Indole reaction
<i>Escherichia coli</i> (25922) (00013*)	50-100	luxuriant	blue-green	+	+
# <i>Klebsiella aerogenes</i> (13048) (00075*)	50-100	luxuriant	blue-green	-	-
M1465 : Cultural characteristics observed after an incubation at 35-37°C for 18-24 hours.					
Organisms (ATCC)	Inoculum (CFU)	Growth	Colour change in medium	Fluorescence (Under UV light)	Indole reaction
<i>Escherichia coli</i> (25922) (00013*)	50-100	luxuriant	blue-green	+	+
# <i>Klebsiella aerogenes</i> (13048) (00075*)	50-100	luxuriant	blue-green	-	-
<i>Klebsiella pneumoniae</i> (13883) (00097*)	50-100	luxuriant	blue-green	-	-
<i>Salmonella Typhimurium</i> (14028) (00031*)	50-100	luxuriant	yellow	-	-

Key: * : Corresponds to WDCM number
 : + : positive reaction, - : negative reaction.
 # : Formerly known as *Enterobacter aerogenes*

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (8, 9).

References

1. Hahn G., and Wittrock E., 1991, Acta Microbiologica Hungarica, 38 (3-4): 265.
2. Manafi M., and Kneifel W., 1989, Zbl. Hygiene and Umweltmedizin, 189:225.
3. Manafi M., 1990, Forum Stadte-Hygiene, 41:181-184.
4. Manafi M., 1991, Ernährung / Nutrition, 15, Nr. 10.
5. Manafi M., Kneifel W., 1991, Acta Microbiologica Hungarica, 33(3-4):293-304
6. Manafi M., Kneifel B., Bascon S., 1991, Microbiol. Rev. 55:335-348.
7. Baird R.B., Eaton A.D., and Rice E.W., (Eds.), 2015, Standard Methods for the Examination of Water and Wastewater, 23rd ed., APHA, Washington, D.C
8. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
9. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1

For Fluorogenic Identification of *E.coli*

HiCrome™ ECD Agar w/ MUG

Recommended for the detection of *Escherichia coli* in water and food samples by using a combination of chromogenic and fluorogenic substrate.

M1488

Composition **

Ingredients	Grams/Litre
Tryptone	20.00
Bile salts mixture	1.50
L-Tryptophan	1.00
Lactose	5.00
Sodium chloride	5.00
Dipotassium hydrogen phosphate	4.00
Potassium dihydrogen phosphate	1.50
Fluorogenic substrate	0.07
Chromogenic substrate	0.10
Agar	15.00

Final pH (at 25°C) 7.0 ± 0.2

** Formula adjusted, standardized to suit performance parameters

Directions

Suspend 53.17 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Mix well and pour into sterile Petri plates.

Principle and Interpretation

HiCrome™ ECD Agar w/ MUG is recommended for rapid detection of *Escherichia coli* by using a combination of chromogenic and fluorogenic substrates. The presence of *Escherichia coli* is indicated by blue coloured colony formation due to cleavage of chromogenic substrate. Fluorogenic substrate permits rapid detection of *Escherichia coli* when medium is observed for fluorescence using UV light (1, 2). Fluorogenic substrate also detects anaerogenic strains, which may not be detected in conventional procedure (1). It is hydrolysed by enzyme β -D-glucuronidase, possessed by *Escherichia coli* to yield a fluorescent end product. The reaction is indicated by a blue fluorescence under UV light. Tryptone provides nitrogenous, carbonaceous compounds, long chain amino acids and other essential nutrients. Lactose is the fermentable carbohydrate. Sodium chloride maintains osmotic equilibrium. The medium has a strong buffering system to control the pH in the presence of fermentive action. The bile salt mixture inhibits gram-positive bacteria especially *Bacillus* species and faecal Streptococci.

Type of specimen

Food samples, Water samples

Specimen Collection and Handling


For food samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (3).

For water samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (4).

After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

Read the label before opening the container. Wear protective gloves/ protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets



M1488 HiCrome ECD Agar w/ MUG

HiCromeVeg™ Freedom from BSE / TSE worries
Single Streak Rapid Differentiation Series

HiCrome™ ECD Agar w/ MUG (M1488) is also available as HiCrome™ ECD HiVeg™ Agar w/ MUG (MV1488) wherein all the animal origin nutrients have been replaced by vegetable base nutrients.

HiCrome™ ECD Agar w/ MUG

Recommended for the detection of *Escherichia coli* in water and food samples by using a combination of chromogenic and fluorogenic substrate.

M1488

Limitations

1. β -glucuronidase especially present in 97% of *Escherichia coli*, however few *E. coli* may be negative.
2. Slight colour variation may be observed depending upon the utilization of the substrate by the organism.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

- Appearance of Powder** : Cream to yellow coloured, homogeneous, free flowing powder.
- Gelling** : Firm, comparable with 1.5% Agar gel.
- Colour and Clarity of prepared medium** : Light amber coloured, clear gel forms in Petri plates.
- Reaction** : Reaction of 5.32% w/v aqueous solution at 25°C. pH : 7.0 ± 0.2
- Cultural Response** : Cultural characteristics observed after an incubation at 44-45°C for 18-24 hours.

Organisms (ATCC)	Inoculum (CFU)	Growth	Recovery	Cololur of colony	Fluorescence under UV	Indole
<i>Escherichia coli</i> (25922) (00013*)	50-100	good	40-50%	bluish green	+	+
<i>Klebsiella pneumoniae</i> (13883) (00097*)	50-100	good	40-50%	colourless	-	-
<i>Staphylococcus aureus subsp aureus</i> (25923) (00034*)	≥10 ³	inhibited	0%	-	-	-
<i>Enterococcus faecalis</i> (29212) (00087*)	≥10 ³	inhibited	0%	-	-	-
<i>Pseudomonas aeruginosa</i> (27853) (00025*)	50-100	good	40-50%	colourless	-	-

Key : * : Corresponds to WDCM number

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (5, 6).

References

1. Feng PCS and Hartman PAS, (1982), Appl. Environ. Microbiol. 43:132.
2. Robinson (1984), Appl. Environ. Microbiol., 48:285.
3. Salfinger Y, and Tortorello M.L. Fifth (Ed.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.
4. Baird R.B., Eaton A.D., and Rice E.W., (Eds.), 2015, Standard Methods for the Examination of Water and Wastewater, 23rd ed., APHA, Washington, D.C
5. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
6. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1

For Fluorogenic Identification of *Pseudomonas*

HiFluoro™ *Pseudomonas* Agar Base

Recommended for selective isolation of *Pseudomonas aeruginosa* from clinical and non-clinical specimens by fluorogenic method.

M1469

Composition **

Ingredients	Grams/Litre
Gelatin peptone	18.00
Magnesium chloride	1.40
Potassium sulphate	10.00
Cetrimide	0.30
Fluorogenic mixture	2.05
Agar	15.00

Final pH (at 25°C) 7.2 ± 0.2

** Formula adjusted, standardized to suit performance parameters

Directions

Suspend 46.75 grams in 1000 ml distilled water containing 10 ml glycerol. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Mix well and pour into sterile Petri plates.

Principle and Interpretation

Pseudomonas aeruginosa (also known as *Pseudomonas pyocyanea*) is a gram-negative, aerobic, rod-shaped bacterium. Like other *Pseudomonads*, *P. aeruginosa* secretes a variety of pigments, including pyocyanin (blue-green), fluorescein (yellow-green and fluorescent), and pyorubin (red-brown). King et al developed *Pseudomonas* Agar P (i.e. King A media) for enhancing pyocyanin and pyorubin production and *Pseudomonas* Agar F (i.e. King B media) for enhancing fluorescein production (1). HiFluoro *Pseudomonas* Agar Base is devised based on the formula described by King et al. (1) except fluorogenic mixture. It is used as the selective medium for the isolation of *P. aeruginosa* from pus, sputum and drains etc.

Cetrimide (Cetyltrimethylammonium bromide) is incorporated in the medium to inhibit bacteria other than *P. aeruginosa*. It acts as a quaternary ammonium compound, cationic detergent that causes nitrogen and phosphorus to be released from bacterial cells other than *P. aeruginosa*. *P. aeruginosa* cleaves the fluorogenic compound to release the fluorogen which produces a visible fluorescence under long wave UV light.

Type of specimen

Clinical samples - Pus, Sputum, Drains ; Water samples

Specimen Collection and Handling


For Clinical samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (2, 3).

For water samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (4).

After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

In Vitro diagnostic use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets



P. aeruginosa

M1469 HiFluoro™ *Pseudomonas* Agar Base

HiCromeVeg™ Freedom from BSE / TSE worries
Single Streak Rapid Differentiation Series

HiFluoro™ *Pseudomonas* Agar Base (M1469) is also available as HiFluoro™ *Pseudomonas* HiVeg™ Agar Base (MV1469) wherein all the animal origin nutrients have been replaced by vegetable based nutrients.

HiFluoro™ *Pseudomonas* Agar Base

Recommended for selective isolation of *Pseudomonas aeruginosa* from clinical and non-clinical specimens by fluorogenic method.

M1469

Limitations

1. Due to variable nutritional requirements, some strains may show poor growth on this medium.
2. Slight colour variation may be observed depending upon strains

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the recommended temperature.

Quality Control

- Appearance of powder** : Cream to yellow coloured, homogeneous, free flowing powder.
- Gelling** : Firm, comparable with 1.5% Agar gel.
- Colour and Clarity of prepared medium** : Light amber coloured, opalescent gel with slight precipitate forms in Petri plates.
- Reaction** : Reaction of 4.67% w/v aqueous solution at 25°C. pH:7.2 ± 0.2
- Cultural Response** : Cultural characteristics observed after an incubation at 35-37°C for 24-48 hours.

Organisms (ATCC)	Inoculum (CFU)	Growth	Recovery	Fluorescence
<i>Pseudomonas aeruginosa</i> (27853) (00025*)	50-100	good-luxuriant	≥50%	+
<i>Escherichia coli</i> (25922) (00013*)	≥10 ³	inhibited	0%	-
<i>Stenotrophomonas maltophilia</i> (13637)	≥10 ³	inhibited	0%	-
<i>Staphylococcus aureus subsp aureus</i> (25923) (00034*)	≥10 ³	inhibited	0%	-

Key : * Corresponds to WDCM number

Storage and Shelf-life

Store between 2-8°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (2, 3).

References

1. King, Ward and Raney, 1954, J. Lab. Clin. Med., 44:301.
2. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
3. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S. and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1
4. Baird R.B., Eaton A.D., and Rice E.W., (Eds.), 2015, Standard Methods for the Examination of Water and Wastewater, 23rd ed., APHA, Washington, D.C

For detection of *Escherichia coli*

HiCrome™ EC Broth w/ RUG

Recommended for detection of *Escherichia coli* in water and food samples by a chromogenic and fluorogenic method

M2073

Composition **

Ingredients	Grams/Litre
Yeast extract	2.00
Acicase#	1.00
Buffers	4.50
Sodium chloride	0.50
Salts	2.45
Chromogenic mixture	0.112

Final pH (at 25°C) 7.0 ± 0.2

** Formula adjusted, standardized to suit performance parameters

Equivalent to Casein acid hydrolysate

Directions

Suspend 10.56 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. DO NOT AUTOCLAVE. Cool to 45-50°C. Dispense into sterile tubes or flasks as desired.

Principle and Interpretation

Escherichia coli is a member of faecal coliform group of bacteria. It is a member of the indigenous faecal flora of warmblooded animals. *E. coli* is considered a specific indicator of faecal contamination and the possible presence of enteric pathogens. *E. coli* can be reliably detected with media that contain a chromogenic or fluorogenic substrate for beta glucuronidase, an enzyme that occurs almost exclusively in *E. coli*.

▲ Resorufin-beta-D-glucuronic acid methyl ester (RUG) is a highly sensitive chromogenic and fluorogenic indicator for *E. coli*. In contrast to MUG, RUG is more specific and does not require fluorescent detection. The released dye Resorufin itself gives intense pink color which can be visually detected. Additional confirmation can be done by observation of fluorescence under uv light.

Yeast extract and Acicase provides carbonaceous, nitrogenous substances, long chain amino acids, vitamins and other essential nutrients. Sodium chloride maintains osmotic equilibrium. The medium has a strong buffering system to control the pH in the medium. Gram-positive bacteria especially *Staphylococcus*, *Bacillus* species and faecal *Streptococcus* are inhibited.

Type of specimen

Food sample; Water samples

Specimen Collection and Handling


For food samples, follow appropriate techniques for sample collection and processing as per guidelines (1).

For water samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (2).

After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

Read the label before opening the container. Wear protective gloves/ protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in




M2073 HiCrome™ EC Broth w/ RUG

1. *E. coli* ATCC 25922
2. *P. aeruginosa* ATCC 27853
3. *S. Typhimurium* ATCC 14028
4. Control

▲ Resorufin-beta-D-glucuronic acid methyl ester (RUG) is a patent of BIOSYNTH

HiCrome™ EC Broth w/ RUG

Recommended for detection of *Escherichia coli* in water and food samples by a chromogenic and fluorogenic method



M2073

individual safety data sheets.

Limitations

1. Slight variation in intensity of colour may be observed depending on the isolates

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

Quality Control

- Appearance of powder** : Pale yellow to orange homogeneous free flowing powder.
- Colour and Clarity of prepared medium** : Yellow coloured clear solution without any precipitate.
- Reaction** : Reaction of 1.05% w/v aqueous solution at 25°C. pH : 7.0±0.2
- Cultural Response** : Cultural characteristics observed after an incubation at 35 - 37°C for 18 - 24 hours.

Organisms (ATCC)	Inoculum (CFU)	Growth	Colour of medium	*Fluorescence (at 366 nm)
<i>Escherichia coli</i> ATCC 25922 (00013*)	50-100	luxuriant	bright pink	positive, throughout the tube
<i>Escherichia coli</i> ATCC 10536	50-100	luxuriant	bright pink	positive, throughout the tube
<i>Citrobacter freundii</i> ATCC 8090	50-100	luxuriant	pale yellow	negative
<i>Salmonella</i> Enteritidis ATCC 10376 (00030*)	50-100	luxuriant	pale yellow	negative
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> ATCC 25923 (00034*)	≥10 ³	inhibited	no colour change	
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> ATCC 6633 (00003*)	≥10 ³	inhibited	no colour change	

Key : * : Corresponds to WDCM number

Storage and Shelf-life

Store dehydrated powder and prepared medium on receipt at 2-8°C. Use before expiry period on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label. Product performance is best if used within stated expiry period

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (3, 4).

References

1. Salfinger Y., and Tortorello M.L. Fifth (Ed.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.
2. Baird R.B., Eaton A.D., and Rice E.W., (Eds.), 2015, Standard Methods for the Examination of Water and Wastewater, 23rd ed., APHA, Washington, D.C.
3. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
4. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

Supplement

Bacillus Selective Supplement

FD324

An antibiotic supplement recommended for the selective isolation of *Bacillus* species.

Formula

(Per vial sufficient for 1000 ml medium)

Polymyxin sulphate	10 mg
Bacitracin	10 mg

Directions

Rehydrate the contents of one vial aseptically with 10 ml sterile distilled water. Mix well and aseptically add it to 1000 ml sterile, molten HiCrome™ Bacillus Agar Base (M1651) / HiCrome™ Bacillus HiCynth™ Agar Base (MCD1651). Mix well and pour into sterile Petri plates / tubes.

Cefoxitin Supplement

FD259

An antimicrobial supplement recommended for the selective isolation of Methicillin Resistant *Staphylococcus aureus* from clinical specimens.

Formula

(Per vial sufficient for 500 ml medium)

Cefoxitin	3.000 mg
-----------	----------

Directions

Rehydrate the contents of 1 vial with 5 ml of sterile distilled water and aseptically add to 500 ml of sterile, molten, cooled (45-50°C) MeReSa Agar Base (M1594). HiCrome™ MeReSa Agar Base (M1674) / HiCrome™ MeReSa HiVeg Agar Base (MV1674) / HiCrome™ MeReSa HiCynth™ Agar Base (MCD1674) / HiCrome™ MeReSa Agar Base, Modified (M1953). This supplement can either be used individually or in combination with MeReSa Selective Supplement (FD229) for more selectivity. Mix well and pour into sterile Petri plates.

Chromogenic Supplement

FD270

Chromogenic Supplement is recommended for the enumeration of faecal coliform by membrane filter technique.

Formula

(Per vial sufficient for 1000 ml medium)

Chromogenic Substrate	0.100g
-----------------------	--------

Directions

Rehydrate the contents of one vial aseptically with 5 ml sterile distilled water. Mix well and aseptically add it to 1000 ml of sterile, molten, cooled (45-50°C) MFC Basal Medium (M1812). Mix well and pour into sterile Petri plates.

Ciprofloxacin Supplement

FD345

Recommended for the selective isolation of *Lactobacillus* species.

Formula

(Per vial sufficient for 1000 ml medium)

Ciprofloxacin	2 mg
---------------	------

Directions

Rehydrate the contents of one vial aseptically with 5 ml of distilled water. Mix well and aseptically add to 1000 ml of sterile molten cooled (45-50°C) HiCrome™ Lacto*Bacillus* Selective Agar Base (M2065). Mix well and pour into sterile Petri plates.

ECC Selective Supplement, Modified

FD344

Recommended for the detection of *E. coli* and other *Enterobacteriaceae* in water samples.

Formula

(Per vial sufficient for 1000 ml medium)

Erythromycin	2 mg
--------------	------

Directions

Rehydrate the contents of one vial aseptically with 5 ml of distilled water. Mix well and aseptically add to 1000 ml of sterile molten cooled (45-50°C) HiCrome™ Coliconfirm Broth Base (M2064). Mix well and pour into sterile Petri plates.

Supplement

Egg Yolk Tellurite Emulsion (50 ml/100 ml per vial)

FD046

Sterile stabilized tellurite emulsion of egg yolk recommended for identification of *Staphylococcus species*.

Formula

	(100ml per vial)	(50ml per vial)
Egg yolk	30.000 ml	15.000 ml
Sterile saline	64.000 ml	32.000 ml
Sterile 3.5% potassium tellurite solution	6.000 ml	3.000 ml

Directions

Warm up the refrigerated Egg Yolk Tellurite Emulsion to room temperature. Shake well to attain uniform emulsion (since on refrigeration emulsion has a tendency to form layers or small lumps). Aseptically add 50 ml in 950 ml of sterile, molten, cooled (45-50°C) Baird Parker Agar Base (M043 / M043S) / Baird Parker HiVeg™ Agar Base (MV043) / Baird Parker Agar Base, Granulated (GM043)/ Baird Parker HiCynth™ Agar Base (MCD043) / Baird Parker Agar Base w/ Sulpha (M1140) / HiCrome™ Aureus Agar Base (M1468). Mix well and pour into sterile Petri plates.

Enterococcus faecium Selective Supplement

FD226

Recommended to differentiate *Enterococcus faecium* from *Enterococcus faecalis*.

Formula

(Per vial sufficient for 500 ml medium)	
Cephalexin	25.000 mg
Aztreonam	37.500 mg

Directions

Rehydrate the contents of 1 vial aseptically with 5 ml sterile distilled water. Mix well and aseptically add to 500 ml of sterile, molten, cooled (45-50°C) Arabinose Agar Base (M1576) / HiCrome™ Enterococcus faecium Agar Base (M1580) / HiCrome™ Enterococcus faecium HiVeg™ Agar Base (MV1580). Mix well and pour into sterile Petri plates.

HiCrome™ Candida Differential Selective Supplement

FD283R

An antibiotic supplement recommended for rapid and direct isolation and identification of *Candida species* from mixed cultures.

Formula

(per vial, sufficient for 500 ml medium)	
Chloramphenicol	250.00 mg

Directions

Rehydrate the contents of 1 vial aseptically with 2 ml of 95% ethanol. Mix well and aseptically add to 500 ml of sterile, molten cooled (45-50°C) HiCrome™ Candida Differential Agar Base (M1297AR). Mix well and pour into sterile Petri plates.

HiCrome™ Candida Selective Supplement

FD192

An antibiotic supplement recommended for the selective isolation of *Candida species* from mixed cultures.

Formula

(Per vial sufficient for 500 ml medium)	
Gentamicin	50.00 mg

Directions

Rehydrate the content of 1 vial with 5 ml of sterile distilled water. Mix well and aseptically add to 500 ml of sterile, molten, cooled (45-50°C) HiCrome™ Candida Differential Agar Base, Modified (M1456A) / HiCrome™ Candida Differential HiVeg™ Agar Base, Modified (MV1456A). Mix well and pour into sterile Petri plates.

Supplement

HiCrome™ ECC Selective Supplement

FD190

An antibiotic supplement recommended for the selective isolation of *Escherichia coli* and coliforms from water and food samples.

Formula

(per vial sufficient for 1000 ml medium)

Cefsulodin	10.00 mg
------------	----------

Directions

Rehydrate the contents of one vial aseptically with 2 ml sterile distilled water. Mix well and aseptically add to 1000 ml of sterile, molten, cooled (45-50°C) HiCrome™ ECC Selective Agar Base (M1294) / HiCrome™ ECC Selective Agar Base, HiVeg™ (MV1294). Mix well and pour into sterile Petri plates.

HiCrome™ EC 0157:H7 Selective Supplement

FD187

Recommended for selective isolation and easy detection of *Escherichia coli* O157:H7 from food samples.

Formula

(per vial sufficient for 1000 ml medium)

Novobiocin	10.00 mg
Potassium tellurite	1.00 mg

Directions

Rehydrate the contents of 1 vial aseptically with 10 ml of sterile distilled water. Mix well and add aseptically to 990 ml sterile, molten, cooled (45-50°C) HiCrome™ EC 0157:H7 Selective Agar Base (M1575) / HiCrome™ EC 0157:H7 Selective HiVeg™ Agar Base (MV1575). Mix well and pour into sterile Petri plates.

HiCrome™ EC 0157: H7 Selective Supplement I

FD230

An antimicrobial supplement recommended for isolation and easy detection of *Escherichia coli* O157:H7 from food and environmental samples.

Formula

(Per vial sufficient for 500 ml medium)

Novobiocin	15.00 mg
Potassium tellurite	1.50mg
Distilled water	5.00 ml

Directions

Warm up refrigerated contents to 45-50°C and aseptically add one vial to 495 ml of sterile, cooled (45-50°C) HiCrome™ Enrichment Broth Base for EC 0157:H7 (M1598). Mix well and dispense in sterile test tubes.

HiCrome™ EC 0157:H7 Selective Supplement, Modified

FD295

A selective supplement recommended for presumptive enumeration of *Escherichia coli* O157:H7 by membrane filtration technique.

Formula

(per vial, sufficient for 1000 ml medium)

Monensin	0.038 gm
Novobiocin	0.0075 gm

Directions

Rehydrate the contents of one vial aseptically with 5 ml sterile distilled water. Mix well and aseptically add it to 1000 ml of sterile, molten cooled (45-50°C) HiCrome™ Modified EC 0157:H7 Selective Agar Base (M1862). Mix well and pour into sterile Petri plates.

HiCrome™ ESBL Agar Supplement

FD278

Recommended for the detection of Extended Spectrum β -Lactamase (ESBL) producing organisms.

Formula

(per vial, sufficient for 500 ml medium)

Ceftazidime	1.50 mg
Cefotaxime	1.50 mg
Ceftriazone	1.00 mg
Aztreonam	1.00 mg
Fluconazole	5.00 mg

Directions

Rehydrate the contents of 1 vial aseptically with 5 ml sterile distilled water. Mix well and aseptically add to 500ml of sterile, molten, cooled (45-50°C) HiCrome™ ESBL Agar (M1829). Mix well and pour into sterile Petri plates.

HiCrome™ KPC Agar Supplement

FD279

Recommended for the detection of Carbapenem resistant gram negative bacteria.

Formula

(Per vial sufficient for 500 ml medium)

Selective Mix	50 mg
---------------	-------

Directions

Rehydrate the contents of 1 vial aseptically with 5 ml sterile distilled water. Mix well and aseptically add to 500ml of sterile, molten, cooled (45-50°C) HiCrome™ KPC Agar Base (M1831). Mix well and pour into sterile Petri plates.

HiCrome™ Listeria Selective Supplement

FD181

An antimicrobial supplement recommended for rapid and direct identification of *Listeria species*.

Formula

(per vial sufficient for 500 ml medium)

Ceftazidime	2.00 mg
Amphotericin B	2.50 mg

Directions

Rehydrate the contents of 1 vial with 5 ml of sterile distilled water. Mix well and aseptically add to 500 ml sterile, molten, cooled (45-50°C) HiCrome™ Listeria Agar Base, Modified (M1417/SM1417) / "HiCrome™ Listeria HiCynth™ Agar Base, Modified (MCD1417) / HiCrome™ Listeria Agar Base (M1417F) / HiCrome™ L.mono Rapid Differential Agar Base (M1924). Mix well and pour into sterile Petri plates.

HiCrome™ Nickels & Leesment Selective Supplement

FD245

An antibiotic supplement used for the selective isolation of *Leuconostoc species*.

Formula

(Per vial sufficient for 500 ml medium)

Vancomycin	100.00 mg
------------	-----------

Directions

Rehydrate the contents of 1 vial aseptically with 5 ml of sterile distilled water. Mix well and aseptically add it to sterile, molten, cooled (45-50°C) HiCrome™ Nickels and Leesment Medium (M1712) / HiCrome™ Nickels and Leesment HiVeg™ Agar (MV1712). Mix well and pour into sterile Petri plates.

Supplement

HiCrome™ Selective Salmonella Agar Supplement

FD274

For the selective isolation and differentiation of *Salmonella species* from coliforms by chromogenic method.

Formula

(Per vial sufficient for 1000 ml medium)

Novobiocin	10.00 mg
Cefsulodin	24.00 mg

Directions

Rehydrate the contents of 1 vial aseptically with 5ml of sterile distilled water. Mix gently to dissolve the contents completely. Aseptically add the rehydrated contents to 1000 ml of sterile, cooled (45-50°C) HiCrome™ Selective Salmonella Agar Base (M1842) / HiCrome™ Selective Salmonella HiCynth™ Agar Base (MCD1842). Mix well and pour into sterile Petri plates.

HiCrome™ Strep B Selective Supplement

FD273

An antibiotic supplement recommended for the selective isolation of Group B Streptococci from clinical samples.

Formula

(per vial, sufficient for 1000 ml medium)

Colistin	10.00 mg
Nalidixic Acid	10.00 mg
Gentamicin	2.00 mg

Directions

Rehydrate the contents of 1 vial aseptically with 5 ml of sterile distilled water. Mix well and aseptically add to 1000 ml of sterile, molten cooled (45-50°C) HiCrome™ Strep B Selective Agar Base (M1840) / HiCrome™ Strep B HiCynth™ Agar Base (MCD1840 / HiCrome™ Strep B Selective Agar Base, Modified (M1966). Mix well and pour into sterile Petri plates.

HiCrome™ VRE Agar Supplement

FD277

HiCrome™ VRE Agar Supplement is recommended for selective isolation of Vancomycin Resistant Enterococci (VRE).

Formula

(Per vial sufficient for 500 ml medium)

Vancomycin	4.00 mg
Fluconazole	5.00 mg

Directions

Rehydrate the contents of 1 vial aseptically with 5 ml sterile distilled water. Mix well and aseptically add to 500ml of sterile, molten, cooled (45-50°C) HiCrome™ VRE Agar Base (M1830) / HiCrome™ VRE Agar Base, Modified (M1925). Mix well and pour into sterile Petri plates.

Klebsiella Selective Supplement

FD225

Recommended for the selective isolation and easy detection of *Klebsiella species* from water and other sources.

Formula

(Per vial sufficient for 500 ml medium)

Carbenicillin	25.000 mg
---------------	-----------

Directions

Rehydrate the contents of 1 vial aseptically with 2 ml of sterile distilled water. Mix well and aseptically add it to 500 ml of sterile, molten, cooled (45-50°C) HiCrome™ Klebsiella Selective Agar Base (M1573) / HiCrome™ Klebsiella Selective HiVeg™ Agar Base (MV1573). Mix well and pour into sterile Petri plates.

Supplement

L. mono Enrichment Supplement I

FD214

For selective differentiation of *Listeria monocytogenes* from other *Listeria* species, as per ISO Committee.

Formula

(Per vial sufficient for 500 ml medium)

L – phosphatidylinositol	1.00 g
Distilled water	25.00 ml

Directions

Thaw the contents of 1 vial of L. mono Enrichment Supplement I at room temperature. Aseptically add the sterile contents to 460 ml of sterile, molten, cooled (45-50°C) L. mono Differential Agar Base (M1540) / L. mono Differential HiVeg™ Agar Base (MV1540) along with sterile rehydrated contents of 1 vial each of L. mono Selective Supplement I (FD212) and L. mono Selective Supplement II (FD213) / or 470ml of HiCrome™ L. mono Rapid Differential Agar Base (M1924). Mix well and pour into sterile Petri plates.

L. mono Enrichment Supplement II

FD227

Recommended for selective differentiation of *Listeria monocytogenes* from other *Listeria* species.

Formula

(Per vial sufficient for 500 ml medium)

L-phosphatidylinositol	0.50 g
Distilled water	15.00 ml

Directions

Thaw the contents of 1 vial of L. mono Enrichment Supplement II at room temperature. Aseptically add the sterile contents of one vial to 470 ml of sterile, molten, cooled (45-50°C) L. mono Confirmatory Agar Base (M1552) / L. mono Confirmatory HiVeg™ Agar Base (MV1552) along with sterile rehydrated contents of one vial each of L. mono Selective Supplement I (FD212) and L. mono Selective Supplement II (FD213). Mix well and pour into sterile Petri plates.

L. mono Selective Supplement I

FD212

A selective supplement recommended by ISO Committee for the isolation of *Listeria species*.

Formula

(Per vial sufficient for 500 ml medium)

Polymyxin B sulphate	38350 IU
----------------------	----------

Directions

Rehydrate the contents of 1 vial aseptically with 10 ml sterile distilled water. Mix well and aseptically add it to 460 ml of sterile, molten, cooled (45-50°C) L. mono Differential Agar Base (M1540) / L. mono Differential HiVeg™ Agar Base (MV1540) along with sterile contents of one vial of L. mono Enrichment Supplement I (FD214) and sterile rehydrated contents of one vial of L. mono Selective Supplement II (FD213) or add in 470 ml of sterile, molten, cooled (45-50°C) L. mono Confirmatory Agar Base (M1552) / L. mono Confirmatory HiVeg™ Agar Base (MV1552) along with sterile contents of one vial of L. mono Enrichment Supplement II (FD227) and rehydrated contents of one vial of L. mono Selective Supplement II (FD213). Mix well and pour into sterile Petri plates.

L. mono Selective Supplement II

FD213

A selective supplement recommended by ISO Committee for the isolation of *Listeria species*.

Formula

(Per vial sufficient for 500 ml medium)

Ceftazidime	10.00 mg
Amphotericin B	5.00 mg
Nalidixic acid, sodium salt	10.00 mg

Directions

Rehydrate the contents of 1 vial aseptically with 2 ml of 0.2 N Sodium hydroxide, further add 3 ml of sterile distilled water. Mix well and aseptically add it to 460 ml of sterile, molten, cooled (45-50°C) L. mono Differential Agar Base (M1540) / L. mono Differential HiVeg™ Agar Base (MV1540) along with sterile contents of 1 vial of L. mono Enrichment Supplement I (FD214) and sterile rehydrated contents of 1 vial of L. mono Selective Supplement I (FD212) or add in 470 ml of sterile, molten, cooled (45-50°C) L. mono confirmatory Agar Base (M1552) / L. mono Confirmatory HiVeg™ Agar Base (MV1552) along with sterile contents of 1 vial of L. mono Enrichment Supplement II (FD227) and sterile rehydrated contents of 1 vial of L. mono selective Supplement I (FD212). Mix well and pour into sterile Petri plates.

Supplement

Lecithin Solution

FD332

Recommended for selective differentiation of *L.monocytogenes* and *Listeria* species.

Formula

(Per vial sufficient for 480 ml medium)

Soya lecithin	20 ml
---------------	-------

Directions

Thaw the content of one vial at room temperature. Mix well and aseptically add it to 480 ml of sterile molten, cooled (45-50°C) HiCrome L.mono Differential Agar Base (M2009) along with sterile contents of 1 vial of Modified L.mono Selective Supplement (FD333). Mix well and pour into sterile Petri plates.

Leeds Acinetobacter Selective Supplement I

FD335

For the selective isolation of MDR *Acinetobacter* species.

Formula

(Per vial sufficient for 1000 ml medium)

Vancomycin	10 mg
Cefsulodin	15 mg
Cefradine	50 mg

Directions

Rehydrate the contents of 1 vial aseptically with 5ml of sterile distilled water. Mix gently to dissolve the contents completely. Aseptically add the rehydrated contents to 1000ml of sterile, molten cooled (45-50°C) Leeds Acinetobacter Agar Base (M1839) / HiCrome™ Acinetobacter Agar Base (M1938). Mix well and pour into sterile Petri plates.

M-CP Selective Supplement - I

FD153

An antibiotic supplement recommended by the Directive of the Council of the European Union 98/83/EC for the selective isolation of *Clostridium perfringens*

Formula

(per vial sufficient for 500 ml medium)

D-cycloserine	200.00 mg
Polymyxin B sulphate	12.50 mg

Directions

Rehydrate the contents of one vial aseptically with 5 ml of sterile distilled water. Mix well and aseptically add it to 485 ml sterile, molten, cooled (45-50°C) M-CP Agar Base (M1354) / M-CP HiVeg™ Agar Base (MV1354) along with one vial of M-CP Selective Supplement II (FD154) / M-CP Selective Supplement II, Modified (FD154A). Mix well and pour into sterile Petri plates.

M-CP Selective Supplement - II

FD154

Filter sterilized 0.5% solution of phenolphthalein diphosphate recommended by the Directive of the Council of the European Union 98/83/EC for the selective isolation of *Clostridium perfringens*.

Formula

(per vial sufficient for 500 ml medium)

Phenolphthalein diphosphate	0.05 g
Distilled water	10.00 ml

Directions

Warm up the refrigerated 0.5% phenolphthalein diphosphate solution to room temperature and add aseptically 10 ml of solution to 485 ml sterile, molten, cooled (45-50°C) M-CP Agar Base (M1354) / M-CP HiVeg™ Agar Base (MV1354) along with rehydrated contents of one vial of M-CP Selective Supplement I (FD153). Mix well and pour into sterile Petri plates.

M-CP Selective Supplement - II , Modified

FD154A

This supplement is recommended for selective isolation of *Clostridium perfringens*.

Formula

(per vial sufficient for 500 ml medium)

Phenolphthalein diphosphate	0.05 g
-----------------------------	--------

Directions

Rehydrate the contents of 1 vial aseptically with 10 ml of sterile distilled water. Mix well. Add aseptically to 485 ml sterile, molten, cooled to 45-50°C M-CP Agar Base (M1354)/ M-CP HiVeg™ Agar Base (MV1354) along with rehydrated contents of one vial of M-CP Selective Supplement I (FD153). Mix well and pour into sterile Petri plates.

Supplement

MDR Acinetobacter Selective Supplement

FD271

This antibiotic supplement is recommended for the selective isolation of MDR strains of *Acinetobacter species*.

Formula

(Per vial sufficient for 500 ml medium)

Ampicillin, sodium salt	5.00 mg
Ceftazidime	5.00 mg

Directions

Rehydrate the contents of 1 vial aseptically with 5ml of sterile distilled water. Mix gently to dissolve the contents completely. Aseptically add the rehydrated contents to 500ml of sterile, molten cooled (45-50°C) Leeds Acinetobacter Agar Base (M1839) / HiCrome™ Acinetobacter Agar Base (M1938). Mix well and pour into sterile Petri plates.

MeReSa Selective Supplement

FD229

An antimicrobial supplement recommended for the selective isolation of Methicillin Resistant *Staphylococcus aureus* from clinical specimens.

Formula

(Per vial sufficient for 500 ml medium)

Methicillin	2.00 mg
-------------	---------

Directions

Rehydrate the contents of 1 vial with 5 ml of sterile distilled water and aseptically add to 500 ml of sterile molten cooled (45-50°C) MeReSa Agar Base (M1594)/ HiCrome™ MeReSa Agar, Base (M1674)/ HiCrome™ MeReSa HiVeg Agar Base (MV1674) / HiCrome™ MRSA Agar Base Modified (M1953). This supplement either be used individually or in combination with (FD259) Cefoxitin Supplement. Mix well and pour into sterile Petri plates.

Modified L.mono Selective Supplement

FD333

Recommended for selective isolation of *L.monocytogenes* and *Listeria species* based on PCPLC activity.

Formula

(Per vial sufficient for 500 ml medium)

Nalidixic acid	0.10 g
Polymyxin B Sulfate	0.000004564 g
Ceftazidime	0.10 g
Amphotericin	0.005 g

Directions

Rehydrate the content of one vial aseptically with 5 ml of sterile distilled water. Mix well and aseptically add it to 480 ml of sterile molten, cooled (45-50°C) HiCrome™ L.mono Differential Agar Base (M2009) along with sterile contents of 1 vial of Lecithin solution(FD332). Mix well and pour into sterile Petri plates.

Monensin Selective Supplement

FD309

Recommended for the selective and differential isolation of coliform bacteria using membrane filtration technique.

Formula

(per vial sufficient for 1000 ml medium)

Monensin	0.038 g
----------	---------

Directions

Rehydrate the contents of one vial aseptically with 2 ml of methanol. Mix well and aseptically add to 1000 ml sterile cooled (45-50°C) HiCrome M-Coliform Differential Agar Base (M1951). Mix well and pour into sterile Petri plates.

Supplement

MUG Supplement

FD092

A fluorogenic substrate recommended for measuring β - glucuronidase activity, for rapid and sensitive identification of *Escherichia coli*.

Formula

(per vial sufficient for 500 ml / 1000 ml medium)
4-Methylumbeliferyl β -D-Glucuronide 50.00 mg
(MUG)

Directions

Rehydrate the contents of one vial aseptically with 5 ml methanol. One vial is sufficient for 500 ml agar media or 1000 ml of broth media. It is also recommended for 1000 ml of M-FC Basal Medium (M1812).

Oxytetra Selective Supplement

FD032

An antibiotic supplement recommended for the selective isolation and cultivation of yeasts and moulds.

Formula

(per vial, sufficient for 500 ml medium)
Oxytetracycline 50.00 mg

Directions

Rehydrate the contents of 1 vial aseptically with 5 ml of sterile distilled water. Mix well and aseptically add to 500 ml of sterile, molten, cooled (45-50°C) Oxytetra Glucose Yeast Extract Agar Base (M639 / M639I) / Oxytetra Glucose Yeast Agar Base, Granulated (GM639 / GM639I) / Oxytetra Glucose Yeast Extract Agar Base with Biotin (M1136) / HiCrome™ OGYE Agar Base (M1467). Mix well and pour into sterile Petri plates.

Polymyxin B Selective Supplement

FD003

An antibiotic supplement recommended for the selective isolation of various microorganisms

Formula

(per vial, sufficient for 500 ml/1000 ml medium)
Polymyxin B sulphate 50000 Units

Directions

Rehydrate the contents of one vial aseptically with 2 ml sterile distilled water. Mix well and aseptically add it to 475 ml of sterile, molten *Bacillus cereus* Agar Base (M833) / *Bacillus cereus* HiVeg™ Agar Base (MV833) / *Bacillus cereus* HiCynth Agar Base (MCD833) or to 450 ml of KG Agar Base (M658) / KG HiVeg™ Agar Base (MV658) / MYP Agar Base (M636/M636S / M636F) / MYP HiCynth™ Agar Base (Phenol Red Egg Yolk Polymyxin HiCynth™ Agar Base) (MCD636) / MYP HiVeg™ Agar Base (MV636) / Modified MYP Agar Base (M1139 / M1139I) / Modified MYP HiVeg™ Agar Base (MV1139) along with 25 ml / 50 ml Egg Yolk Emulsion (FD045) to make a total volume of 500 ml or to 500 ml of SDS Agar (M1155) / SDS HiVeg™ Agar (MV1155) / Salt Polymyxin Broth Base (M821/M821I) / Salt Polymyxin HiVeg™ Broth Base (MV821) / HiCrome™ Staph Agar Base, Modified (M1837) / Soyabean Casein Digest Medium Base (M011F) or to 1000 ml of HiCrome™ *Bacillus* Agar (M1651). Mix well and pour into sterile Petri plates / tubes.

Supplement

Potassium Tellurite 1%

FD052

Recommended for the selective isolation of *Staphylococci* and *Corynebacteria*.

Formula

(to achieve 1% solution dilute the contents in 8.9 ml sterile distilled water)

Potassium tellurite concentrate	1.1 ml
---------------------------------	--------

Directions

Warm up the refrigerated contents of one vial to room temperature. Add aseptically 8.9 ml sterile distilled water, mix well and add in sterile, molten, cooled (45-50°C) Baird Parker Agar Base (M043B / MM043 / MU043 / ME043)/ Vogel Johnson Agar Base w/o Tellurite (M023 / MM023 / MU023) / Vogel- Johnson Agar Base w/ 1.5% Agar (M023F) / Vogel Johnson HiVeg™ Agar Base w/o Tellurite (MV023) / Vogel Johnson HiCynth™ Agar Base w/o Tellurite (MCD023) / Mycoplasma Broth Base w/ CV (M268) / Mycoplasma HiVeg™ Broth Base w/ CV (MV268) / TPEY Agar Base (M402)/ TPEY HiVeg™ Agar Base (MV402)/ Tellurite Glycine Agar Base (M448) / Cholera Medium Base (M558) / Cholera HiVeg™ Medium Base (MV558) / Giolitti-Cantoni Broth Base (M584) / Dextrose Proteose Peptone Agar Base (M734) / Dextrose Proteose Peptone HiVeg™ Agar Base (MV734) / Cystine Tellurite Agar Base (M881) / Diphtheria Virulence Agar Base (M882) / Diphtheria Virulence HiVeg™ Agar Base (MV882) / Tryptone Tellurite Agar Base (M1056) / Baird Staphylococcus Enrichment Broth Base (M1091) / Baird Staphylococcus Enrichment Broth Base, Granulated (GM1091) / Tellurite Blood Agar Base (M1260) / Mitis Salivarius Agar Base (M259)/ Mitis Salivarius HiVeg™ Agar Base (MV259) / Monsur Medium Base (M474) as desired. Mix well and dispense in sterile Petri plates or tubes.

Tellurite - Cefixime Supplement

FD147

A selective supplement recommended by ISO Committee for the isolation of *Escherichia coli* 0157:H7.

Formula

(Per vial sufficient for 500 ml medium)

Potassium tellurite	1.250 mg
Cefixime	0.025 mg

Directions


Rehydrate the contents of 1 vial aseptically with 5 ml of sterile distilled water. Mix gently to dissolve the contents completely. Aseptically add the contents to 495 ml of sterile, molten, cooled (45-50°C) MacConkey Sorbitol Agar Base (M2981) / HiCrome™ MacConkey Sorbitol Agar Base (M1340). Mix well and pour into sterile Petri plates.

DT001 HiDtect™ UTI Identification Disc

For rapid detection and confirmation of microorganisms mainly causing urinary tract infections, For eg. *E.coli*, *Proteus*, *Klebsiella*, *Pseudomonas*, *Staphylococcus aureus* and *Enterococcus* species.

- Appearance** : White coloured sterile identification disc.
Cultural Response : Identification observed within 1-4 hours on replication and incubation at 35-37°C, when disc is placed on an 18 hour old grown culture plate of any general media.

Organisms (ATCC)	Colour of Colony
<i>Escherichia coli</i> (25922) (00013*)	Pink-purple
<i>Staphylococcus aureus</i> (25923) (00034*)	colourless to green
<i>Pseudomonas aeruginosa</i> (27853) (00025*)	colourless (greenish pigment is observed)
<i>Enterococcus faecalis</i> (29212) (00087*)	blue - blue green (small)
<i>Klebsiella pneumoniae</i> (13883) (00097*)	blue to purple mucoid
<i>Proteus mirabilis</i> (12453)	light brown




Key : * corresponds to WDCM numbers

DT002 HiDtect™ Salmonella Identification Disc

For rapid detection of *Salmonella* species from coliforms

- Appearance** : White coloured sterile identification disc.
Cultural Response : Identification observed within 1-4 hours on replication and incubation at 35-37°C, when disc is placed on an 18 hour old grown culture plate of any general media.

Organisms (ATCC)	Colour of Colony
<i>Escherichia coli</i> (25922) (00013*)	blue - greenish blue
<i>Salmonella</i> Typhimurium (14028) (00031*)	light purple
<i>Salmonella</i> Enteritidis (13076) (00030*)	light purple
<i>Klebsiella pneumoniae</i> (13883) (00097*)	blue to purple mucoid



Key : * corresponds to WDCM numbers


DT003 HiDtect™ Pseudomonas Identification Disc

For rapid detection of *Pseudomonas aeruginosa* from clinical and nonclinical specimen

- Appearance** : White coloured sterile identification disc.
Cultural Response : Identification observed within 1-4 hours on replication and incubation at 35-37°C, when disc is placed on an 18 hour old grown culture plate of any general media.

Organisms (ATCC)	Colour of Colony	Fluorescence
<i>Escherichia coli</i> (25922) (00013*)	colourless	Negative
<i>Pseudomonas aeruginosa</i> (27853)	colourless	Positive
<i>Enterococcus faecalis</i> (29212) (00087*)	colourless	Negative
<i>Klebsiella pneumoniae</i> (13883) (00097*)	colourless, mucoid	Negative

Flourescence under uv (365 nm)



Key : * corresponds to WDCM numbers


DT005 HiDtect™ Universal Microbial Limit Test Disc

For detection of pathogenic microorganisms such as *E.coli*, *S. aureus*, *P.aeruginosa*, and *Salmonella* species from pharmaceutical preparations, raw materials and cosmetic samples etc.

Appearance : White coloured sterile identification disc.

Cultural Response : Identification observed within 1-4 hours on replication and incubation at 35-37°C, when disc is placed on an 18 hour old grown culture plate of any general media.

Organisms (ATCC)	Colour of Colony
<i>Pseudomonas aeruginosa</i> (9027) (00026*)	Pink-purple colourless (greenish pigment is observed)
<i>Escherichia coli</i> (8739) (00012*)	Pink-purple
<i>Staphylococcus aureus subsp aureus</i> (6538) (00032*)	green to bluish green
<i>Salmonella</i> Typhimurium (14028) (00031*)	colourless
<i>Salmonella</i> Abony (NCTC 6017) (00029*)	colourless
<i>Proteus mirabilis</i> (12453)	light brown



Key : * corresponds to WDCM numbers

DT008 HiDtect™ Dual Confirmation of E. coli Identification Disc

For rapid detection and confirmation of *Escherichia coli* in water and food samples based on chromogenic and fluorogenic confirmation.


Appearance : White coloured sterile identification disc.

Cultural Response : Identification observed within 1-4 hours on replication and incubation at 35-37°C, when disc is placed on an 18 hour old grown culture plate of any general media.

Organisms (ATCC)	Colour of Colony	Fluorescence under uv light
<i>Pseudomonas aeruginosa</i> (27853) (00025*)	colourless	negative
<i>Escherichia coli</i> (25922) (00013*)	blue	positive
<i>Staphylococcus aureus subsp aureus</i> (25923) (00034*)	colourless	negative
<i>Salmonella</i> Typhimurium (14028) (00031*)	colourless	negative

Key : * corresponds to WDCM numbers

Fluorescence under uv (365 nm)



DT010 HiDtect™ Universal Food pathogen Identification Disc


For rapid detection of food pathogens such as *E.coli*, *E.coli* O157:H7, *Staphylococcus aureus*, *Salmonella*, *Listeria* and *Shigella* species etc. from various food, dairy, fish and meat products.

Appearance : Light pink coloured sterile identification disc.

Cultural Response : Identification observed within 1-4 hours after replication and incubation at 35-37°C, when disc is placed on an 18 hour old grown culture plate of any general media.

Organisms (ATCC)	Colour of Colony
<i>Escherichia coli</i> (25922) (00013*)	purple
<i>Staphylococcus aureus subsp aureus</i> (25923) (00034*)	colourless - green
<i>Salmonella</i> Typhimurium (14028) (00031*)	colourless
<i>Listeria monocytogenes</i> (19111)	blue - green
<i>Escherichia coli</i> O157:H7 (NCTC 12900)	purple - pink
<i>Shigella flexneri</i> (12022)	colourless

-  *Listeria monocytogenes* (ATCC 19111)
Blue coloured
-  *Escherichia coli* (ATCC 25922)
Purple coloured
-  *Escherichia coli* O157:H7 (NCTC 12900)
Pink purple coloured
-  *Pseudomonas aeruginosa* (ATCC 27853)
colourless, greenish pigment may be observed



Key : * corresponds to WDCM numbers

DT011 HiDtect™ Bacillus Identification Disc


For rapid detection and differentiation between various species of *Bacillus* such as *B.subtilis*, *B.cereus*, *B.thuringiensis*, from food, meat, fish, cosmetic and pharmaceutical preparations.

Appearance : Pale pink coloured sterile identification disc.

Cultural Response : Identification observed within 1-4 hours after incubation at 35-37°C, when disc is placed on an 18 hour old grown culture plate of any general media.

Organisms (ATCC)	Colour of Colony
<i>Bacillus cereus</i> (10876)	light blue
<i>Bacillus subtilis</i> subsp. spizizenii (6633) (00003*)	yellowish green to green
<i>Bacillus thuringiensis</i> (10792)	light blue

-  *Bacillus subtilis* (ATCC 6633)
green coloured




Key : * corresponds to WDCM numbers

DT012 HiDtect™ Total Coliform Identification Disc

For qualitative detection of coliforms from water, pharmaceutical preparations, dairy and food products.

Appearance : Pink coloured sterile identification disc.
Cultural Response : Identification observed within 1-4 hours after replication and incubation at 35-37°C, when disc is placed on an 18 hour old grown culture plate of any general media.

Organisms (ATCC)	Colour of Colony
<i>Escherichia coli</i> (25922) (00013*)	dark blue
<i>Enterobacter cloacae</i> (23355)	reddish to purple
<i>Citrobacter freundii</i> (8090)	reddish to purple
<i>Klebsiella pneumoniae</i> (13883) (00097*)	pink to purple




Key : * corresponds to WDCM numbers

DT013 HiDtect™ Differential Coli-E. coli Identification Disc

For or rapid detection of *E.coli*, *Klebsiella*, *Pseudomonas* and *Salmonella* species in food, environmental and water samples.

Appearance : White coloured sterile identification disc.
Cultural Response : Identification observed within 1-4 hours after replication and incubation at 35-37°C, when disc is placed on an 18 hour old grown culture plate of any general media.

Organisms (ATCC)	Colour of Colony
<i>Escherichia coli</i> (25922) (00013*)	dark blue to violet
<i>Klebsiella pneumoniae</i> (13883) (00097*)	light pink
<i>Pseudomonas aeruginosa</i> (27853) (00025*)	colourless
<i>Salmonella</i> Enteritidis (13076)	colourless




Key : * corresponds to WDCM numbers

DT015 HiDtect™ Universal Enviro Identification Disc

For rapid detection of *Pseudomonas*, *Enterococcus*, *E.coli* and *Salmonella* species etc. from food environmental samples, samples of clinical origin such as nosocomial samples.

Appearance : Light pink coloured sterile identification disc.
Cultural Response : Identification observed within 1-4 hours on replication and incubation at 35-37°C, when disc is placed on an 18 hour old grown culture plate of any general media.

Organisms (ATCC)	Colour of Colony
<i>Escherichia coli</i> (25922) (00013*)	Pink-purple
<i>Staphylococcus aureus subsp aureus</i> (25923) (00034*)	colourless - green
<i>Pseudomonas aeruginosa</i> (27853) (00025*)	colourless (greenish pigment is observed)
<i>Enterococcus faecalis</i> (29212) (00087*)	blue - blue green (small)
<i>Salmonella</i> Typhimurium (14028) (00031*)	colourless




Key : * corresponds to WDCM numbers

HiCrome™
Single Streak Rapid Differentiation Series

Chromogenic identification





Antimicrobial susceptibility



Features :

***Rapid & Reliable Method
with Dual advantage***





K092 Hi*E.coli* Test Kit

Compartment Bag Test replaces MPN Technique

Recommended for easy detection of *E. coli* from water samples.

Features

Easy detection and enumeration of *E. coli* using Compartment Bags



- Economical, user friendly
- Self equipped test, only water sample to be tested is required
- No consumption of power or laboratory equipment
- Incubation at Ambient temperature

Simple Steps to Perform

- Step 1 : Collection of water sample
- Step 2 : Prepare for the test, collect all accessories
- Step 3 : Add Hi*E.coli* test bud to the water sample.
Mix well to release the contents of bud.
- Step 4 : Open the compartment bag. Transfer the sample into compartment bag.
Water levels filled should be at equal in levels each compartment of the bag.
- Step 5 : Seal the bag with clip provided.
- Step 6 : Incubate at 24-30°C for 24-30 hours or 35-44.5°C for 20-24 hours.
- Step 7 : Check for colour change from yellow to blue-green colour.
Yellow indicates negative and blue-green indicates positive test.
Correlate for outcomes of positive and negative results with MPN chart.
- Step 8 : Decontaminate the contents by opening the bag and adding the chlorine tablet.



Bottle for sample collection



HiFast™ Food Pathogen Detection Kit

Food Testing Becomes Easier Confirmation in 4-6 hours*

HiFast™ Food Pathogen Detection Kit, K097S/M/L

Food safety is a major focus of food microbiology. It is important to be able to detect microorganisms in food, in particular pathogenic microorganisms which include *Escherichia coli*, *Klebsiella pneumoniae*, *Salmonella*, *Escherichia coli* O157: H7, *Listeria* and *Clostridium*. So this test was designed to identify these pathogens rapidly as the conventional method of food testing takes longer time.



Salient Features :

- Simple, Fast and Convenient
- Ready to use
- Covers Major food pathogens

Kit contents :

- Stomacher bag for sample processing {available in three sizes : small (390 ml capacity), medium (710 ml capacity) and large (1.63 L capacity) {3 Packings available}}
- Enrichment medium for sample processing
- Differential Food pathogen Testing medium
- Aureus Confirmation test
- Listeria Confirmation test
- Clostridium Confirmation test



1. Control
2. *Escherichia coli*
3. *Salmonella*
4. *Escherichia coli* O157
5. *Vibrio*



Aureus Confirmation test



Listeria Confirmation test



Clostridium Confirmation test

* After Preincubation



Dehydrated Culture Media &
Ready Prepared Media
Comprehensive range of media formulations,
standard media and customer specified
media for diagnostic purpose.



MICROBIOLOGY



"All that the CELLS need."



ANIMAL CELL CULTURE



"All the flowers and fruits of tomorrow are in the
SEEDS of today." Nurture them with supreme
quality of Plant Tissue Culture media and
Plant Tissue Culture tested chemicals.
Media available in Powder and Liquid Form.



PLANT TISSUE CULTURE



"Somewhere, something incredible
known as **GENES** waiting to be acknowledged."
Unzipping Genes in Molecular Biology and
Diagnostics with perfect choice.



MOLECULAR BIOLOGY



HiMedia Laboratories Pvt. Ltd.

www.himedialabs.com

- CORPORATE OFFICE -

A-516, Swastik Disha Business Park, Via Vadhani Indl Est, LBS Marg,
Mumbai - 400 086, India.
Tel : +91-22-6147 1919 / 2500 3747 | Fax : +91-22-6147 1920 / 2500 5764
Email : info@himedialabs.com


- OVERSEAS OFFICES -

USA & Canada

HiMedia Laboratories LLC, 107 W Dorothys Way, Lincoln University, West Chester,
Pennsylvania 19352, USA.
Tel : +1-484-734-4401 | Fax : +1-484-734-4402
Email : infous@himedialabs.com

Europe

HiMedia Laboratories GmbH, Marie-Curie-Str. 3, 64683,
Einhausen, Germany.
Tel : +49 6251 989 24 26 | Fax : +49 6251 989 24 27
Email : infoeu@himedialabs.com



HIMEDIA®

For Life is Precious