



Esculin Mannitol Agar

M2072

Intended use

Recommended as a selective and differential media for the isolation of *Staphylococci* and *Enterococcus* based on mannitol fermentation and esculin hydrolysis.

Composition**

Ingredients	Gms / Litre
Peptone	23.000
Mannitol	10.000
Sodium chloride	5.000
Corn starch	1.000
Esculin	1.000
Ferric ammonium citrate	0.500
Phenol red	0.025
Nalidixic acid	0.015
Colistin sulphate	0.010
Agar	15.000
Final pH (at 25°C)	7.3±0.2

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 55.55 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Mix well and pour into sterile Petri plates.

Principle And Interpretation

Staphylococci are widespread in nature, although they are mainly found on the skin, skin glands and mucous membranes of mammals and birds. The coagulase-positive species i.e. *Staphylococcus aureus* is well documented as a human opportunistic pathogen and also considered potent pathogen from the point of view of food hygiene. The ability to clot plasma continues to be the most widely used and accepted criterion for the identification of pathogenic *staphylococci* associated with acute infections (1).

Group D Streptococci possess the group D lipoteichoic acid antigen in their cell walls. Former Group D species, which are predominant normal inhabitants of the human gastrointestinal tract, are termed as faecal Streptococci or Enterococci (2). The unique ability of Enterococci to split esculin was reported by Meyer and Schonfeld (3). The use of esculin hydrolysis in identification of Enterococci was first cited by Rochaix (4).

Peptone supplies nitrogenous and carbonaceous compounds, long chain amino acids, other essential growth factors and trace nutrients to the growing bacteria. Sodium chloride maintains osmotic balance. Mannitol is the fermentable carbohydrate, fermentation of which leads to acid production, detected by phenol red indicator. Corn starch helps in neutralizing the toxic compounds. Esculin in the medium is hydrolyzed to esculetin and dextrose. Esculetin reacts with ferric ammonium citrate to form a dark brown or black complex, visualized as a zone of black precipitate around the colonies. Esculin hydrolysis is shown by *Enterococcus* species. Nalidixic acid and Colistin sulphate helps in inhibiting gram negative bacteria. *S.aureus* ferment mannitol and produce yellow coloured colonies surrounded by yellow zones. Coagulase-negative strains of *S.aureus* are usually mannitol non-fermenters and therefore produce pink to red colonies surrounded by red-purple zones. Presumptive coagulase-positive yellow colonies of *S.aureus* should be confirmed by performing the coagulase test [tube or slide](1). *Enterococcus* species hydrolyse esculin and hence black precipitate is observed around the colonies. This medium helps in simultaneous differentiation of *Enterococcus* and *Staphylococcus* based on mannitol fermentation and esculin hydrolysis.

Type of specimen

Food and dairy samples ; Water samples

Specimen Collection and Handling

For food and dairy samples, follow appropriate techniques for sample collection and processing as per guidelines (5,6,7).

For water samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards.(8)

After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions :

In Vitro diagnostic Use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/ face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets

Limitations :

Due to nutritional variations, some strains may show poor growth

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

Quality Control

Appearance

Light yellow to pink homogeneous free flowing powder

Gelling

Firm,comparable with 1.5% Agar gel

Colour and Clarity of prepared medium

Red coloured clear to slightly opalescent gel forms in Petri plates

Reaction

Reaction of 5.55% w/v aqueous solution at 25°C. pH : 7.3±0.2

pH

7.10-7.50

Cultural Response

Cultural characteristics observed after an incubation at 35-37°C for 18-72 hours. Recovery rate is considered as 100% for bacteria growth on Soybean Casein Digest Agar.

Organism	Inoculum (CFU)	Growth	Recovery	Mannitol fermentation	Esculin hydrolysis
<i>Staphylococcus aureus</i> ATCC 6538 (00032*)	50 -100	luxuriant	≥50 %	yellow/white colonies surrounded by yellow zone	
<i>Staphylococcus aureus</i> ATCC 25923 (00034*)	50 -100	luxuriant	≥50 %	yellow/white colonies surrounded by yellow zone	
<i>Enterococcus faecalis</i> ATCC 29212 (00087*)	50-100	luxuriant	≥50%		positive reaction,blackening of medium around the colony
<i>Escherichia coli</i> ATCC 25922 (00013*)	≥10 ³	inhibited	0%		
<i>Enterobacter aerogenes</i> ATCC 13048 (00175*)	≥10 ³	inhibited	0%		

Key:- (*) Corresponding WDCM numbers

Storage and Shelf Life

Store between 10-30°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition. Seal the container tightly after use. Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (9,10).

Reference

1. Murray P. R., Baron J. H., Pfaller M. A., Tenover J. C. and Tenover F. C., (Eds.), 2003, Manual of Clinical Microbiology, 8th Ed., American Society for Microbiology, Washington, D.C. ,
2. Koneman E. W., Allen S. D., Janda W. M., Schreckenberger P. C., Winn W. C. Jr., 1992, Colour Atlas and Textbook of Diagnostic Microbiology, 4th Ed., J. B. Lippincott Company
3. Meyer and Schonfeld, 1926, Zentralbl. Bakteriol, Parasitenk. Infektionskr. Hyg. Abt. Orig. 99:402.
4. Rochaix, 1924, Comt. Rend. Soc. Biol., 90:771.
5. American Public Health Association, Standard Methods for the Examination of Dairy Products, 1978, 14th Ed., Washington D.C.
6. Salfinger Y., and Tortorello M.L. Fifth (Ed.), 2015, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.
7. Wehr H. M. and Frank J. H., 2004, Standard Methods for the Microbiological Examination of Dairy Products, 17th Ed., APHA Inc., Washington, D.C.
8. Baird R.B., Eaton A.D., and Rice E.W., (Eds.), 2015, Standard Methods for the Examination of Water and Wastewater, 23rd ed., APHA, Washington, D.C.
9. Isenberg, H.D. Clinical Microbiology Procedures Handbook, 2nd Edition.
10. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

Revision : 00/ 2018

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.