



# Technical Data

## Sodium Bicarbonate Solution (20 ml per vial)

FD080

A selective solution used in media is recommended for the selective isolation of *Enterococcus* species.

### Composition

Per vial sufficient for 980 ml medium

### Ingredients

Sodium bicarbonate (10% w / v)

### Concentration

20ml

### Directions:

Warm up the refrigerated solution to room temperature. Add content of 1 vial aseptically with 980 ml sterile, molten, cooled (45-50°C) GTC Agar Base [M922](#) along with 2 vials of GTC Supplement [FD044](#) . Mix well and pour into sterile petri plates.

### Storage and Shelf Life

Store at 2-8°C. Use before the expiry date on the label.

Revision : 01 / 2015

### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.