



# Technical Data

## Mycoplasma Supplement - P

FD182

A selective supplement for the isolation of *Mycoplasma pneumoniae* when used with Mycoplasma Agar Base / Mycoplasma HiVeg™ Agar Base.

### Composition

*Ingredients	Concentration
Horse serum	6ml
Yeast extract	3ml
Thallos acetate	0.008g
Penicillin	12000IU
Glucose	0.300g
Phenol red	0.0012g
Methylene blue	0.0003g
Mycoplasma Broth Base w/o CV (M267)	0.145g

### Directions:

Aseptically add 20 ml sterile distilled water to the vial and mix gently. Separately prepare Mycoplasma Agar Base [M266](#) / Mycoplasma HiVeg™ Agar Base [MV266](#) without supplement addition. Heat to boiling to dissolve the medium completely. Add 1 ml of the medium to each 10 small bottles. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. After sterilization allow the medium to set or solidify. Then aseptically add 2 ml of the reconstituted Mycoplasma Supplement-P [FD182](#) to each bottle containing agar. This will form a biphasic medium.

### Storage and Shelf Life

Store at 2-8°C. Use before the expiry date on the label.

Revision : 2 / 2016

\* Not For Medicinal Use

### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.