

Technical Data

Meropenem Supplement

An antimicrobial supplement recommended for the selective enrichment of Enterococci.

Composition

Per vial sufficient for 500 ml medium

*Ingredients

Meropenem

Directions:

 Rehydrate the contents of 1 vial with 5 ml of sterile distilled water and aseptically add to 1000 ml of sterile molten, cooled

 (45-50°C)
 , Vancomycin Resistant Enterococci (VRE) Agar Base
 M1763
 along with 2 vials of Vancomycin

 Supplement
 FD261
 or aseptically add 2 vials to 1000 ml of Vancomycin Resistant Enterococci (VRE) Broth Base
 M1762

 . Mix well and dispense as desired.
 Mix well and dispense as desired.
 M1762

Concentration

1mg

Storage and Shelf Life

Store at 2 - 8°C. Use before expiry date on the label.

* Not For Medicinal Use

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia[™] publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia[™] Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. A-516,Swastik Disha Business Park,Via Vadhani Ind. Est., LBS Marg, Mumbai-400086, India. Customer care No.: 022-6147 1919 Email: techhelp@himedialabs.com Website: www.himedialabs.com

FD262