



Technical Data

MDR Acinetobacter Selective Supplement

FD271

This antibiotic supplement is recommended for selective isolation of MDR strains of Acinetobacter.

Composition

Per vial sufficient for 500 ml medium

*Ingredients	Concentration
Ampicillin, sodium salt	5mg
Ceftazidime	5mg

Directions:

Rehydrate the contents of 1 vial aseptically with 5ml of sterile distilled water. Mix gently to dissolve the contents completely. Aseptically add the rehydrated contents to 500ml of sterile, molten, cooled (45-50°C) Leeds Acinetobacter Agar Base [M1839](#) / 500ml of HiCrome™ Acinetobacter Agar Base [M1938](#) . Mix well and pour into sterile Petri plates.

Storage and Shelf Life

Store at 2 - 8°C. Use before expiry date on the label.

Revision : 1 / 2012

* Not For Medicinal Use

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.