



Technical Data

NAD Supplement

FD327

For preparation of Mueller Hinton Agar F medium as per EUCAST guidelines

Composition

Per vial sufficient for 5000 ml medium

Ingredients

Beta NAD

Distilled water

Concentration

100 mg

5ml

Directions:

Thaw the contents of NAD Supplement at room temperature. Aseptically add 1ml of NAD solution to , molten, cooled (45-50°C) Mueller Hinton Agar [M173](#) / [M1084](#) along with 5% defibrinated horse blood . Mix well and pour into sterile petri plates.

Storage and Shelf Life

On receipt product should be stored at 2-8°C. For prolonged storage store at -20°C. Use before the expiry date on the label.

Revision : 00 / 2014

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.