



Actidione Agar Base without Actidione

M058

Actidione Agar Base without Actidione is used for the enumeration and detection of bacteria in specimens containing large number of yeasts and moulds.

Composition**

Ingredients	Gms / Litre
Casein enzymic hydrolysate	5.000
Yeast extract	4.000
Dextrose	50.000
Monopotassium phosphate	0.550
Potassium chloride	0.425
Calcium chloride	0.125
Magnesium sulphate	0.125
Ferric chloride	0.0025
Manganese sulphate	0.0025
Bromo cresol green	0.022
Agar	15.000
Final pH (at 25°C)	5.5±0.2

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 75.25 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45°C and aseptically add 10 mg cycloheximide to the medium. Mix well before pouring into sterile Petriplates.

Warning : Cycloheximide (Actidione) is very toxic. Avoid skin contact or aerosol formation and inhalation.

Principle And Interpretation

Actidione Agar was formulated by Green and Gray (1) which may be used for microbiological investigation during brewing and baking. Actidione® (Cycloheximide) at a concentration of 0.001% permits the growth of bacteria and inhibits the growth of most yeasts and moulds except dermatophytes. This medium may be used for the estimation of bacterial contamination of pitching yeast. Addition of penicillin or streptomycin may be used for selective isolation of dermatophytes.

Casein enzymic hydrolysate acts as source of nitrogen while yeast extract is a rich reservoir of vitamins. Dextrose in high amount along with mineral salts at acidic pH favours sugar fermentation.

Quality Control

Appearance

Light yellow to light green homogeneous free flowing powder

Gelling

Firm, comparable with 1.5% Agar gel

Colour and Clarity of prepared medium

Greenish blue clear to slightly opalescent gel forms in Petri plates

Reaction

Reaction of 7.53% w/v aqueous solution at 25°C. pH : 5.5±0.2

pH

5.30-5.70

Cultural Response

M058: Cultural characteristics observed after an incubation at 25-30°C for 40-48 hours.

Organism	Inoculum (CFU)	Growth (Plain)	Recovery	Growth (w/ Actidione)	Recovery (w / Actidione)
----------	----------------	----------------	----------	-----------------------	--------------------------

Cultural Response

<i>Escherichia coli</i> ATCC 25922	50-100	luxuriant	$\geq 70\%$	good-luxuriant	$\geq 50\%$
<i>Lactobacillus fermentum</i> ATCC 9338	50-100	luxuriant	$\geq 70\%$	good-luxuriant	$\geq 50\%$
<i>Proteus mirabilis</i> ATCC 25933	50-100	luxuriant	$\geq 70\%$	good-luxuriant	$\geq 50\%$
<i>Saccharomyces cerevisiae</i> ATCC 9763	50-100	luxuriant	$\geq 70\%$	inhibited	0%
<i>Saccharomyces uvarum</i> ATCC 9080	50-100	luxuriant	$\geq 70\%$	inhibited	0%

Storage and Shelf Life

Store below 30°C in tightly closed container and the prepared medium at 2 - 8°C. Use before expiry date on label.

Reference

1.Green, S.R. and Gray, P.P., 1950, Wallerstein Lab. Communication, 13:357.

Revision : 2 / 2015

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.