



## Lactic Bacteria Differential Broth

M1086

Lactic Bacteria Differential Broth is used for differentiation of homofermentative and heterofermentative lactic acid bacteria.

### Composition\*\*

Ingredients	Gms / Litre
Casein enzymic hydrolysate	10.000
Papaic digest of soyabean meal	1.500
Casein acid hydrolysate	3.000
Yeast extract	1.000
Fructose	2.500
Monopotassium phosphate	2.500
Bromocresol green	0.055
Final pH ( at 25°C)	7.0±0.2

\*\*Formula adjusted, standardized to suit performance parameters

### Directions

Suspend 20.5 grams in 1000 ml distilled water. Add 1 gram of polysorbate 80. Heat if necessary to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Mix well and dispense as desired.

### Principle And Interpretation

Lactic Bacteria Differential Broth is formulated as per McDonald et al (1) for differentiation of homofermentative Lactobacilli and heterofermentative Streptococci. Lactobacilli and Streptococci are used as starter cultures in food and dairy industry. Streptococci grow first and produce metabolites, lowering redox potential which enables Lactobacilli to grow. Lactobacilli synthesize products which stimulate growth of Streptococci.

Medium constituents like casein acid hydrolysates, papaic digest of soyabean meal and yeast extract supply all the necessary nutrients for the growth of lactic bacteria. Fructose is the fermentable carbohydrate in the medium. Bromo cresol green is the pH indicator.

Heterofermentative lactic acid bacteria produce CO<sub>2</sub>, lactic acid, acetic acid, ethanol and mannitol. Homofermentative lactic acid bacteria produce only lactic acid. Homofermentative lactic acid bacteria produce lactic acid from fructose and is indicated by yellow colour formation. Heterofermentative lactic acid bacteria induce lesser acidification and thus vary in the colour formation by the indicator in the medium

### Quality Control

#### Appearance

Light yellow to bluish grey homogeneous free flowing powder

#### Colour and Clarity of prepared medium

Blue coloured clear solution in tubes

#### Reaction

Reaction of 2.05% w/v aqueous solution at 25°C. pH : 7.0±0.2

#### pH

6.80-7.20

#### Cultural Response

M1086: Cultural characteristics observed after an incubation at 35-37°C for 18-48 hours

Organism	Inoculum (CFU)	Growth	Colour of medium
----------	----------------	--------	------------------

#### Cultural Response

<i>Lactobacillus casei</i> ATCC 9595	50-100	luxuriant	green
--------------------------------------	--------	-----------	-------

---

<i>Lactobacillus plantarum</i> ATCC 8014	50-100	luxuriant	green
<i>Streptococcus cremoris</i> ATCC 19257	50-100	luxuriant (incubated at 30°C)	blue
<i>Streptococcus thermophilus</i> ATCC 14485	50-100	luxuriant (incubated at 45°C)	bluish-green

### Storage and Shelf Life

Store at 30°C in tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label.

### Reference

1. McDonald L.C., McFecters R.F., Daeschel M.A. and Fleming H.P., 1987, Appl. Environ. Microbiol., 53:1382.

Revision : 2 / 2015

#### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.