



## CPC Agar Base

M1241

CPC Agar Base is used for the cultivation and identification of *Vibrio* species from foods.

### Composition\*\*

Ingredients	Gms / Litre
Peptone	10.000
Meat extract B #	5.000
Cellobiose	15.000
Sodium chloride	20.000
Bromothymol blue	0.040
Cresol red	0.040
Agar	15.000
Final pH ( at 25°C)	7.6±0.2

\*\*Formula adjusted, standardized to suit performance parameters

# Equivalent to Beef extract

### Directions

Suspend 32.54 grams in 500 ml of distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C and aseptically add the rehydrated contents of 1 vial of CPC Supplement (FD110). Mix well and pour into sterile Petri plates

### Principle And Interpretation

*Vibrio* species are natural inhabitants of brackish and salt water. Human disease is associated with ingestion of contaminated water or consumption of contaminated seafood. Wound and systemic infections develop following contact with contaminated water (1). CPC (Cellobiose, Polymyxin and Colistin) Agar Base formulated as per APHA (2) is recommended for the cultivation and identification of *Vibrio* species from foods. CPC Agar is a selective and differential agar medium, designed to differentiate *Vibrio vulnificus* from other *Vibrios* (2). *Vibrio cholerae* strains except *V.cholerae* 01-classical biotype grow on CPC Agar while most *Vibrio parahaemolyticus* strains do not grow on CPC Agar.

CPC Agar contains meat extract B and peptone, which supply the essential nitrogenous, carbonaceous compounds, long chain amino peptides, vitamins and other growth nutrients to *Vibrios*. Cellobiose is fermented by some *Vibrios* producing acid and is indicated by the pH indicator bromothymol blue, which turns yellow at acidic pH. Cresol red is the pH indicator of alkaline range, which turns red at alkaline pH. Alkaline pH of the medium enhances the recovery of *Vibrios*.

Blend approximately 25 grams of food sample with 225 ml Alkaline Peptone Water (M618). Transfer a loopful from the surface growth of either Alkaline Peptone Water (M618) or Gelatin Phosphate Salt Broth to the surface of the dried plates of CPC Agar. Streak in a manner that will yield isolated colonies. Incubate CPC Agar at 40 - 42°C for 18 to 24 hours.

Typical colonies of *V.cholerae* on CPC Agar are small, smooth, opaque and green to purple in colour as CPC Agar contains two pH indicators viz. bromothymol blue and cresol red. A purple background will also develop in the CPC Agar upon extended incubation.

### Quality Control

#### Appearance

Light yellow to light brown homogeneous free flowing powder

#### Gelling

Firm, comparable with 1.5% Agar gel

#### Colour and Clarity of prepared medium

Olive-green to light brown coloured, clear to slightly opalescent gel forms in Petri plates

**Reaction**

Reaction of 6.5% w/v aqueous solution at 25°C. pH : 7.6±0.2

**pH**

7.40-7.80

**Cultural Response**

M1241: Cultural characteristics observed after an incubation at 40±2°C for 18-24 hours.

Organism	Inoculum (CFU)	Growth	Recovery	Colour of colony
<b>Cultural Response</b> <i>Vibrio cholerae</i> ATCC 15748	50-100	good - luxuriant	≥50%	green-purple
<i>Vibrio parahaemolyticus</i> ATCC 17802	≥10 <sup>3</sup>	inhibited	0%	
<i>Vibrio vulnificus</i>	50-100	good - luxuriant	≥50%	yellow

**Storage and Shelf Life**

Store below 30°C in tightly closed container and prepared medium below 2-8°C. Use before expiry period on the label.

**Reference**

1. Murray P. R., Baron J. H., Pfaller M. A., Tenover J. C. and Tenover F. C., (Eds.), 2003, Manual of Clinical Microbiology, 8th Ed., American Society for Microbiology, Washington, D.C.
2. Vanderzant C. and Splittstoesser D. F., (Eds), 1992, Compendium of Methods for the Microbiological Examination of Foods, 3rd Ed., APHA, Washington DC.

Revision : 02 / 2016

**Disclaimer :**

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.