



# YT Broth (2X YT Broth)

YT Broth is used for the cultivation of recombinant strains of Escherichia coli .

Composition**	
Ingredients	Gms / Litre
Casein enzymic hydrolysate	16.000
Yeast extract	10.000
Sodium chloride	5.000
Final pH ( at 25°C)	$7.0\pm0.2$
**Formula adjusted, standardized to suit performance parameters	

## **Directions**

Suspend 31.0 grams in 1000 ml distilled water. Heat if necessary to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Mix well and dispense as desired.

## **Principle And Interpretation**

YT Broth is recommended for use in the cultivation of recombinant strains of *Escherichia coli* (1, 2, 4). It is also used in culture of *E. coli* strains for propagation of M 13 bacteriophages (2-4).

These media contain casein enzymic hydrolysate and yeast extract, which supply nitrogenous compounds, vitamin B complex and other essential nutrients and co-factors necessary for the luxuriant growth of recombinant *E. coli* and allows the bacteria to recover from the stress of transformation and grow well. Sodium chloride helps in maintaining isotonic conditions in the medium.

## **Quality Control**

#### Appearance

Light yellow to beige homogeneous free flowing powder

#### **Color and Clarity of Prepared Medium**

Light amber coloured clear solution in tubes.

#### Reaction

Reaction of 3.1% w/v aqueous solution at 25°C. pH : 7.0±0.2

## pH

6.80-7.20

#### **Cultural Response**

M1251: Cultural characteristics observed after an incubation at 35-37°C for 18-24 hours.

Organism	Inoculum (CFU)	Growth
Cultural Response		
Escherichia coli ATCC	50-100	good-luxuriant
23724		
Escherichia coli ATCC	50-100	good-luxuriant
53868		

## **Storage and Shelf Life**

Store below 30°C in tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label.

#### Reference

1. Miller H., 1987, Meth. Enzymol; 152, 145.

2. Ausubel F. M., Brent R., Kingston R. E., Moore B. D., Seidman J. G., Smith J. A. and Strohl K., 1994, Current Protocols in Molecular Biology, Vol. I, Current Protocols, New York, N.Y.

3. Davis L. G., Dibner M. D., Battey J. F., 1986, Basic Methods in Molecular Biology, Elsevier, New York, N.Y.

4. Sambrook J., Fritsch E. E. and Maniatis T., 1989, Molecular Cloning: A Laboratory Manual, 2nd Ed., Cold Spring Harbor Laboratory, Cold Spring Harbor, N.Y.

Revision : 02 / 2015

#### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia<sup>™</sup> publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia<sup>™</sup> Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. A-516, Swastik Disha Business Park, Via Vadhani Ind. Est., LBS Marg, Mumbai-400086, India. Customer care No.: 022-6147 1919 Email: techhelp@himedialabs.com