



## Acetobacter Broth (Mannitol)

M1718

Acetobacter Broth (Mannitol) is recommended as a cultivation media for mannitol positive *Acetobacter* species.

### Composition\*\*

Ingredients	Gms / Litre
Peptic digest of animal tissue	3.000
Yeast extract	5.000
Mannitol	25.000
Final pH ( at 25°C)	7.4±0.2

\*\*Formula adjusted, standardized to suit performance parameters

### Directions

Suspend 33.0 grams in 1000 ml distilled water. Heat if necessary to dissolve the medium completely. Dispense in test tubes and sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes.

### Principle And Interpretation

Acetic acid bacteria are found in fruits with high carbohydrate concentration, which is selective for yeasts, that produce ethanol. This ethanol forms the substrate for acetic acid bacteria and may oxidize ethanol to acetic acid (1). Various synthetic and maintenance media for *Acetobacter* cultures have been cited (2). A typical maintenance medium is Acetobacter Broth (2) Acetobacter Broth is formulated as per Manual of Microbiological Methods (3) and used for the maintenance of *Acetobacter* species utilizing mannitol (4).

Peptic digest of animal tissue, yeast extract in the medium provides nitrogen, vitamins and minerals necessary to support bacterial growth. Mannitol acts as energy source. Calcium carbonate acts as a buffer.

### Quality Control

#### Appearance

Cream to yellow homogeneous free flowing powder

#### Colour and Clarity of prepared medium

Yellow coloured, clear to slightly opalescent solution in tubes

#### Reaction

Reaction of 3.3% w/v aqueous solution at 25°C. pH : 7.4±0.2

#### pH

7.20-7.60

#### Cultural Response

Organism	Inoculum (CFU)	Growth
<b>Cultural Response</b> <i>Acetobacter Hansenii</i> ATCC 35959	50-100	luxuriant
<i>Acetobacter pasteurianus</i> ATCC 6033	50-100	luxuriant

### Storage and Shelf Life

Store below 30°C in tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label.

### Reference

1. Vanderzant C., Splittstoesser D. F., (Eds.), 1992, Compendium of Methods for the Microbiological Examination of Foods, 3rd Ed., APHA, Washington, D. C.
2. Asai, 1968, Univ. of Tokyo Press, Tokyo, Japan and Univ. Park Press, Baltimore, MD.
3. Manual of Microbiological Methods, 1957, Society of American Bacteriologists, McGraw-Hill Book Company, New York.

4. Catalogue of Bacteria and Bacteriophages, 1992, 18th Ed., American Type Culture Collection, Rockville, MD.

Revision : 2 / 2015

**Disclaimer :**

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.