



Anaerobic Agar without Dextrose and Eh Indicator

M229

Anaerobic Agar without Dextrose and Eh Indicator is recommended for the isolation and identification of anaerobic pathogens and for studies of haemolytic activity of Clostridia, Streptococci and other anaerobic organisms.

Composition**

Ingredients	Gms / Litre
Casein enzymic hydrolysate	20.000
Sodium chloride	5.000
Sodium thioglycollate	2.000
Sodium formaldehyde sulphonylate	1.000
Agar	15.000
Final pH (at 25°C)	7.2±0.2

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 43 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Mix well and pour into sterile Petri plates.

Principle And Interpretation

Anaerobic Agar without dextrose and methylene blue (Eh indicator) is used for studies of haemolytic activity of Clostridia, Streptococci and other anaerobes (1). For isolation or cultivation of the highly fermentative butyric types, 1% dextrose may be added prior to sterilization. These media contain sodium thioglycollate and sodium formaldehyde sulphonylate which provides adequate anaerobiosis. Casein enzymic hydrolysate provides essential nutrients while sodium chloride maintains osmotic equilibrium. For isolation or cultivation of the highly fermentative butyric types, 1% dextrose may be added prior to sterilization.

For haemolytic tests anaerobic blood agar plates may be prepared in one of the following ways;

- 1) Sterile blood in about 0.7 ml amount and small inoculum may be mixed with 25-50 ml of cooled agar and mixture is poured into the Petri plate filling it up to 3/4. After solidification the lid is replaced with Brewer Anaerobic Petri plate cover,
- 2) An ordinary sterile Blood Agar plate (made from Blood Agar Base or Tryptone Soya Agar) may be streaked with a culture. Melted and cooled Anaerobic Agar without Dextrose is then poured over the Blood Agar to provide the proper depth. After solidification the lid is replaced with anaerobic Petri plate cover.

The anaerobic cover should not rest on the Petri plate bottom: its inner ridge should seal the agar, and the medium within the ridge should not touch the cover at any point. The medium should be cherry red in colour after addition of blood.

Quality Control

Appearance

Cream to yellow homogeneous free flowing powder

Gelling

Firm, comparable with 1.5% Agar gel

Colour and Clarity of prepared medium

Yellow to light amber coloured clear to slightly opalescent gel forms in Petri plates

Reaction

Reaction of 4.3% w/v aqueous solution at 25°C. pH : 7.2±0.2

pH

7.00-7.40

Cultural Response

M229: Cultural characteristics observed under anaerobic condition after an incubation at 35-37°C for 18-48 hours.

Organism	Inoculum (CFU)	Growth	Recovery
Cultural Response			
<i>Clostridium butyricum</i> ATCC 13732	50-100	good-luxuriant	>=50%
<i>Clostridium perfringens</i> ATCC 12919	50-100	good-luxuriant	>=50%
<i>Clostridium sporogenes</i> ATCC 11437	50-100	good-luxuriant	>=50%

Storage and Shelf Life

Store below 30°C in tightly closed container and the prepared medium at 2 - 8°C. Use before expiry date on the label.

Reference

1. Vera J., 1942, J. Bact., 44:497.

Revision : 1 / 2011



Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.