

Technical Data

Legionella Selective Supplement III

An antimicrobial supplement recommended for the selective isolation of *Legionella* species from mixed cultures.

Composition

Per vial sufficient for 500 ml medium

*Ingredients

Vancomycin Glycine Polymyxin B sulphate

Directions:

Rehydrate the contents of one vial aseptically with 5 ml of 50% ethanol. Mix well and aseptically add to 500 ml of sterile, molten, cooled (45-50°C) Buffered Charcoal Yeast Extract Agar Base <u>M813</u> / Modified Buffered Charcoal Agar Base <u>M892</u> / Modified Buffered Charcoal HiVegTM Agar Base <u>MV892</u> alongwith rehydrated contents of Legionella Selective Supplement IV <u>FD040</u> and Legionella Supplement <u>FD041</u>. Mix well and pour into sterile petri plates.

Concentration

2.500mg 1.500mg

50000Unit

Storage and Shelf Life

Store at 2 - 8°C. Use before expiry date on the label.

Revision : 1 / 2012

* Not For Medicinal Use

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia[™] publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia[™] Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. A-516,Swastik Disha Business Park,Via Vadhani Ind. Est., LBS Marg, Mumbai-400086, India. Customer care No.: 022-6147 1919 Email: techhelp@himedialabs.com Website: www.himedialabs.com

FD038