



Technical Data

Esculin (0.5gms per vial)

FD050

This supplement is recommended for detection of group D streptococci by means of esculin hydrolysis.

Composition

Per vial sufficient for medium

Ingredients

Esculin

Concentration

0.500g

Directions:

Rehydrate the content of one vial in 2 ml distilled water. Mix well and add it to 500 ml of sterile, molten, cooled (45-50°C) Bile Esculin Agar Base [M340](#) / Bile Esculin HiVeg™ Agar Base [MV340](#) . Mix well and sterilize the medium by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Mix well and pour into sterile petri plates.

Storage and Shelf Life

Store between 10-30°C. Use before expiry date on the label.

Revision : 1 / 2012

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.