



# Technical Data

## Mueller Tellurite Serum (25 ml per vial)

FD100

This supplement is recommended for isolation, cultivation and differentiation of *Corynebacterium diphtheria* .

### Composition

Per vial sufficient for 975 ml medium

#### Ingredients

Potassium tellurite  
Calcium D pantothenate  
Horse serum  
Sodium lactate solution  
Ethanol

#### Concentration

100mg  
0.050mg  
12.500ml  
10ml  
2.500ml

### Directions:

Warm up the refrigerated solution to room temperature. Add the contents of 1 vial to 975 ml of sterile, molten, cooled (45-50°C)

Mueller Tellurite Agar Base [M1202](#) . Mix well and pour into sterile petri plates.

### Storage and Shelf Life

On receipt product should be stored at -20°C. Use before the expiry date on the label.

Revision : 1 / 2012

#### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.