



# Technical Data

## \*\*L. mono Enrichment Supplement II

FD227

Recommended for selective differentiation of *Listeria monocytogenes* from other Listeria species.

### Composition

Per vial sufficient for 500 ml medium

#### Ingredients

L-phosphatidylinositol  
Distilled water

#### Concentration

0.500g  
15ml

### Directions:

Thaw the contents of 1 vial of L. mono Enrichment Supplement II at room temperature. Aseptically add the sterile contents of one vial to 470 ml of sterile, molten, cooled (45-50°C) L. mono Confirmatory Agar Base [M1552](#) / L. mono Confirmatory HiVeg™ Agar Base [MV1552](#) along with sterile rehydrated contents of one vial each of L. mono Selective Supplement I [FD212](#) and L. mono Selective Supplement II [FD213](#). Mix well and pour into sterile petri plates.

### Storage and Shelf Life

On receipt product should be stored at -20°C. Use before the expiry date on the label.

Revision : 1 / 2012

#### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.