



Technical Data

Bifido Selective Supplement B

FD251

A selective supplement recommended for the isolation of Bifidobacterium species by colony count technique from milk products.

Composition

Per vial sufficient for 500 ml medium

Ingredients

Glacial acetic acid (1%)

Concentration

5ml

Directions:

Aseptically add the contents of 1 vial in 500 ml sterile, molten cooled (45-50°C) Bifidobacteria Selective Count Agar Base (BSC Propionate Agar Base) [M1734](#) along with 1 vial of Bifido Selective Supplement A [FD250](#) . Mix well and pour into sterile Petri plates.

Note : After addition of FD251 pH of the medium will shift to the acidic side which does not affect the performance of the medium.

Storage and Shelf Life

Store at 2-8°C. Use before expiry date on label.

Revision : 1 / 2012

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.