



# Technical Data

## Salmonella Selective Enrichment Supplement

FD275

Selective supplement recommended for Selective isolation and differentiation of Salmonella species from coliforms by chromogenic method.

### Composition

Per vial sufficient for 1000 ml medium

#### \*Ingredients

Novobiocin

#### Concentration

12mg

### Directions:

Rehydrate the contents of 1 vial aseptically with 5ml of sterile distilled water. Mix gently to dissolve the contents completely. Aseptically add the rehydrated contents of one vial to 1000ml of sterile, cooled (45-50°C) of Salmonella Selective Enrichment Broth Base [M1843](#) . Mix well and dispense as desired.

### Storage and Shelf Life

Store at 2 - 8°C. Use before expiry date on the label.

\* Not For Medicinal Use

#### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.