



Hayflick Supplement

FD300

A selective supplement recommended for detection of *Mycoplasma* species.

Composition

(per vial, sufficient for 500 ml medium)

Ingredients

Penicillin 19770 IU

Directions

Rehydrate the contents of one vial aseptically with 5ml sterile distilled water. Mix well and aseptically add to 416 ml sterile, cooled (45-50⁰C) Hayflick Broth Base (ME1885) or molten Hayflick Agar Base, ME1886 cooled to 45-50⁰C. Then add 79 ml of unheated Horse serum (RM1239). Mix well and dispense as desired.

Unused rehydrated contents can be stored at -20°C for future use.

Disclaimer :

User must ensure of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ Publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.