



# Technical Data

## Modified L.mono Selective Supplement

FD333

Modified L.mono Selective Supplement is recommended for selective isolation of *L.monocytogenes* and *Listeria* species based on PCPLC activity.

### Composition

Per vial sufficient for 500 ml medium

Ingredients	Concentration/vial
Nalidixic acid	0.10gm
Polymyxin B Sulfate	0.000004564gm
Ceftazidime	0.10gm
Amphotericin	0.005gm

### Directions:

Rehydrate the content of one vial aseptically with 5 ml of sterile distilled water. Mix well and aseptically add it to 480 ml of sterile molten, cooled (45-50°C) HiCrome L.mono Differential Agar Base (M2009) along with sterile contents of 1 vial of Lecithin solution(FD332). Mix well and pour into sterile Petri plates.

### Storage and Shelf Life

Store at 2-8°C. After rehydration, it is advisable to use immediately within one hour of preparation. Use before the expiry date on the label.

Revision : 00 / 2016

#### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.