



Aspergillus Differentiation Medium Base

M1127

Aspergillus Differentiation Medium Base is used for detection of aflatoxin producing *Aspergillus* species from food samples.

Composition**

Ingredients	Gms / Litre
Peptic digest of animal tissue	10.000
Yeast extract	20.000
Ferric ammonium citrate	0.500
Dichloran	0.002
Agar	15.000
Final pH (at 25°C)	6.3±0.2

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 22.75 grams in 500 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C and aseptically add sterile rehydrated contents of 1 vial of Chloramphenicol Selective Supplement (FD033). Mix well and pour into sterile Petri plates.

Principle And Interpretation

Aspergilli are hyaline moulds that commonly cause opportunistic infections in humans. Allergic bronchopulmonary disease is a manifestation of hypersensitivity to fungal spores or products, a common manifestation of *Aspergillus* species (particularly *A. flavus*) (1). Aspergillus Differentiation Medium Base formulated by Pitt et al (2) is a modification of the medium formulated by Bothast and Fennel (3). *Aspergillus flavus* develops intense yellow orange colour at the base of the colonies, which is a differential characteristic of this species. This pigmentation helps in differentiating *A. flavus* from other *Aspergillus* species (3-5). Assante et al (6) showed that the orange yellow coloration was due to the reaction of ferric ions (from ferric ammonium citrate) with aspergillilic acid or neoaspergillilic acid forming a colored complex.

A mixture of chloramphenicol and dichloran restricts the spreading of moulds. It also inhibits bacterial growth and helps in the identification of fungi. Peptic digest of animal tissue and yeast extract serve as sources of nitrogen, amino acids and B complex vitamins. Ferric ammonium citrate aids in the production of yellow orange pigment characteristic of *A. flavus*. *A. parasiticus*, associated with aspergillosis also produces a yellow orange pigment similar to the one produced by *A. flavus* (7).

Quality Control

Appearance

Cream to yellow homogeneous free flowing powder

Gelling

Firm, comparable with 1.5% Agar gel

Colour and Clarity of prepared medium

Medium amber coloured clear to slightly opalescent gel forms in Petri plates

Reaction

Reaction of 4.55% w/v aqueous solution at 25°C. pH : 6.3±0.2

pH

6.10-6.50

Cultural Response

M1127: Cultural characteristics observed with added 1 vial of Chloramphenicol Selective Supplement (FD033) after an incubation at 25-30°C for 48-72 hours.

Organism	Inoculum (CFU)	Growth	Colour of colony
----------	-------------------	--------	---------------------

Cultural Response

* <i>Aspergillus brasiliensis</i> ATCC 9642	50-100	good-luxuriant	pale yellow colour on the reverse side of colonies with black heads on the top of the colonies
<i>Aspergillus flavus</i> ATCC 22547	50-100	good-luxuriant	yellowish orange colour on the reverse side of colonies
<i>Aspergillus parasiticus</i> ATCC 28285	50-100	good-luxuriant	yellowish orange colour on the reverse side of colonies

*Key: Formerly known as *Aspergillus niger*

Storage and Shelf Life

Store below 30°C in tightly closed container and the prepared medium at 2 - 8°C. Use before expiry date on label.

Reference

1. Koneman E. W., (Ed.), Mycology, Colour Atlas and Textbook of Diagnostic Microbiology, 4th Ed, 1992, J. B. Lippincott Company.
2. Pitt J., Hocking D., and Glenn D. R., 1983
3. Bothast and Fennel, 1974, Mycologia. 66:365.
4. Haley and Callaway, 1978, Laboratory methods in medical mycology, 4th Ed., Center for Disease Control, Atlanta, Ga.
5. McGinnis, 1980, Laboratory Handbook of Medical Mycology, Academic Press, New York, N.Y.
6. Assante G. et al., 1981, J. Ag. Food Chem., 29:785
7. Murray P. R., Baron E. J., Jorgensen J. H., Pfaller M. A., Tenover F. C., Tenover J. C., (Eds.), 8th Ed., 2003, Manual of Clinical Microbiology, ASM, Washington, D.C.

Revision : 2 / 2015

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.