



Bryant and Burkey Agar

M1690

Bryant and Burkey Agar is used for detecting and enumerating spores of lactate fermenting *Clostridium* in milk and dairy products.

Composition**

Ingredients	Gms / Litre
Casein enzymic hydrolysate	15.000
Yeast extract	5.000
Beef extract	7.500
Sodium acetate	5.000
L-Cystine hydrochloride	0.600
Resazurin	0.0025
Sodium thioglycollate	0.200
Agar	0.750
Final pH (at 25°C)	5.9±0.2

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 34.05 grams in 1000 ml distilled water containing 3 grams of Sodium lactate. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes.

Note : If more than the upper one-third of the medium has acquired a pink colour, the medium may be restored once by heating in a water bath or in free flowing steam until the pink colour disappears.

Principle And Interpretation

Bryant and Burkey Agar is used to enumerate the spores of gasogenic clostridia that are responsible for swelling of cheese in dairy industry (1, 2). In normal conditions of use, the medium allows the growth of other microorganisms also which are not directly related to cheese alteration, e.g. *Clostridium sporogenes* or *Clostridium butyricum*. Germination and growth of butyric acid bacteria (BAB) causes severe defects in cheese with silage being the main source of BAB spores in cheese milk (4). Clostridia spores are heat resistant and therefore can contaminate cheese brines. The gas produced by the growth of clostridia swells the cheese and is responsible for defect known as butyric swelling, resulting in bad taste. The main species causing this butyric swelling defect is *Clostridium tyrobutyrium*.

Bryant and Burkey Agar is the modification of Bryant and Burkey Medium with addition of little amount of agar that makes the medium thicker to be used in greater volumes. Addition of sodium thioglycollate maintains anaerobic conditions and also the amount of lactate is reduced so that the density of the medium will be retained easily for the production of gas bubbles.

Resazurin is a redox indicator and monitors the oxygen level turning from pink under aerobic conditions to colourless under anaerobic conditions. The nutrient composition of the basal medium, particularly of casein enzymic hydrolysate, yeast extract, beef extract and L-cystine HCl help in rapid growth of *Clostridium* species. Sodium acetate promotes spore germination (3) and improves the selectivity of the medium. Sodium lactate is fermented under anaerobic conditions by *Clostridium tyrobutyricum*. A pink colour indicates the presence of oxygen in the media which turns colourless on boiling.

Quality Control

Appearance

Cream to yellow homogeneous free flowing powder

Colour and Clarity of prepared medium

Light amber coloured, clear to slightly opalescent solution with upper 10% or less medium pink on standing.

Reaction

Reaction of 3.4% w/v aqueous solution at 25°C. pH : 5.9±0.2

pH

5.70-6.10

Cultural Response

M1690: Cultural characteristics observed under anaerobic condition, after an incubation at 35-37°C for 6 days.

Organism	Inoculum (CFU)	Growth	Gas production
Cultural Response <i>Clostridium tyrobutyricum</i> ATCC 25755	50-100	good	positive reaction

Storage and Shelf Life

Store below 30°C in tightly closed container and the prepared medium at 2 - 8°C. Use before expiry date on the label.

Reference

1. Bryant M. P. and Burkey L. A., 1956, J. Bacteriol., 71: 43-46.
2. Bergeres J. L. and Sivela S., 1989, Detection and enumeration of Clostridial spores related to cheese quality-Classical and new methods, FIL-IDF Bull. 51:18-23.
3. Touraille C. and Bergere J. L., 1974, Biochimie, 56: 404-422.
4. Vissers M. M. M., Drichuis F., Giffel T., John P. D., Lankveld J. M. G., 2007, J. Dairy Sci., 90:928

Revision : 2 / 2015

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.