



Eugonic Broth

M429

Eugonic broth is recommended for the cultivation of fastidious microorganisms like *Haemophilus*, *Nesseria*, *Pasteurella*, *Brucella* and *Lactobacillus* species.

Composition**

Ingredients	Gms / Litre
Casein enzymic hydrolysate	15.000
Papaic digest of soyabean meal	5.000
Dextrose	5.000
Sodium chloride	4.000
Sodium sulphite	0.200
L-Cystine	0.200
Final pH (at 25°C)	7.0±0.2

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 29.4 grams in 1000 ml distilled water. Heat if necessary with frequent stirring to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 40-45°C and add 5 -10% v/v sterile defibrinated blood if desired. The blood may be chocolated by heating.

Principle And Interpretation

Eugonic Broth was developed by Pelczar and Vera (1) for cultivation of fastidious organisms like *Brucella*. This medium can also be used to grow *Mycobacteria* and various pathogenic fungi including *Nocardia*, *Histoplasma* and *Blastomyces*, when enriched with blood. Niven used this media for detection of spoilage of meats (2). Eugonic Broth was developed to obtain eugonic (luxuriant) growth of fastidious microorganisms like *Brucella* which are otherwise difficult to cultivate (5). The unenriched medium supports rapid growth of lactobacilli associated with cured meat products, dairy products and other foods. APHA recommends Eugonic Broth, which is also used in germinating anaerobic spores pasteurized at 104°C (3, 4). Organisms like *Bordetella* and *Neisseria* proliferate in Eugonic Broth because large amount of sulfur and carbon sources have been added in the formulation. Therefore Eugonic Broth is recommended for the direct isolation of *Bordetella pertussis* and *Neisseria meningitides* from the test materials such as throat mucus, blood, cerebrospinal fluid, pleural fluid and other specimens. For the isolation of *Bacillus pumilus*, Eugonic Broth can be supplemented with 0.1% starch, prior to sterilization (5)

Casein enzymic hydrolysate and papaic digest of soyabean meal provide the nitrogen, vitamins and amino acids, which supports the growth of fastidious microbial species. The high concentration of dextrose is the energy source for rapid growth of bacteria. L-Cystine and sodium sulphite are added to stimulate growth. Sodium chloride maintains the osmotic balance of the media. The high carbohydrate content along with high sulfur (cystine) content improves growth with chromogenicity (4).

Quality Control

Appearance

Cream to yellow homogeneous free flowing powder

Colour and Clarity of prepared medium

Yellow coloured, clear solution in tubes

Reaction

Reaction of 2.94% w/v aqueous solution at 25°C. pH : 7.0±0.2

pH

6.80-7.20

Cultural Response

M429: Cultural characteristics observed with added 5-10% v/v sterile defibrinated blood after an incubation at 35-37°C for 48 hours (fungal cultures incubated at 25-30°C).

Organism	Inoculum (CFU)	Growth
<i>Bacillus pumilus</i> ATCC 14884	50-100	good (with 0.1% starch)
<i>Brucella abortus</i> ATCC 43155	50-100	good (under 3-5% CO ₂)
<i>Candida albicans</i> ATCC 26790	50-100	good
<i>Lactobacillus fermentum</i> ATCC 9338	50-100	good
<i>Neisseria meningitidis</i> ATCC 13090	50-100	good
<i>Streptococcus pneumoniae</i> ATCC 6303	50-100	luxuriant (under 3-5% CO ₂)
<i>Streptococcus pyogenes</i> ATCC 19615	50-100	luxuriant (under 3-5% CO ₂)

Storage and Shelf Life

Store below 30°C in tightly closed container and the prepared medium at 2 - 8°C. Use before expiry date on the label.

Reference

1. Pelczar and Vera J., 1949, Milk Plant Monthly 38:30
2. Niven C. F., Castellani A. G., and Allanson V., 1949, J. Bacteriol., 58:633.
3. Downes F. P. and Ito K., (Ed.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 4th Ed., American Public Health Association, Washington, D.C.
4. Frank H. A., 1955, J. Bacteriol., 70:269.
5. MacFaddin J. F., 1985, Media for Isolation-Cultivation-Identification-Maintenance of Medical Bacteria, Vol. 1, Williams & Wilkins, Baltimore, Md.

Revision : 1 / 2011



Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.