



Lactic Agar

M599

Lactic Agar is recommended for enumeration and identification of lactic Streptococci and Lactobacilli by pour plate technique.

Composition**

Ingredients	Gms / Litre
Casein enzymic hydrolysate	20.000
Yeast extract	5.000
Gelatin	2.500
Dextrose	5.000
Lactose	5.000
Sucrose	5.000
Sodium chloride	4.000
Sodium acetate	1.500
Ascorbic acid	0.500
Agar	15.000

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 63.5 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Mix well and pour into sterile Petri plates.

Principle And Interpretation

Lactic Agar was developed by Elliker et al (1) and recommended by APHA (2) for cultivation of lactic bacteria to promote the colony development of Lactobacilli and lactic Streptococci. Samples are analyzed by pour plate technique. Lactic acid bacteria are fastidious in nature and hence Lactic Agar is designed to satisfy their growth requirement. Lactic acid bacteria survive at low pH, but are very sensitive to other adverse conditions.

Samples to be examined for enumeration of viable lactic acid bacteria should not be frozen prior to analysis. Many of the lactic acid bacteria are easily killed or injured by freezing. For dilution of products it is best to use sterile 0.1% Peptone Water (M028) as the diluent since it protects bacteria during the dilution process (3, 4).

Casein enzymic hydrolysate and yeast extract provide amino acids, other nitrogenous nutrients, vitamin B complex etc. Dextrose, lactose and sucrose are the fermentable carbohydrates. Ascorbic acid provides vitamin C required by lactic acid bacteria. Sodium chloride maintains the osmotic equilibrium of the medium. Sodium acetate inhibits contaminating bacteria and restricts the swarming of lactic acid bacteria. Upon incubation, the colonies are examined for gram staining and catalase production. Gram-positive, catalase-negative cocci or rods are tentatively considered to be lactic acid bacteria (2).

Quality Control

Appearance

Cream to yellow homogeneous free flowing powder

Gelling

Firm, comparable with 1.5% Agar gel

Colour and Clarity of prepared medium

Yellow coloured clear to slightly opalescent gel forms in Petri plates

Cultural Response

M599: Cultural characteristics observed after an incubation at 35-37°C for 18-48 hours.

Organism	Inoculum (CFU)	Growth	Recovery
----------	-------------------	--------	----------

Cultural Response

<i>Lactobacillus bulgaricus</i> ATCC 11842	50-100	good-luxuriant	>=50%
<i>Lactobacillus casei</i> ATCC 9595	50-100	good-luxuriant	>=50%
<i>Lactobacillus lactis</i> ATCC 8000	50-100	good-luxuriant	>=50%
<i>Streptococcus cremoris</i> ATCC 19257	50-100	good-luxuriant	>=50%
<i>Streptococcus thermophilus</i> ATCC 14485	50-100	good-luxuriant	>=50%

Storage and Shelf Life

Store below 30°C in tightly closed container and the prepared medium at 2 - 8°C. Use before expiry date on the label.

Reference

1. Elliker P. R., Anderson A. W. and Hanesson G., 1956, J. Dairy Science, 39:1611.
2. Downes F. P. and Ito K., (Eds.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 4th Ed., APHA, Washington, D.C.
3. Hartman P. A., and Huntsberger D. V., 1961, Appl. Microbiol., 9-324. Jayne-Williams D. J., 1963, J. Appl. Bacteriol., 26:398

Revision : 2 / 2015

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.