



## Liquoid Broth

M817

Liquoid Broth is used for screening of blood specimens from suspected bacteremic cases.

### Composition\*\*

| Ingredients                   | Gms / Litre |
|-------------------------------|-------------|
| Calf brain, infusion from     | 200.000     |
| Beef heart, infusion from     | 250.000     |
| Proteose peptone              | 10.000      |
| Sodium chloride               | 5.000       |
| Disodium phosphate            | 2.500       |
| Dextrose                      | 2.000       |
| Sodium polyanethol sulphonate | 0.500       |
| Final pH ( at 25°C)           | 7.4±0.2     |

\*\*Formula adjusted, standardized to suit performance parameters

### Directions

Suspend 37.5 grams in 1000 ml distilled water. Heat if necessary to ensure complete solution. Dispense into bottles or tubes and sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. If desired, 1 gm/litre agar can be added to encourage growth of anaerobic organisms. For best results, use the medium on the day it is prepared otherwise boil or steam it to remove dissolved oxygen just before use.

### Principle And Interpretation

In most bacteriemic conditions in man, the organisms are not numerous. Therefore for their demonstration by blood culture, relatively large amount of blood e.g. 5-10 ml should be used as inoculum. As the bloods natural bactericidal or bacteriostatic action may interfere with the growth of any bacteria present, diluting the blood with medium should annul this effect. The technology of blood culture was revised by Gould and Duerden (1). Upto 10 ml or more blood may be added to 100 ml of broth without a detectable antibacterial effect. The antibacterial effect may be further prevented by incorporation of substances such as sodium polyanethol sulphonate (SPS). Liquoid Broth is used for the culturing of blood specimens from suspected bacteremia cases (2). Liquoid (Sodium polyanethol sulphonate) is a good anticoagulant. Moreover it is not inhibitory and has the added advantage of annulling the natural bactericidal action of blood (3).

The medium is composed of rich ingredients for blood culture. Beef heart and calf brain infusion and proteose peptone provide the necessary carbonaceous and nitrogenous nutrients, vitamins and growth factors to the organisms. Dextrose is the carbon source and sodium chloride maintains the osmotic equilibrium of the medium. It is advisable to seed more than one medium for blood culture. One of each set of bottles should be incubated in an atmosphere of air with 10% CO<sub>2</sub>. It is essential to loosen the caps of bottles during incubation. Growth may produce a generalized turbidity; make subculture from all bottles to solid media.

### Quality Control

#### Appearance

Cream to yellow homogeneous free flowing powder

#### Colour and Clarity of prepared medium

Light amber coloured, clear solution without any precipitate

#### Reaction

Reaction of 3.75% w/v aqueous solution at 25°C. pH : 7.4±0.2

#### pH

7.20-7.60

#### Cultural Response

M817: Cultural characteristics observed after an incubation at 35-37°C for 24-48 hours.

---

| Organism                                 | Inoculum (CFU) | Growth         |
|--|----------------|----------------|
| <i>Escherichia coli</i> ATCC 25922       | 50-100         | luxuriant      |
| <i>Salmonella Typhi</i> ATCC 6539        | 50-100         | luxuriant      |
| <i>Staphylococcus aureus</i> ATCC 25923  | 50-100         | luxuriant      |
| <i>Streptococcus pyogenes</i> ATCC 19615 | 50-100         | good-luxuriant |

### Storage and Shelf Life

Store below 30°C in tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label.

### Reference

1. Gould J. C., Duerden B. I., 1983 (Ed.), J. Clin. Pathol., 36: 963-977
2. Collee J. G., Fraser A. G., Marmion B. P., Simmons A., (Eds.), Mackie and McCartney, Practical Medical Microbiology, 1996, 14th Edition, Churchill Livingstone
3. Von Haebler T., Miles A. A., The Journal of Pathology and Bacteriology, Vol. 46, Issue 2, Pages 245- 252.

Revision : 1 / 2011



#### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.