



A C Broth

M875

A C Broth is recommended for cultivation of common aerobes and sterility testing of solutions and biological products without mercurial preservatives.

Composition**

Ingredients	Gms / Litre
Proteose peptone	20.000
Beef extract	3.000
Yeast extract	3.000
Malt extract	3.000
Dextrose	5.000
Ascorbic acid	0.200
Final pH (at 25°C)	7.2±0.2

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 34.2 grams in 1000 ml of distilled water. Heat if necessary to dissolve the medium completely. Distribute in tubes or bottles to give the desired depth and sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes.

If the medium is not used on same day, it is advisable to drive off dissolved gases by boiling or steaming in the autoclave and cool without agitation.

Principle And Interpretation

AC Broth support an early and luxuriant growth of aerobic, anaerobic and microaerophilic microorganisms . Many pathogenic and saprophytic aerobes can also be cultivated using AC Broth (1). This medium can also be used for sterility testing of solutions and biological products not containing mercurial preservatives. Some of the media containing sodium thioglycollate exhibit toxicity for some organisms. This toxicity is not seen in the case of AC Broth as reported by Christensen (2) and Malin and Finn (3). Earlier studies performed have reported the usefulness of using this medium for the cultivation of a wide variety of organisms (4, 5).

Proteose peptone, beef extract, yeast extract and malt extract serve as the carbon and nitrogen sources in addition to being a source of vitamins and cofactors. Dextrose serves as the fermentable carbohydrate and source of energy. Ascorbic acid in the media helps to improve the clarity of the medium.

Quality Control

Appearance

Cream to yellow homogeneous free flowing powder

Colour and Clarity of prepared medium

Medium amber coloured clear to slightly opalescent solution.

Reaction

Reaction of 3.42% w/v aqueous solution at 25°C. pH : 7.2±0.2

pH

7.00-7.40

Cultural Response

M875: Cultural characteristics observed after an incubation at 35-37°C for 18-24 hours (*Clostridium* species incubated anaerobically).

Organism	Inoculum (CFU)	Growth
Cultural Response <i>Clostridium perfringens</i> ATCC 12919	50-100	luxuriant

<i>Escherichia coli</i> ATCC 25922	50-100	luxuriant
<i>Neisseria meningitidis</i> ATCC 13090	50-100	luxuriant
<i>Staphylococcus aureus</i> ATCC 25923	50-100	luxuriant
<i>Streptococcus mitis</i> ATCC 9811	50-100	luxuriant
<i>Streptococcus pneumoniae</i> ATCC 6303	50-100	luxuriant

Storage and Shelf Life

Store below 30°C in tightly closed container and the prepared medium at 2-8°C. Use before expiry date on label.

Reference

1. MacFaddin J. F., 1985, Media for Isolation-Cultivation-Identification-Maintenance of Medical Bacteria, Vol. I. Williams & Wilkins, Baltimore, Md.
2. Christensen, 1944, Paper read at New York Meeting, American Public Health Association.
3. Malin and Finn, 1951, J. Bacteriol., 62:349.
4. Reed and Orr, 1943, J. Bacteriol., 45:309.
5. Schneiter, Dunn and Caminita, 1945, Public Health Rep., 60:789.

Revision : 2 / 2015

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.