



Campylobacter Enrichment Broth Base (Preston Enrichment Broth Base)

M899

Campylobacter Enrichment Broth Base is recommended for selective enrichment and cultivation of *Campylobacter* species.

Composition**

Ingredients	Gms / Litre
Peptic digest of animal tissue	10.000
Beef extract	10.000
Sodium chloride	5.000
Final pH (at 25°C)	7.5±0.2

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 12.5 grams in 470 ml distilled water. Heat if necessary to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to room temperature and aseptically add sterile 25 ml lysed horse blood and reconstituted contents of 1 vial of Campylobacter Supplement IV (Preston Selective Supplement) (FD042). Mix well and dispense as desired.

Principle And Interpretation

Balton and Robertson (1) described this as a selective medium for the cultivation of *Campylobacter* species. It is recommended by APHA (2) for enrichment of thermotolerant *Campylobacter* species from foods. Preston Enrichment Broth has a rich basal medium to aid resuscitation of sublethally damaged *Campylobacter*. Preliminary incubation of the medium complete with antibiotics for 4 hours at 37°C was recommended to aid resuscitation of injured organisms followed by 42°C for 18- 48 hours (3).

Peptic digest of animal tissue and beef extract in the medium provide nitrogen, vitamins and minerals necessary to support bacterial growth. Sodium chloride provides essential ions.

Preston Selective Supplement (FD042) contains antibacterial and antifungal agents. Polymyxin B is active only against gram-negative bacteria and *Proteus* species are sometimes resistant. Trimethoprim usually inhibits *Proteus* species as well as other gram-negative bacteria. Rifampicin is also active against gram-negative organisms. Cycloheximide acts as antifungal agent.

Direct plating without enrichment is adequate for fresh faecal samples, fecal contents or intestinal specimens as high numbers of the organisms may be anticipated. For food sample s enrichment is required. Humphrey (1989) suggested that pre-enrichment at 37°C should be continued for 4 hours and that addition of all antibiotics should be delayed until the 4 hours pre- enrichment had been completed. Enrichment medium with rifampicin was recommended in parallel with similar plating medium.

The *Campylobacter* species grow well in microaerobic conditions i.e. in 5% O₂ at 42°C in about 48 hours. Addition of about 4 drops of glycerol to a filter paper kept within the jar/container will hamper confluent and swarming growth of *Campylobacter* (3).

Quality Control

Appearance

Cream to yellow homogeneous free flowing powder

Colour and Clarity

Basal medium :Light yellow coloured clear solution. After addition of sterile lysed horse blood : Cherry red coloured opaque solution in tubes

Reaction

Reaction of 2.5% w/v aqueous solution at 25°C. pH : 7.5±0.2

pH

7.30-7.70

Cultural Response

Cultural characteristics observed with added 25ml sterile lysed horse blood and Campylobacter Supplement IV (Preston Selective Supplement),(FD042), after an incubation at 42°C for 48 hours(5% O₂ + 10% CO₂ + 85% N₂).

Cultural Response**Organism****Growth****Cultural Response**

Bacillus cereus ATCC 10876 inhibited

Campylobacter coli ATCC 33559 good-luxuriant

Campylobacter jejuni ATCC 29428 good-luxuriant

Campylobacter lari ATCC 35221 good-luxuriant

Escherichia coli ATCC 25922 inhibited

Proteus mirabilis ATCC 25933 inhibited

Staphylococcus aureus ATCC 25923 inhibited

Storage and Shelf Life

Store below 30°C in tightly closed container and the prepared medium at 2 - 8°C. Use before expiry date on the label.

Reference

1. Balton F.J. and Robertson L., 1982, J. Clin. Pathol., 35:462.
2. Vanderzant C. and Splittstoesser D. (Eds.), 1992, Compendium of Methods for the Microbiological Examination of Foods, 3rd ed., APHA, Washington, D.C.
3. Humphrey T. J., 1989, J. Appl. Bacteriol. 66, 119-126

Revision : 1 / 2011

**Disclaimer :**

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.