



Islams Medium Base for Group B Streptococci

M998

Islams Medium Base for Group B Streptococci is used for identification and cultivation of group B Streptococci from clinical specimens.

Composition**

Ingredients	Gms / Litre
Proteose peptone	23.000
Starch, soluble	5.000
Monosodium phosphate	1.482
Disodium phosphate	5.749
Agar	10.000
Final pH (at 25°C)	7.4±0.2

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 45.23 grams in 950 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 115°C for 10 minutes. Cool to 50°C and aseptically add 50 ml sterile inactivated horse serum (RM1239), (inactivated by heating at 56°C for 30 minutes). Mix well and pour into sterile Petri plates.

Principle And Interpretation

Islam (1) formulated this medium to exploit the ability of most Group B Streptococci to produce orange/red-pigmented colonies when incubated under anaerobic conditions. Lancefield first noted carotenoid pigmentation, characteristic of group B Streptococci when incubated under anaerobic conditions (2). This medium also supports growth of other genital bacteria that cause neonatal infection (1) such as anaerobic *Streptococcus*, *Bacteroides* and *Clostridium* species.

Proteose peptone provides the necessary nutrients for the growth of Group B Streptococci. Disodium and monosodium phosphates provide buffering to the medium.

Pigmentation can be enhanced by adding trimethoprim / sulphonamides (3). No inhibition of growth occurs and the pigmentation is seen clearly over a radius of 10-20 mm. The medium must have the correct pH to ensure good pigmentation but some strains of Group B Streptococci do not produce pigmented colonies (4). Other organisms that can grow on this medium do not produce the characteristic orange-red pigment. Inoculate the specimen swab onto the surface of Islams Medium. If desired, apply a disc containing 300 or 500µg of sulphafurazole onto an area of the plate where growth can be expected to be moderately profuse. Incubate the plates anaerobically at 35°C for 24 to 48 hours.

Quality Control

Appearance

Cream to yellow homogeneous free flowing powder

Gelling

Firm, comparable with 1.0% Agar gel.

Colour and Clarity of prepared medium

Yellow coloured, clear to slightly opalescent gel forms in Petri plates

Reaction

Reaction of 4.5% w/v aqueous solution at 25°C. pH : 7.4±0.2

pH

7.20-7.60

Cultural Response

M998: Cultural characteristics observed with added sterile inactivated horse serum (RM1239), after an incubation at 35-37°C for 24-48 hours.

Organism	Inoculum (CFU)	Growth	Recovery	Pigmentation
<i>Bacteroides fragilis</i> ATCC 25285	50-100	fair-good	40-50%	no pigmentation
<i>Enterococcus faecalis</i> ATCC 29212	50-100	luxuriant	>=70%	no pigmentation
<i>Streptococcus agalactiae</i> ATCC 13813	50-100	luxuriant	>=70%	orange/red

Storage and Shelf Life

Store below 30°C in tightly closed container and prepared medium at 2-8°C. Use before expiry period on the label.

Reference

1. Islam A. K. M. S., 1927, Lancet, i: 256 (letter).
2. Merrit K. and Jacobs N. J., 1978, J. Clin. Microbiol., 8:105.
3. de al Rosa M., Villareal R., Vega D., Miranda C. and Martinezbrocal A., 1983, J. Clin. Microbiol., 18:779.
4. Islam A. K. M. S., 1981, J. Clin. Pathol., 34:78.

Revision : 2 / 2015

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.