

Technical Data

Legionella BMPA Selective Supplement

FD144

An antibiotic supplement recommended for selective isolation of Legionella species.

Composition

Per vial sufficient for 450 ml medium

*Ingredients Concentration
Cefamandole 0.400mg
Polymyxin B 8000IU
Anisomycin 8mg

Directions:

Rehydrate the contents of 1 vial aseptically with 5 ml of 50% ethanol. Mix well and aseptically add to 450 ml of sterile, molten, cooled (45-50°C) Legionella Agar Base M809A along with 50 ml of sterile Legionella Growth Supplement (BCYE) FD142 . Mix well and pour into sterile petri plates.

Storage and Shelf Life

Store at 2-8°C.Use before expiry date on the label.

Revision: 1 / 2012

* Not For Medicinal Use

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia[™] publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia[™] Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.