

Technical Data

NYC Supplement

An antibiotic supplement recommended for the isolation and cultivation of pathogenic *Neisseria* species.

Composition

Per vial sufficient for 500 ml medium

*Ingredients

Colistin methane sulphonate Trimethoprim Vancomycin Nystatin

Directions:

Rehydrate the contents of one vial aseptically with 10 ml of sterile distilled water. Mix well and aseptically add it to 320 ml of sterile, molten, cooled (45-50°C) NYC Agar Base M1348 along with 100 ml of sedimented horse blood cells, 60 ml of citrated horse plasma and rehydrated contents of 1 vial of Yeast Autolysate Supplement FD027. Mix well and pour into sterile petri plates.

Concentration

3.750mg

1.500mg

6250Units

1mg

Storage and Shelf Life

Store at 2-8°C. Use before the expiry date on the label.

Revision : 1 / 2012

* Not For Medicinal Use

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia[™] publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia[™] Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. A-516,Swastik Disha Business Park,Via Vadhani Ind. Est., LBS Marg, Mumbai-400086, India. Customer care No.: 022-6147 1919 Email: techhelp@himedialabs.com Website: www.himedialabs.com

FD150