

Technical Data

Lactic Supplement, Modified

FD155

A selective supplement recommended for the isolation of lactic acid bacteria in beer and brewing process.

Composition

Per vial sufficient for 500 ml medium

*Ingredients Concentration
Polysorbate 80 5ml
Amphotericin B 2.500mg
Phenyl ethanol 1.500g

Directions:

Aseptically add the contents of 1 vial to 500 ml sterile, molten, cooled (45-50°C) Lactic Acid Bacteria Selective Agar Base M1072 / Lactic Acid Bacteria Selective Broth Base M1384 . Mix well and dispense as desired.

Storage and Shelf Life

Store at 2-8°C.Use before expiry date on the label.

Revision: 02 / 2016

* Not For Medicinal Use

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia[™] publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia[™] Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.