



# **Doyle's Antibiotic Supplement, Modified**

FD159

An antibiotic supplement recommended for the selective isolation of Campylobacter species.

# Composition

Per vial sufficient for 500 ml medium

*Ingredients	Concentration
Vancomycin	7.500mg
Trimethoprim	2.500mg
Polymyxin B sulphate	10000IU
Amphotericin B	5mg

### **Directions:**

Rehydrate the contents of 1 vial aseptically with 5 ml of 50% ethanol. Mix well and aseptically add to 460 ml of sterile, molten, cooled (45-50°C) Doyles Enrichment Broth Base  $\underline{M916}$  / Doyles Enrichment HiVeg<sup>TM</sup> Broth Base  $\underline{MV916}$  along with 35 ml sterile lysed horse blood. Mix well and dispense as desired.

# **Storage and Shelf Life**

Store at 2 - 8°C. Use before expiry date on the label

Revision : 02 / 2016

\* Not For Medicinal Use

#### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia<sup>™</sup> publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia<sup>™</sup> Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. A-516,Swastik Disha Business Park,Via Vadhani Ind. Est., LBS Marg, Mumbai-400086, India. Customer care No.: 022-6147 1919 Email: techhelp@himedialabs.com Website: www.himedialabs.com