



Bushnell Haas Broth

M350

Bushnell Haas Broth is recommended for the examination of fuels for microbial contamination and for studying microbial hydrocarbon deterioration.

Composition**

Ingredients	Gms / Litre
Magnesium sulphate	0.200
Calcium chloride	0.020
Monopotassium phosphate	1.000
Dipotassium phosphate	1.000
Ammonium nitrate	1.000
Ferric chloride	0.050
Final pH (at 25°C)	7.0±0.2

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 3.27 grams in 1000 ml distilled water. Heat if necessary to dissolve the medium completely. Dispense into tubes or flasks as desired. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. A white precipitate prior to sterilization becoming yellow to orange after sterilization is normal.

Principle And Interpretation

Bushnell Haas Broth is prepared as per the formula of Bushnell and Haas (1) and recommended for the microbiological examination of fuels by the SIM Committee on microbiological deteriorations of fuels (2). These media contain all nutrients except carbon source, necessary for the growth of bacteria. Only those bacteria that are able to decompose hydrocarbon will grow in these media. Specific carbon source i.e. hydrocarbon can be added to this medium and their utilization by different microorganisms can be studied.

These bacteria can decompose a variety of hydrocarbons like kerosene, mineral oil, paraffin wax and gasoline. For liquid hydrocarbon the hydrocarbon is layered on the surface of inoculated agar. For testing volatile hydrocarbons such as gasoline the Petri-plates containing the medium are inverted and the hydrocarbon is poured into the lid.

Magnesium sulphate, calcium chloride and ferric chloride provide trace elements. Ammonium nitrate is a nitrogen source while monopotassium phosphate and potassium phosphate buffers the medium.

Quality Control

Appearance

White to cream homogeneous free flowing powder

Colour and Clarity of prepared medium

Colourless clear to slightly opalescent solution with slight brown precipitate in tubes.

Reaction

Reaction of 0.33% w/v aqueous solution at 25°C. pH : 7.0±0.2

pH

6.80-7.20

Cultural Response

M350: Cultural characteristics observed after an incubation at 25-30°C within 1 week.

Organism	Inoculum (CFU)	Growth (Plain)	Growth w/ minerals
----------	-------------------	-------------------	-----------------------

Cultural Response

<i>Pseudomonas aeruginosa</i> ATCC 27853	50-100	poor	good-luxuriant
<i>Pseudomonas aeruginosa</i> ATCC 9027	50-100	poor	good-luxuriant

Storage and Shelf Life

Store below 30°C in tightly closed container and the prepared medium at 2 - 8°C. Use before expiry date on the label.

Reference

1. Bushnell and Haas, 1941, J. Bacteriol., 41:653.
2. Allred, DeGray, Edwards, Hedrick, Klemme, Rogers, Wulf and Hodge, 1963, Proposed Procedures for Microbiological Examination of Fuels, SIM Special Publications, No. 1. Merck, Sharp & Dohme Research Laboratories, Rahway, N.J.

Revision : 02 / 2016



Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.