



Thiobacillus Broth

M789

Thiobacillus Broth is used for cultivation of *Thiobacillus* species

Composition**

Ingredients	Gms / Litre
Ammonium sulphate	0.400
Monopotassium phosphate	4.000
Calcium chloride	0.250
Ferrous sulphate	0.010
Magnesium sulphate	0.500
Sodium thiosulphate	5.000

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 10.16 grams in 1000 ml distilled water. Heat if necessary to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Mix well and dispense as desired.

Principle And Interpretation

The genus *Thiobacillus* is also known under the name of *Acidithiobacillus*. *Thiobacillus* are obligate autotrophic organisms, as they require organic carbon both as an electron and carbon source. Thiobacilli produce high quantity of sulphuric acid as a byproduct during oxidation of thiosulphates, sulphur and related inorganic sulphur-containing compounds to generate metabolic energy. *Thiobacillus*, by its production of sulphuric acid is involved in the destruction of concrete sewers and the acid corrosion of metals (2).

Thiobacillus Broth is a modification of formulation described by Starkey (1) and is used for the isolation and maintenance of *Thiobacillus* species.

The medium contains three inorganic sulphates and a thiosulphate. Phosphate serves as a buffer while sodium chloride maintains the osmotic balance of the medium.

Samples are inoculated into Thiobacillus Broth. After incubation at 25-30°C for about 7 days or more, turbidity or sulphur precipitation on the surface of the liquid or against the walls of the flasks, indicates growth of bacteria. Isolation is subsequently done on Thiobacillus Agar.

Quality Control

Appearance

White to cream homogeneous free flowing powder

Colour and Clarity of prepared medium

Colourless clear solution in tubes.

Cultural Response

M789: Cultural characteristics observed after an incubation at 25-30°C for upto 7 days.

Organism	Growth
<i>Thiobacillus thioparus</i> ATCC 8158	luxuriant
<i>Thiobacillus thiooxidans</i> ATCC 8085	luxuriant

Storage and Shelf Life

Store below 30°C in tightly closed container and the prepared medium at 2 - 8°C. Use before expiry date on the label.

Reference

1. Starkey R. L., 1935, Science, 39:197.

2.Eaton A. D., Clesceri L. S. and Greenberg A. E., (Ed.), 1995, Standard Methods for the Examination of water and Wastewater, 19th Ed., American Public Health Association, Washington, D.C.

Revision : 03 / 2015

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory,diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.