

# **Technical Data**

GTC Supplement FD044

An antibiotic supplement recommended for selective isolation of Enterococci.

### **Composition**

Per vial sufficient for 500 ml medium

\*Ingredients Concentration
Gentamicin 1.250mg

#### **Directions:**

Rehydrate the content of 1vial aseptically with 2 ml sterile distilled water and aseptically add to 500 ml sterile, molten, cooled (45-50°C) GTC Agar Base M922 alongwith 10 ml of sterile 10% Sodium bicarbonate solution FD080 just before use. Mix well and pour into sterile petriplates.

## **Storage and Shelf Life**

Store at 2-8°C. Use before expiry period on the label.

Revision: 1 / 2012

\* Not For Medicinal Use

#### Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia<sup>TM</sup> publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia<sup>TM</sup> Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.