

Technical Data

GBS Supplement FD054

An antibiotic supplement recommended for the selective isolation and presumptive identification of group B Streptococci.

Composition

Per vial sufficient for 450 ml medium

*Ingredients Concentration
Metronidazole 5mg
Gentamicin 1mg

Directions:

Rehydrate the contents of 1 vial aseptically with 2 ml of sterile distilled water. Mix well and add it to 475 ml sterile, molten, cooled (45-50°C) GBS Medium Base M1073 along with 25 ml inactivated Horse Serum RM1239. Mix well.

Storage and Shelf Life

Store at 2-8°C. Use before expiry period on the label.

Revision: 1/2012

* Not For Medicinal Use

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia[™] publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia[™] Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.