



Technical Data

Campylobacter Selective Supplement

FD077

An antibiotic supplement recommended for the selective isolation of *Campylobacter* species.

Composition

Per vial sufficient for 1000 ml medium

Ingredients

Rifampicin

Cefsulodin

Polymyxin B Sulphate

Concentration

25mg

6.250mg

20000Unit

Directions:

Rehydrate the contents of 1 vial aseptically with 10 ml of 50% ethanol. Mix well and aseptically add to 1000 ml sterile, molten, cooled (45-50°C) Wesley Broth Base [M1152](#) / Wesley HiVeg™ Broth Base [MV1152](#) along with 6.25 ml of cooled alkaline Hematin Solution. Mix well and dispense as desired.

Storage and Shelf Life

Store at 2-8°C. Use before the expiry date on the label.

Revision : 1 / 2012

* Not For Medicinal Use

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.