



Technical Data

Park and Sanders Selective Supplement I

FD104

An antibiotic supplement recommended for the selective isolation and enumeration of thermotolerant *Campylobacter* species .

Composition

Per vial sufficient for 950 ml medium

*Ingredients	Concentration
Vancomycin	10mg
Trimethoprim	10mg

Directions:

Rehydrate the contents of 1 vial aseptically with 5 ml of 50% ethanol. Mix well and aseptically add to 940 ml of sterile, cooled Park and Sanders Enrichment Broth Base [M1185](#) / Park and Sanders Broth Base [M1185I](#) / Park and Sanders Enrichment HiVeg™ Broth Base [MV1185](#) alongwith 50 ml sterile defibrinated blood and Park and Sanders Selective Supplement II [FD105](#) . Mix well and dispense into sterile tubes.

Storage and Shelf Life

Store at 2-8°C. Use before the expiry date on the label.

Revision : 1 / 2012

* Not For Medicinal Use

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.