

# **Technical Data**

## Listeria Selective Supplement II, Modified

FD163

An antimicrobial supplement recommended for the selective cultivation of Listeria species.

#### **Composition**

Per vial sufficient for 1000 ml medium

*Ingredients	Concentration
Acriflavin hydrochloride (Trypaflavin)	15mg
Nalidixic acid	40mg
Amphotericin B	5mg

#### **Directions:**

Rehydrate the contents of 1 vial aseptically with 5 ml of sterile distilled water. Mix well and aseptically add to 1000 ml sterile, cooled (45-50°C) Listeria Selective Broth Base M889 / Listeria Selective HiVeg<sup>TM</sup> Broth Base MV889 or to 1000 ml sterile, molten, cooled (45-50°C) Listeria Selective Agar Base M1474 . Mix well and dispense as desired.

### **Storage and Shelf Life**

Store at 2-8°C. Use before expiry date on the label.

Revision: 02/2016

\* Not For Medicinal Use

#### Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia<sup>TM</sup> publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia<sup>TM</sup> Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.