



Technical Data

Campylobacter Selective Supplement, Abeyta

FD294

An antibiotic Supplement recommended by FDA BAM for the selective isolation of *Campylobacter* species

Composition**

Formula : (One vial,sufficient for 1000ml medium)

| Ingredients | mg/vial |
|---------------------|---------|
| Sodium cefoperazone | 16.00 |
| Rifampicin | 5.00 |
| Amphotericin B | 1.00 |

Directions

Rehydrate the contents of 1 vial aseptically with 1 ml of 0.2 N Sodium hydroxide, further add 4 ml of sterile distilled water. Mix well and aseptically add to 470 ml of sterile, molten, cooled (45-50°C) Heart Infusion Yeast extract Agar Base along with 25 ml of sterile lysed horse blood and reconstituted contents of one vial of Campylobacter Growth Supplement. Mix well and pour into sterile Petri plates.

Storage and Shelf Life

Store at 2-8°C. Use before expiry date on the label.

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.